# English

#### A. DEVICE DESCRIPTION

Arthrex remote controls are used to operate only Arthrex arthroscopy pumps. Refer to the pump User Guide for the proper remote control for use with that device.

#### B. REFERENCES

These instructions were developed using the guidance given in the following standards:

- ANSI/AAMI ST79: Comprehensive guide to steam sterilization and sterility assurance in healthcare facilities.
- ISO 17664: Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices
- ISO 17665-1: Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices
- AMI TIR30: A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices
- AMI TIR34: Water for the reprocessing of medical devices
- Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling Guidance for Industry and Food and Drug Administration Staff*

#### C. LIMITATIONS ON REPROCESSING

Repeated processing does affect the useful life of these devices. End of useful life (which is typically 1 year from the start of use) is normally determined by wear due to repeated processing, wear due to use, and/or damage.

#### D. VALIDATION

The recommended cleaning and sterilization methods in this DFU have been validated in compliance with federal and international guidance/standards. In accordance with ISO 17665, the “overkill” approach was used for sterilization validation, and demonstrates a sterility assurance level (SAL) of 10<sup>4</sup>. Cleaning and sterilizing equipment and materials vary in performance characteristics; therefore, it is the responsibility of the facility/end user to perform the appropriate validation testing for any use beyond recommended performance characteristics. Distilled (critical) water was utilized as the terminal rinse water quality to ensure that residuals will not interfere with subsequent processing steps and to prevent spotting or discoloration.

#### E. CAUTIONS

- Never immerse the device in water or cleaning solution.
- Never clean the device with bleach, chlorine-based detergents, or any products containing sodium hydroxide (such as alkaline cleaning agents). These fluids will degrade the device and potentially damage functionality.
- Never clean the device in an ultrasonic cleaner or an automated washer/disinfector.
- To avoid damaging the device, do not impact or subject to blunt force. Rough handling may cause damage.
- Do not use the device for any purpose other than its intended use. Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
- Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

#### F. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

#### G. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbols glossary](http://www.arthrex.com/symbols glossary).

#### H. CLEANING

Devices must be adequately cleaned and sterilized prior to use or re-use. All devices are to be cleaned and sterilized prior to each application; this is required as well for the first use after delivery of the unsterile devices. Effective cleaning is an indispensable requirement for effective sterilization of the devices.

- POINT-OF-USE PREPARATION, CONTAINMENT, AND TRANSPORTATION**

It is recommended that devices are reprocessed within a maximum of 2 hours of use. At point of use, soiled devices must be removed from trays and moistened to prevent debris from drying before transportation to the reprocessing area for cleaning procedures. Avoid soaking devices in solutions; instead, apply low-linting cloths or towels dampened in water or neutral-pH detergent. Devices should be contained and transported in a closed, puncture-proof device to ensure safety. Do not clean soiled devices while in cases or trays.

#### II. DETERGENT SELECTION

Consider the following points during selection of the cleaning agent:

- Compatibility of the cleaning agent with the devices. Arthrex recommends the use of neutral pH or enzymatic cleaning agents or quaternary-based cleaning/disinfecting agents. Alkaline cleaning agents should be avoided since these agents will degrade the device’s exterior material composition and potentially damage the handpiece’s device’s functionality. Follow the instructions of the detergent manufacturer regarding use concentration and temperature. Use only freshly prepared cleaning solutions as well as critical (e.g. RO/DI or distilled) water for final rinse and filtered medical grade (instrument) air for drying.

#### III. MANUAL CLEANING

- Prepare a suitable cleaning solution (reference “Detergent Selection” section).
- Thoroughly scrub the device with a clean, soft-bristled brush that is dampened with the cleaning solution until all visible soil has been removed.
- Alternatively, if using a cleaning/disinfecting towelette or a low-linting cloth dipped in cleaning solution, thoroughly wipe down all surfaces until all visible soil has been removed. The use of additional towelettes or cloths may be used to ensure the surface is thoroughly cleaned.
- Rinse thoroughly with critical water.
- Check devices for visible soil. Repeat cleaning if soil is visible, and re-inspect.
- Dry devices thoroughly utilizing filtered medical grade air or a soft, clean, and low-linting cloth.

#### I. INSPECTION AND MAINTENANCE

- Arthrex devices are precision medical instruments and must be used with care.
- Inspect the device for damage prior to use and at all stages of handling thereafter.
- If damage is detected, consult your Arthrex representative for guidance.

#### J. STERILIZATION

##### I. STERILE PACKAGING

**Singly:** Single devices should be packed as to ensure that the pack is large enough to contain the device without stressing the seals. Packaging should be completed utilizing a pouch or wrap which conforms to the recommended specifications for steam sterilization as outlined below. If a wrap is utilized, it should be completed following AAMI double-wrap or equivalent guidelines with an appropriate wrap (cleared by the FDA or the local governing body).

**Sets:** Where appropriate, cleaned, disinfected and inspected devices should be placed into trays/cases as provided or in general-purpose sterilization trays. The total weight of trays/cases should not exceed 11.4 kg/25 lbs. (other local limits below 11.4 kg/25 lbs. may apply). Trays/cases should be double wrapped following AAMI or equivalent guidelines with an appropriate wrap (cleared by the FDA or the local governing body).

##### II. STERILIZATION

Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in the table below. Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer’s instructions.

RECOMMENDED STEAM STERILIZATION PARAMETERS				
Cycle Type	Minimum Exposure Temperature	Minimum Exposure Time	Minimum Drying Time <sup>1</sup>	Minimum Cooling Time <sup>2</sup>
<b>Preva-cuum Cycle</b>	132°C (270°F)	4 Minutes	30 Minutes	30 minutes
<b>Gravity-Displacement Cycle</b>	132°C (270°F)	15 Minutes	30 Minutes	30 minutes
<b>Gravity-Displacement Cycle</b>	121°C (250°F)	30 minutes	30 Minutes	30 minutes

<sup>1</sup>Drying times vary according to load size and should be increased for larger loads.

<sup>2</sup>Cooling times vary according to the sterilizer used device design, temperature and humidity of ambient environment, and type of packaging used. Cooling process should comply with ANSI/AAMI ST79.

**CAUTION: Ensure that the connector cap is secured before sterilization. NEVER allow the console plug pins to have any contact with liquids. Remove dust or moisture, if necessary, with dry compressed air. ONLY dry connectors may be plugged into the console.**

**CAUTION: NEVER immerse the cable or handpiece in liquid to cool, as this can damage the cable, handpiece, connectors, and/or seals.**

##### III. SPECIAL PRECAUTION - TRANSMISSIBLE SPONGIFORM ENCEPHALOPATHY AGENTS

It is outside the scope of this document to describe in detail the precautions that should be taken for Transmissible Spongiform Encephalopathy (TSE) Agents. The agents for transmission of Creutzfeldt-Jakob disease (CJD) are believed to be resistant to normal processes of disinfection and sterilization and therefore the normal processing methods of decontamination and sterilization as described above may not be appropriate where CJD transmission is a risk. In general, the tissues that come into contact with orthopedic surgical instruments are those of low TSE infectivity. However, particular precautions should be taken when handling instruments that have been used on known, suspected, or at-risk patients. Refer to ANSI/AAMI ST79 for further information.

#### K. STORAGE

Non-sterile devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions. It is the responsibility of the end-user to ensure devices,

once sterilized, are stored in such a way as to maintain the sterility of the device until use. Sterile, packaged devices should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, and temperature/humidity extremes. Sterile device packages should be carefully examined prior to opening to ensure that package integrity has not been compromised. Maintenance of sterile package integrity is generally event related. If a sterile wrap is torn, shows any evidence of tampering, or has been exposed to moisture, the device or set must be cleaned, repackaged, and sterilized.

- STORAGE CONDITIONS**
  - Temperature -30° to +70°C (-22° to 158°F)
  - Relative Humidity 10% to 90%, non-condensing
  - Barometric pressure 500 hPa (7.25 PSI) to 1060 hPa (15.37 PSI)

- OPERATING CONDITIONS**
  - Temperature 10° to 40 °C (50° to 104 °F)
  - Relative Humidity 20% to 75%, non-condensing
  - Barometric pressure 700 hPa (10.15 PSI) to 1060 hPa (15.37 PSI)

#### L. INFORMATION

- In CE accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE accepting Countrries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

#### M. DIRECTIONS FOR USE

- Only use with Arthrex arthroscopic pumps.
- Remove the connector cover.
- Rotate the connector until the red dot is pointed up.
- Insert the connector from the device into the appropriate receptacle of the console. The remote connects and locks to the console to prevent accidental separation during use.
- When connected, the console will recognize the remote and display an icon on the screen when the console power is on.
- Disconnect the remote by pulling on the connector, not the cable.

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