

ALLOGRAFT PACKAGE INSERT

DONATED HUMAN TISSUE

THIS ALLOGRAFT IS SUPPLIED STERILE

This human tissue allograft is processed and manufactured by CellRight Technologies. All tissue was retrieved, processed, stored and distributed for use in accordance with the standards of the American Association of Tissue Banks (AATB), FDA requirements for Human Cellular and Tissue Based Products (HCT/Ps 21 CFR Part 1271), and applicable State regulations. The Donor has been determined to be eligible based on the results of screening and testing. Screening includes a review of medical and social history, available hospital records, infectious disease screening, autopsy report (if performed), and physical exam. The Donor has been tested and was found negative or non-reactive for:

- Human Immunodeficiency Virus Types 1 and 2 Antibody (anti-HIV-1/anti-HIV-2)
- Hepatitis B Surface Antigen (HBsAg)
- Hepatitis B Core Antibody - Total (anti-HBc)
- Hepatitis C Virus Antibody (anti-HCV)
- Human Immunodeficiency Virus 1 and Hepatitis C Virus Nucleic Acid Test (HIV 1/HCV NAT)
- Syphilis Rapid Plasma Reagin or Treponemal Specific Assay

Additional tests, including but not limited to HTLV I/II and HBV Nucleic Acid Testing, may have been performed and were found to be acceptable for transplantation. U.S. Food and Drug Administration (FDA) licensed, approved, or cleared donor screening test kits are used when available. Communicable disease testing has been performed by a laboratory registered with the FDA to perform donor testing in accordance with the Clinical Laboratory Improvement Amendments of 1988 (CLIA) and 42 CFR part 493, or that has met equivalent requirements as determined by the Centers for Medicare and Medicaid Services (CMS). A list of additional communicable disease test(s) performed will be provided upon request.

The CellRight Technologies Medical Director has determined this donor tissue to be suitable for transplantation. The testing and medical release records are maintained by CellRight Technologies. The names and addresses of the testing laboratories, the interpretation of all required infectious disease tests, and a listing of the documents reviewed as part of the relevant medical records are kept on file at CellRight Technologies and are available upon request.

Tissue has been sterilized, using Cobalt 60, to a SAL of 10^{-6} (Sterility Assurance Level). Allografts are processed using some or all of the following agents: physiological buffers, acids, alcohols, surfactants, hydrogen peroxide, Gentamicin Sulfate, Vancomycin HCl, Amphotericin B, Polymyxin B, and/or Ciprofloxacin and traces may remain.

WARNINGS AND PRECAUTIONS

- Intended for use in one patient, on a single occasion only.
- Do not use if package integrity has been compromised. Once the user breaks the seal on the inner-most pouch, the tissue grafts must be transplanted or discarded.
- Tissue may not be sterilized or re-sterilized by your facility.
- This tissue is intended for use by qualified healthcare specialists such as physicians, dentists, or podiatrist.
- Although this tissue has been tested and screened for human pathogens, and processed under aseptic conditions, human derived tissue may still transmit infectious agents.
- It is the responsibility of the Tissue Dispensing Service, Tissue Distribution Intermediary, and/or End-User clinician to maintain tissue intended for transplantation in appropriate storage conditions prior to further Distribution or transplant.

STORAGE

DEMINERALIZED BONE MATRIX (DBM) GRAFTS – Maintain tissue at room temperature (15°C - 30°C).

FREEZE-DRIED or DEHYDRATED TISSUE – Maintain tissue at room temperature (15°C - 30°C).

FROZEN TISSUE – Maintain tissue at -40°C or colder. Short term storage of up to 6 months is acceptable if tissue is maintained at -20°C to -39°C.

TISSUE PREPARATION

BEFORE USE – Examine Allograft Packaging – Do Not Use This Allograft If:

1. Any of the package elements appear to be missing, tampered with or damaged.
2. The product label or identifying bar code is severely damaged, illegible or missing.
3. The expiration date shown on the package label has passed.
4. Frozen allograft has not been stored according to storage temperature requirements or the allograft has been prematurely thawed.

If any of the above conditions exist or are suspected, this allograft should NOT be used.

PREPARATION OF DBM FOR USE

1. Ready to use – DBM is provided as a paste, putty, gel or crunch product that is ready for implantation. It does not require thawing or rehydration.
2. Opening Peel Packages: peel outer package down and aseptically deliver inner peel pouch to the sterile field or sterile team member.
3. DBM is contained inside the inner pouch in a jar, syringe, or other storage container.
4. Remove container of DBM from the Inner peel pouch.
 - a. Jar – Unscrew the top. Remove DBM from the jar. Mold into desired shape and press into defect.
 - b. Syringe – Remove protective cap from syringe tip or remove the syringe end cap completely, apply pressure to the plunger to extrude the DBM. Mold into desired shape and press into defect.
5. Irrigation resistant once molded and pressed into the defect.
6. For best results. The DBM must fill the defect and contact as much viable bone as possible.

PREPARATION OF FREEZE-DRIED or DEHYDRATED ALLOGRAFT TISSUE FOR USE

1. Opening Peel Packages: peel outer package down and aseptically deliver inner peel pouch to the sterile field or sterile team member.
2. Remove tissue from Inner peel pouch.
3. Tissue may be maintained within the inner pouch in a jar, syringe, or other storage container.
 - a. Jar – Unscrew the top. Rehydrate tissue in jar or transfer tissue to a basin for rehydration.
 - b. Syringe – Rehydrate tissue in syringe or transfer to a basin for rehydration.

**PREPARATION OF FREEZE-DRIED
or DEHYDRATED ALLOGRAFT TISSUE FOR USE**

Continued from previous page

4. Rehydrate the tissue, when applicable.
 - a. Final determination of allograft reconstitution should be made by the physician prior to use. Rehydrate using a sterile isotonic solution or solution of physicians' choice.
 - b. Recommendation – Non-weight bearing osseous grafts and soft tissue should be reconstituted for a minimum of 30 minutes.
 - c. Recommendation – Weight bearing grafts should be reconstituted approximately 1 hour.
 - d. Recommendation – Grafts that are to be manipulated by drilling or cutting or require force to insert may require a longer period of reconstitution prior to manipulation to reduce the chance of fracturing.
5. Tissue should be used as soon as possible after reconstitution. If tissue is to be stored for longer than 2 hours after reconstitution, it should be refrigerated at 1°C to 10°C in an aseptic container for no longer than 6 hours.
6. **IMPORTANT!** Peel away and remove all internal packaging materials, if present, from the graft (ie, gauze or mesh) prior to implantation.

RETURNS

With prior approval, unused, unopened tissue may be returned to CellRight Technologies provided Arthrex personnel have authorized the return and issued a return authorization number. The responsible individual at your facility must obtain a Tissue Return Authorization Form from CellRight Technologies, complete the required information and provide a signature declaring the unopened tissue has been continuously stored according to instructions and that proper transportation has been utilized to ensure tissue integrity during the return. This form must be completed for credit to be issued. Allograft tissue products may not be returned to Arthrex.

Contact Customer Service at Arthrex by email or phone.

Email: returns@arthrex.com

Phone: 1-800-934-4404

ADVERSE OUTCOMES

Adverse outcomes potentially attributable to this tissue must be reported promptly to Arthrex, Inc. at 1-800-933-7001 ext. 78718 or email at: Complaints@Arthrex.com

TISSUE TRACKING

Complete the enclosed Allograft Tracking Form and mail to Arthrex, Inc. US Federal Regulations (21 CFR 1271.290(b)) and Joint Commission Standards (TS.03.02.01, EP 7) require proper tracking of this tissue. It is the responsibility of the end-user to provide this information, which enables Arthrex, Inc. to maintain records for the purpose of tracing the tissue post-transplant.

Manufactured for:



www.Arthrex.com

**1370 Creekside Blvd.
Naples, FL 34108**

Processed and Manufactured by:



**1808 Universal City Blvd.
Universal City, TX 78148
210-659-9353
Fax: 210-659-9556**

CellRight Technologies holds:

AATB Accreditation No. 00212
US FDA Registration No. 3009234552
Canadian Registration No. 100228
California Tissue Bank ID No. CNC80949
Florida License No. 212
Maryland Tissue Bank No. TB1898
New York State Tissue Bank ID No. CPI73TP141TS145



www.Arthrex.com