English

A. DEVICE DESCRIPTION

The Arthrex Low Profile Screws are headed and self-tapping. They are available for full or partially threaded, solid or cannulated. The screw family ranges from 1.0 mm to 6.7 mm in diameter and from 6 mm to 120 mm in length (in 1, 2 or 5 mm increments).

B. INDICATIONS

The **Arthrex Mini GFS (Comprehensive Fixation System) Screw** (1.0-2.4 mm solid) is intended for use in selective trauma, revision, and general surgery of the hand, wrist, and other small bones. When used with a plate, the screws may be used with the Arthrex Comprehensive Fixation System Plates (1.4-2.4 mm).

The **Arthrex Mini GFS Screws (2.0-2.4 mm solid)** are intended for fixation of fractures, osteotomies, nonunions, replantations, and general surgery of the hand, wrist, and ankle, particularly in osteopenic bone in the hand, wrist, foot, and ankle. When used with a plate, the screws may be used with the Arthrex Comprehensive Fixation System Plates (2.0-2.4 mm).

The **Arthrex Low Profile Screws (2.0-3.0 mm solid)** are intended to be used as stand-alone bone screws, or in a plate-screw system for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, and wrist.

The **Arthrex Low Profile Screws (2.7 mm screws only)** Use with Intracortical Nails (Except in EU).

The **Arthrex Low Profile Screws (2.0-3.0 mm cannulated)** are intended to be used as stand-alone bone screws for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, and wrist.

The **Arthrex Low Profile Screws (3.5 mm and larger, solid)** are intended to be used as stand-alone bone screws, or in a plate-screw system for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, femur and fibula. When used with a plate, the screws may be used with the Arthrex Low Profile, Small Fragment Plates, Humeral Fracture Plates, and Distal Radius Plates. 2.7 mm screws only: Use with Intracortical Nails (Except in EU).

The **Arthrex Low Profile Screws (3.5 mm and larger, cannulated)** are intended to be used as stand-alone bone screws for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, femur and fibula.

The **Arthrex Vented Cannulated Locking Screws (4.0 mm, cannulated)** are intended to be used in a plate-screw system for internal bone fixation for bone fractures in the humerus. The screws are used with the Arthrex Humeral Fracture Plates. Additional augmentation, as needed, may be obtained through the screw.

These screws are for EU distribution only.

The **Arthrex Blunt Tip Screws** are intended to be used as stand-alone bone screws for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus, femur and fibula. When used in conjunction with FiberGrip™, they can be used to treat patella fractures (Except in EU).

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply, osteomyelitis and previous infections, which may require healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are potentially immature. The use of this immature device and the placement of hardware or implants must not bridge, distort or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.

E. WARNINGS

- Caution: Federal Law restricts this device to sale by or on the order of a physician.
- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be reused.
- All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other non-supportive stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of devices, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- Devices that have been implanted for a long period of time may require the use of screw removal instrumentation.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeons should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management and avoid re-fracture.
- Biobazard waste, such as explanted devices, needles and contaminated surgical equipment, should be disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex, Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex Low Profile Screws are MR Conditional. A patient with the following can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5 Tesla and 3-Tesla, only
 - Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
 - Maximum MR system reported, while body averaged specific absorption rate (SAR) of 1W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the Arthrex Low Profile Screws are exposed from the scan head field is not recommended. Scan times of 15-minutes of continuous scanning at a temperature rise of up to 6°C

In non-clinical testing, the image artifact caused by the Arthrex Low Profile Screws can extend up to approximately 120 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate screw size based on the specific technique, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Use the appropriately sized drill bit for the screw.
- Damage to the driver may cause screw result from failure to seat the driver fully into the screw or to align the driver properly with the screw.
- QuickFit™** Screws Only: Bending the insertion post to remove it from the screw head is not recommended. Screws should be inserted by hand and not with powered equipment.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolglossary.

I. VALIDATION

The recommended cleaning, disinfection, and sterilization methods in this DFI have been validated in compliance with federal and international guidance/standards. In accordance with ISO 17665, the "overkill" approach was used for sterilization validation, and demonstrates a sterility assurance level (SAL) of 10⁻⁶. Cleaning, disinfecting, and sterilizing equipment and materials vary in performance characteristics. Therefore, it is the responsibility of the facility-end user to perform the appropriate validation testing for any use beyond recommended performance characteristics.

In accordance with ISO 17664 and ANSI/AAMI TR30, limit values and a means for monitoring chemical residues following cleaning have been established for the product. In assessing the level of contamination and the cleaning and disinfection process or the machine (automated) cleaning and disinfection process, a clinically relevant method was utilized for testing the safety of residuals as part of the validation protocol. Deionized (critical) water was utilized as the terminal rinse water quality to ensure that residuals will not interfere with subsequent processing steps. Repeat cleaning (soil is visible and re-inspect, otherwise, proceed to sterilization stage).

IV. MANUAL CLEANING AND DISINFECTION
Following preliminary cleaning, the instructions for Manual Cleaning and Disinfection may be followed as an alternative cleaning method to Machine (Automated) Cleaning and Thermal Disinfection if an automated procedure is not available.

After preliminary cleaning is complete, repeat steps 1-5, including in the Preliminary Cleaning section with this DFI, provided cleaning and disinfection may be followed as an alternative cleaning method to Machine (Automated) Cleaning and Thermal Disinfection if an automated procedure is not available.

Use references to those single use devices that have come into contact with blood, bone, tissue, or other body fluids. Any unused single-use devices that has been exposed to blood, bone, tissue, or body fluids **must not be reprocessed and must be discarded**.

The instructions in this DFI were developed using the guidance given in the following standards:

- ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities"
- Repeted processing has minimal effect on these devices.
- Soak the devices for the given soaking time (provided by the distal manufacturer) in disinfecting solution so that the devices are sufficiently cleaned. Make sure that there is no contact between the devices. Ensure that the device is in the open position during soaking. Activate movable parts at least five times during disinfection, as applicable.
- Remove the devices from the disinfectant solution and rinse per disinfectant manufacturer's instructions.
- Drives thoroughly utilizing filtered medical grade air or a soft, low-linting cloth. Proceed to Sterilization section.

J. CLEANING AND DISINFECTION

Certain Arthrex devices that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. All devices are to be cleaned, disinfected, and sterilized prior to each application; this is required as well for the first use after delivery of the sterile devices. An effective cleaning and disinfection is an indispensable requirement for an effective sterilization of the devices. Single-use devices **must** be cleaned separately from solid devices.

If possible, the Machine (Automated) procedure should be used for cleaning and disinfection of the instruments. The Manual Cleaning procedure should only be used if an automated procedure is not available; in this case, the significantly lower efficiency and reproducibility of the manual procedure should be considered. The Preliminary Cleaning steps are to be performed in both cases. Manual cleaning may require onsite validation by the healthcare facility and appropriate procedures should be in place to avoid human factor variability.

L. DETERGENT SELECTION

Consider the following points during selection of the cleaning detergent:

- Suitability of the cleaning agent for ultrasonic cleaning (no foam development).
- Compatibility of the cleaning agent with the instruments. Arthrex recommends the use of neutral pH or enzymatic cleaning agents. Alkaline agents may be used to clean devices in countries where required by law or local ordinance, or where prior diseases such as Transmissible Spongiform Encephalopathy (TSE) or Creutzfeldt-Jakob disease (CJD) are a concern. **Caution: Low or alkali solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as PEP (Fluorinatedethylenepropylene), ABS (Acrylonitrile Butadiene Styrene), Ultem™, Lexan™, and Cycloc™™. Non-neutral pH cleaning chemicals are utilized, care should be taken to ensure appropriate rinsing, as validated by the end-user facility, and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.**
- Follow the instructions of the detergent manufacturer regarding use concentration and temperature for either manual or automated cleaning. Use only freshly prepared solutions as well as only purified/highly purified water at least for final rinse, and a soft, low-linting cloth and/or filtered medical grade air for drying, respectively.

NOTE: No assembly/disassembly of these devices is required unless stated on the labeling, directions for use, or literature assemblations (LA) pertaining to cleaning, disinfection, and sterilization. Devices that require disassembly are to be disassembled prior to cleaning.

- Remove excess soil from devices, especially in areas such as joints and crevices, by cleaning the surfaces with a sponge or brush under cold running water or with a non-shedding disposable wipe for a minimum of 30 seconds.
- Rinse the devices at least 1 minute under running utility water (temperature < 35 °C/95 °F). Special attention should be given to lumens, joints, crevices, and other hard-to-reach areas. Immerse the device in cleaning solution inside an ultrasonic bath. While immersed in solution, brush the devices for 1 minute using a soft-bristled brush. Special attention should be given to lumens, joints, crevices, and other hard-to-reach areas. Lumens should be brushed with appropriate diameter and length bristle sonics for the particular lumen. Activate movable parts at least (5) times during soaking, as applicable.
- After brushing, rinse on ultrasonic power and soak and sonicate for 10 minutes at a minimum of 40-2.5 kHz. Ensure devices are in the open position and that lumens have complete contact with cleaning solution during soaking.
- Remove the devices from the cleaning solution and rinse at least 1 minute with ultra water. Thoroughly and aggressively rise lumens, joints, crevices, and other hard-to-reach areas.
- After the completion of preliminary cleaning, the end user has the option to perform either Manual Cleaning and Disinfection Or Machine (Automated) Cleaning and Thermal Disinfection (preferred).

III. MACHINE (AUTOMATED) CLEANING AND THERMAL DISINFECTION

Considerations for the selection of the washer-disinfector:

- Capable of providing an approved program for thermal disinfection (appropriate exposure time and temperature according to A2 concept)
- Final rinse completed with purified (critical, e.g. RO or DI) water, and utilizes only filtered air for drying.
- After preliminary cleaning is complete, load the devices in the washer-disinfector such that all design features of the device are accessible to cleaning and that design features that might retain liquid can drain (for example, hinges should be open and cannulations/holes positioned to drain).
- Use appropriate cleaning agents; a neutralization step should be utilized as appropriate.
- Run an automated wash cycle with fundamentally approved efficiency of the washer-disinfector (for example, CE marking according to EN ISO 15883 or FDA approved/clearance/registration). The following minimum recommended automatic wash cycle parameters were utilized by Arthrex during the validation of these instructions.

RECOMMENDED WASHING CYCLE PARAMETERS			
Phase	Recirculation Time	Temperature	Detergent
Pre-Wash	3 Minutes	Gold Water	N/A
Cleaning Wash	10 Minutes	Follow detergent manufacturer's recommendation	Enzymatic or alkaline detergent
Neutralization Rinse (optional)	2 Minutes	Follow detergent manufacturer's recommendation	Neutralizing agent (as needed)
Rinse	3 Minutes	Gold Water	N/A
Thermal Disinfection Rinse	5 Minutes	90°C (194°F)	N/A
Drying	Minimum 6 Minutes or until visibly dry	Minimum 100°C (212°F)	N/A

- Remove the devices from the washer-disinfector following the completion of the program and check devices for visible soil. Repeat cleaning (soil is visible and re-inspect, otherwise, proceed to sterilization stage).

IV. MANUAL CLEANING AND DISINFECTION

Following preliminary cleaning, the instructions for Manual Cleaning and Disinfection may be followed as an alternative cleaning method to Machine (Automated) Cleaning and Thermal Disinfection if an automated procedure is not available.

- After preliminary cleaning is complete, repeat steps 1-5, including in the Preliminary Cleaning section with this DFI, provided cleaning and disinfection may be followed as an alternative cleaning method to Machine (Automated) Cleaning and Thermal Disinfection if an automated procedure is not available.
- Use references to those single use devices that have come into contact with blood, bone, tissue, or other body fluids. Any unused single-use devices that has been exposed to blood, bone, tissue, or body fluids **must not be reprocessed and must be discarded**.
- The instructions in this DFI were developed using the guidance given in the following standards:
 - ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities"
 - Repeted processing has minimal effect on these devices.
 - Soak the devices for the given soaking time (provided by the distal manufacturer) in disinfecting solution so that the devices are sufficiently cleaned. Make sure that there is no contact between the devices. Ensure that the device is in the open position during soaking. Activate movable parts at least five times during disinfection, as applicable.
 - Remove the devices from the disinfectant solution and rinse per disinfectant manufacturer's instructions.
 - Drives thoroughly utilizing filtered medical grade air or a soft, low-linting cloth. Proceed to Sterilization section.

K. STERILIZATION

This device may be provided either sterile or non-sterile. Check the package labeling for more information. For devices that are not provided in a terminally sterilized configuration, sterilization is to be performed following cleaning, disinfection, and sterile packaging prior to use, and may be re-sterilized (if unused) following cleaning, disinfection, and sterile packaging prior to use.

Devices that are provided in a terminally sterilized configuration should never be re-sterilized under any conditions.

Certain Arthrex instruments and devices that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

L. STERILE PACKAGING

Singly: Single devices should be packed to ensure that the pack is large enough to contain the device without stressing the device. Packaging should be completed utilizing a pouch or wrap which conforms to the recommended specifications for steam sterilization as outlined below. If a wrap is utilized, it should be completed following AAMI double-wrap or equivalent guidelines with an appropriate wrap. An appropriate wrap is one that, for example, is cleared by the FDA or the local governing body or is an approved reusable rigid sterilization container. Aesculap SteriContainer™ rigid containers with perforated bottoms and lids are approved for use with Arthrex, Inc. sets. Areas designated for specific devices shall contain only devices intended for those areas. Devices should not be stacked or placed in close contact. Only Arthrex devices should be included in the trays or cases. These validated reprocessing instructions are not applicable to trays or cases that include devices not intended to be used with Arthrex trays or cases.

Sets: Where appropriate, cleaned, disinfected and inspected devices should be placed into trays/cases as provided or in general-purpose sterilion trays. The total weight of trays/cases should not exceed 11.4kg/25 lbs. (other local limits follow 11.4kg/25 lbs. may apply). Trays/cases should be double wrapped following AAMI or equivalent guidelines with an appropriate wrap. An appropriate wrap is one that, for example, is cleared by the FDA or the local governing body at the point of use. Sets may also be placed into an approved reusable rigid sterilization container. Aesculap SteriContainer™ rigid containers with perforated bottoms and lids are approved for use with Arthrex, Inc. sets.

Areas designated for specific devices shall contain only devices intended for those areas. Devices should not be stacked or placed in close contact. Only Arthrex devices should be included in the trays or cases. These validated reprocessing instructions are not applicable to trays or cases that include devices not intended to be used with Arthrex trays or cases.

II. STERILIZATION

Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in the table below. Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions.

RECOMMENDED STEAM STERILIZATION PARAMETERS				
Cycle Type	Minimum Exposure Temperature	Minimum Exposure Time	Minimum Drying Time ¹	Minimum Cooling Time ²
US Pre-vacuum Cycle	132°C (270°F)	4 Minutes	30 Minutes	30 Minutes
UK Pre-vacuum Cycle	134°C (273°F)	3 Minutes	30 Minutes	30 Minutes
Pre-vacuum Cycle ³ (Prion Cycle)	134°C (273°F)	18 Minutes	30 Minutes	30 Minutes

¹Drying times vary according to load size and should be increased for larger loads.

²Cooling times vary according to the sterilizer used, device design, temperature and humidity of ambient environment, and type of packaging used. Cooling process should comply with ANSI/AAMI ST79.

³Reprocessing parameters recommended by the World Health Organization (WHO) where there is concern regarding TSE/GJD contamination.

L. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

This device is made of titanium or stainless steel.

M. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under recommended conditions.

It is the responsibility of the end-user to ensure devices, once sterilized, are stored in such a way as to maintain the sterility of the device until use. Sterile, packaged devices should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, and temperature/humidity extremes. Sterile device packages should be carefully examined prior to opening to ensure that package integrity has not been compromised. Maintenance of sterile package integrity is generally event related. If the wrap is torn, shows an evidence of tampering, or has been exposed to moisture, the device or set must be cleaned, repackaged, and sterilized.

It is the responsibility of the end-user to ensure devices, once sterilized, are stored in such a way as to maintain the sterility of the device until use. Sterile, packaged devices should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, and temperature/humidity extremes. Sterile device packages should be carefully examined prior to opening to ensure that package integrity has not been compromised. Maintenance of sterile package integrity is generally event related. If the wrap is torn, shows an evidence of tampering, or has been exposed to moisture, the device or set must be cleaned, repackaged, and sterilized.

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Nota: Si el producto no está etiquetado como apto para su uso con los dispositivos de Arthrex, Inc.