English

A. DEVICE DESCRIPTION

The Arthrex Syndesmosis TightRope® XP Buttress Plate Implant System is comprised of one 2 Hole Buttress Plate, two metal buttons and #5 suture. The metal buttons and #5 suture are preloaded onto a disposable inserter. The ancillary 2 Hole Buttress Plate is a metal plate with two holes used with the Syndesmosis TightRope XP devices. The devices are packeaged with various ancillary instruments to aid in insertion.

B. INDICATIONS

The Arthrex Syndesmosis TightRope XP Buttress Plate Implant System is intended as an adjunct in fracture repair involving metaphyseal and periarticular small bone fragments where screws are not indicated, and as an adjunct in external and intramedullary fixation systems involving plates and rods, with fracture braces and casting.

Specifically, the Arthrex Syndesmosis TightRope XP Buttress Plate Implant System is intended to provide fixation during the healing process following a syndesmosis trauma, such as but not limited to, fixation of syndesmosis (syndesmosis disruptions) in connection with Weber B and C ankle fractures.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeris other than those indicated.

D. ADVERSE EFFECTS

- Osteomyelitis surrounding the Syndesmosis TightRope XP implant.
- Rediastasis resulting from a failure of the implant insertion technique.
- Infections, both deep and superficial.
- Foreign body reactions.
- Silicone sensitivity, though very rare, has been reported.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be re-used.
- Do not re-sterilize this device.
- All metallic implant devices used for this surgical procedure should have the same metallurgical composition, including the plates that might be used in conjunction with the Arthrex Syndesmosis TightRope XP Buttress Plate Implant System. Therefore, different material implants are provided.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device and bone.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device, the patient leaflet (www.arthrex.com/patientleaflets) and the patient implant card should be given to the patient. Your surgeon will guide you in deciding what particular treatment is best for you and explain the benefits, risks, and contraindications associated with the treatment.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion,

with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.

- Once open, discard any unused suture.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex Syndesmosis TightRope XP Buttress Plates and Suture Buttons manufactured from Titanium or Stainless Steel are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15-minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the Arthrex Syndesmosis TightRope XP Buttress Plates and Suture Buttons are expected to produce a maximum temperature rise of ≤2.5 °C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Arthrex Syndesmosis TightRope XP Buttress Plates and Suture Buttons can extend up to approximately 40 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

2. MR Safe

Artbrex Suture manufactured from UHMWPE is MR safe.

G. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

Buttons/Plates: The buttons and plates are available in two different metals in order to match the metallurgic composition of implants that may be used in the surgical procedure. The materials available are titanium and stainless steel.

Suture: Ultra High Molecular Weight Polyethylene (UHMWPE). Suture meets or exceeds U.S.P. and European standards for non-absorbable surgical sutures.

K. STORAGE CONDITIONS

Devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

- In CE accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE accepting Countries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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