

---

# BoneSync™ Cement: A Biomechanical Analysis of Rotator Cuff Repair Augmentation in Open-Cell Foam Block

Arthrex Orthopedic Research

---

## Background

---

Poor bone quality induces several challenges for arthroscopic rotator cuff repairs. Poor-quality bone of the proximal humerus can increase the chances of anchor pullout and require revisional surgery.<sup>1</sup> Calcium phosphate bone cement, such as BoneSync cement, can be useful for augmenting poor-quality bone.<sup>2</sup>

## Materials and Methods

---

BoneSync cement has a set time of 6-10 minutes, and hardening occurs 24 hours postimplantation.<sup>2</sup> Several timepoints (0 minutes, 5 minutes, 10 minutes, and 24 hours) were selected to evaluate the fixation strength of Biocomposite 4.75 mm SwiveLock® anchors augmented with BoneSync cement throughout the hardening process. Timepoints were defined as the amount of time from the start of mixing until tensile testing was conducted.

## Sample Preparations

---


Open-cell foam block was selected to mimic poor bone quality. A hole was punched into the 7.5 pcf open-cell foam block using a SwiveLock punch.

FiberTape® suture was threaded through the eyelet of a SwiveLock anchor to form a loop. For control group samples (no BoneSync cement), SwiveLock anchors were immediately inserted into the punched hole with the loop held around a 0.5-inch dowel rod.

For all other samples, BoneSync cement was mixed with 0.8 cc phosphate-buffered saline (PBS) according to the packaging instructions. The mixture was injected into the punched hole using a delivery cannula, and then the SwiveLock anchor was inserted directly following injection (Figure 1).

Each sample was closed in a jar of PBS and placed in a 37 °C (body temperature) environmental chamber for a set amount of time (5 minutes, 10 minutes, or 24 hours). Zero-minute timepoint samples were not immersed in PBS.

Figure 1. Prepared BoneSync cement sample



## Tensile Testing Method

---

Once removed from the PBS, the suture loop was attached to a hook and the foam block was secured under a box fixture. A preload of ~0.1 N was applied. Following the preload, tensile testing was performed at a rate of 5 mm/s until failure. The maximum load (up to 3 mm of displacement), displacement at maximum load, and failure mode were recorded.


## Results

---

The maximum load was defined as the maximum load achieved at or before 3 mm displacement. Displacement at 3 mm was used to represent clinical failure for this application.

Figure 2 shows the mean maximum load at each time point. Significant differences were found between the control group and 10-minute ( $P = .009$ ) and 24-hour ( $P < .001$ ) timepoints. The failure mode for all samples was anchor pullout (Figure 3).


**Figure 2.** Maximum load at or before 3 mm of displacement for each timepoint<sup>3</sup>



**Table 1.** Maximum Load Data (Mean ± SD)<sup>3</sup>

| Timepoint | Maximum Load (N) |
|-----------|------------------|
| Control   | 6.67 ± 2.27      |
| 0 min     | 15.4 ± 2.02      |
| 5 min     | 17.1 ± 5.61      |
| 10 min    | 22.2 ± 1.38      |
| 24 hr     | 74.2 ± 14        |

**Figure 3.** Failure mode for all samples was anchor pullout



## Conclusions

The presented porous sawbones model with BoneSync™ augmentation demonstrates significantly higher ultimate loads at the 10-minute and 24-hour timepoints compared to the models without BoneSync augmentation. The time-zero analysis provides evidence that BoneSync augmentation in poor-quality bone may be more resistant to pullout forces after the specified setting time.

### References

1. Kirchoff C, Braunstein V, Milz S, et al. Assessment of bone quality within the tuberosities of the osteoporotic humeral head: relevance for anchor positioning in rotator cuff repair. *Am J Sports Med.* 2010;38(3):564-569. doi:10.1177/0363546509354989
2. Miyamoto Y, Ishikawa K, Takechi M, et al. Tissue response to fast-setting calcium phosphate cement in bone. *J Biomed Mater Res.* 1997;37(4):457-64. doi:10.1002/(sici)1097-4636(19971215)37:4<457::aid-jbm3>3.0.co;2-k
3. Arthrex, Inc. Data on file (APT-06573). Naples, FL; 2024.