English

A. DEVICE DESCRIPTION

The Arthrex Corkscrew™, PushLock®, and SwiveLock® suture anchors consist of anchors with an integral or separate eyelet. The PushLock Tenodesis anchor is a two-piece push-in anchor system with either a forked or closed eyelet. They are pre-loaded on a handled inserter. Suture, with or without needles, and a suture threader may also be provided.

The Arthrex Implant System, Tenodesis Screw Eyelet is a kit comprised of attachable eyelets, FiberTape™ suture and suture threading device. This Implant System is used in conjunction with Arthrex Tenodesis Screws and a Tenodesis driver for use as a SwiveLock (Tenodesis) suture anchor.

B. INDICATIONS

The Arthrex Corkscrew, PushLock, and SwiveLock suture anchors are intended to be used for sutures (soft tissue) fixation to bone in the foot, ankle, knee, hand, wrist, elbow, shoulder, and hip (devices with FiberWire™ only):

I. IN NON CE ACCEPTING COUNTRIES:

- Corkscrew:**
 - Shoulder:** Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tendosis, Capsular Shift and Capsulolabral Reconstruction, Subscapularis Tendon Tears
 - Elbow:** Biceps Tendon Reattachment, Ulnar/Radial Collateral Ligament Reconstruction, and Lateral Epicondylitis repair (**except** Micro and Mini Corkscrew FT suture anchors)
 - Foot/Ankle:** Achilles Tendon Repair, Lateral Stabilization, Medial Stabilization, Achilles Tendn Repair, Mid-foot Reconstruction, Hallux Valgus Reconstruction, Metatarsal Ligament Repair, Metatarsal Tendon Repair, and Digital tendon transfers (with **Micro and Mini Corkscrew FT suture anchors only**)
 - Hand/Wrist:** Scapholunate Ligament Reconstruction and Ulnar/Radial Collateral Ligament Reconstruction, Digital Tendon Transfers with **Micro and Mini Corkscrew suture anchors**, Carpometaarpal Joint Arthroplasty (Basal Thumb Joint Arthroplasty) with **Micro and Mini Corkscrew FT suture anchors only**.
 - Knee:** Anterior Cruciate Ligament Repair, Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, and Iliotibial Band Tenodesis, Quadriceps Tendon Repair (Except in Canada) and Meniscal Root Repair (Except in Canada). Secondary or adjunct fixation for ACL/PCL reconstruction or repair (Except in Canada).
 - Hip:** Gluteus Medius Repair (5.5 - 6.5 mm PEEK and Titanium Corkscrew suture anchors only), Acetabular labral repair (**except** Micro and Mini Corkscrew FT suture anchors) and Proximal Hamstring Repair (Except in Canada)
- PushLock:**
 - Shoulder:** Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Symptomatic Acromioclavicular Joint Instability Biceps Tenodesis, Capsular Shift, and Capsulolabral Reconstruction
 - Elbow:** Biceps Tendon Reattachment, Ulnar/Radial Collateral Ligament Reconstruction, and Lateral Epicondylitis repair (**except** PushLock Tenodesis anchor, 2.4 mm and 2.5 mm PushLocks)
 - Foot/Ankle:** Achilles Tendon Repair, Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Mid-Foot Reconstruction, Hallux Valgus Reconstruction, Metatarsal Ligament Repair, Metatarsal Tendon Repair (**except** PushLock Tenodesis anchor), and Digital tendon transfers.
 - Hand/Wrist:** Carpometaacarpal Joint Arthroplasty (Basal Thumb Joint Arthroplasty) with PushLock Tenodesis anchor, Scapholunate Ligament Reconstruction and Ulnar/Radial Collateral Ligament Reconstruction, Digital Tendon Transfers with **2.5 mm PushLock**.
 - Knee:** Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, and Iliotibial Band Tenodesis, Quadriceps Tendon Repair (Except in Canada) and Meniscal Root Repair (Except in Canada). Secondary or adjunct fixation for ACL/PCL reconstruction or repair (Except in Canada).
 - Hip:** Acetabular labral repair (**except** PushLock Tenodesis anchor), Acetabular labral repair (**except** PushLock Tenodesis anchor), Proximal Hamstring Repair (Except in Canada), Gluteus Medius repair (Except in Canada)
- SwiveLock:**
 - Shoulder:** Rotator Cuff Repair, Biceps Tendosis, Capsulolabral Reconstruction, Subscapularis Tendon Tears, Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Capsular Shift, and Capsulolabral Reconstruction
 - Elbow:** Biceps Tendon Reattachment, Ulnar/Radial Collateral Ligament Reconstruction, and Lateral Epicondylitis repair (**except** DX SwiveLock SL, and Nano suture anchors only)
 - Foot/Ankle:** Achilles Tendon Repair, Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Mid-Foot Reconstruction, Hallux Valgus Reconstruction, Metatarsal Ligament Repair, Metatarsal Tendon Repair (**except** PushLock Tenodesis anchor), and Digital tendon transfers.
 - Hand/Wrist:** Scapholunate Ligament Reconstruction and Ulnar/Radial Collateral Ligament Reconstruction, Digital Tendon Transfers with **2.5 mm PushLock, DX SwiveLock SL, and Nano suture anchors only**.
 - Knee:** Anterior Cruciate Ligament Repair (4.75- 5.5 SwiveLock Only), Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, and Iliotibial Band Tenodesis, Quadriceps Tendon Repair (4.75 SwiveLock C Only) (Except in Canada) and Meniscal Root Repair (4.75 SwiveLock C Only) (Except in Canada). Secondary or adjunct fixation for ACL/PCL reconstruction or repair (Except in Canada).
 - Hip:** Acetabular labral repair (**except** PushLock Tenodesis anchor), Acetabular labral repair (**except** PushLock Tenodesis anchor), Proximal Hamstring Repair (Except in Canada), Gluteus Medius repair (Except in Canada)

labral repair **ONLY**.

- SwiveLock:**
 - Shoulder:** Rotator Cuff Repair, Biceps Tendosis, Capsulolabral Reconstruction, Subscapularis Tendon Tears, Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Capsular Shift, and Capsulolabral Reconstruction
 - Elbow:** Biceps Tendon Reattachment, Ulnar/Radial Collateral Ligament Reconstruction, and Lateral Epicondylitis repair (**except** DX SwiveLock SL, and Nano suture anchors only)
 - Foot/Ankle:** Achilles Tendon Repair, Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Mid-Foot Reconstruction, Hallux Valgus Reconstruction, Metatarsal Ligament Repair, Metatarsal Tendon Repair, and Digital tendon transfers (with **DX SwiveLock SL, and Nano suture anchors only**).
 - Hand/Wrist:** Scapholunate Ligament Reconstruction and Ulnar/Radial Collateral Ligament Reconstruction, Digital Tendon Transfers with **2.5 mm PushLock, DX SwiveLock SL, and Nano suture anchors only**.
 - Knee:** Anterior Cruciate Ligament Repair (4.75- 5.5 SwiveLock Only), Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, and Iliotibial Band Tenodesis, Quadriceps Tendon Repair (4.75 SwiveLock C Only) (Except in Canada) and Meniscal Root Repair (4.75 SwiveLock C Only) (Except in Canada). Secondary or adjunct fixation for ACL/PCL reconstruction or repair (Except in Canada).
 - Hip:** Acetabular labral repair (**except** PushLock Tenodesis anchor), Acetabular labral repair (**except** PushLock Tenodesis anchor), Proximal Hamstring Repair (Except in Canada), Gluteus Medius repair (Except in Canada)
- SwiveLock:**
 - Shoulder:** Rotator Cuff Repair, Biceps Tendosis, Capsulolabral Reconstruction, Subscapularis Tendon Tears, Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Capsular Shift, and Capsulolabral Reconstruction
 - Elbow:** Biceps Tendon Reattachment, Ulnar/Radial Collateral Ligament Reconstruction, and Lateral Epicondylitis repair (**except** DX SwiveLock SL, and Nano suture anchors only)
 - Foot/Ankle:** Achilles Tendon Repair, Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Mid-Foot Reconstruction, Hallux Valgus Reconstruction, Metatarsal Ligament Repair, Metatarsal Tendon Repair, and Digital tendon transfers (with **DX SwiveLock SL, and Nano suture anchors only**).
 - Hand/Wrist:** Scapholunate Ligament Reconstruction and Ulnar/Radial Collateral Ligament Reconstruction, Digital Tendon Transfers with **2.5 mm PushLock, DX SwiveLock SL, and Nano suture anchors only**.
 - Knee:** Anterior Cruciate Ligament Repair (4.75- 5.5 SwiveLock Only), Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, and Iliotibial Band Tenodesis, Quadriceps Tendon Repair (4.75 SwiveLock C Only) (Except in Canada) and Meniscal Root Repair (4.75 SwiveLock C Only) (Except in Canada). Secondary or adjunct fixation for ACL/PCL reconstruction or repair (Except in Canada).
 - Hip:** Acetabular labral repair (**except** DX SwiveLock SL, and Nano suture anchors), Gluteus Medius Repair (4.75- 5.5mm PEEK and Titanium SwiveLock suture anchors only), and Proximal Hamstring Repair (4.75-5.5 PEEK SwiveLock Only) (Except in Canada).

II. IN CE ACCEPTING COUNTRIES:

- Corkscrew:**
 - Shoulder:** Rotator Cuff Repair, SLAP Lesion Repair, Biceps Tenodesis, Capsulolabral Reconstruction, Subscapularis Tendon Tears
 - Foot/Ankle:** Achilles Tendon Repair
 - Hip:** Gluteus Medius Repair
- PushLock:**
 - Shoulder:** Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Symptomatic Acromioclavicular Joint Instability
 - Hand/Wrist:** Carpometaacarpal Joint Arthroplasty (Basal Thumb Joint Arthroplasty) with PushLock Tenodesis anchor, DX SwiveLock SL, Nano, Micro and Mini Corkscrew FT suture anchors only.
 - Foot/Ankle:** Achilles Tendon Repair
 - Hip:** Acetabular labral repair (**except** PushLock Tenodesis anchor)
- SwiveLock:**
 - Shoulder:** Rotator Cuff Repair, Biceps Tenodesis, Capsulolabral Reconstruction, Subscapularis Tendon Tears,
 - Knee:** Anterior Cruciate Ligament Repair (4.75- 5.5 SwiveLock Only)
 - Foot/Ankle:** Lateral Stabilization (for SwiveLock), Achilles Tendon Repair (for SpeedBridge)
 - Hip:** Gluteus Medius Repair

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Bioabsorbable only:** Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.

- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate. Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Bioabsorbable only:** Allergic-like reactions to PLA materials (PLLA, PDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
- Titanium only:** Shoulder dislocation/subluxation.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional. An internal fixation device must never be re-used.
- Do not re-sterilize this device.
- Titanium only:** All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
- Bioabsorbable only:** Attempting implantation into hard, dense bone and/or drilling/punching smaller diameter holes than recommended may cause failure (breakage) of the implant during insertion.
- Bioabsorbable Corkscrew suture anchor only:** The Arthrex 6.5 mm anchor should be used in soft bone only.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device, the patient leaflet (www.arthrex.com/patientleaflets) and the patient implant card that is given to the patient. Your surgeon will guide you in deciding what particular treatment is best for you and explain the benefits, risks, and contraindications associated with the treatment.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Bioabsorbable only:** Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the metal (titanium and stainless steel) Corkscrew, PushLock, and SwiveLock suture anchors are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 min of scanning in the Normal Operating Mode of operation for the MR system.

Under the scan conditions defined, the Corkscrew, PushLock, and SwiveLock suture anchors are expected to produce a maximum temperature rise of up to 1.8°C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Corkscrew, PushLock, and SwiveLock suture anchors can extend up to approximately 17 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

2. MR Safe

The Corkscrew, PushLock, and SwiveLock suture anchor devices manufactured from only polyetheretherketone (PEEK) poly (L-Lactide, PLLA), poly (L-Lactide-co-D-acid, PLDLA), and/or Poly (L-Lactide acid, PLLA) and tri-calcium phosphate (TCP) are MR safe. Arthrex Suture manufactured from UHMWPE and polyester with or without silicone elastomer coating, cyanoacrylate, and nylon are MR safe.

G. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate suture anchor size based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Make sure to use the recommended drill bit or punch to create the bone socket.
- PushLock and SwiveLock suture anchors only:** During anchor insertion, ensure that the angle of anchor insertion is coaxial to that of the previously prepared bone socket.
- PushLock and SwiveLock suture anchors only:** Insert the driver into the bone socket until the anchor body makes contact with the bone. Preview and adjust suture tension, if necessary. Tension will not increase during final advancement of the anchor body.
- PushLock and SwiveLock suture anchors only:** Ensure that the anchor body is in full contact with the bone before advancing the anchor body into the prepared bone socket.
- Self-punching suture anchors only:** Insertion in very hard bone may require pre-punching a bone socket to avoid damage to the implant.
- Self-Punching PushLock and SwiveLock suture anchors only:** Ensure that the angle of anchor insertion is perpendicular to the bone.
- Implant System, Tenodesis Screw Eyelet and PushLock Tenodesis anchor only:** Under insertion of the device may leave the proximal end of the implant protruding beyond the cortical bone, which could potentially cause soft tissue irritation and/or pain post-operatively.
- Implant System, Tenodesis Screw Eyelet only:** Ensure that the eyelet post is securely engaged into the end of the Tenodesis driver and held in place by tensioning the FiberTape towards the driver handle before inserting into the prepared bone socket.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolslossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Obtain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI S179 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

These devices consist of one or more components. Each component is manufactured from either/or titanium alloy, polyetheretherketone (PEEK), Poly (L-Lactide acid, PLLA), Poly (L-Lactide-co-D-acid, PLDLA), or Poly (L-Lactide acid, PLLA) and tri-calcium phosphate (TCP).

Suture (if supplied): See package label for size and type of suture provided with device.

The FiberWire, FiberWire CL, TigerWire®, TigerWire CL, TigerTape, and TigerTape™ sutures are made of Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester. The FiberWire and TigerWire sutures include FiberStick, TigerStick, FiberChain, FiberLink, TigerLink, TigerTail, FiberLoop and TigerLoop sutures. The FiberLink, TigerLink, FiberLoop and TigerLoop sutures may also be braided from SutureTape sutures. Additional materials may include silicone elastomer coating (except suture with the suffix -Tape), cyanoacrylate, and may include nylon.

SutureTape is made from Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester. Additional materials may include nylon and/or cyanoacrylate.

The sutures supplied meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements). The suture does may include: D&C Blue No. 6, D&C Green No.6, Logwood Black and FD&C Blue No. 2. Suture strands that are dyed with Logwood black are made of nylon. Other non-absorbable suture is made from polyester, polytetrafluoroethylene (PTFE) coating, and cyanoacrylate.

K. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.

Non-bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

- In CE Accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE Accepting Countries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE Accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- For Russia:** registered names of the medical devices: “Screw”, “Anchor”

- INFORMATION**
 - In CE Accepting Countries:** Procedures carried out using these devices may be used on the general population.
 - In CE Accepting Countries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
 - In CE Accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
 - For Russia:** registered names of the medical devices: “Screw”, “Anchor”

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