English
<p><b>Symbols glossary can be found at <a href="http://www.arthrex.com/symbolsglossary">www.arthrex.com/symbolsglossary</a>.</b></p>

#### A. DEVICE DESCRIPTION

The Unvers™ II Total Shoulder System consists of a stem for fixation within the humerus, a spherical head for replacing the humeral head, and a trunnion to connect the head to the stem. The stems and heads are available in a variety of sizes. The Unvers™ Apex Total Shoulder System consists of a stem for fixation within the humerus, a spherical head for replacing the humeral head, and a trunnion to connect the head to the stem. The stems and heads are available in a variety of sizes. The Unvers Apex stem is identical to the existing Unvers II stem, except for the overall length. The stem length has been shortened to optimize bone conservation. The stem and trunnions in the Unvers II and Unvers Apex Systems are comprised of titanium. The heads are available in cobalt alloy.

The Unvers VaultLock®, Keeled and Pegged glenoids are comprised of UHMWPE and available in a variety of sizes.

#### B. INDICATIONS

The Arthrex Unvers II Total Shoulder System and Unvers Apex Total Shoulder System are indicated in replacement(s) when conditions including severe pain or significant disability resulting from degenerative, rheumatoid, traumatic disease, or injury of the glenohumeral joint; non-union humeral head fractures of long duration; irreducible 2- and 4-part proximal humeral fractures; avascular necrosis of the humeral head, or other difficult clinical management problems where arthrodesis or resection arthroplasty is not acceptable.

The polyethylene glenoid components are intended for cemented fixation in the joint and must only be used with appropriate bone cement.

The Arthrex Titanium Humeral Head is indicated for patients with suspected cobalt alloy sensitivity. The wear properties of Titanium and Titanium alloys are inferior to that of cobalt alloy. A Titanium humeral head is not recommended for patients who lack suspected material sensitivity to cobalt alloy.

#### C. CONTRAINDICATIONS

- Insufficient quantities or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period, including severe neuro-arthropathy.
- Do not use for surgeries other than those indicated.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.

#### D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Allergies and other reactions to device materials.
- Loosening of the implant as a result of changed conditions in load transfer, respectively, fatigue wear and breakage of the cement bed and / or tissue reaction to the implant. Loosening is frequently a consequence of one or several of the above listed risk conditions, but can also be caused by inadequate anchoring technique (see below).
- Dislocation, subluxation, or inadequate scope of movement as a result of failure to achieve optimum positioning of the implant.
- Bone fractures as a result of one-sided overload or weakened bone substance.
- Temporary or permanent nerve damage as a result of pressure or hematoma.
- Cardiovascular diseases including venous thrombosis, pulmonary embolism, and cardiac arrest.
- Wound hematoma and delayed wound healing.
- Tissue reactions caused by allergic reactions to the implanted material, particularly metal, or caused by accumulations of wear particles or cement particles.

#### E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- Failure to achieve the appropriate torque requirements when tightening locking screws may result in the premature loosening of the device.
- The Unvers Apex Total Shoulder System stems (lengths 55-65 mm) are not recommended for fractures of the proximal humerus.
- Postoperatively, until healing is complete the fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the implant.
- Detailed instructions on the use and limitations of this device, the patient leaflet ([www.arthrex.com/patientleaflets](http://www.arthrex.com/patientleaflets)) and the patient implant card should be given to the patient. Your surgeon will guide you in deciding what particular treatment is best for you and explain the benefits, risks, and contraindications associated with the treatment.

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device.
- The following operative situations may cause premature loosening and complications:
  - Extreme weakening of the bone structure in preparing the bone bed;
  - Unsuitable selection of the implant size;
  - Inadequate cleaning of the bone bed prior to implantation; and,
  - Excessive use of force in placing or fastening the implant, provoking splintering fractures, or causing the bone to tear.
- An internal fixation device must never be re-used.
- Do not re-sterilize this device.
- The appropriate Arthrex delivery system is required for proper insertion of the implant.
- Only Arthrex delivery systems, instruments, and trial prostheses should be used for the implantation procedure.
- Endoprostheses may not be altered mechanically, or changed in any other way.
- Do not implant any parts that have been scratched or damaged.
- An artificial joint is subject to wear and/or can loosen over a period of time. Wear and loosening may make it necessary to re-operate on an artificial joint.
- An infection in an artificial joint may lead to implant removal.
- This device should only be used in conjunction with other implants designed specifically for use with this system.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

#### F. MRI SAFETY INFORMATION

##### 1. MR Conditional

*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Unvers II Total Shoulder System and Unvers Apex Total Shoulder System are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla
- Maximum spatial gradient magnetic field of 3000 Gauss/cm (30 T/m)
- Maximum MR system reported, with body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined above, the Unvers II Total Shoulder System and Unvers Apex Total Shoulder System are expected to produce a maximum temperature rise less than 2.87°C after 15-minutes of continuous scanning.

*In non-clinical testing, the image artifact caused by the device extends approximately 60 mm from the Unvers II Total Shoulder System and Unvers Apex Total Shoulder System when imaged with a gradient echo pulse sequence and a 3-Tesla MRI system.*

#### G. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

#### H. FACTORS AND RISKS IMPACTING THE SAFETY AND SERVICE LIFE OF THE IMPLANT

- Patient weight. An overweight patient may present additional risk.
- Extreme stress or strain resulting from work or sport related activity.
- Patients with increased risk of fractures due to repeated strain or trauma, or medical conditions that increase the patient’s risk for trauma, including falls.
- Osteoporosis or osteomalacia.
- Exposure to infectious diseases with possible manifestation in the joints.
- Deformation of the operative site, which can prevent or impede anchoring of the implant.
- Tumors that weaken the supportive structure.
- Allergic reactions to implant materials.
- Thrombosis and pulmonary infarction caused during implant site preparation and subsequent treatment.

#### I. PACKAGING AND LABELING

- Arthrex implants should be accepted only if the factory packaging and labeling arrive intact.
- All prosthesis parts and components should be stored unopened in the original packaging, and should be kept in protective packaging until immediately prior to use.
- Contact Customer Service if the package has been opened or altered.

#### J. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

#### K. MATERIAL SPECIFICATIONS

The device consists of a stem body manufactured of titanium, a metal head support of titanium alloy, and a head manufactured of cobalt-chromium (Co-Cr) alloy. Humeral heads are also offered in Titanium Alloy 6L-4V Eli. The glenoid component is composed of ultra-high molecular weight polyethylene (UHMWPE).

**Worldwide except Canada:** The head is also available as TiNbN coated cobalt-chromium (Co-Cr) alloy.

#### L. STORAGE CONDITIONS

Joint prostheses must be treated with great care prior to implantation. Scratches or nicks on the joint surface of the prosthesis can result in excessive warp and complications.

- Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.
- Do not remove the device from the package until immediately before use.

#### M. INFORMATION

- In CE Accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE Accepting Countries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE Accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

**1. MR Conditional**  
*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Unvers II Total Shoulder System and Unvers Apex Total Shoulder System are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:*

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