# English

#### A. DEVICE DESCRIPTION

These devices consist of the Bio-FASTak<sup>®</sup>, FASTak<sup>™</sup>, SutureTak<sup>™</sup>, SutureTak<sup>®</sup> and FiberTak<sup>®</sup> Suture Anchors.

These suture anchor families consist of rigid and threaded implants and all-suture implants (i.e. FiberTak). FiberTak suture anchors are comprised of a polyester sheath with single or multiple sutures assembled through the sheath. All implants are provided preloaded on an inserter and may include eyelets which are preloaded with various types of suture, with or without needles.

#### B. INDICATIONS

The Bio-FASTak, FASTak, SutureTak and FiberTak Suture Anchors are intended to be used for suture or tissue fixation in the foot, ankle, knee, hand, wrist, elbow, shoulder, and hip (devices with FiberWire<sup>®</sup> or #2 polyblend Suture).

The Bio-FASTak, FASTak, SutureTak and FiberTak Suture Anchors are intended for the following indications:

**Shoulder:** Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Capsular Shift, and Capsulolabral Reconstruction

**Foot/Ankle:** Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Mid-foot Reconstruction, Hallux Valgus Reconstruction (Except in EU), Metatarsal Ligament Repair (Except in EU), and Digital Tendon Transfers (**with SutureTak and FiberTak Suture Anchors only, except in EU**)

**Knee (Except in EU):** Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Iliotibial Band Tenodesis, and Joint Capsule Closure (**with Knotless SutureTak only**)

**Hand/Wrist (Except in EU):** Scapholunate Ligament Reconstruction, Digital Tendon Transfers, and Carpometacarpal joint arthroplasty (basal thumb joint arthroplasty) (**with Knotted FiberTak and Micro BioComposite SutureTak Suture Anchors only, except in EU**)

**Elbow:** Biceps Tendon Reattachment, Ulnar/Radial Collateral Ligament Reconstruction (Except in EU), and Lateral Epicondylitis Repair (**with Knotless SutureTak only, except in EU**)

**Hip:** acetabular labral repair, Capsular repair (**excludes FiberTak Suture Anchors**). Exclusions to hip indications include the 2.4 mm FASTak Suture Anchors.

The 2.0 mm Mini Hip SutureTak Suture Anchor is indicated for acetabular labral repair **ONLY**.

#### C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

#### D. ADVERSE EFFECTS

- Infection, both deep and superficial.
- Foreign body reactions.
- Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.

- FASTak Suture Anchors only:** Shoulder dislocation/subluxation.

#### E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be re-used.
- Do not re-sterilize this device.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.

#### F. MRI SAFETY INFORMATION

**1. MR Conditional**  
*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the metal (titanium) FASTakSuture Anchors are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, FASTak Suture Anchors are expected to produce a maximum temperature rise of up to 1.8<sup>o</sup> C after 15-minutes of continuous scanning.

*In non-clinical testing, the image artifact caused by the FASTak Suture Anchors can extend up to approximately 17 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*

- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device, the patient leaflet ([www.arthrex.com/patientleaflets](http://www.arthrex.com/patientleaflets)) and the patient implant card should be given to the patient. Your surgeon will guide you in deciding what particular treatment is best for you and explain the benefits, risks, and contraindications associated with the treatment..
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

#### F. MRI SAFETY INFORMATION

##### 1. MR Conditional

*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the metal (titanium) FASTakSuture Anchors are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, FASTak Suture Anchors are expected to produce a maximum temperature rise of up to 1.8<sup>o</sup> C after 15-minutes of continuous scanning.

*In non-clinical testing, the image artifact caused by the FASTak Suture Anchors can extend up to approximately 17 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*

##### 2. MR Safe

*The Bio-FASTak, SutureTak and FiberTak Suture Anchors devices manufactured from only polyester, polyetheretherketone (PEEK) poly (L-Lactide, PLLA), Poly (L-Lactide-co-D acid, PLDLA), and/or Poly (L-Lactide acid, PLLA) and tri-calcium phospbate (TCP) are MR safe.*

#### G. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate suture-anchor type and size based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- FASTak Suture Anchors only:** Stopping and restarting the drill may result in excess torque being applied to the implant, which may cause the device to fail.
- FASTak Suture Anchors only:** It is important to stop drilling as soon as the chuck contacts the guide or grasper. Failure to do so may result in damage to the suture and/or implant.
- Insert the anchor with the same orientation as that of the prepared drilled bone hole to avoid damage to the anchor.
- SutureTak and FiberTak Suture Anchors only:** Anchors Loaded on Flexible Drivers: The drill guide tip must remain in contact with the bone surface during the drilling and implant impact steps of the procedure. Failure to do so may result in difficulty seating the implant to its intended depth.

#### H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbolsglossary](http://www.arthrex.com/symbolsglossary).

#### I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

#### J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

**Anchor:** The device is made of titanium, PEEK (polyetheretherketone), Poly (L-Lactide) or Poly (L-Lactide-co-D, L-Lactide). FiberTak is made with a polyester sheath.

Biocomposite anchors are made of Poly (L-Lactide) or Poly (L-Lactide-co-D, L-Lactide) and tricalcium phosphate (TCP).

**Suture:** See package label for size and type of suture provided with device.

FiberWire, FiberWire CL, TigerWire<sup>®</sup>, TigerWire CL, FiberTape<sup>®</sup>, and TigerTape<sup>™</sup> sutures are made of Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester. Additional materials may include silicone elastomer coating (except suture with the suffix -Tape), cyanoacrylate, and may include nylon.

LabralTape<sup>™</sup> suture is made from Ultra High Molecular Weight Polyethylene (UHMWPE) and may include polyester or nylon.

SutureTape suture is made from Ultra High Molecular Weight Polyethylene (UHMWPE), and polyester. Additional materials may include nylon and/or cyanoacrylate.

The #2 Suture and #5 Suture are a polyblend made from Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester.

##### or

The suture is made from Polyester, polytetrafluoroethylene (PTFE) coating, cyanoacrylate.

The sutures supplied with the suture anchor meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements). The suture dyes may include: D&C Blue No. 6, D&C Green No. 6, and Logwood Black. Suture strands that are dyed black are made of nylon.

#### K. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32<sup>o</sup> C / 90<sup>o</sup> F, and should not be used after the expiration date.

Non-Bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

##### L. INFORMATION

- In EU Only:** Procedures carried out using these devices may be used on the general population.
- In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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