# English

**A. DEVICE DESCRIPTION**

The Arthrex Corkscrew®, PushLock®, and SwiveLock® suture anchors consist of anchors with an integral or separate eyelet. The PushLock Tenodesis anchor is a two-piece push-in anchor system with either a forked or closed eyelet. They are pre-loaded on a handled inserter. Suture, with or without needles, and a suture threader may also be provided.

The Arthrex Implant System, Tenodesis Screw Eyelet is a kit comprised of attachable eyelets, FiberTape® suture and a suture threading device. This Implant System is used in conjunction with Arthrex Tenodesis Screws and a Tenodesis driver for use as a SwiveLock (Tenodesis) suture anchor.

**B. INDICATIONS**

The Arthrex Corkscrew, PushLock, and SwiveLock suture anchors are intended to be used for suture (soft tissue) fixation to bone in the foot, ankle, knee, hand, wrist, elbow, shoulder, and hip (devices with FiberWire® only): Acetabular labral repair (**except** DX SwiveLock SL, 2.4 mm and 2.5 mm PushLock, Nano, Micro and Mini Corkscrew FT suture anchors). The Arthrex PushLock Tenodesis anchor is intended to provide soft tissue reattachment (i.e. fixation of ligament and tendon graft tissue) in surgeries of the shoulder, elbow, knee, foot/ankle, and hand/wrist.

**Shoulder:** Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Capsular Shift, and Capsulolabral Reconstruction

**Foot/Ankle:** Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Mid-foot Reconstruction, Hallux Valgus Reconstruction (Except in EU), Metatarsal Ligament Repair (Except in EU), Metatarsal Tendon Repair (**except** PushLock Tenodesis anchor, except in EU), and Digital tendon transfers (with **DX SwiveLock SL, Nano, Micro and Mini Corkscrew FT suture anchors only**, (Except in EU)

**Knee:** Anterior Cruciate Ligament Repair (4.75- 5.5 **SwiveLock Only**), Medial Collateral Ligament Repair (Except in EU), Lateral Collateral Ligament Repair (Except in EU), and Iliotibial Band Tenodesis (Except in EU)

**Hand/Wrist:** Scapholunate Ligament Reconstruction (Except in EU) and Ulnar/Radial Collateral Ligament Reconstruction (Except in EU). Digital Tendon Transfers with **2.5 mm PushLock, DX SwiveLock SL, Nano, Micro and Mini Corkscrew suture anchors only**, (Except in EU). Carpometacarpal Joint Arthroplasty (Basal Thumb Joint Arthroplasty) with **PushLock Tenodesis anchor, DX SwiveLock SL, Nano, Micro and Mini Corkscrew FT suture anchors only**.

**Elbow:** Biceps Tendon Reattachment, Ulnar/Radial Collateral Ligament Reconstruction (Except in EU), and Lateral Epicondylitis repair (**except** PushLock Tenodesis anchor, 2.4 mm and 2.5 mm PushLocks, DX SwiveLock SL, Nano, Micro and Mini Corkscrew FT suture anchors, except in EU)

**Hip:** Acetabular labral repair (**except** PushLock Tenodesis anchor, DX SwiveLock SL, Nano, Micro and Mini Corkscrew FT suture anchors)

The 2.4 mm Hip PushLock suture anchor is indicated for acetabular labral repair **ONLY**.

**Gluteus Medius (Except in EU):** 4.75 – 5.5 mm PEEK and Titanium SwiveLock suture anchors and 5.5 - 6.5 mm PEEK and Titanium Corkscrew suture anchors

**C. CONTRAINDICATIONS**

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Bioabsorbable only:** Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

**D. ADVERSE EFFECTS**

- Infections, both deep and superficial.
- Foreign body reactions.
- Bioabsorbable only:** Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
- Titanium only:** Shoulder dislocation/subluxation.

**E. WARNINGS**

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be reused.
- Do not re-sterilize this device.
- Titanium only:** All metallic implant devices used for this surgical procedure should have the same metallurgical composition.

- Bioabsorbable only:** Attempting implantation into hard, dense bone and/or drilling/punching smaller diameter holes than recommended may cause failure (breakage) of the implant during insertion.

- Bioabsorbable Corkscrew suture anchor only:** The Arthrex 6.5 mm anchor should be used in soft bone only. Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.

- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.

- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.

- Detailed instructions on the use and limitations of this device, the patient leaflet and the patient implant card should be given to the patient.

- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.

- Bioabsorbable only:** Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.

- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

**F. MRI SAFETY INFORMATION**

**1. MR Conditional**

*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the metal (titanium and stainless steel) Corkscrew, PushLock, and SwiveLock suture anchors are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the Corkscrew,PushLock, and SwiveLock suture anchors are expected to produce a maximum temperature rise of up to 1.8° C after 15-minutes of continuous scanning.

*In non-clinical testing, the image artifact caused by the Corkscrew, PushLock, and SwiveLock suture anchors can extend up to approximately 17 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*

**2. MR Safe**

*The Corkscrew, PushLock, and SwiveLock suture anchor devices manufactured from only polyetheretherketone (PEEK) poly (L-Lactide, PLLA), Poly (L-Lactide-co-D acid, PLDLA), and/or Poly (L-Lactide acid, PLLA) and tri-calcium phosphate (TCP) are MR safe.*

**G. PRECAUTIONS**

- Surgeons must apply their professional judgment when determining the appropriate suture anchor size based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Make sure to use the recommended drill bit or punch to create the bone socket.
- PushLock and SwiveLock suture anchors only:** During anchor insertion, ensure that the angle of anchor insertion is coaxial to that of the previously prepared bone socket.
- PushLock and SwiveLock suture anchors only:** Insert the driver into the bone socket until the anchor body makes contact with the bone. Preview and adjust suture tension, if necessary. Tension will not increase during final advancement of the anchor body.
- PushLock and SwiveLock suture anchors only:** Ensure that the anchor body is in full contact with the bone before advancing the anchor body into the prepared bone socket.
- Self-punching suture anchors only:** Insertion in very hard bone may require pre-punching a bone socket to avoid damage to the implant.
- Self-Punching PushLock and SwiveLock suture anchors only:** Ensure that the angle of anchor insertion is perpendicular to the bone.
- Implant System, Tenodesis Screw Eyelet and PushLock Tenodesis anchor only:** Under insertion of the device may leave the proximal end of the implant protruding beyond the cortical bone, which could potentially cause soft tissue irritation and/or pain post-operatively.
- Implant System, Tenodesis Screw Eyelet only:** Ensure that the eyelet post is securely engaged into the end of the Tenodesis driver and held in place by tensioning the FiberTape towards the driver handle before inserting into the prepared bone socket.

**H. PACKAGING AND LABELING**

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.

- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbolsglossary](http://www.arthrex.com/symbolsglossary).

**I. STERILIZATION**

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

**J. MATERIAL SPECIFICATIONS**

Refer to the package label for the materials.

These devices consist of either one or two components. Each component is manufactured from either/or titanium alloy, polyetheretherketone (PEEK), Poly (L-Lactide acid, PLLA), Poly (L-Lactide-co-D acid, PLDLA), or Poly (L-Lactide acid, PLLA) and tri-calcium phosphate (TCP).

**Suture (if supplied):** See package label for size and type of suture provided with device.

The FiberWire, FiberWire CL, TigerWire®, TigerWire CL, FiberTape, and TigerTape™ sutures are made of Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester. Additional materials may include silicone elastomer coating (except suture with the suffix -Tape), cyanoacrylate, and may include nylon.

SutureTape is made from Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester. Additional materials may include nylon and/or cyanoacrylate.

The sutures supplied meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements). The suture dyes may include: D&C Blue No. 6, D&C Green No. 6, and Logwood Black. Suture strands that are dyed black are made of nylon. Other non-absorbable suture is made from polyester, polytetrafluoroethylene (PTFE) coating, and cyanoacrylate.

**K. STORAGE CONDITIONS**

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.

Non-bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

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