English

A. REFERENCES

These instructions were developed using the guidance given in the following standards:

- ANSI AAMI S797, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities"
- ISO 17664: Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices
- ISO 17665-1: Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices
- ISO 10993-5: Biological evaluation of medical devices – Part 5: Tests for in vivo cytotoxicity
- AAMI TR30:2011: A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices

B. DEVICE DESCRIPTION AND INFORMATION

This device may be a reusable or a single-use non-sterile instrument. Or, it may be a reusable sterile or non-sterile instrument that may be attached to power. Check the package labeling.

Users of this device are encouraged to contact their Arthrex representatives if, in their professional judgment, they require a more comprehensive surgical technique. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations.

C. LIMITATIONS ON REPROCESSING

For reusable devices: Repeated processing has minimal effect on these instruments. End of life is normally determined by wear and damage due to use.

A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that may include, but are not limited to cross-infection, breakage resulting in irremovable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.

D. VALIDATION

The recommended cleaning, disinfection, and sterilization methods in this IFU have been validated in compliance with federal and international guidance/standards. In accordance with ISO 17665, the half cycle "overkill" approach was used for sterilization validation, and demonstrates a sterility assurance level (SAL) of 10⁻⁶. Cleaning, disinfecting, and sterilizing equipment and materials vary in performance characteristics. Therefore, it is the responsibility of the facility/end user to ensure equipment cleared by the FDA or the local governing body and perform operations in accordance with directions from this IFU.

In accordance with DIN EN ISO 17664 limit values and a means for monitoring chemical residues following cleaning shall be established for the product. In assessing the level of cleaning residues following the manual cleaning and disinfection process or the machine (automatic) cleaning and disinfection process, it is recommended that cytotoxicity testing be utilized as a clinically relevant method for testing the safety of residuals. Testing for cytotoxicity should be completed and scored in accordance with ISO 10993-5. Biological evaluation of medical devices – Part 5: Tests for in vivo cytotoxicity.

The water quality for the cleaning and disinfection process has been validated to ensure that residuals will not interfere with subsequent processing steps. Testing for cytotoxicity as the primary means of assessing the safety of the detergent and/or disinfectant levels after the cleaning and washing processes was carried out in the validation protocol in accordance with AAMI TR30:2011. Deionized water is the validated water quality type used for manual and automated cleaning.

E. CONTAINMENT AND TRANSPORTATION

It is recommended that instruments are reprocessed as soon as is reasonably practical following use. Any visible soil must be removed at point-of-use to prevent the drying of soil. Instrument cases and trays are considered reusable devices. Trays should be inspected for visible soil and must be cleaned prior to use. They can be cleaned manually or in an automatic washer using a detergent.

F. CLEANING

I. PRELIMINARY CLEANING

When properly formed, cleaning, disinfection and/or sterilization do not compromise the use and mechanical performance of these instruments. These instruments are used with or on patients who may harbor unrecognized and unrecognized infections. To prevent the spread of infection, all reusable instruments must be thoroughly cleaned, disinfected, and sterilized after each patient.

Note: No assembly/disassembly of these instruments is required unless stated on the labeling, directions for use, or literature assembly instructions (LAI) pertaining to cleaning, disinfection, and sterilization.

- Devices that require disassembly should be disassembled prior to cleaning.
- Remove dried-on soil from devices, especially in areas such as joints and crevices, prior to washing. Clean the surfaces with a sponge or brush under cold running water until visible contamination has been removed. Pre-cleaning may also be completed utilizing an ultrasonic bath. If cleaning ultrasonically, place the device in the ultrasonic unit for a minimum of 10 minutes and process according to the ultrasonic unit's directions.
- Check the instruments for visible soil. Repeat the preliminary cleaning if soil is visible, then re-inspect.
- After the completion of preliminary cleaning, the end user has the option to perform either Manual Cleaning (Section F II) and Disinfection or Machine (Automatic) Cleaning and Thermal Disinfection (Section F III).

II. MANUAL CLEANING AND DISINFECTION

After the preliminary cleaning, the instructions for Manual Cleaning and Disinfection may be used as an alternative cleaning method to the Machine (Automatic) Cleaning and Thermal Disinfection.

CAUTION: Manual cleaning is not an appropriate method for cleaning hand instruments containing non-exposed actuating components, e.g. wishbone handle or standard handgrip thersers, punches, etc. The Machine (Automatic) Cleaning and Thermal Disinfection procedure should be followed for these instruments.

- Immerse the instrument in an enzymatic or alkaline cleaning detergent solution. Cleaning solutions may include, but are not limited to: ENZOL® enzymatic, neodisher® Mediclean forte, and ThermoSep® alk alean. **CAUTION:** Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinated ethylenepropylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™. If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing, as validated by the end-user facility, and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device. Disinfection solutions should always be mixed to the manufacturer's specification for concentration and cleaning should be conducted at ambient temperature unless otherwise stated in the cleaning solution manufacturer's instructions.
- Scrub the instrument with a soft bristled brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the instrument. Immerse the instrument in detergent, agitate and allow it to soak for at least one minute.
- Rinse the instrument thoroughly with cold distilled water for at least one minute after the cleaning process.
- Immerse instruments in disinfection solution for a minimum of 20 minutes. Suitable disinfection solutions may include, but are not limited to: CIDEX®/ WAWCID®/0-1, GigaSep®/0, Kohrsolid®/0, and GigaSep®/0. Use the supplier's instructions for preparing the solution. **CAUTION:** Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP

(Fluorinated ethylenepropylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™. If non-neutral pH disinfection chemistries are utilized, care should be taken to ensure appropriate rinsing, as validated by the end-user facility, and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device. Disinfection solutions should always be mixed to the manufacturer's specification for concentration and cleaning should be conducted at ambient temperature unless otherwise stated in the cleaning solution manufacturer's instructions.

- After disinfection, the instruments should be rinsed for at least one minute with cold distilled water or deionized sterile water.
- Dry instruments thoroughly utilizing compressed air, wipes, or an oven.
- Check the instruments for visible soil. Repeat the cleaning if soil is visible, then re-inspect.

III. MACHINE (AUTOMATIC) CLEANING, AND THERMAL DISINFECTION

After the preliminary cleaning, the instructions for Machine (Automatic) Cleaning and Thermal Disinfection may be used as an alternative cleaning method to the Manual Cleaning and Disinfection.

- Load the instruments in the washer so that all design features of the device are accessible to cleaning and any that may retain liquid can drain (things should be open and cam/latches/holes positioned to drain).
- Run the automatic wash cycle. The minimum cycle parameters are listed below:

| Phase | Recirculation Time | Temperature | Detergent |
|-----------------------------------|--------------------|------------------------|-----------------------------|
| Cold Pre-Wash | 3 Minutes | 20 ± 5°C (68 ± 9°F) | N/A |
| Cleaning Wash | 10 Minutes | 65.5 ± 5°C (150 ± 9°F) | Enzymatic or alkaline agent |
| Rinse 1 | 3 Minutes | 50 ± 5°C (122 ± 9°F) | N/A |
| Rinse 2 | 3 Minutes | 50 ± 5°C (122 ± 9°F) | N/A |
| Thermal Disinfection Rinse | 5 Minutes | 90 ± 5°C (194 ± 9°F) | N/A |
| Drying | 6 Minutes | 115°C (239°F) | N/A |

- Automatic wash cleaning solutions may include, but are not limited to: neodisher® Mediclean forte, ThermoSep® alk alean, Prolystica® Ultra Concentrate Enzymatic Cleaner, and Proline® Nyl Instrument Lubricant and Antifog. **Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinated ethylenepropylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™.** If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing, as validated by the end-user facility, and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.

- Check the instruments for visible soil. Repeat the cleaning if soil is visible, then re-inspect.
- Automatic wash cleaning solutions may include, but are not limited to: neodisher® Mediclean forte, ThermoSep® alk alean, Prolystica® Ultra Concentrate Enzymatic Cleaner, and Proline® Nyl Instrument Lubricant and Antifog. **Low acid or high alkaline solutions are not recommended, as they corrode metal parts and anodized aluminum and compromise polymer plastics, such as FEP (Fluorinated ethylenepropylene), ABS (Acrylonitrile Butadiene Styrene), Utem™, Lexan™, and Cyoclear™.** If non-neutral pH cleaning chemistries are utilized, care should be taken to ensure appropriate rinsing, as validated by the end-user facility, and neutralization steps are taken so as to not negatively impact the fit, finish, or function of the device.

- The water quality for the cleaning and disinfection process has been validated to ensure that residuals will not interfere with subsequent processing steps. Testing for cytotoxicity as the primary means of assessing the safety of the detergent and/or disinfectant levels after the cleaning and washing processes was carried out in the validation protocol in accordance with AAMI TR30:2011. Deionized water is the validated water quality type used for manual and automated cleaning.
- Devices with cutting functions or sharp points become dull with continuous use. This condition does not indicate a device defect. This condition indicates normal wear. Dull devices may require replacement, if they no longer perform as designed. Inspection prior to use should include visually verifying the cutting ability and sharpness of these tools.
- If damage is detected, do not use the device prior to consulting the manufacturer for guidance.
- Dry the instruments thoroughly and lubricate all moving parts with a water-soluble instrument lubricant, prior to sterilization. Acceptable lubricants may include, but are not limited to: Sieris Hinge-Free® Instrument Lubricant and neodisher® IP Spray. Apply lubricants in accordance with manufacturer's instructions.

G. INSPECTION AND MAINTENANCE

- Arthrex non-sterile instruments are precision mechanical instruments and must be used and handled prior to use and at all stages of handling.
- Devices with cutting functions or sharp points become dull with continuous use. This condition does not indicate a device defect. This condition indicates normal wear. Dull devices may require replacement, if they no longer perform as designed. Inspection prior to use should include visually verifying the cutting ability and sharpness of these tools.
- If damage is detected, do not use the device prior to consulting the manufacturer for guidance.
- Dry the instruments thoroughly and lubricate all moving parts with a water-soluble instrument lubricant, prior to sterilization. Acceptable lubricants may include, but are not limited to: Sieris Hinge-Free® Instrument Lubricant and neodisher® IP Spray. Apply lubricants in accordance with manufacturer's instructions.

H. PACKAGING

Singley: Single device should be packed to ensure that the pack is large enough to contain the instrument without stressing the seals. Packaging should be completed using a pouch or wrap, which conforms to the recommended specifications for steam sterilization as outlined above. If a wrap is used, it should be completed following AAMI or other applicable standards. The wrap should be wrapped with an appropriate wrap. An appropriate wrap is one that, for example, is cleared by the FDA or the local governing body at the point of use.

- Devices that require disassembly should be disassembled prior to cleaning.
- Remove dried-on soil from devices, especially in areas such as joints and crevices, prior to washing. Clean the surfaces with a sponge or brush under cold running water until visible contamination has been removed. Pre-cleaning may also be completed utilizing an ultrasonic bath. If cleaning ultrasonically, place the device in the ultrasonic unit for a minimum of 10 minutes and process according to the ultrasonic unit's directions.
- Check the instruments for visible soil. Repeat the preliminary cleaning if soil is visible, then re-inspect.
- After the completion of preliminary cleaning, the end user has the option to perform either Manual Cleaning (Section F II) and Disinfection or Machine (Automatic) Cleaning and Thermal Disinfection (Section F III).

I. STERILIZATION

This device may be provided either sterile or non-sterile. Check the package labeling for more information. For instruments that are not provided in a terminally sterilized configuration, sterilization must be performed following cleaning, disinfection, and sterile packaging.

Certain Arthrex devices that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or reuse. (See above instructions for cleaning.) Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions.

Follow your country-specific guidelines, standards, and requirements.

| STERILIZATION PARAMETERS: FOR THE USA ONLY: | | | |
|---|----------------------|---------------|------------------|
| | Exposure Temperature | Exposure Time | Drying Time |
| Gravity-Displacement Steam Sterilization Cycle | 121°C (250°F) | 30 Minutes | 15 to 30 Minutes |
| | 132°C (270°F) | 15 Minutes | 15 to 30 Minutes |
| | 135°C (275°F) | 10 Minutes | 30 Minutes |
| Pre-vacuum Cycle | 132°C (270°F) | 4 Minutes | 20 to 30 Minutes |
| | 135°C (275°F) | 3 Minutes | 16 Minutes |

| STERILIZATION PARAMETERS: FOR OUTSIDE THE USA ONLY: | | | |
|---|-------------------------------|---------------|------------------|
| | Exposure Temperature | Exposure Time | Drying Time |
| Gravity-Displacement Steam Sterilization Cycle | 132°C – 135°C (270°F – 275°F) | 18 Minutes | 15 to 30 Minutes |
| | 121°C (250°F) | 30 Minutes | 15 to 30 Minutes |
| | 132°C – 135°C (270°F – 275°F) | 4 Minutes | 20 to 30 Minutes |

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L. CAUTIONS

- Users of this device are encouraged to contact their Arthrex representatives if, in their professional judgment, they require a more comprehensive surgical technique or more information. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations.
- To avoid damaging the instruments, do not impact or subject to blunt force any instruments that are designed to be turned or screwed in. When two devices are intended to be threaded together, ensure that they are fully engaged prior to use.
- Do not use Arthrex instruments for any purpose other than their intended use. Manipulating soft tissue or bone with an instrument not intended for that use may result in damage to the instrument.
- Instruments with adjustable components must be handled with care. Overtightening or rough handling of the instrument may damage the locking mechanism. Locking mechanisms with internal polymer components may become weakened after repeated autoclaving.
- Do not use an instrument that is intended to be used with a specific implant on another implant.
- Flexion of the joint with the instrument in position in the joint may result in bending or breakage of the instrument.

- Instruments with adjustable components must be handled with care. Overtightening or rough handling of the instrument may damage the locking mechanism. Locking mechanisms with internal polymer components may become weakened after repeated autoclaving.
- Do not use an instrument that is intended to be used with a specific implant on another implant.
- Flexion of the joint with the instrument in position in the joint may result in bending or breakage of the instrument.

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