English

Symbols glossary can be found at
www.arthrex.com/symbolsglossary

A. DEVICE DESCRIPTION

The TightRope® devices are comprised of various configurations of one or two metal buttons, a metal or bioabsorbable anchor, and suture. Single buttons (clavicle and coracoid) and washers may be provided in individual packages without sutures. These devices may be preloaded onto an inserter, or may be packaged with various ancillary instruments to aid in insertion.

The Syndesmosis TightRope XP Buttress Plate Implant System is comprised of one two-hole buttress plate, four metal buttons and #5 suture. The metal buttons and #5 suture are preloaded onto a disposable inserter. The ancillary two-hole syndesmosis buttress plate is a metal plate with two holes used with the Syndesmosis TightRope XP devices. The devices are packaged with various ancillary instruments to aid in insertion.

The ACL/PCL/BTB/RT/ABS TightRope products are comprised of suture with or without a button, wedge or inserter. TightRope ABS implants are used with TightRope ABS buttons (standard, round and oblong buttons) and TightRope button extender.

The ancillary buttress plates are metal plates with two or four holes, depending on design, for TightRope or screw fixation.

B. INDICATIONS

The TightRope devices are intended as adjuncts in fracture repair involving metatarsal and periarticular small bone fragments where screws are not indicated, and as adjuncts in external and intramedullary fixation systems involving plates and rods, with fracture braces and casting.

Specifically, they are intended to provide fixation during the healing process following:

- Mini TightRope, Mini TightRope FT and CMC Mini TightRope: hallux valgus reconstruction (correction) by providing for the reduction of 1st metatarsal-2nd metatarsal intermetatarsal angle
- Mini TightRope, Mini TightRope FT and CMC Mini TightRope: tarsometatarsal (TMT) injury, such as fixation of foot soft tissue separations due to a Lisfranc injury (midfoot reconstruction)
- Mini TightRope, Mini TightRope FT and CMC Mini TightRope: syndesmotic trauma, such as fixation of dorsal distal radioulnar ligament (DRUL) disruptions
- Syndesmosis TightRope: syndesmotic trauma, such as fixation of syndesmosis (syndesmosis disruptions) in connection with Weber B and C ankle fractures
- Acromioclavicular (AC) TightRope, twin-tail TightRope, and Dog Bone™ button used with FiberTape™: syndesmotic trauma, such as acromioclavicular separation due to coracoclavicular ligament disruption
- Arthrex recommends the use of #5 FiberWire® or equivalent with the clavicle and coracoid buttons and washers.

The Arthrex Mini TightRope, Mini TightRope FT and CMC Mini TightRope, are used for fixation of bone-to-bone or soft-tissue-to-bone, are intended as fixation posts, distribution bridges, or for distributing suture tension over areas of ligament or tendon repair. Specifically, the Arthrex Mini TightRope, the Mini TightRope FT and the CMC Mini TightRope are indicated for carpal metacarpal (CMC) joint arthroplasty as an adjunct in the healing process of the reconstruction of the ligament at the base of the thumb metacarpal by providing stabilization between the base of the first and second metacarpal when the trapezium has been excised due to osteoarthritis. The Mini TightRope, Mini TightRope FT and CMC Mini TightRope are also indicated for use as adjuncts in the suspension of the thumb metacarpal during the healing process of hematoma distraction arthroplasty by providing stabilization at the base of the first and second metacarpal when the trapezium has been excised due to osteoarthritis (except AR-S917DS Mini TightRope FT).

The Arthrex Syndesmosis TightRope XP Buttress Plate Implant System is intended as an adjunct in fracture repair involving metapsal and periarticular small bone fragments where screws are not indicated, and as an adjunct in external and intramedullary fixation systems involving plates and rods, with fracture braces and casting.

Specifically, the Arthrex Syndesmosis TightRope XP Buttress Plate Implant System is intended to provide fixation during the healing process following a syndesmotic trauma, such as fixation of syndesmosis (syndesmosis disruptions) in connection with Weber B and C ankle fractures.

The Arthrex ACL/PCL/ BTB/RT/ABS TightRope and TightRope II are used for fixation of bone to bone or soft tissue to bone, and are intended as fixation posts, a distribution bridge, or for distributing suture tension over areas of ligament or tendon repair. Specifically, Arthrex will be offering these for ACL/PCL repair and reconstruction.

- Arthrex recommends the use of #2 or #5 FiberWire or equivalent with the suture button.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign body reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, distort or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Syndesmosis TightRope XP Syndesmosis TightRope, Mini TightRope, Mini TightRope FT and CMC Mini TightRope Only:** Osteomyelitis surrounding the TightRope.
- Syndesmosis TightRope XP, Syndesmosis TightRope, Mini TightRope, Mini TightRope FT and CMC Mini TightRope Only:** Rediastasis resulting from a failure of the implant insertion technique.
- Syndesmosis TightRope, Mini TightRope, Mini TightRope FT and CMC Mini TightRope Only:** FiberWire polyethylene wear-related painful neuropathic osteoarthritis.
- Infections, both deep and superficial.
- Foreign body reactions.
- Bio only** will be allergic-like reactions to PLA materials (PLLA, PLDLA) are being reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
- Silicone sensitivity, though very rare, has been reported.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- The ACL/PCL/ BTB/RT/ABS TightRope devices are not intended to be used as a ligament replacement.
- Do not add additional suture to the ACL TightRope. The extra suture may impede passage of the device through the femur.
- Do not re-sterilize this device.

- All metallic implant devices used for this surgical procedure should have the same metallurgical composition, including the plates that might be used in conjunction with the Arthrex TightRope syndesmosis repair kit implants. Therefore, different material implants are provided.
- AC TightRope repair kit, twin-tail TightRope, Dog Bone button with FiberTape:** Do not use this device as the sole means of reconstructing a chronic acromioclavicular joint dislocation.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device and bone.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device, the patient leaflet (www.arthrex.com/patientleaflets) and the patient implant card should be given to the patient. Your surgeon will guide you in deciding what particular treatment is best for you and explain the benefits, risks, and contraindications associated with the treatment.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
- BIOABSORBABLE:** Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.
- Once open, discard any unused suture.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-company representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the buttress plates and suture buttons manufactured from titanium or stainless steel are MR Conditional. A patient with this device can be scanned safely in an MR system using the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla
 - Maximum spatial field gradient 3000 Gauss/cm (30 T/m)
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)
- Under the scan conditions defined above, the buttress plates and suture buttons are expected to produce a maximum temperature rise of less than 2.5 °C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the buttress plates and suture buttons, can extend approximately 40 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

2. MR Safe

Arthrex suture manufactured from UHMWPE and polyester with or without silicone elastomer coating, cyanacrylate, and nylon are MR safe.

G. PRECAUTIONS

- ACL/PCL/ BTB/RT/ABS TightRope only:** Excessive force on the shortening suture strands may break the strands and impair ability to fully seat the implant. No additional force on the shortening suture strands is required when the graft/wedge construct reaches the desired position in the femoral socket and the graft stability is verified by pulling distally on the graft.
- Lead the ACL/PCL TightRope button over the unspliced, thinner portion of the ACL/PCL TightRope suture to assist assembly.** Once assembled slide the ACL/PCL TightRope button down to the thicker, spliced portion of the ACL/PCL TightRope suture to help prevent disassembly.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

Anchor: The anchor is made of titanium or poly (L-lactic acid) PLLA. **Buttons/Plates:** The buttons and plates are available in two different metals in order to match the metallurgic composition of implants that may be used in the surgical procedure. The materials available are titanium and stainless steel.

- Needles (if applicable):** Stainless steel (may be coated in silicone).
- Suture:** Ultra-high molecular weight polyethylene (UHMWPE), or ultra-high molecular weight polyethylene and polyester.
- FiberWire and TigerWire®:** a braided suture made of UHMWPE and polyester over a core of UHMWPE yarns. Some braided polybled sutures or UHMWPE sutures may be corless. Additional materials to the Arthrex sutures include silicone elastomer coating, and may include cyanocrylate and/or nylon. The coating acts as a lubricant for suture sliding, knot tying, and ease of passing suture through tissue. FiberWire CL and other corless polybled suture may be uncoated.
- FiberTape, TigerTape™, LabraTape™ and SutureTape™** are flat braided constructs. FiberTape and TigerTape suture is composed of UHM-WPE and polyester yarns over a core of FiberWire or TigerWire. LabralTape and SutureTape are made of UHMWPE but may include polyester yarns and/or nylon. SallTape™ is a woven construct made of UHMWPE and polyester with tapered suture tails made of UHMWPE. Additional materials to the Arthrex tape constructs may include cyanocrylate and/or nylon.

The FiberWire and TigerWire sutures include FiberStick™, TigerStick®, FiberChain™, FiberLink™, TigerLink™, FiberLoop™, TigerLoop™, FiberTape loop, TigerTape loop, and FiberRing™ sutures and are assembled from FiberWire, FiberTape, or SutureTape sutures.

The sutures meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements). The diameters of the sutures typically span the diameter requirements of two USP sizes. For that reason, Arthrex tests the sutures to conform to the greater USP size specification requirement. SutureTape sutures meet or exceed USP standards for size 2 or size 2-0 non-absorbable surgical sutures, except for an oversize in diameter for the round tails and width of flat tape.

Additional materials to the sutures include silicone elastomer coating (except suture with the suffix -Tape), and may include cyanocrylate and/ or nylon. The coating acts as a lubricant for suture sliding, knot tying, and ease of passing suture through tissue. Dyes may include: D&C Blue No. 6, FD&C Blue No. 2, D&C Green No. 6, Chromium-cobalt-aluminum oxide (Co-Cr-Al oxide) blue, D&C Black No. 4 and Logwood Black.

#2 Suture: Ultra-high molecular weight polyethylene (UHMWPE) and polyester.

Wedge (if present): Polyetheretherketone (PEEK).

K. STORAGE CONDITIONS
Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.

Non-bioabsorbable devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

- In CE Accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE Accepting Countries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE Accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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