English

A. DEVICE DESCRIPTION

The Univers Revers™ Shoulder Prosthesis System has an articular design that is inverted, compared to traditional total shoulder prosthesis. The system is comprised of two main components: the Arthrex Univers Revers Shoulder Prosthesis and the Universal Glenoid™ Shoulder Prosthesis. The Arthrex Univers Revers Shoulder Prosthesis is a humeral stem and epiphysis or humeral cup, a spacer, and humeral cup liner. The humeral stem and epiphysis are available uncoated or coated.

The Universal Glenoid Shoulder Prosthesis consists of a glenoid baseplate, a glenosphere, and screws.

The Univers Revers Modular Glenoid System consists of a monoblock baseplate or a modular baseplate; both baseplates are available with either a central screw or central post. The baseplate is designed to be used cementless with peripheral screws and a glenosphere. The glenoid system is designed to be used as the glenoid side of the existing Univers Revers Shoulder Prosthesis System.

The Univers Revers Cuff Arthropathy System is designed to be used with the existing well-fixed Revers stem, or to convert an existing reversed shoulder prosthesis to a hemi anatomic configuration. The Revers Cuff Arthropathy Humeral Head is designed with a larger area of articulation to allow for articulation with the acromion in patients with gross rotator cuff deficiency.

B. INDICATIONS

The Univers Revers Shoulder Prosthesis System and the Univers Revers Modular Glenoid System are indicated for use in a grossly rotator cuff deficient glenohumeral joint with severe arthropathy or a previously failed joint replacement with a gross rotator cuff deficiency. The patient’s joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device.

The Univers Revers Shoulder Prosthesis and the Univers Revers Modular Glenoid System are indicated for primary, fracture, or revision total shoulder replacement for the relief of pain and significant disability due to gross rotator cuff deficiency. (Note: The Univers Revers Shoulder Prosthesis and the Univers Revers Modular Glenoid System are not indicated for fracture in the EU.)

(Humeral) Stems are intended for cementless applications with Arthrex humeral Suture Cups. The Universal Glenoid baseplate is CaP coated and is intended for cementless use with the addition of screws for fixation. The Univers Revers Modular Glenoid System are porous coated and are intended for cementless use with the addition of screws for fixation.

The **Arthrex Univers Revers Cuff Arthropathy Humeral Heads and Adapters** are indicated for:

- salvage of a failed reverse total shoulder, with an irreparable rotator cuff tear and a well-fixed humeral stem, to an anatomic hemi-shoulder replacement; or
- conversion of a primary reverse total shoulder, for the relief of pain secondary to severe rotator cuff arthropathy and an irreparable rotator cuff tear, to anatomic hemi-shoulder replacement when insufficient glenoid bone stock is encountered intraoperatively after the humeral stem has been implanted.

The patient’s joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery.
- Do not use with other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Allergies and other reactions to device materials.
- Loosening of the implant as a result of changed conditions in load transfer, respectively, fatigue wear and breakage of the cement bed and/ or tissue reaction to the implant. Loosening is frequently a consequence of one or several of the above listed risk conditions, but can also be caused by inadequate anchoring technique (see below).
- Dislocation, subluxation, or inadequate scope of movement as a result of failure to achieve optimum positioning of the implant.
- Bone fractures as a result of one-sided overload or weakened bone substance.
- Temporary or permanent nerve damage as a result of pressure or hematoma.
- Cardiovascular diseases including venous thrombosis, pulmonary embolism, and cardiac arrest.
- Wound hematoma and delayed wound healing.

E. WARNINGS

- 6mm offset Humeral Liners must **not** be used in combination with Humeral Spacers. Humeral Spacers **should only** be used with 3mm offset Humeral Liners.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Detailed instructions on the use and limitations of the device should be given to the patient.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- Removal of the device should be performed using standard surgical practices for device removal.
- Preoperative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device. The following operative situations may cause premature loosening and complications:
 - Extreme weakening of the bone structure in preparing the bone bed;
 - Unsuitable selection of the implant size;
 - Inadequate cleaning of the bone bed prior to implantation; and,
 - Excessive use of force in placing or fastening the implant, provoking splintering fractures, or causing the bone to tear.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user. Do not re-sterilize this device.

- The appropriate Arthrex delivery system is required for proper insertion of the implant.
- Only Arthrex delivery systems, instruments, and trial prostheses should be used for the implantation procedure.
- Endoprostheses may not be processed mechanically or changed in any other way.
- Do not implant any parts that have been scratched or damaged. An artificial joint is subject to wear and/or can loosen over a period of time. Wear and loosening may make it necessary to re-operate on an artificial joint.
- An infection in an artificial joint may lead to implant removal.
- This device should only be used in conjunction with other implants designed specifically for use with this system.
- CaP coated device – Fluid contact other than patient’s blood should be avoided to achieve the best ongrowth results.
- Failure to properly align and completely seal the components together can lead to disassociation. Proper technique must be followed to ensure there is no bony or soft tissue interference between modular components. All screws must be adequately tightened to ensure they are recessed to prevent a mechanical interference between modular components.
- Thoroughly clean and dry tapers, prior to attachment of modular components to avoid crevice corrosion and improper seating. Glenosphere forceps are required to verify integrity of the Morse Taper connection between Glenosphere and Baseplate.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex Univers Revers is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system
- Under the scan conditions defined, the Arthrex Univers Revers is expected to produce a maximum temperature rise of 3 °C after 15-minutes of continuous scanning.

a. Artifact Information

In non-clinical testing, the image artifact caused by the Arthrex Univers Revers extends approximately 60 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print format, video format, and electronically. Or, contact your Arthrex representative for an onsite demonstration. The Arthrex website also provides detailed surgical technique information and demonstrations.
- Patient weight. An overweight patient may present additional risk.
- Extreme stress or strain resulting from work or sport related activity. Patients with increased risk of fractures due to repeated strain or trauma, or medical conditions that increase the patient’s risk for trauma, including falls.
- Osteoporosis or osteomalacia.
- Exposure to infectious diseases with possible manifestation in the joints.
- Deformation of the operative site, which can prevent or impede anchoring of the implant.
- Tumors that weaken the supportive structure.
- Allergic reactions to implant materials.
- Thromboses and pulmonary infarction caused during implant site preparation and subsequent / treatment.
- If a revision to the Revers prosthesis is performed and the glenosphere is removed from the baseplate for any reason, a new glenosphere should be implanted.
- If a revision to the Revers prosthesis is performed and the baseplate is not removed, the baseplate Morse Taper socket should be inspected for damage. If damage is identified or suspected, the baseplate should be removed and replaced.

H. PACKAGING AND LABELING

- Arthrex implants should be accepted only if the factory packaging and labeling arrive intact.
- All prosthesis parts and components should be stored unopened in the original packaging, and should be kept in protective packaging until immediately prior to use.
- Contact customer service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at **www.arthrex.com/symbols glossary**.

I. STERILIZATION

This device is provided sterile. Refer to the package label for the sterilization method.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79, “Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities”.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials. This device is made of:

- The humeral device consists of a stem, cup, screw, and spacer manufactured of titanium. The stem body may be partly coated with a calcium phosphate (CaP) coating.
- The humeral inlay component is composed of ultra high molecular weight polyethylene (UHMWPE).
- The glenoid device consists of a glenosphere manufactured of cobalt-chromium. The Univers Revers glenoid baseplate consists of the baseplate, bushings, a central screw and cancellous screws – all manufactured of titanium. The Univers Revers glenoid baseplate is partially coated with calcium phosphate (CaP). The Univers Revers Modular Glenoid System consists of a baseplate, central screw or post and cancellous bone screws – all manufactured of titanium. The Univers Revers Modular Glenoid baseplate is partially coated with CP (commercially pure) titanium. **Worldwide except US and Canada:** The glenosphere is also available as TiNbN coated cobalt-chromium (Co-Cr).

K. STORAGE CONDITIONS

Joint prosthesis must be handled carefully prior to implantation. Scratches or nicks on the joint surface of the prostheses can result in excessive wear and complications.

- Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.
- Do not remove the device from the package until immediately before use.

L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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