


# Apollo<sup>RF</sup>® i90 Probe

Enhanced Ablation. Optimized Aspiration.



**Arthrex**® 

## Apollo<sup>RF</sup> i90 Probe (AR-9831)

The Apollo<sup>RF</sup> i90 device is the first probe with a torpedo-shaped electrode for optimized edge control. Its performance-enhanced features include an improved flow rate and flow path,<sup>1</sup> together with the most efficient rate of tissue removal when compared to the leading market competitors.<sup>2</sup>


### Features and Benefits

- Torpedo-shaped electrode\* promotes easier access and maneuverability into and within the joint
- Optimized flow rate / flow path\* designed to efficiently remove bubbles for better visual clarity and reduce tissue clogging,<sup>3</sup> while maintaining temperature with increases of  $<1^{\circ}\text{C}$ <sup>4</sup>
- Enhanced ablation for the most efficient removal rate of tissue on the market<sup>2</sup>
- Edge control for a defined plasma edge that enables greater precision and gentle ablation distally and along the sides of the electrode
- Shorter working length of 140 mm and ability to be used with a 5 mm cannula

### Ordering Information

Apollo<sup>RF</sup> i90 probe

AR-9831



The Apollo<sup>RF</sup> i90 probe prepping the tuberosity with distal end ablation while working through the lateral portal.

\*Patent pending

# Application Examples

## Rotator Cuff Repair

- › Release/remove subdeltoid bursas and/or adhesions before placing the lateral row of anchors
- › Use the lateral portal to mark and clear ideal locations for anchors with distal end ablation

## Subacromial Decompression

- › Market leader in efficient resection,<sup>2</sup> which is optimal for removing bursal tissue and adipose off the acromion

## Bankart Preparation

- › Precisely debride degenerative tissue of the labrum with the side or distal end of the probe

## Adhesive Capsulitis

- › Dissect the middle glenohumeral ligament (MGHL) and capsule from the undersurface of the subscapularis
- › Maintain precise division of the capsule while preserving the adjacent axillary nerve

## Chondroplasty

- › Target treatment of chondral lesions for smoothing and stabilization

Bankart Preparation



Stabilization of Chondral Lesions



## References

1. Arthrex, Inc. Data on file (PLM85639). Naples, FL; 2021.
2. Arthrex, Inc. LA1-000176-en-US\_C. Naples, FL; 2024.
3. Arthrex, Inc. Data on file (AR-9831 clog performance design verification). Naples, FL; 2021.
4. Arthrex, Inc. LA1-000202-en-US\_A. Naples, FL; 2024.

This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.



Arthrex manufacturer, authorized  
representative, and importer  
information (Arthrex eIFUs)



US patent  
information