



Autograft Harvesting and Augmentation Recommendations

Blades and Biologics Quick Reference Chart

	Technique	Blade	Setting	Outcome
	AutoCart™ Technique	4 mm Bone Cutter	3000 RPM, AGG mode	Harvest medium-sized cartilage particulate while maintaining high cell viability (>80%) ¹
		4 mm Torpedo™ Shaver	3000 RPM, AGG mode	
	Microfracture	PowerPick™ XL Attachment	6200 RPM, FWD direction	Patient-reported outcomes were significantly improved for those treated with the PowerPick attachment vs conventional chondral picks ²
	BioACL™ Technique	5.5 mm Bone Cutter	6200 RPM, FWD direction	Autograft bone has shown to be the gold standard for use in orthopedic procedures as it inherently includes osteogenic, osteoinductive, and osteoconductive properties ³
	Bone Collection for BTB Backfill			
	Rotator Cuff Repair Augmentation With Subacromial Bursa	4 mm Excalibur	1750 RPM, EFF mode	The 4 mm Excalibur shaver blade maintains 87% CTP cell viability after subacromial bursal tissue harvest ⁴
	Preparing Tuberosity for RCR	4 mm PoweRasp™ Attachment	5000 RPM, FWD direction	Creation of a bleeding bone surface to accelerate tendon-to-bone healing



References

1. Arthrex, Inc. Data on file (APT 03989). Naples, FL; 2019.
2. Beletsky A, et al. *Arthrosc Sports Med Rehabil.* 2021;3(3):e629-e638. doi:10.1016/j.asmr.2020.10.006
3. Lavender C. *Arthrosc Tech.* 2019;8(9):e1019-e1023. doi:10.1016/j.eats.2019.05.013
4. Wellington IJ, et al. *J Clin Med.* 2022;11(5):1272. doi:10.3390/jcm11051272