English

A. DEVICE DESCRIPTION
The Univrs Revers™ Shoulder Prosthesis System has an articular design that is inverted, compared to traditional total shoulder prosthesis. The system is comprised of two main components; the Arthrex Univrs Revers Shoulder Prosthesis and the Universal Glenoid™ Shoulder Prosthesis. The Arthrex Univrs Revers Shoulder Prosthesis is a humeral stem and epiphysis or humeral cup, a spacer, and humeral cup liner. The humeral stem and epiphysis are available uncoated or coated.

The Universal Glenoid Shoulder Prosthesis consists of a glenoid baseplate, a glenosphere, and screws.

The Univrs Revers™ Modular Glenoid System consists of a monoblock baseplate or a modular baseplate; both baseplates are available with either a central screw or central post. Except in CE Accepting Countries, modular baseplates are also available with wedge augments (full or half wedge). The augmented modular baseplates are used with central posts. The baseplates are designed to be used cementless with peripheral screws and a glenosphere. The glenoid system is designed to be used as the glenoid side of the existing Univrs Revers Shoulder Prosthesis System.

The Univrs Revers™ Cuff Arthropathy System is designed to be used with the existing well-fixed Univrs Revers™ stem, or to convert an existing reversed shoulder prosthesis to a hemi anatomic configuration. The Univrs Revers Cuff Arthropathy (CA) Humeral Head is designed with a larger area of articulation to allow for articulation with the acromion in patients with gross rotator cuff deficiency.

B. INDICATIONS

The Univrs Revers Shoulder Prosthesis System and the Univrs Revers Modular Glenoid System are indicated for use in a grossly rotator cuff deficient gleno-humeral joint with severe arthropathy or a previously failed joint replacement with a gross rotator cuff deficiency. The patient’s joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device.

The Univrs Revers Shoulder Prosthesis and the Univrs Revers Modular Glenoid System are indicated for primary, fracture, or revision total shoulder replacement for the relief of pain and significant disability due to gross rotator cuff deficiency. (Note: The Univrs Revers Shoulder Prosthesis and the Univrs Revers Modular Glenoid System are not indicated for fracture in CE Accepting Countries.)

(Humeral) Stems are intended for cementless applications with Arthrex humeral Suture Cups. The Universal Glenoid baseplate is CaP coated and is intended for cementless use with the addition of screws for fixation. The Univrs Revers Modular Glenoid System are porous coated and are intended for cementless use with the addition of screws for fixation.

The *Arthrex Univrs Revers CA Humeral Head and Adapters* are indicated for:

- salvage of a failed reverse total shoulder, with an irreparable rotator cuff tear and a well-fixed humeral stem, to an anatomic hemi-shoulder replacement; or
- conversion of a primary reverse total shoulder, for the relief of pain secondary to severe rotator cuff arthropathy and an irreparable rotator cuff tear, to anatomic hemi-shoulder replacement when insufficient glenoid bone stock is encountered intraoperatively after the humeral stem has been implanted

The patient’s joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity: Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Injury to nerves, blood vessels, muscles, tendons.
- Malunion/nonunion of the bone.
- Necrosis of the humeral head.
- Risk of chronic pain.
- Allergies and other reactions to device materials.
- Loosening of the implant as a result of changed conditions in load transfer, respectively, fatigue wear and breakage of the cement bed and/ or tissue reaction to the implant. Loosening is frequently a consequence of one or several of the above listed risk conditions, but can also be caused by inadequate anchoring technique (see below).
- Dislocation, subluxation, or inadequate scope of movement as a result of failure to achieve optimum positioning of the implant.
- Bone fractures as a result of one-sided overload or weakened bone substance.
- Temporary or permanent nerve damage as a result of pressure or hematoma.
- Cardiovascular diseases including venous thrombosis, pulmonary embolism, and cardiac arrest.
- Wound hematoma and delayed wound healing.
- Tissue reactions caused by allergic reactions to the implanted material, particularly metal, or caused by accumulations of wear particles or cement particles.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional. 6mm offset Humeral Liners must not be used in combination with Humeral Spacers. Humeral Spacers should only be used with 3mm offset Humeral Liners.
- Failure to achieve the appropriate torque requirements when tightening locking screws may result in the premature loosening of the device.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Detailed instructions on the use and limitations of this device, the patient leaflet (www.arthrex.com/patientleaflets) and the patient implant card should be given to the patient. Your surgeon will guide you in deciding what particular treatment is best for you and explain the benefits, risks, and contraindications associated with the treatment.
- Any action to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- Removal of the device should be performed using standard surgical practices for device removal.

- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device. The following operative situations may cause premature loosening and complications:
 - Extreme weakening of the bone structure in preparing the bone bed;
 - Unsuitable selection of the implant size;
 - Inadequate cleaning of the bone bed prior to implantation; and,
 - Excessive use of force in placing or fastening the implant, provoking splintering fractures, or causing the bone to tear.

- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Do not re-sterilize this device.
- The operation is to be planned based on the pre-operative X-rays.
- The specific Arthrex implantation instruments are to be used both for preparing the bone bed and for adjusting and inserting the joint prosthesis. The appropriate Arthrex delivery system is required for proper insertion of the implant.
- Only Arthrex delivery systems, instruments, and trial prostheses should be used for the implantation procedure.
- Endoprosthesis may not be processed mechanically or changed in any other way.
- Do not implant any parts that have been scratched or damaged. An artificial joint is subject to wear and/or can loosen over a period of time. Wear and loosening may make it necessary to re-operate on an artificial joint.
- An infection in an artificial joint may lead to implant removal. This device should only be used in conjunction with other implants designed specifically for use with this system.
- Proper anchoring is of decisive importance for firm, permanent positioning of the prosthesis.
- CaP coated device – Fluid contact other than patient’s blood should be avoided to achieve the best on grossly results.
- In the case of joint endoprosthesis intended for cemented anchoring, the surgeon must comply with the instructions and recommendations of the cement manufacturer when it comes to preparation and cementing techniques. Failure to properly align and completely seat the components together can lead to disassociation. Proper technique must be followed to ensure there is no bony or soft tissue interference between modular components. All screws must be adequately tightened to ensure they are recessed to prevent a mechanical interference between modular components.
- Thoroughly clean and dry tapers, prior to attachment of modular components to avoid crevice corrosion and improper seating. Glenosphere forceps are required to verify integrity of the Morse Taper connection between Glenosphere and Baseplate.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.
- For oblique augmented baseplates peripheral screw lengths of 24mm or longer must be used if the screw passes through the augment. For oblique half augmented baseplates, any peripheral screw length can be utilized if the screw passes through the non-augmented portion of the baseplate.

- CAP coating device – Fluid contact other than patient’s blood should be avoided to achieve the best on grossly results.
- In the case of joint endoprosthesis intended for cemented anchoring, the surgeon must comply with the instructions and recommendations of the cement manufacturer when it comes to preparation and cementing techniques. Failure to properly align and completely seat the components together can lead to disassociation. Proper technique must be followed to ensure there is no bony or soft tissue interference between modular components. All screws must be adequately tightened to ensure they are recessed to prevent a mechanical interference between modular components.
- Thoroughly clean and dry tapers, prior to attachment of modular components to avoid crevice corrosion and improper seating. Glenosphere forceps are required to verify integrity of the Morse Taper connection between Glenosphere and Baseplate.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.
- For oblique augmented baseplates peripheral screw lengths of 24mm or longer must be used if the screw passes through the augment. For oblique half augmented baseplates, any peripheral screw length can be utilized if the screw passes through the non-augmented portion of the baseplate.

F. FACTORS AND RISKS IMPACTING THE SAFETY AND SERVICE LIFE OF THE IMPLANT

- Patient weight. An overweight patient may present additional risk.
- Extreme stress or strain resulting from work or sport related activity.
- Patients with an increased risk of fractures due to repeated strain or trauma, or medical conditions that increase the patient’s risk for trauma, including falls.
- Osteoporosis or osteomalacia.
- Exposure to infectious diseases with possible manifestation in the joints.
- Deformation of the operative site, which can prevent or impede anchoring of the implant.
- Tumors that weaken the supportive structure.
- Allergic reactions to implant materials.
- Thromboses and pulmonary infarction caused during implant site preparation and subsequent treatment.

G. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex Univrs Revers Shoulder System is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whose body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the Arthrex Univrs Revers Shoulder System is expected to produce a maximum temperature rise of 3 °C after 15-minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Arthrex Univrs Revers Shoulder System extends approximately 60 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

H. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print format, video format, and electronically. Or, contact your Arthrex representative for an onsite demonstration. The Arthrex website also provides detailed surgical technique information and demonstrations.
- Patient weight. An overweight patient may present additional risk.
- Extreme stress or strain resulting from work or sport related activity.
- Patients with increased risk of fractures due to repeated strain or trauma, or medical conditions that increase the patient’s risk for trauma, including falls.
- Osteoporosis or osteomalacia.
- Exposure to infectious diseases with possible manifestation in the joints.
- Deformation of the operative site, which can prevent or impede anchoring of the implant.
- Tumors that weaken the supportive structure.
- Allergic reactions to implant materials.
- Thromboses and pulmonary infarction caused during implant site preparation and subsequent / treatment.
- If a revision to the Univrs Revers prosthesis is performed and the glenosphere is removed from the baseplate for any reason, a new glenosphere should be implanted.
- If a revision to the Univrs Revers prosthesis is performed and the baseplate is not removed, the baseplate Morse Taper socket should be inspected for damage. If damage is identified or suspected, the baseplate should be removed and replaced.

I. PACKAGING AND LABELING

- Arthrex implants should be accepted only if the factory packaging and labeling arrive intact.
- All pre-assembled parts and components should be stored unopened in the original packaging, and should be kept in protective packaging until immediately prior to use.
- Contact customer service if the package has been opened or altered.

- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

J. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

K. MATERIAL SPECIFICATIONS

Refer to the package label for the materials. This device is made of:

- The humeral device consists of a stem, cup, screw, and spacer manufactured of titanium. The stem body may be partly coated with a calcium phosphate (CaP) coating.
- The humeral alloy component is composed of ultra-high molecular weight polyethylene (UHMWPE).
- The glenoid device consists of a glenosphere manufactured of cobalt-chromium. The Univrs Revers glenoid baseplate consists of the baseplate, bushings, a central screw and cancellous screws – all manufactured of titanium. The Univrs Revers glenoid baseplate is partially coated with calcium phosphate (CaP). The Univrs Revers Modular Glenoid System consists of a baseplate, central screw or post and cancellous bone screws – all manufactured of titanium. The Univrs Revers Modular Glenoid baseplates and modular posts are partially coated with CP (commercially pure) titanium. The glenosphere is also manufactured of Ti6N coated cobalt-chromium (Co-Cr) but this device is not available in all countries.

L. STORAGE CONDITIONS

Joint prosthesis must be handled carefully prior to implantation. Scratches or nicks on the joint surface of the prostheses can result in excessive wear and complications.

- Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.
- Do not remove the device from the package until immediately before use.

M. INFORMATION

- In CE Accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE Accepting Countries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE Accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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