

3DAnatomy™

TTC ArthroLab™ Model



TTC ArthroLab™ Model

Developed by 3DAnatomy™, the tibiotocalcaneal (TTC) arthrodesis ArthroLab model, is a high-fidelity anatomical model engineered to replicate the structural and mechanical properties of human bone, enabling precise, repeatable execution of the hindfoot procedure.

- › Simulates realistic cortical and cancellous response under instrumentation.
- › Purpose-built for the DualCompression hindfoot fusion nail procedure.
- › Engineered for consistent and repeatable performance.
- › Supports procedural training, technique refinement, and education.
- › For use with 10.5 mm × 210 mm or 10.5 mm × 180 mm hindfoot nails



Assembly Instructions

Step 7: Insert the pin into the pin hole, securing the tibia to the TRIMANO puck.



Step 6: Insert the Proximal Dowel of the tibia into the TTC TRIMANO puck.



Step 5: Slide the proximal tibia onto the compression tab cylinders from the posterior side of the assembly.



Step 4: Slide the compression tab into the cylindrical slots, making sure the side labeled DIST is facing down (the PROX label should face up).




Step 3: Place the distal tibia onto the K-wires, then gently push the wires downward until they seat fully against the built-in stop in the tibia.




Note: The compression tab must be removed before disassembling the parts.

Step 2: Secure the talus onto the K-wires from the proximal end, ensuring the smooth cartilage surfaces of both the talus and calcaneus are facing forward (anterior).



Step 1: Insert the two nitinol K-wires through the distal calcaneus (medial or lateral hole) until it exits proximally.



Tips for a Successful Lab

- › Follow the Setup Guide to assemble the model prior to the session.
- › Confirm orientation using molded labels (eg, PROX/DIST, anterior/posterior).
- › If resistance is felt, realign and re-seat components—do not force parts together.
- › Stage required instruments and perform a quick dry run before attendees arrive.
- › Secure the model using a compatible holder (eg, TRIMANO®, or vise clamp).
- › Components should be stored in a cool, dry environment and protected from direct light or sunlight when not in use.

Note: Ensure joint spaces are collapsed prior to the Hindfoot Nail insertion, to ensure maximum nitinol compression.

Ordering Information

TTC Full Kit

AMS-3DC-04-005-01K

- › Calcaneus
- › Talus
- › Distal Tibia
- › Proximal Tibia
- › TRIMANO® Puck
- › Compression Pull-Out Tab
- › Carabiner
- › Assembly Wires

TTC Disposable Kit

AMS-3DC-04-005-02K

- › Calcaneus
- › Distal Tibia
- › Compression Pull-Out Tab

Exclusively available for Arthrex employees and agencies through the business portal or CSSynergy@arthrex.com



Learn more about
3DAnatomy models