


Shoulder and Elbow

Next Generation in Repair and Reconstruction



Arthrex® 

Helping Surgeons Treat Their Patients Better[®]

Since its inception, Arthrex has been committed to one mission: Helping Surgeons Treat Their Patients Better. We are strategically focused on constant product innovation through scientific research, surgeon collaboration, and medical education to make less invasive surgical procedures simple, safer, and more reproducible. Each year, we develop more than 1,000 new innovative products and procedures to advance minimally invasive orthopedics worldwide.

Arthrex has always remained a privately held company, which allows for the rapid evaluation of new technologies and ideas and the freedom to develop products and techniques that truly make a difference. Our experienced team of dedicated professionals represents a shared passion and commitment to delivering uncompromising quality to the health care providers who use our products and the millions of patients whose lives we impact.

The medical significance of our contributions serves as our primary benchmark of success and will continue into the future as the legacy of Arthrex.


Table of Contents

Rotator Cuff Repair	05
SLAP and Bankart Repair	17
FiberWire® Suture	27
Cannulas.....	35
Instruments.....	41
Biceps and Pectoralis Repair	57
Elbow Arthroscopy	67
AC Reconstruction Systems.....	71
Glenoid Bone Loss	77
Fracture Management	81
Cardiothoracic.....	91
Limb Positioning Devices.....	95
Anchor Charts.....	101

Rotator Cuff Repair

SwiveLock® C Suture Anchor	06
SwiveLock SP Suture Anchor	07
2.6 FiberTak® Soft Anchors	08
SpeedBridge™ Implant Systems	09
CuffMend™ Rotator Cuff Augmentation System	10
Autograft Tissue Compression System	11
PushLock® Suture Anchor	11
Superior Capsular Reconstruction (SCR).....	12
Biologic Tuberoplasty Technique	12
PASTA Bridge	13
Knotless Corkscrew® Suture Anchor.....	13
Corkscrew FT Suture Anchors	14
Corkscrew Suture Anchor.....	15

SwiveLock® C Suture Anchor



This closed-eyelet, twist-in knotless anchor can be used as the lateral row of the SutureBridge™ repair. Additionally, it can be combined with FiberTape® suture for SpeedFix™ and SpeedBridge™ repairs. Knotless versions incorporate a knotless tensionable mechanism in the eyelet.

Biocomposite

Product Description	Item Number
3.9 mm × 17.9 mm (a)	AR-2326BCC
Knotless , 4.75 mm × 19.1 mm w/ #2 suture (blue) (c)	AR-2324KBCC
Knotless , 4.75 mm × 19.1 mm w/ FiberTape Loop (blue) and #2 suture (white/black)	AR-2324KBCCT
Knotless , 4.75 mm × 19.1 mm w/ TigerTape™ Loop (white/black) and #2 suture (blue)	AR-2324KBCCTT
4.75 mm × 19.1 mm	AR-2324BCC
4.75 mm × 22 mm, w/ two #2 TigerWire® CL suture (1 white/blue, 1 white/black) (e)	AR-2324BCC-2
4.75 mm × 22 mm, w/ two #2 TigerTail® sutures (1 white/black, 1 blue/black)	AR-2324BCT-2
4.75 mm × 19.1 mm, w/ blue FiberTape Loop (d)	AR-2324BCCT
4.75 mm × 19.1 mm, w/ white/black TigerTape Loop	AR-2324BCCTT
Knotless , 5.5 mm × 19.1 mm w/ #2 suture (blue)	AR-2323KBCC
5.5 mm × 22 mm, w/ two #2 TigerTail sutures (1 white/black, 1 blue/black)	AR-2323BCT-2
5.5 mm × 19.1 mm (f)	AR-2323BCC

PEEK

Product Description	Item Number
3.9 mm × 17.9 mm (b)	AR-2326PSLC
Knotless , 4.75 mm × 19.1 mm w/ FiberTape Loop (blue) and #2 suture (white/black)	AR-2324KPCT
Knotless , 4.75 mm × 19.1 mm w/ TigerTape Loop (white/black) and #2 suture (blue)	AR-2324KPCTT
Knotless , 4.75 mm × 19.1 mm w/ #2 suture (blue)	AR-2324KPSLC
4.75 mm × 19.1 mm	AR-2324PSLC
4.75 mm × 22 mm, w/ two #2 TigerTail sutures (1 white/black, 1 blue/black)	AR-2324PSLC-2
4.75 mm × 19.1 mm, w/ blue FiberTape Loop	AR-2324PCT
4.75 mm × 19.1 mm, w/ white/black TigerTape Loop	AR-2324PCTT
Knotless , 5.5 mm × 19.1 mm w/ #2 suture (blue)	AR-2323KPSLC
5.5 mm × 22 mm, w/ two #2 TigerTail sutures (1 white/black, 1 blue/black)	AR-2323PSLC-2
5.5 mm × 19.1 mm (g)	AR-2323PSLC

Bio



Product Description	Item Number
4.75 mm × 19.1 mm	AR-2324BSLC

Implant System

Product Description	Item Number
SwiveLock Implant System	AR-2324BCCS
4.75 mm BioComposite SwiveLock anchor preloaded w/ blue FiberTape Loop, HD Scorpion needle	

Product Description	Item Number
Required Instruments	
Punch for 5.5 mm Corkscrew® FT and 4.75 mm and 5.5 mm SwiveLock anchors	AR-1927PB
Punch for 3.9 mm SwiveLock anchor	AR-2326P
Optional Instruments	
Punch/Tap for 5.5 mm Corkscrew FT and 5.5 mm SwiveLock anchors	AR-1927CTB
Punch/Tap for 4.75 mm SwiveLock anchors	AR-2324PTB
Disposable Punch for 5.5 mm Corkscrew FT and 4.75 mm and 5.5 mm SwiveLock anchors	AR-1927PBS
Spade Tip Drill for 5.5 mm Corkscrew FT and 4.75 mm and 5.5 mm SwiveLock anchors	AR-1927D
Punch/Tap for 3.9 mm SwiveLock anchor	AR-2326PTB
Spade Tip Drill for 3.9 mm SwiveLock anchor	AR-2326D

SwiveLock® SP Suture Anchor



Eliminate the need for prepunching a bone socket with the SwiveLock SP (self-punching) anchor. These anchors have a variety of options, including PEEK or titanium tip with a biocomposite, PEEK, or titanium anchor body. The self-punching design can help save valuable OR time, while increasing the precision of the final construct. Combine the SwiveLock SP anchor with the 2.6 FiberTak® RC for a fully self-punching FiberTak SpeedBridge™ knotless rotator cuff repair.

Biocomposite

Product Description	Item Number
Knotless , 4.75 mm × 24.5 mm w/ #2 suture (blue), SP PEEK eyelet (a)	AR-2324KBCSP
4.75 mm × 24.5 mm w/ 1.3 mm SutureTape (white/blue), SP PEEK eyelet	AR-2324BCSP
4.75 mm × 24.5 mm, SP titanium eyelet (b)	AR-2324BCM
Knotless , 5.5 mm × 24.5 mm w/ #2 suture (blue), SP PEEK eyelet	AR-2323KBCSP
5.5 mm × 24.5 mm w/ 1.3 mm SutureTape (white/blue), SP PEEK eyelet	AR-2323BCSP
5.5 mm × 24.5 mm, SP titanium eyelet	AR-2323BCM

PEEK


Product Description	Item Number
Knotless , 4.75 mm × 24.5 mm w/ #2 suture (blue), SP PEEK eyelet	AR-2324KPSP
4.75 mm × 24.5 mm w/ 1.3 mm SutureTape (white/blue), SP PEEK eyelet	AR-2324PSP
4.75 mm × 24.5 mm, SP titanium eyelet (c)	AR-2324PSLM
Knotless , 5.5 mm × 24.5 mm w/ #2 suture (blue), SP PEEK eyelet	AR-2323KPSP
5.5 mm × 24.5 mm w/ 1.3 mm SutureTape (white/blue), SP PEEK eyelet	AR-2323PSP
5.5 mm × 24.5 mm, SP titanium eyelet	AR-2323PSLM

Titanium

Product Description	Item Number
4.75 mm × 19.1 mm, SP titanium eyelet (d)	AR-2324SLM

Product Description	Item Number
Recommended FiberTape® Suture Options for SwiveLock C and SP Anchors	
36 in (blue) tape, w/ each end tapered to #2 FiberWire® suture, 54 in	AR-7237
7 in (blue) tape, w/ each end tapered to #2 FiberWire suture, 30 in (e)	AR-7237-7
7 in (blue) tape, w/ each end tapered to #2 FiberWire suture, collagen-coated	AR-7237-7B
7 in (white/black) tape, w/ each end tapered to #2 TigerWire® suture, 30 in (f)	AR-7237-7T
17 in (blue) tape, w/ each end tapered to #2 FiberWire suture w/ tapered needle, 26.5 mm, ½ circle	AR-7237-17N
FiberTape Suture Instruments	
Mini FiberTape Suture Retriever w/ self-release (SR) handle (g)	AR-12974SR
Mini FiberTape Suture Retriever w/ nonratcheting (NR) handle	AR-12974NR
Mini FiberTape Suture Retriever w/ WishBone™ handle	AR-12974W
FiberTape Suture Retriever w/ SR handle	AR-13974SR
FiberTape Suture Retriever w/ NR handle	AR-13974NR
FiberTape Suture Retriever w/ WishBone handle	AR-13974W
Mini KingFisher Retriever/Grasper, 3.4 mm w/ SR handle	AR-13972SR
Mini KingFisher Retriever/Grasper, 3.4 mm w/ NR handle	AR-13972NR
Mini KingFisher Retriever/Grasper, 3.4 mm w/ WishBone handle	AR-13972W
KingFisher® FiberTape Suture Retriever/Grasper w/ SR handle	AR-13971SR
KingFisher FiberTape Suture Retriever/Grasper w/ NR handle	AR-13971NR
KingFisher FiberTape Suture Retriever/Grasper w/ WishBone handle (h)	AR-13971W
Penetrator™ FiberTape Suture Retriever, 15° up curved	AR-2167-3
Penetrator FiberTape Suture Retriever, straight	AR-2167ST-3
Penetrator FiberTape Suture Retriever w/ WishBone handle, 15° up curved	AR-2167W-3
Penetrator FiberTape Suture Retriever w/ WishBone handle, straight (i)	AR-2167STW-3
FiberTape Suture Cutter	AR-13250
FiberTape Suture Cutter w/ WishBone handle	AR-13250W

2.6 FiberTak® Soft Anchors



2.6 FiberTak RC soft anchors are specifically designed for surgeons who prefer an all-suture anchor for the medial row of a double-row, bridging rotator cuff repair. They are combined with knotless SwiveLock® C anchors for the lateral row. The unique 2.6 FiberTak RC anchor has a 1.7 FiberTape® suture fixed to the sheath so that it cannot slide, making it easy to tension the lateral row for a totally knotless construct. The 2.6 FiberTak RC anchor is available with a tensionable knotless mechanism or a sliding 1.3 mm SutureTape that can be incorporated into the repair, or discarded, for ultimate flexibility.

Knotless 2.6 FiberTak soft anchors—which offer versatility for rotator cuff repairs—can be used in multiple configurations (simple stitches, mattress stitches, or interconnected). These anchors are ideal for repairs such as upper border subscapularis, PASTA, or remplissage. Double- and triple-loaded 2.6 FiberTak options feature the latest SutureTape technology for knotted repairs for single- or double-row repairs.

The self-punching inserter for select FiberTak anchors allows for direct insertion into bone without creating a prepunched or predrilled socket for anchor insertion. The self-punching FiberTak inserter may be used with a guide and can be placed directly into bone or used with a transtendon insertion technique.

2.6 FiberTak RC Soft Anchor

Product Description	Item Number
With 1.7 mm TigerTape™ Loop (white/black) and 1.3 mm SutureTape (white/blue), self-punching, qty 5 (a)	AR-3652SP
With 1.7 mm FiberTape® Loop (blue) and 1.3 mm SutureTape (white/black), self-punching, qty 5	AR-3652TSP
With 1.7 mm FiberTape Loop (black/blue) and 1.3 mm SutureTape (blue), self-punching, qty 5	AR-3652TTSP
Knotless with 1.7 mm TigerTape Loop (white/black) and #2 suture (blue), self-punching, qty 5	AR-3653SP
Knotless with 1.7 mm FiberTape Loop (blue) and #2 suture (white/black), self-punching, qty 5	AR-3653TSP
Knotless with 1.7 mm FiberTape Loop (black/blue) and #2 suture (blue/black), self-punching, qty 5 (b)	AR-3653TTSP


2.6 FiberTak Soft Anchor

Product Description	Item Number
Knotless w/ two #2 knotless repair sutures (white/blue & white/black), self-punching, qty. 5 (c)	AR-3642SP
Knotless w/ #5 suture, self-punching, qty. 5	AR-3641SP
Knotless w/ #5 suture, qty. 5	AR-3641
With Two 1.3 mm SutureTapes (white/blue, white/black), self-punching, qty. 5	AR-3632SP
With Three 1.3 mm SutureTapes (white/blue, black/white, blue), self-punching, qty. 5 (d)	AR-3633SP

Product Description	Item Number
Disposables for 2.6 FiberTak anchors	
Disposable Kit Includes angled circumferential teeth spear, sharp obturator, and drill	AR-3650DS
Instruments for 2.6 FiberTak anchors	
Angled Spear w/ Circumferential Teeth	AR-3655
Blunt Obturator	AR-3658B
Sharp Obturator	AR-3658T
Circumferential Teeth Spear w/ Sharp Obturator	AR-1941CT
Fishmouth Spear w/ Sharp Obturator	AR-1941DGF
Drill, 2.6 mm	AR-3657
2.6 mm Shaver Drill	AR-3657SD
Punch	AR-3656

Literature


Product Description	Item Number
FiberTak® SpeedBridge™ Rotator Cuff Repair Surgical Technique	LT1-000204-en-US
FiberTak SpeedFix™ Rotator Cuff Repair Surgical Technique	LT1-000296-en-US



SpeedBridge™ Implant Systems



(a)



Tensionable Knotless
SpeedBridge Repair



SpeedBridge Repair



SwiveLock® anchor

SpeedBridge implant systems are intended to facilitate a more convenient and efficient SpeedBridge repair. Implant systems are packaged in a single “peel-pack” blister and contain all implants and sutures required to perform a typical SpeedBridge repair.


Product Description	Item Number
Knotless BioComposite SpeedBridge Implant System	AR-2600SBS-10
Two 4.75 mm Knotless SwiveLock anchors w/ 1 preloaded FiberTape® Loop (1 blue, 1 white/black) for medial row, two 4.75 mm Knotless SwiveLock anchors w/ #2 FiberWire® suture (blue) for lateral row, disposable punch, and HD Scorpion needle	
Knotless BioComposite SpeedBridge Implant System	AR-2600SBS-11
Two 4.75 mm Knotless SwiveLock anchors w/ 1 preloaded FiberTape Loop (1 blue, 1 white/black) for medial row, two 5.5 mm Knotless SwiveLock anchors w/ #2 FiberWire suture (blue) for lateral row, disposable punch, and HD Scorpion needle	
BioComposite SpeedBridge Implant System	AR-2600SBS-8
Two 4.75 mm SwiveLock C anchors w/ 1 preloaded FiberTape Loop (1 blue, 1 white/black) for medial row, two 4.75 mm SwiveLock C anchors for lateral row, disposable punch, and HD Scorpion needle	
BioComposite SpeedBridge Implant System	AR-2600SBS-9
Two 4.75 mm SwiveLock C anchors w/ 1 preloaded FiberTape Loop (1 blue, 1 white/black) for medial row, two 5.5 mm SwiveLock C anchors for lateral row, disposable punch, and HD Scorpion needle	
BioComposite SpeedBridge Implant System	AR-2600SBS-4
Two 4.75 mm SwiveLock C anchors w/ 1 preloaded FiberTape Loop (1 blue, 1 white/black) for medial row, two 4.75 mm SwiveLock C anchors for lateral row, disposable punch	
BioComposite SpeedBridge Implant System	AR-2600SBS-5
Two 4.75 mm SwiveLock C anchors w/ 1 preloaded FiberTape Loop (1 blue, 1 white/black) for medial row, two 5.5 mm SwiveLock SP anchors for lateral row, disposable punch	
PEEK SpeedBridge Implant System	AR-2600SBS-6
Two 4.75 mm SwiveLock C anchors w/ 1 preloaded FiberTape Loop (1 blue, 1 white/black) for medial row, two 4.75 mm SwiveLock C anchors for lateral row, disposable punch	
PEEK SpeedBridge Implant System	AR-2600SBS-7
Two 4.75 mm SwiveLock C anchors w/ 1 preloaded FiberTape Loop (1 blue, 1 white/black) for medial row, two 5.5 mm SwiveLock SP anchors for lateral row, disposable punch	

Product Description	Item Number
FiberTak SpeedBridge Implant System (a)	AR-2600SBS-7
Three 2.6 FiberTak RC Self-Punching anchors w/1.7 FiberTape Loop and 1.3 mm SutureTape (1 white/black, white/blue, 1 blue, white/black, and 1 black/blue, blue), two 4.75 Knotless BioComposite SwiveLock SP anchors, HD Scorpion needle, and MegaLoader	
Knotless FiberTak SpeedBridge Implant System	AR-2600FSB-2
Three 2.6 Knotless FiberTak RC Self-Punching anchors w/ 1.7 mm FiberTape Loop and #2 suture (1 white/black, blue, 1 blue, white/black, and 1 black/blue, blue/black), two 4.75 Knotless BioComposite SwiveLock SP anchors, HD Scorpion needle, and MegaLoader	
FiberTak SpeedBridge Implant System	AR-2600FSB-3
Three 2.6 FiberTak RC Self-Punching anchors w/1.7 FiberTape Loop and 1.3 mm SutureTape (1 white/black, white/blue, 1 blue, white/black, and 1 black/blue, blue), two 5.5 Knotless BioComposite SwiveLock SP anchors, HD Scorpion needle, and MegaLoader	
Knotless FiberTak SpeedBridge Implant System	AR-2600FSB-4
Three 2.6 Knotless FiberTak RC Self-Punching anchors w/ 1.7 mm FiberTape Loop and #2 suture (1 white/black, blue, 1 blue, white/black, and 1 black/blue, blue/black), two 5.5 Knotless BioComposite SwiveLock SP anchors, HD Scorpion needle, and MegaLoader	
FiberTak SpeedBridge Implant System	AR-2600FSB-5
Two 2.6 FiberTak RC Self-Punching anchors (1 blue, white/black and 1 black/blue, blue), two 4.75 Knotless BioComposite SwiveLock SP anchor, HD Scorpion needle, and MegaLoader	
FiberTak SpeedBridge Implant System	AR-2600FSB-6
Three 2.6 Knotless FiberTak RC Self-Punching anchors w/ 1.7 mm FiberTape Loop and #2 suture (1 white/black, white/blue, 1 blue, white/black, and 1 black/blue, blue), two 4.75 BioComposite SwiveLock C anchors, HD Scorpion needle, and MegaLoader	

Literature

Product Description	Item Number
SpeedBridge™ Knotless Rotator Cuff Repair Surgical Technique	LT1-000170-en-US
FiberTak® SpeedBridge™ Rotator Cuff Repair Surgical Technique	LT1-000204-en-US

CuffMend™ Rotator Cuff Augmentation System



The CuffMend™ system provides an efficient, simplified approach for augmentation of partial- and full-thickness rotator cuff tears. The system employs the benefits of an ArthroFLEX® dermal allograft, which provides mechanical strength and supports healing.^{1,2} The graft is introduced using the graft spreader and attached medially with the FiberStitch™ RC implants for soft-tissue fixation to the rotator cuff tendon. This allows the augmentation to span the complete tendon length. PushLock® anchors are used, spanning the graft over the tuberosity footprint, for lateral bony fixation. Scientific literature supports the use of acellular dermal grafts as a safe and effective solution for rotator cuff augmentation.³

Implants

Product Description	Item Number
CuffMend Implant System	AR-19042
Two FiberStitch 1.5 curved RC Implants, two 3.5 mm BC PushLock anchors w/ SP PEEK eyelet, a 4 Layer ATCS Plate w/ three #0 FiberWire sutures w/ needle and two 0.9 mm FiberLink SutureTapes w/ needle, one each MegaLoader, graft spreader, knot pusher/suture cutter w/ portal skid	
CuffMend Rotator Cuff Augmentation Implant System	AR-19043
Two FiberStitch 1.5 curved RC Implants, two 3.5 mm BC PushLock w/ SP PEEK eyelet, one each MegaLoader, graft spreader, and knot pusher/suture cutter w/portal skid	
FiberStitch 1.5 RC Implant, straight w/ two polyester implants and 2-0 FiberWire®	AR-19032S
FiberStitch 1.5 RC Implant, curved w/ two polyester implants and 2-0 FiberWire	AR-19032C
FiberStitch RC Implant, simple, straight	AR-19033S
FiberStitch RC Implant, simple, curved	AR-19033C
BioComposite PushLock anchor, 3.5 mm × 19.5 mm, self-punching PEEK eyelet	AR-1926BCSP
PEEK PushLock anchor, 3.5 mm x 19.5 mm, self-punching PEEK eyelet	AR-1926PSSP
FiberLink™ SutureTape, 0.9 mm, w/ loop (white/blue)	AR-7559
TigerLink™ SutureTape, 0.9 mm, w/ loop (white/black)	AR-7559T
0 FiberWire® Suture, 38 in (blue)	AR-7254

References

- Ely EE, et al. *Orthopedics*. 2014;37(9):608-614. doi:10.3928/01477447-20140825-05
- Smith MJ, et al. *J Shoulder Elbow Surg*. 2020;29(8):1573-1583. doi:10.1016/j.jse.2019.11.028
- Gilot GJ, et al. *Arthroscopy*. 2015;31(8):1459-1465. doi:10.1016/j.arthro.2015.02.032

Graft Options

Product Description	Item Number
ArthroFlex® Options	
ArthroFlex Dermal Allograft, 25 mm × 30 mm × 2 mm, w/ Matracell® Decellularized Dermis	AFLEX202
ArthroFlex Dermal Allograft, 20 mm × 25 mm × 1.0 mm, w/ Matracell® Decellularized Dermis	AFLEX402
ArthroFlex Dermal Allograft, 25 mm × 30 mm × 1.0 mm, w/ Matracell® Decellularized Dermis	AFLEX403
ArthroFlex Dermal Allograft, 20 mm × 30 mm × 3.0 mm, Matracell Decellularized Dermis	AFLEX352
SpeedFLEX® Options	
SpeedFLEX Implant, 20 mm × 25 mm × 1 mm, pre-sutured graft	ABS-4402
SpeedFLEX Implant, 25 mm × 30 mm × 1 mm, pre-sutured graft	ABS-4403
SpeedFLEX Implant, 20 mm × 25 mm × 2 mm, pre-sutured graft	ABS-2202
SpeedFLEX Implant, 25 mm × 30 mm × 2 mm, pre-sutured graft	ABS-2203
SuturePatch Options	
SuturePatch Rectangle, 60 × 30 mm, 5/pack	AR-7290
SuturePatch Trapezoid, 60 × 35 mm, 5/pack	AR-7291
SuturePatch Trapezoid, 45 × 35 mm, 5/pack	AR-7292

Instruments


Product Description	Item Number
Graft Spreader	AR-19007GS
PassPort Button™ Cannula, 10 mm ID × 2 cm	AR-6592-10-20
PassPort Button Cannula, 10 mm ID × 3 cm	AR-6592-10-30
PassPort Button Cannula, 10 mm ID × 4 cm	AR-6592-10-40
PassPort Button Cannula, 10 mm ID × 5 cm	AR-6592-10-50
10 mm PassPort Inserter	AR-6592-10PI
FastPass Scorpion™ SL suture passer	AR-13999MF
HD Scorpion Needle w/ MegaLoader	AR-13999HDN
FiberWire Scissor	AR-11796
Knot Pusher / Suture Cutter with Portal Skid	AR-5845

Literature

Product Description	Item Number
CuffMend™ Rotator Cuff Augmentation System Surgical Technique	LT1-000201-en-US

ArthroFLEX and MatrACELL are trademarks of LifeNet Health.


Autograft Tissue Compression System



The Autograft Tissue Compression System is designed to create a graft using autologous tissue, such as biceps tendon. The graft can be used in various procedures, including rotator cuff augmentation and subscapularis repair.

Product Description	Item Number
Graft Preparation	
ATCS Instrument Set Includes press and case	AR-19081S
Baseplate Taper Assembly Press (b)	AR-9622
ATCS Instrument Case	AR-19081C
ATCS Compression Plates, single use, 5/pack	AR-19083
ATCS Compression Plates, 4 layer, single-use, 5/pack (a)	AR-19082
CuffMend Rotator Cuff Augmentation Repair	
Graft Spreader	AR-19007GS
0 FiberLoop® Suture, 20 in (blue), 76 mm straight needle w/ 7 mm loop	AR-7253
FiberLink™ SutureTape, 0.9 mm w/ loop (white/blue)	AR-7559
TigerLink™ SutureTape, 0.9 mm w/ loop (white/black)	AR-7559T
Subscap Repair	
ATCS Compression Plates, Long, single use, 5/pack	AR-19087
FiberTape® Tendon Compression Bridge	AR-7297-2
Two FiberTape tendon compression bridge assemblies, two 1.3 mm SutureTapes with needles, one free needle, and a drill bit	
FiberTape Cerclage Disposable Tensioner	AR-7820

PushLock® Suture Anchor



The unique knotless design allows tension to be adjusted intraoperatively while helping to ensure precise tissue reduction and fixation.

The self-punching anchor was developed to help speed completion of a SutureBridge™ repair. The titanium and PEEK tips minimize the need to prepare a bone socket for the lateral row, where soft tissue can sometimes obscure the view.


Product Description	Item Number
Biocomposite, 3.5 mm × 19.5 mm	AR-1926BC
Biocomposite, 3.5 mm × 19.5 mm, self-punching PEEK eyelet (e)	AR-1926BCSP
Biocomposite, 4.5 mm × 24 mm (a)	AR-1922BC
Biocomposite, 4.5 mm × 28 mm, self-punching* (c)	AR-1922BCM
PEEK, 3.5 mm × 19.5 mm	AR-1926PS
PEEK, 3.5 mm × 19.5 mm, self-punching PEEK eyelet	AR-1926PSSP
PEEK, 4.5 mm × 24 mm (b)	AR-1922PS
PEEK, 4.5 mm × 28 mm, self-punching* (d)	AR-1922PSM
Bio, 3.5 mm × 19.5 mm	AR-1926B
Required Instruments	
Punch for 3.5 mm PushLock anchor	AR-1926P
Punch for 4.5 mm PushLock and 4.5 mm Corkscrew® FT anchors	AR-1922P
Optional Instruments	
Disposable Punch for 3.5 mm PushLock anchor	AR-1926PBS
Disposable Punch for 4.5 mm PushLock and 4.5 mm Corkscrew FT anchors	AR-1922PBS

*Instrumentation not required for PushLock SP anchor

Literature

Product Description	Item Number
SutureBridge™ Double Row Rotator Cuff Repair Technique	LT1-0515-EN

Superior Capsular Reconstruction (SCR)



The SCR technique reinforces the superior capsule by attaching ArthroFLEX® dermal allograft from the glenoid to the greater tuberosity, keeping the humeral head centered on the glenoid.

Product Description	Item Number
SCR Knotless Implant System	AR-2600SCR
BioComposite SwiveLock® C Implant System w/ SCORPION-multifire needle Two 4.75 mm SwiveLock C anchors w/ 1 preloaded FiberTape® Loop (1 blue, 1 white/black) for medial row Two 4.75 mm SwiveLock C anchors for lateral row Disposable punch SCORPION-multifire needle Two 3.9 mm PEEK Knotless Corkscrew® w/ #5 suture Fishmouth Spear, Drill, Portal Dilator, 17-ga Spinal Needle w/ Stylet, Guide Pin, PassPort Cannula Divider	
Other	
3 mm ArthroFlex® Dermal Allograft w/ 40 mm × 70 mm Matrancell® Decellularized Dermis	AFLEX301
SCR Guide	AR-16950SR
Back Grasper w/ SR handle	AR-12531SR
PassPort Button™ Cannula, 12 mm ID × 3 cm	AR-6592-12-30
PassPort Button Cannula, 12 mm ID × 4 cm	AR-6592-12-40
PassPort Button Cannula, 12 mm ID × 5 cm	AR-6592-12-50
12 mm PassPort Inserter	AR-6592-12PI
PassPort Divider, 12 mm	AR-6592-12D
SutureTape, 1.3 mm, w/ tapered needle (white/blue)	AR-7500

ArthroFLEX and MatrACELL are trademarks of LifeNet Health.

Literature

Product Description	Item Number
Superior Capsular Reconstruction for Massive Irreparable Rotator Cuff Tears Surgical Technique	LT1-00020-EN

Biologic TuberoPlasty Technique



Biologic tuberoPlasty harnesses the power of ArthroFlex HDA to provide a biologic cushion between acromion and tuberosity, preventing bone on bone contact in an irreparable rotator cuff tear. Self punching tensionable knotless technology allows the ArthroFlex graft to be placed in a quick, effective, and reproducible technique preventing shoulder impingement syndrome, often associated with shoulder pain.

Product Description	Item Number
Implants, Medial Row Options	
Knotless 2.6 FiberTak® SP anchor w/ #5 suture, self-punching	AR-3641SP
Knotless BioComposite SwiveLock anchor, 4.75 mm × 19.1 mm w/ #2 suture (blue)	AR-2324KBCC
Knotless PEEK SwiveLock anchor, 4.75 mm × 19.1 mm w/ #2 suture (blue)	AR-2324KPSLC
Implants, Lateral Row Options	
BioComposite SwiveLock® SP, 4.75 mm × 24.5 mm w/ 1.3 mm SutureTape (white/blue), self-punching PEEK eyelet	AR-2324BCSP
PEEK SwiveLock SP anchor, 4.75 mm × 24.5 mm w/ 1.3 mm SutureTape, self-punching PEEK eyelet	AR-2324PSP
Bone Preparation	
PowerRasp™ Instrument, 5.5 mm × 13 cm	AR-8550PR
PowerPick™ Instrument, 45°, 6 mm drill depth	AR-8150PX-45
Other	
Arthroscopic Measurement Probe, 60°, 220 mm	AR-4070-01
SCR Guide	AR-16950SR
Back Grasper w/ SR handle	AR-12531SR
PassPort Button™ Cannula, 12 mm ID × 3 cm	AR-6592-12-30
PassPort Button Cannula, 12 mm ID × 4 cm	AR-6592-12-40
PassPort Button Cannula, 12 mm ID × 5 cm	AR-6592-12-50
12 mm PassPort Inserter	AR-6592-12PI
PassPort Divider, 12 mm	AR-6592-12D
FiberLink™ SutureTape, 1.3 mm, w/ loop (white/blue)	AR-7535
TigerLink™ SutureTape, 1.3 mm, w/ loop (white/black)	AR-7535T


Biologic TuberoPlasty Technique (Cont.)

Product Description	Item Number
Graft Options	
ArthroFlex® 20 mm x 30 mm x 3 mm, w/ Matracell® Decellularized Dermis	AFLEX352
ArthroFlex 25 mm x 35 mm x 4 mm, w/ Matracell® Decellularized Dermis	AFLEX600
ArthroFlex 20 mm x 30 mm x 4 mm, w/ Matracell® Decellularized Dermis	AFLEX602
ArthroFlex 20 mm x 30 mm x 5 mm, w/ Matracell® Decellularized Dermis	AFLEX652

Literature

Product Description	Item Number
Biologic TuberoPlasty Surgical Technique	LT1-000279-en-US


PASTA Bridge



Use this technique for transtendon repair of partial, articular-sided supraspinatus tendon avulsions, providing a secure bridging repair intended to enhance footprint compression and maximize contact between tendon and bone.

Product Description	Item Number
PASTA Bridge Repair Implant System	AR-1934PRK
Two 2.4 BioComposite SutureTak® anchors	
4.75 mm BioComposite SwiveLock® anchor	
Disposable spear	
Disposable punch	
17-ga spinal needle	
1.1 mm nitinol guidewire	
Dilator	

Knotless Corkscrew® Suture Anchor



The partially threaded Knotless Corkscrew anchor provides knotless soft-tissue fixation and has a tensionable suture-locking mechanism. Suture is passed and shuttled into the locking mechanism, allowing users to control the tension of knotless repairs. These anchors can be used for PASTA, remplissage, upper subscapularis, and SCR repairs.

3.9 mm Knotless Corkscrew Anchor

Product Description	Item Number
Biocomposite, 3.9 mm x 11.2 mm, w/ #5 suture (a)	AR-1941BC
PEEK, 3.9 mm x 11.2 mm, w/ #5 suture	AR-1941PSV
Instrumentation	
Drill, sterile	AR-1941D
Slotted Spear w/ Dilator	AR-1941DG
Fishmouth Spear w/ Sharp Obturator	AR-1941DGF
Circumferential Teeth Spear w/ Sharp Obturator	AR-1941CT
Punch	AR-1941P
Punch/Tap	AR-1941PTB
Disposables	
Disposable Kit Includes fishmouth spear, sharp obturator, and punch	AR-1941DS
TransTendon Percutaneous Instrument Kit Includes 17-ga spinal needle, 1.5 mm nitinol wire, dilator, fishmouth spear, and punch for 3.9 mm Knotless Corkscrew anchor	AR-1941PK

4.75 mm Double-Loaded Tensionable Knotless Corkscrew Anchor

Product Description	Item Number
BioComposite, 4.75 mm x 14.7 mm w/ two #2 knotless repair sutures (b)	AR-1943BC
PEEK, 4.75 mm x 14.7 mm w/ two #2 knotless repair sutures	AR-1943PS
Instrumentation	
Punch	AR-1927KPB-475

Literature

Product Description	Item Number
Rotator Cuff Repair with Double-Loaded Tensionable Knotless Anchors	LB1-000560-en-US
Knotless Corkscrew/PASTA Bridge Surgical Technique	LT1-00094-EN

Corkscrew® FT Suture Anchors



A fully threaded (FT) design over the entire length of the suture anchor provides cortical fixation. When inserting into hard bone, the internal drive mechanism substantially increases resistance to stripping. The recessed FiberWire® suture eyelet is self-aligning and minimizes suture abrasion during knot tying. The Corkscrew FT suture anchors are available in biocomposite (β -TCP and PLLA), PEEK, titanium, and bio materials.

Biocomposite

Product Description	Item Number
4.5 mm × 14 mm, w/ two #2 FiberWire sutures (c)	AR-1927BCF-45
4.5 mm × 14 mm, w/ two #2 TigerTail® sutures	AR-1927BCFT-45
4.5 mm × 14 mm, w/ two #2 FiberWire sutures and 4 needles	AR-1927BCNF-45
4.75 mm × 14 mm w/ two 1.3 mm SutureTape (white/black, white/blue)	AR-1927BCT-475
4.75 mm × 14 mm, w/ three #2 TigerTail sutures (d)	AR-1927BCF-475
5.5 mm × 14.7 mm, w/ two 1.3 mm SutureTape (white/black, white/blue) (a)	AR-1927BCT
5.5 mm × 14.7 mm w/ three 1.3 mm SutureTapes (b)	AR-1927BCT-3
5.5 mm × 14.7 mm, w/ two #2 FiberWire sutures (e)	AR-1927BCF
5.5 mm × 14.7 mm, w/ two #2 TigerTail sutures	AR-1927BCFT
5.5 mm × 14.7 mm, w/ two #2 FiberWire suture and 4 needles	AR-1927BCNF
5.5 mm × 14.7 mm, w/ three #2 FiberWire sutures	AR-1927BCF-3
6.5 mm × 14.7 mm, w/ two #2 FiberWire sutures (f)	AR-1927BCF-65

PEEK

Product Description	Item Number
4.5 mm × 14 mm, w/ two #2 FiberWire sutures (g)	AR-1927PSF-45
4.5 mm × 14 mm, w/ two #2 FiberWire suture and 4 needles	AR-1927PNF-45
4.75 mm × 14 mm, w/ two 1.3 mm SutureTape (white/black, white/blue)	AR-1927PST-475
4.75 × 14 mm, w/ three #2 TigerTail sutures	AR-1927PSF-475
5.5 mm × 14.7 mm w/ two 1.3 mm SutureTapes and 4 needles	AR-1927PNT
5.5 mm × 14.7 mm w/ two 1.3 mm SutureTapes	AR-1927PST
5.5 mm × 14.7 mm w/ three 1.3 mm SutureTapes	AR-1927PST-3
5.5 mm × 14.7 mm, w/ two #2 FiberWire sutures (h)	AR-1927PSF
5.5 mm × 14.7 mm, w/ two #2 FiberWire sutures and 4 needles	AR-1927PNF
5.5 mm × 14.7 mm, w/ three #2 FiberWire sutures	AR-1927PSF-3
6.5 mm × 14.7 mm, w/ two #2 FiberWire sutures (i)	AR-1927PSF-65

Titanium (No Instrumentation Required)

Product Description	Item Number
4.5 mm × 14 mm, w/ two #2 FiberWire sutures (j)	AR-1928SF-45
5.5 mm × 16.3 mm w/ two 1.3 mm SutureTapes	AR-1928ST
5.5 mm × 16.3 mm w/ three 1.3 mm SutureTapes	AR-1928ST-3
5.5 mm × 16.3 mm w/ two 1.3 mm SutureTapes and 4 needles	AR-1928STN
5.5 mm × 16.3 mm, w/ two #2 FiberWire sutures (k)	AR-1928SF-2
5.5 mm × 16.3 mm, w/ two #2 FiberWire sutures and 4 needles	AR-1928SNF-2
5.5 mm × 16.3 mm, w/ two #2 TigerTail sutures	AR-1928SFT-2
5.5 mm × 16.3 mm, w/ three #2 FiberWire sutures	AR-1928SF-3
6.5 mm × 16.3 mm, w/ three #2 FiberWire sutures (l)	AR-1929SF-3


Bio

Product Description	Item Number
5.5 mm × 14.7 mm, w/ two #2 FiberWire sutures (m)	AR-1927BF


Corkscrew FT Instrumentation

Product Description	Item Number
4.5 mm Instrumentation	
Punch for 4.5 mm Corkscrew FT anchors (n)	AR-1922P
Punch/Tap for 4.5 mm Corkscrew FT anchor	AR-1927PTB-45
Disposable Punch for 4.5 mm Corkscrew FT and 4.5 mm PushLock® anchors	AR-1922PBS
Punch w/ Cortical Tap for 4.5 mm Corkscrew FT anchor	AR-1927PTB-45-2
4.75 mm Instrumentation	
Punch/Tap for 4.75 mm Corkscrew FT anchor (required) (o)	AR-1927PTB-475
5.5 mm Instrumentation	
Punch for 5.5 mm Corkscrew FT and 4.75 mm and 5.5 mm SwiveLock® anchors (p)	AR-1927PB
Punch/Tap for 5.5 mm Corkscrew FT and 5.5 mm SwiveLock anchors (q)	AR-1927CTB
Disposable Punch for 5.5 mm Corkscrew FT and 4.75 mm and 5.5 mm SwiveLock anchors	AR-1927PBS
Spade-Tip Drill for 5.5 mm Corkscrew FT anchor (In hard bone, a drill may be used to prevent bone cracks)	AR-1927D

Corkscrew® FT Suture Anchors (Cont.)



Corkscrew Suture Anchor



The titanium Corkscrew anchor design incorporates a cancellous thread with a very small core diameter to maximize pullout strength in cancellous or osteopenic bone.

Product Description	Item Number
3.5 mm × 12 mm, w/ one #2 FiberWire® suture (a)	AR-1915SF
5 mm × 15.5 mm, w/ two #2 FiberWire sutures (2 suture eyelets) (b)	AR-1902SF
5 mm × 15.5 mm, w/ two #2 FiberWire suture (c)	AR-1920SF
5 mm × 15.5 mm, w/ needles and two #2 FiberWire sutures	AR-1920NSF
6.5 mm × 15.5 mm, w/ two #2 FiberWire sutures (d)	AR-1925SF

Note: Instrumentation not required for titanium Corkscrew anchors

SLAP and Bankart Repair

Knotless 1.8 FiberTak® Soft Anchor.....	18
FiberTak Soft Anchor.....	19
PushLock® Suture Anchor	20
Labral SwiveLock® Suture Anchor	21
SutureTak® Suture Anchor.....	22
3 mm Knotless SutureTak Suture Anchor.....	23
FASTak™ Suture Anchor	23
Spears	24
Offset Guide	24
Cannulated Guide.....	25

Knotless 1.8 FiberTak® Soft Anchor



This tensionable knotless suture anchor combines the benefits of soft anchors with knotless soft-tissue fixation. Curved and straight guide delivery systems are available with 1.8 mm drills.

Product Description	Item Number
Knotless 1.8 FiberTak Soft Anchor w/ #2 suture, qty. 5	AR-3636
Knotless 1.8 FiberTak Implant System (a)	AR-3637
Includes: three 1.8 Knotless FiberTak anchors (blue, blue/black, white/blue), Disposable Curved Instrument Kit (1.8 mm drill, tapered curved spear and obturator), and Low-Profile Percutaneous Instrument Kit (needle with stylet, K-wire, cannula, and obturator)	
Disposables	
Curved Disposable Kit, blunt Includes 1.8 mm flexible drill w/ hub, curved circumferential teeth spear, and blunt obturator	AR-3638DC
Tapered Curved Disposable Kit, sharp Includes 1.8 mm flexible drill w/ hub, tapered curved circumferential teeth spear, and sharp obturator	AR-3610DC-3
Straight Disposable Kit Includes 1.8 mm rigid drill, straight circumferential teeth spear, and blunt obturator	AR-3638DS
Percutaneous Insertion Kit Includes 1.8 mm drill, fishmouth spear, dilator, needle with stylet, and guide pin	AR-3610PK-3
Low-Profile Percutaneous Instrument Kit Includes 5 mm × 7 cm low-profile cannula, obturator, needle with stylet, and k-wire	AR-6548PK
1.8 mm Flexible Drill, w/ hub and sharp obturator, sterile	AR-3610ND-2
1.8 mm Rigid Drill, sterile	AR-3600D-2
1.8 mm Flexible ShaverDrill Device, w/ sharp obturator, sterile	AR-3610NSD-2
1.8 mm Flexible Fluted Drill, w/ hub and sharp obturator, sterile	AR-3610NDF-2
1.9 mm Flexible Drill, w/ hub and sharp obturator, hard bone, sterile	AR-3610ND-4
1.9 mm Rigid Drill, hard bone, sterile	AR-3600D-4
1.9 mm Flexible Fluted Drill, w/ hub and sharp obturator, sterile	AR-3610NDF-4
Reusables	
Curved Circumferential Teeth Spear, w/ flexible sharp obturator, reusable	AR-3610CTC
Tight Curved Circumferential Teeth Spear, w/ flexible sharp obturator, reusable	AR-3610CTC-2

Fishmouth Spear, reusable	AR-3610F
Circumferential Teeth Spear, reusable	AR-3610CT
Slotted Spear, reusable	AR-3610ST
Pivot Guide, reusable	AR-3610R

Short Depth Options

Product Description	Item Number
Disposables	
Knotless FiberTak Curved Spear Disposable Kit, Short Includes: 1.8 mm flexible drill w/ hub, curved circumferential teeth spear, short w/ flexible sharp obturator	AR-3610DCS
Knotless FiberTak Tight Curved Spear Disposable Kit, Short Includes: 1.8 mm flexible drill w/ hub, tight curved circumferential teeth spear, short w/ flexible sharp obturator	AR-3610DTCS
Reusables	
Curved Circumferential Teeth Spear, short, w/ flexible sharp obturator, reusable	AR-3610CTCS
Tight Curved Circumferential Teeth Spear, short, w/ flexible sharp obturator, reusable	AR-3610CTCS-2
Fishmouth Spear, short, reusable	AR-3610FS
Circumferential Teeth Spear, short, reusable	AR-3610CTS
Slotted Spear, short, reusable	AR-3610STS
Sharp Obturator, short, reusable	AR-3610RTS

2.6 FiberTak Soft Anchor


Product Description	Item Number
Knotless 2.6 FiberTak Soft Anchor, w/ #5 suture, qty. 5 (may be used as bail out anchor for glenoid)	AR-3641

2.6 FiberTak instruments listed on page 8

Literature

Product Description	Item Number
Knotless 1.8 FiberTak® Suture Anchor Glenoid Labrum Surgical Technique	LT1-000002-en-US

FiberTak® Soft Anchor



FiberTak soft anchors are available single- and double-loaded each with multiple innovative suture options. These anchors are complemented with both curved- and straight-guide delivery systems and small drills to minimize bone removal.


Product Description	Item Number
With #2 coreless suture (a)	AR-3600
With two #2 coreless sutures (b)	AR-3600-2
With 1.3 mm SutureTape (c)	AR-3602
With two 1.3 mm SutureTape (d)	AR-3602-2
Disposables	
Disposables Kit , straight spear Includes circumferential teeth spear, blunt obturator, and 1.6 mm drill	AR-3600D
Curved Disposables Kit , curved spear Includes curved circumferential teeth spear, blunt obturator, and 1.6 mm flexible drill with hub	AR-3600DC
Tapered Curved Disposables Kit , curved spear Includes tapered curved circumferential teeth spear, sharp obturator, and 1.6 mm flexible drill w/hub	AR-3610DC-1
Disposables Kit , ShaverDrill™ accessory, flexible Includes curved circumferential teeth spear, blunt obturator, and 1.6 mm flexible ShaverDrill™ accessory	AR-3600NDSR-1
Percutaneous Insertion Kit w/ spear, dilator, needle with stylet, guide pin, and 1.6 mm drill	AR-3610PK-1
Low-Profile Percutaneous Instrument Kit Includes 5 mm × 7 cm low-profile cannula, obturator, Needle with stylet, and k-wire	AR-6548PK
1.6 mm Drill, sterile	AR-3600D-1
1.6 mm Flexible Drill w/ hub and sharp obturator, sterile	AR-3610ND-1
1.6 mm Flexible Fluted Drill, w/ hub and sharp obturator, sterile	AR-3610NDF-1
1.7 mm Drill, sterile	AR-3600D-3
1.7 mm Flexible Drill w/ hub and sharp obturator, sterile	AR-3610ND-3
1.7 mm Flexible Fluted Drill, w/ hub and sharp obturator, sterile	AR-3610NDF-3

Product Description	Item Number
Reusables	
Curved Circumferential Teeth Spear w/ flexible sharp obturator, reusable	AR-3610CTC
Tight Curved Circumferential Teeth Spear w/ flexible sharp obturator, reusable	AR-3610CTC-2
Fishmouth Spear, reusable	AR-3610F
Circumferential Teeth Spear, reusable	AR-3610CT
Slotted Spear, reusable	AR-3610ST
Pivot Guide, reusable	AR-3610R
Short Depth Reusable Options	
Curved Circumferential Teeth Spear, short, w/ flexible sharp obturator, reusable	AR-3610CTCS
Tight Curved Circumferential Teeth Spear, short, w/ flexible sharp obturator, reusable	AR-3610CTCS-2
Fishmouth Spear, short, reusable	AR-3610FS
Circumferential Teeth Spear, short, reusable	AR-3610CTS
Slotted Spear, short, reusable	AR-3610STS
Sharp Obturator, short, reusable	AR-3610RTS

Literature

Product Description	Item Number
FiberTak® Anchor Glenoid Labrum Repair Surgical Technique	LT1-00021-EN

PushLock® Suture Anchor



PushLock suture anchors are designed for simple and secure arthroscopic glenohumeral joint instability repair. This knotless technique saves valuable time and eliminates the possibility of knot impingement. The “suture-first” technique allows suture to be passed independently through a desired amount of tissue prior to anchor implantation. Multiple suture configurations, including simple, cinch, and mattress stitches, can be accomplished with FiberWire®, FiberLink™, FiberStick™, and LabralTape™ sutures using the Labral Scorpion™ or QuickPass™ SutureLasso™ suture retriever. Tissue tension can be visualized and adjusted, if necessary, prior to final anchor implantation.

2.4 mm PushLock Anchor

Product Description	Item Number
Biocomposite, 2.4 mm × 11.3 mm (a)	AR-2922BC
PEEK, 2.4 mm × 11.3 mm (b)	AR-2922PS
Note: Use with FiberLink SutureTape [AR-7535]	
Required	
Spear and Trocar- and Blunt-Tip Obturators, for 2.4 mm and 2.9 mm PushLock anchors	AR-1949
Drill for 2.4 mm PushLock anchor, sterile	AR-2922D-24-2
Drill for 2.4 mm PushLock anchor, hard bone, sterile	AR-2922D-24-3
Disposable	
Percutaneous Insertion Kit	AR-2922PK
Includes spear, dilator, needle w/ stylet, guide pin, drill, and cannula	
ShaverDrill™ Accessory	AR-2922DSR-24
ShaverDrill Accessory, hard bone	AR-2922DTSR-24

Optional Instruments for 2.4 mm and 2.9 mm PushLock Anchors

Product Description	Item Number
Offset Guide for 2.4 mm and 2.9 mm PushLock anchors	AR-1934R
Spear w/ Circumferential Teeth and Sharp Obturator for 2.4 mm and 2.9 mm PushLock anchors	AR-1946
Mini Suture Cutter, 3.4 mm, straight Note: May be used w/ metal cannula [AR-1923MCS]	AR-13255
Mini Suture Cutter, 3.4 mm, straight w/ WishBone™ handle Note: May be used w/ metal cannula [AR-1923MCS]	AR-13255W
Mini FiberTape® Retriever w/ SR handle	AR-12974SR
Mini FiberTape Retriever w/ NR handle	AR-12974NR
Mini FiberTape Retriever w/ WishBone handle	AR-12974W

2.9 mm PushLock Anchor

Product Description	Item Number
Biocomposite, 2.9 mm × 12.5 mm (short) (c)	AR-2923BC
PEEK, 2.9 mm × 12.5 mm (short) (d)	AR-2923PS
Implant System	
Short 2.9 mm BioComposite PushLock Implant System 2.9 mm short BioComposite PushLock anchor FiberLink SutureTape 2.9 mm short PushLock Percutaneous Kit (Includes spear, dilator, 17-ga spinal needle, guide pin, drill, and cannula)	AR-2923BCS
Required Instruments	
Spear, Sharp and Blunt-Tip Obturators, for 2.4 mm and 2.9 mm PushLock anchors	AR-1949
Drill for short 2.9 mm PushLock anchor	AR-2923D
Drill for short 2.9 mm PushLock anchor, sterile	AR-2923D-ST
Drill for short 2.9 mm PushLock anchor, hard bone	AR-2923DT
Drill for short 2.9 mm PushLock anchor, hard bone, sterile	AR-2923DT-ST
ShaverDrill Accessory for short anchor	AR-2923DSR



PushLock® Suture Anchor (Cont.)

Product Description	Item Number
Disposables	
Disposables Kit for short anchor, w/ metal spear, trocar, and drill	AR-2923DS
Percutaneous Insertion Kit Includes disposable 17-ga spinal needle, 1.1 mm nitinol wire, portal dilator, spear, drill, and percutaneous metal cannula	AR-1923PK
ShaverDrill™ Disposables Kit for short anchor	AR-2923DSSK
ShaverDrill Disposables Kit for short anchor, hard bone, 5/pack	AR-2923DTSK
ShaverDrill for short anchor, hard bone, 5/pack	AR-2923DTSR
Optional Instruments	
Metal Cannula Set for 2.9 mm PushLock anchor	AR-1923MCS
Disposable Silicone Dam for AR-1923MCS, 10/box	AR-1923MC-03
Spade-Tip Drill for 2.9 mm PushLock anchor	AR-1923SD


3.5 mm PushLock Anchor

Product Description	Item Number
Biocomposite, 3.5 mm × 19.5 mm (e)	AR-1926BC
PEEK, 3.5 mm × 19.5 mm (f)	AR-1926PS
Bio, 3.5 mm × 19.5 mm (g)	AR-1926B
Required Instruments	
Spear and Sharp Obturator for 3.5 mm PushLock and SwiveLock anchors	AR-1907
Drill for 3.5 mm PushLock anchor	AR-1912
Disposables	
Disposables Kit for 3.5 mm PushLock anchor, w/ metal spear and drill	AR-1926DS
Disposables Kit for 3.5 mm PushLock anchor, w/ offset guide and drill	AR-1926DS-2
Optional Instruments	
Offset Guide for 3.5 mm PushLock and SwiveLock anchors	AR-1909R
Spear w/ Circumferential Teeth and Sharp Obturator for 3.5 mm PushLock and SwiveLock anchors	AR-1906
Spade-Tip Drill for 3.5 mm PushLock anchor	AR-1911

Literature

Product Description	Item Number
PushLock® Knotless Instability Repair	LB1-0520-EN

Labral SwiveLock® Suture Anchor




Screw-in, knotless SwiveLock anchors have been miniaturized to treat shoulder instability. Used like a labral PushLock anchor, this anchor offers all of the same knotless advantages (low-profile, quick, "suture-first" technique with multiple stitch configuration options) for surgeons who prefer a threaded anchor.

Product Description	Item Number
Biocomposite, 3.5 mm × 15.8 mm (a)	AR-2325BCC
PEEK, 3.5 mm × 15.8 mm (b)	AR-2325PSLC
Required Instruments	
Spear and Sharp Obturator for 3.5 mm PushLock® and SwiveLock anchors	AR-1907
Drill for 3.5 mm SwiveLock anchor	AR-2325D
Drill for 3.5 mm SwiveLock anchor, sterile	AR-2325D-ST
Optional Instruments	
Offset Guide for 3.5 mm PushLock and SwiveLock anchors	AR-1909R
Spear w/ Circumferential Teeth and Sharp Obturator for 3.5 mm PushLock and SwiveLock anchors	AR-1906

Recommended Suture for Use With PushLock and SwiveLock Anchors

Product Description	Item Number
SutureTape Link, 1.3 mm, w/ loop (white/blue)	AR-7535
SutureTape Link, 1.3 mm, w/ loop (white/black)	AR-7535T
LabralTape™ Suture, 1.5 mm, 36 in tape (white)	AR-7276
LabralTape Suture, 1.5 mm, 36 in tape (white/black)	AR-7276T
#2 FiberLink® Suture, w/ closed loop, 26 in (blue)	AR-7235
#2 TigerLink® Suture, w/ closed loop, 26 in (white/black)	AR-7235T
FiberLink SutureTape, 1.7 mm w/ loop (white/blue)	AR-7538
TigerLink SutureTape, 1.7 mm w/ loop (white/black)	AR-7538T
FiberStick® Suture, #2 FiberWire® suture, 50 in (blue), one end stiffened, 12 in	AR-7209
TigerStick® Suture, #2 TigerWire® suture, 50 in (white/black), one end stiffened, 12 in	AR-7209T
#1 FiberWire Suture for 2.4 mm PushLock anchor, 38 in (blue)	AR-7216

SutureTak® Suture Anchor



SutureTak anchors are available in a variety of sizes, materials, and suture options. Biocomposite anchor options contain a molded-in suture to reinforce the strength of the anchor and suture eyelet. The nonabsorbable PEEK SutureTak anchor has a material eyelet that provides superior abrasion resistance due to PEEK's low coefficient of friction. Anchor insertion and delivery are made simple by drilling through a spear and inserting the anchor through the same spear into the drilled hole.

2.4 mm SutureTak

Product Description	Item Number
Biocomposite, 2.4 mm × 12 mm, w/ two #2 FiberWire sutures (a)	AR-1934BCF-24-2
Required Instruments	
Instrumentation Set Includes	AR-1934-24S
Spear and Sharp and Blunt Obturators	AR-1948
Step Drill	AR-1934D-24
SutureTak Instrumentation Case	AR-1934-24C
Disposables	
Percutaneous Insertion Kit	AR-1934PI
Includes disposable 17-ga spinal needle, 1.1 mm nitinol wire, portal dilator, spear, and drills	
Disposables Kit	AR-1934-24DS
Includes AR-1945S and AR-1934D-24-1	
Optional	
Drill, hard bone	AR-1934D-24-1
Drill, very hard bone	AR-1934D-24-2
Offset Guide	AR-1948R
Spear w/ Circumferential Teeth and Sharp Obturator	AR-1948CT

Literature

Product Description	Item Number
SutureTak® Anchor Bankart and SLAP Repair Surgical Technique	LT1-0507-EN

3 mm SutureTak Anchor

Product Description	Item Number
Biocomposite, 3 mm × 14.5 mm w/ 1.3 mm SutureTape (white/blue) (b)	AR-1934BCT
Biocomposite, 3 mm × 14.5 mm, w/ two 1.3 mm SutureTape (blue, black/white) (c)	AR-1934BCT-2
Biocomposite, 3 mm × 14.5 mm, w/ two #2 FiberWire suture	AR-1934BCF-2
PEEK, 3 mm × 14.5 mm, w/ two 1.3 mm SutureTape (blue, black/white) (d)	AR-1934PST-2
PEEK, 3 mm × 12.2 mm, w/ two #2 FiberWire suture	AR-1934PS-2
Required Instruments	
Step Drill	AR-1250LT
Instrumentation Set Includes	AR-1934S
Spear and Sharp and Blunt Obturators	AR-1949
SutureTak Instrumentation Case	AR-1934C
Disposable Kits	
Disposables Kit	AR-1934DS-2
Includes disposable spear and AR-1250LT	
Percutaneous Insertion Kit	AR-1934PI-30
Includes disposable 17-ga spinal needle, 1.1 mm nitinol wire, portal dilator, spear, and drill	
Optional Instruments	
Offset Guide	AR-1934R
Spear w/ Circumferential Teeth and Sharp Obturator	AR-1946
Disposable Spear and Sharp Obturator	AR-1949S
Disposable Offset Guide	AR-1934GS
Disposable ShaverDrill™ Accessory, 5/pack	AR-1934DSSK
ShaverDrill Accessory	AR-1250LTSR
Spade-Tip Drill	AR-1257
Punch for 3 mm SutureTak	AR-1934P
Needle for Portal Dilator	AR-6521
For Open Procedures	
Short Spear	AR-1326G
Short Spade-Tip Drill	AR-1256

3 mm Knotless SutureTak Suture Anchor



The knotless SutureTak anchor locking mechanism uses a simple and reproducible technique, allowing controlled tensioning of a knotless soft-tissue fixation repair under direct visualization.

Product Description	Item Number
Biocomposite, 3 mm × 12.7 mm w/ #2 suture	AR-1936BC
PEEK, 3 mm × 12.7 mm, w/ #2 suture	AR-1936PS
Biocomposite, 3 mm × 12.7 mm, w/ #2 suture	AR-1938BC
PEEK, 3 mm × 12.7 mm, w/ #2 suture	AR-1938PS
Required Instruments	
Spear and Sharp and Blunt Obturators	AR-1949
2.6 mm Drill, hard bone, sterile	AR-1938D
Disposable Kits	
Disposables Kit Includes spear, drill, and cutter	AR-1938DS
Percutaneous Insertion Kit Includes disposable 17-ga spinal needle, 1.1 mm nitinol wire, portal dilator, spear, drill, PEEK cannula, and cutter	AR-1938PK
Optional Instruments	
Suture Tensioner and Cutter	AR-1938TC
Drill, very hard bone, sterile	AR-1938DL
ShaverDrill™ Accessory, hard bone	AR-1938DSR
Disposable ShaverDrill Kit, hard bone, 5/pack	AR-1938DSSK
For Open Procedures	
Biocomposite, 3 mm × 12.7 mm, w/ #2 suture	AR-1936BCS
PEEK, 3 mm × 12.7 mm, w/ #2 suture	AR-1936PSS
Biocomposite, 3 mm × 12.7 mm, w/ #2 suture	AR-1938BCS
PEEK, 3 mm × 12.7 mm, w/ #2 suture	AR-1938PSS
Disposable Kit Includes step drill and spear	AR-1938DSS

Literature

Product Description	Item Number
Knotless SutureTak® anchor Instability Repair Surgical Technique	LT1-0511-EN



FASTak™ Suture Anchor



Titanium FASTak suture anchors come preloaded with #2 FiberWire® suture. In most cases, the implant can be manually inserted using the handled version. A drill is available for use in hard bone if desired.

Product Description	Item Number
2.8 mm × 11.7 mm, w/ #2 FiberWire suture and handle	AR-1324HF
2.8 mm × 11.7 mm, w/ #2 FiberWire suture	AR-1324SF
Drill for 2.8 mm FASTak II anchor	AR-1324D

Spears

Spears with removable trocar and V-shaped tips allow precise control and accurate anatomical placement of the suture anchors on the glenoid rim, in one simple step. The small-diameter shaft with trocar facilitates percutaneous anchor placement without using a cannula for SLAP repairs or a direct approach through the subscapularis.

Product Description	Item Number
Spear w/ sharp and blunt obturators for 2.8 mm FASTak™ II, 3 mm SutureTak®, and 2.4 mm and 2.9 mm PushLock® anchors	AR-1949
Disposable Spear w/ sharp obturator for 2.8 mm FASTak II, 3 mm SutureTak, 2.4 mm and 2.9 mm PushLock anchors, single use	AR-1949S
Circumferential Teeth Spear w/ sharp obturator for 2.8 mm FASTak II, 3 mm SutureTak, and 2.4 mm and 2.9 mm PushLock anchors	AR-1946
Spear w/ sharp and blunt obturators for FiberTak®, 2.4 mm FASTak, and 2.4 mm SutureTak anchors	AR-1948
Circumferential Teeth Spear w/ sharp obturator for FiberTak, 2.4 mm FASTak, and 2.4 mm SutureTak anchors	AR-1948CT
Fishmouth Spear for FiberTak and Knotless FiberTak anchors	AR-3610F
Circumferential Teeth Spear for FiberTak and Knotless FiberTak anchors	AR-3610CT
Curved Circumferential Teeth Spear w/ flexible sharp obturator for FiberTak and Knotless FiberTak anchors	AR-3610CTC
Tight Curved Circumferential Teeth Spear w/ flexible sharp obturator for FiberTak and Knotless FiberTak anchors (a)	AR-3610CTC-2
Pivot Guide w/ internal curved shaft for FiberTak and Knotless FiberTak anchors, reusable	AR-3610R
Slotted Spear for FiberTak and Knotless FiberTak anchors (b)	AR-3610ST
Curved Circumferential Teeth Spear, short, w/ flexible sharp obturator for FiberTak and Knotless FiberTak anchors, reusable	AR-3610CTCS
Tight Curved Circumferential Teeth Spear, short, w/ flexible obturator for FiberTak and Knotless FiberTak anchors, reusable	AR-3610CTCS-2
Fishmouth Spear, short, for FiberTak and Knotless FiberTak anchors, reusable	AR-3610FS
Circumferential Teeth Spear, short, for FiberTak and Knotless FiberTak, reusable	AR-3610CTS
Short Trocar, for FiberTak and Knotless FiberTak, reusable	AR-3610RTS
Short Slotted Spear, for FiberTak and Knotless FiberTak, reusable	AR-3610STS
Spear w/ sharp obturator for 3.5 mm PushLock and 3.5 mm SwiveLock anchors (c)	AR-1907
Circumferential Teeth Spear w/ sharp obturator for 3.5 mm PushLock and 3.5 mm SwiveLock® anchors (d)	AR-1906

Slotted Spear w/ dilator for 3.9 mm Knotless Corkscrew anchor (e)	AR-1941DG
Fishmouth Spear w/ sharp obturator for 3.9 mm Knotless Corkscrew and 2.6 FiberTak anchors	AR-1941DGF
Circumferential Teeth Spear w/ sharp obturator for 3.9 mm Knotless Corkscrew® and 2.6 FiberTak anchors	AR-1941CT
Angled Circumferential Teeth Spear for 2.6 FiberTak anchors (f)	AR-3655

Offset Guide


The offset guides allow surgeons to accurately position suture anchor implants on the glenoid for anatomic reconstruction of glenohumeral ligament structures during Bankart and SLAP repairs. All offset guides enable surgeons to easily reproduce a 1.5 mm medial offset position relative to the glenoid rim to complete an anatomic reconstruction of the labral tissue.

Product Description	Item Number
Offset Guide for 3.5 mm PushLock® and 3.5 mm SwiveLock® anchors (g)	AR-1909R
Offset Guide for 2.8 mm FASTak™ II, 3 mm SutureTak®, and 2.4 mm and 2.9 mm PushLock® anchors	AR-1934R
Offset Guide for FiberTak®, 2.4 mm FASTak, and 2.4 mm SutureTak anchors	AR-1948R
Disposable Offset Guide for 2.8 mm FASTak II, 3 mm SutureTak, and 2.4 mm and 2.9 mm PushLock anchors	AR-1934GS

Cannulated Guide

Cannulated guides are designed for inserting anchors during open or arthroscopic Bankart repairs. These guides have a larger wall thickness than the spears to increase instrument strength; the guides also include a larger dovetail tip.


Product Description	Item Number
Cannulated Guide for 2.8 mm FASTak II, 3 mm SutureTak, and 2.4 mm and 2.9 mm PushLock anchors	AR-1317



FiberWire® Suture

SutureTape Suture	28
SutureTag	28
FiberWire® and TigerWire® Sutures.....	29
Collagen-Coated FiberWire Suture.....	29
FiberTape® Suture.....	30
LabralTape™ Suture.....	30
FiberLink™ and TigerLink™ Suture.....	31
FiberLoop® Suture	31
FiberStick™ and TigerStick® Suture.....	32
FiberSnare® Suture.....	32
FiberButton.....	33
FiberWire® Suture Cart.....	33
FingerShield™ Finger Guard	33
Free Needles.....	33

SutureTape Suture




SutureTape is a flat-braided suture that is similar in strength to #2 FiberWire® suture, but is 21% more resistant to tissue pull-through.¹ It offers improved handling and knot security and smaller comparable knot stacks.

Product Description	Item Number
SutureTape, 1.3 mm, 40 in (white/blue), w/ tapered ends and tapered needle, 26.5 mm, ½ circle (b)	AR-7500
SutureTape, 1.3 mm, 40 in (white/black) w/ tapered ends and tapered needle, 26.5 mm, 1/2 circle	AR-7500T
SutureTape, X-pattern, 1.3 mm, two each 40 in strands (1 white/blue, 1 white/black) (a)	AR-7501
SutureTape, X-pattern, 1.3 mm, two 40 in strands (1 blue/black, 1 black/white) w/ tapered needle, 36.6 mm, ½ circle	AR-7501T
SutureTape, 1.3 mm, 40 in (white/blue), w/ reverse cutting needle, 36.6 mm, ½ circle	AR-7502
SutureTape, 1.3 mm, 40 in (blue), w/ tapered needle, 36.6 mm, ½ circle (d)	AR-7506
SutureTape, 1.3 mm, 40 in (black/white), w/ tapered needle, 36.6 mm, ½ circle (e)	AR-7506T
FiberStick SutureTape, 1.3 mm (blue/blue), one end stiffened, 12 in	AR-7509
FiberStick SutureTape, 1.3 mm (black/white), one end stiffened, 12 in	AR-7509T
SutureTape, 1.7 mm, 38 in (white/blue), w/ conventional cutting needle, 48 mm, ½ circle	AR-7511
SutureTape, 1.7 mm, 38 in (white/black) w/ conventional cutting needle, 48 mm 1/2 circle	AR-7511T
SutureTape, X-pattern, 1.7 mm, two 42 in strands (1 white/blue, 1 white/black) w/ tapered needle, 36.6 mm, ½ circle	AR-7518
SutureTape, 0.9 mm, 38 in (white/blue)	AR-7521
FiberLoop® SutureTape, 1.3 mm, 20 in loop (white/blue) w/ 76 mm straight needle	AR-7534
TigerLoop™ SutureTape, 1.3 mm, 20 in loop (white/black) w/ 76 mm straight needle (c)	AR-7534T
FiberLink™ SutureTape, 1.3 mm, w/ loop (white/blue)	AR-7535
TigerLink™ SutureTape, 1.3 mm, w/ loop (white/black)	AR-7535T
FiberLink SutureTape, 1.7 mm w/ loop (white/blue)	AR-7538
TigerLink SutureTape, 1.7 mm w/ loop (white/black)	AR-7538T

Reference

1. Arthrex, Inc. LA1-00038-EN_C. Naples, FL; 2021.


SutureTag



The SutureTag looping stitch creates a running ripstop to reinforce the tendon with minimal passes through the tissue. This novel stitch is lower profile, faster to complete, and stronger than a typical Krackow repair stitch. The SutureTag splits into two #2 sutures that can be easily incorporated into the FiberTak® button repair technique.

Product Description	Item Number
SutureTag w/ double 14 in #2 tails and 36.6 mm ½ circle tapered needle (white/black)	AR-7274T
SutureTag w/ double 14 in #2 tails and 36.6 mm ½ circle tapered needle (white/blue/black)	AR-7274TT

FiberWire® and TigerWire® Sutures



FiberWire suture is constructed of a multistranded long chain ultra-high-molecular weight polyethylene (UHMWPE) core with a braided jacket of polyester and UHMWPE that gives FiberWire suture superior strength,¹ soft feel, and abrasion resistance² that is unequalled in orthopedic surgery. Suture breakage during knot tying is virtually eliminated, especially critical during arthroscopic procedures.³

Product Description	Item Number
#5 FiberWire Suture, 38 in (blue)	AR-7210
#5 FiberWire Suture, 38 in (blue) w/ conventional cutting needle, 48 mm, ½ circle	AR-7211
#5 FiberWire Suture, 38 in (1 blue, 1 white, 1 white/black) w/ conventional cutting needle, 48 mm, ½ circle	AR-7212
#5 FiberWire Suture, 38 in (blue) w/ conventional cutting needle, 48 mm, ½ circle, and tapered needle, 26.5 mm, ½ circle	AR-7213
#5 FiberWire Suture, 38 in (white)	AR-7215
#2 FiberWire Suture, 38 in (blue) w/ tapered needle, 26.5 mm, ½ circle	AR-7200
#2 FiberWire Suture, 38 in (blue) w/ tapered needle, 26.5 mm, ½ circle, collagen-coated	AR-7200B
#2 FiberWire Suture, 38 in (blue) w/ reverse cutting needle, 36.6 mm, ½ circle	AR-7202
#2 FiberWire Suture, 38 in (blue) w/ two tapered needles, 26.5 mm ½ circle	AR-7205
#2 TigerWire Suture, 38 in (white/black) w/ two tapered needles, 26.5 mm, ½ circle	AR-7205T
#2 FiberWire Suture, 38 in (blue) w/ tapered needle, 36.6 mm, ½ circle	AR-7206
#2 FiberWire Suture, 38 in (1 blue, 1 white/black) w/ tapered needle, 26.5 mm, ½ circle	AR-7208
#2 FiberWire Suture, 38 in (1 blue, 1 white/black) w/ reverse cutting needle, 36.6 mm, ½ circle	AR-7217
#2 FiberWire Suture, 38 in (blue) w/ diamond-point needle, straight	AR-7246
#2 FiberWire Suture, 38 in (blue) w/ two diamond-point needles, straight	AR-7246-02
#2 FiberWire Suture, 38 in (blue)	AR-7233
#2 TigerWire Suture, 38 in (white/black)	AR-7203
#2 FiberWire Suture, 38 in (1 blue, 1 white/black)	AR-7201
#2 FiberWire Suture, 38 in (1 white, 1 blue/black)	AR-7240

Product Description	Item Number
#1 FiberWire Suture, 38 in (blue)	AR-7216
0 FiberWire Suture, 38 in (blue)	AR-7254
0 FiberWire Suture, 38 in (white)	AR-7255
0 FiberWire Suture, 38 in (blue) w/ tapered needle, 22.2 mm, ½ circle	AR-7250
0 FiberWire Suture, 38 in (blue) w/ diamond point needle, 22.2 mm, ½ circle	AR-7251
FiberWire suture Kit	AR-7219

Includes four #5 FiberWire (blue), four #5 FiberWire (white), four #5 FiberWire (white/black), and six #2 FiberWire (blue) sutures w/ tapered needle and 3 free needles

Collagen-Coated FiberWire Suture


Collagen-coated FiberWire suture is coated with type I bovine collagen, offering a softer coating and improved handling characteristics. These sutures are ideal for knotted and knotless repairs.

Product Description	Item Number
#2 FiberWire Suture, 38 in w/ tapered needle, 26.5 mm, ½ circle, collagen-coated	AR-7200B
#2 FiberLoop® Suture w/ straight needle, 20 in (blue), 76 mm needle w/ 7 mm loop, collagen-coated	AR-7234B
FiberTape® Suture, 2 mm, 7 in (blue) tape w/ each end tapered to #2 FiberWire suture, 30 in, collagen-coated	AR-7237-7B

References

1. Wüst DM, Meyer DC, Favre P, Gerber C. Mechanical and handling properties of braided polyblend polyethylene sutures in comparison to braided polyester and monofilament polydioxanone sutures. *Arthroscopy*. 2006;22(11):1146-1153. doi:10.1016/j.arthro.2006.06.013.
2. Lo IK, Burkhart SS, Athanasiou K. Abrasion resistance of two types of nonabsorbable braided suture. *Arthroscopy*. 2004;20(4):407-413. doi:10.1016/j.arthro.2004.01.022.
3. Lo IK, Burkhart SS, Chan KC, Athanasiou K. Arthroscopic knots: determining the optimal balance of loop security and knot security. *Arthroscopy*. 2004;20(5):489-502. doi:10.1016/j.arthro.2004.03.005.

FiberTape® Suture



Ultra-high-strength 2 mm and 1.7 mm FiberTape sutures have a similar long-chain polyethylene structure as FiberWire® suture. In addition to high-demand applications, like AC joint reconstruction, the broad footprint of FiberTape suture is ideal for repairs where tissue pull-through may be a concern.

Product Description	Item Number
FiberTape Suture, 2 mm, 36 in (blue), each end tapered to #2 FiberWire suture	AR-7237
FiberTape Suture, 2 mm, 7 in (blue), each end tapered to #2 FiberWire suture (d)	AR-7237-7
FiberTape Suture, 2 mm, 7 in (blue), each end tapered to #2 FiberWire suture, collagen-coated	AR-7237-7B
FiberTape Suture, 2 mm, 17 in (blue), each end tapered to #2 FiberWire suture w/ tapered needle 26.5 mm, ½ circle	AR-7237-17N
TigerTape™ Suture, 2 mm, 7 in (white/black), each end tapered to #2 FiberWire suture (b)	AR-7237-7T
FiberTape Suture, 1.7 mm, 36 in (blue), each end tapered to 0 FiberWire suture	AR-7273
TigerTape Suture, 1.7 mm, 36 in (white/black), each end tapered to 0 FiberWire suture	AR-7273T
FiberTape Suture, 1.7 mm, 8 in working length, 30 in overall (blue), each end tapered to 0 FiberWire suture (c)	AR-7273-8
TigerTape Suture, 1.7 mm, 8 in working length, 30 in overall (white/black), each end tapered to 0 FiberWire suture (a)	AR-7273-8T
FiberTape Suture, 1.7 mm, 17 in working length, 34 in overall (blue), each end tapered to 0 FiberWire w/ tapered needle 36.6 mm, ½ circle	AR-7273-17N
TigerTape Suture, 1.7 mm, 17 in working length, 34 in overall (white/black), each end tapered to 0 FiberWire w/ tapered needle 36.6 mm, ½ circle	AR-7273-17TN
FiberTape Loop, 2 mm, 8 in (blue), tapered to #2 FiberWire suture	AR-7275
TigerTape Loop, 2 mm, 8 in (white/black), tapered to #2 TigerWire® suture	AR-7275T

LabralTape™ Suture




This smooth, low-profile 1.5 mm specialty tape made from polyethylene, provides a wider footprint, and is 37% more resistant to tissue pull-through compared to a #2 suture.¹ LabralTape suture is ideal for the restoration of labral height in instability repairs.

Product Description	Item Number
LabralTape Suture, 1.5 mm, 36 in (white)	AR-7276
LabralTape Suture, 1.5 mm, 36 in (white/black)	AR-7276T

Reference

1. Arthrex, Inc. Data on file (LA1-0209-EN_A). Naples, FL; 2012.

FiberLink™ and TigerLink™ Suture



FiberLink suture is constructed from blue #2 suture and transitions from a single strand to an extended loop. This design is intended to allow 2 single strands of suture within a “cinch” stitch to gain purchase around tissue, while only managing one strand through the PushLock® anchor.

Product Description	Item Number
FiberLink Suture, #2 FiberWire® suture w/ loop (blue) (a)	AR-7235
TigerLink Suture, #2 TigerWire® suture w/ loop (white/black)	AR-7235T
FiberLink SutureTape, 1.3 mm, w/ loop (white/blue) (b)	AR-7535
TigerLink SutureTape, 1.3 mm, w/ loop (white/black)	AR-7535T
FiberLink SutureTape, 1.7 mm w/ loop (white/blue)	AR-7538
TigerLink SutureTape, 1.7 mm w/ loop (white/black)	AR-7538T


FiberLoop® Suture



FiberLoop suture is a continuous loop of FiberWire® suture with multiple needle options. The needles are easy to handle and move freely on the suture to recenter themselves after passing through tissue, facilitating even tension.

Product Description	Item Number
#2 FiberLoop® Suture, 20 in (blue), 76 mm straight needle w/ 7 mm loop (a)	AR-7234
#2 FiberLoop Suture, 20 in (blue), 76 mm straight needle w/ 7 mm loop, collagen-coated	AR-7234B
#2 FiberLoop Suture, 20 in (white), 76 mm straight needle w/ 7 mm loop	AR-7234-01
#2 TigerLoop™ Suture, 20 in (white/green), 76 mm straight needle w/ 7 mm loop	AR-7234T
#2 FiberLoop Suture, 20 in (blue) w/ tapered needle, 26.5 mm, ½ circle	AR-7234C
FiberLoop SutureTape, 1.3 mm, 20 in loop (white/blue) w/ 76 mm straight needle	AR-7534
TigerLoop SutureTape, 1.3 mm, 20 in loop (white/black) w/ 76 mm straight needle	AR-7534T
#2 FiberLoop Suture w/ FiberTag® suture (blue), w/ 20 in FiberTape suture w/ 65 mm straight needle	AR-7266
#2 FiberLoop Suture, 20 in (blue), 65 mm straight needle w/ 7 mm loop	AR-7284
0 FiberLoop Suture, 20 in (blue), 76 mm straight needle w/ 7 mm loop	AR-7253
0 TigerLoop Suture, 20 in (white/green), 76 mm straight needle w/ 7 mm loop	AR-7253T


FiberStick™ and TigerStick® Suture



The 12-in stiff “waxed” section of FiberStick suture allows convenient and easy advancement through most cannulated instruments, alleviating the need for a monofilament suture or wire suture shuttle. FiberStick sutures come with a thin plastic tube that helps protect the stiffened suture until use.

Product Description	Item Number
FiberStick Suture, #2 FiberWire® suture, 50 in (blue), one end stiffened, 12 in	AR-7209
TigerStick Suture, #2 TigerWire® suture, 50 in (white/black), one end stiffened, 12 in	AR-7209T
FiberStick SutureTape , 1.3 mm, (blue/blue), one end stiffened 12"	AR-7509
FiberStick SutureTape , 1.3 mm (black/white), one end stiffened 12"	AR-7509T
2-0 FiberStick Suture, 2-0 FiberWire suture, 50 in (blue), one end stiffened, 12 in	AR-7222

FiberSnare® Suture



FiberSnare suture provides an easy one-step approach to creating a FiberWire® loop on the tip of the tenodesis driver. FiberSnare suture can also be used as a suture shuttle for passage of traction sutures through bone tunnels.

Product Description	Item Number
#2 FiberSnare Suture, #2 FiberWire Suture, 26 in (green) w/ ¾ in closed loop, 12 in stiffened	AR-7209SN

FiberButton



FiberButton is a all-polyester embroidered button designed for superior strength and versatility in soft tissue repairs. It comes in two sizes to accommodate a range of surgical needs. Ideal for tissue approximation, FiberButton functions as a washer to distribute load effectively in both knotted and knotless constructs, reducing stress concentration and enhancing repair integrity. Its durable, low-profile design ensures reliable performance while maintaining simplicity and ease of use.

Product Description	Item Number
FiberButton, oval, 10 mm × 6 mm	AR-7285
FiberButton, oval, 14 mm × 8 mm	AR-7286

FingerShield™ Finger Guard



Suture slides easily over the finger guard, a woven polyester sleeve with an embedded radiopaque blue marker. The finger guard will stand up to repetitive hand tying during a case without constraining the fingers.

Product Description	Item Number
FingerShield Finger Guard, 2/pk	AR-7199

Free Needles

Product Description	Item Number
Reverse Cutting Needle w/ Nitinol loop, C-13, 36.6 mm, ½ circle	AR-7280
Tapered Needle w/ Nitinol loop, T-5, 26.5 mm, ½ circle	AR-7281
Spring Eye Needle, 60 mm	AR-7283

FiberWire® Suture Cart




A great way to store and organize sutures in one convenient location, the cart also offers additional storage for miscellaneous items. Wheels make it portable and easy to move around the facility.

Product Description	Item Number
FiberWire Suture Cart, cart only	AR-7200CART

Cannulas

PassPort Button™ Cannula	36
Gemini™ SR8 Cannula	37
Expanula™ Cannula	37
Shoehorn Cannula	37
Triple-Dam Twist-In™ Cannulas	38
Twist-In Cannulas	38
Crystal Cannula® Cannulas	39
Low-Profile 5 mm Cannula	39
Low-Profile Percutaneous Cannula Kit.....	39

PassPort Button™ Cannula



PassPort Button cannulas help maximize visibility and maneuverability inside and outside of the arthroscopic workspace. The double-dam, one-piece, molded design has low-profile flanges that seat flush to the skin and soft tissue. The flanges are designed to create a stable portal that allows the insertion and removal of instruments while preventing cannula loss. The PassPort Button cannula with auxiliary lumen adds a convenient side channel for use with NanoScope™ instruments, or suture management. The Passport Divider also aids in suture management by creating divided channels in the 12 mm cannula.


They are easily introduced using the PassPort inserter or by grasping the inside flange with a curved hemostat and inserting into the incision. PassPort Button cannulas can be used for the shoulder, knee, hip, and elbow. Each cannula is supplied with a 5 mm spacer to ensure accurate length.



Product Description	Item Number
6 mm ID × 2 cm	AR-6592-06-20
6 mm ID × 3 cm	AR-6592-06-30
6 mm ID × 4 cm	AR-6592-06-40
6 mm ID × 5 cm	AR-6592-06-50
8 mm ID × 2 cm	AR-6592-08-20
8 mm ID × 3 cm	AR-6592-08-30
8 mm ID × 3 cm, w/ auxiliary lumen, 1/box (c)	AR-6592-08-30A
8 mm ID × 4 cm	AR-6592-08-40
8 mm ID × 4 cm, w/ auxiliary lumen, 1/box	AR-6592-08-40A
8 mm ID × 5 cm	AR-6592-08-50
8 mm ID × 5 cm, w/ auxiliary lumen, 1/box	AR-6592-08-50A
8 mm ID × 6 cm	AR-6592-08-60
10 mm ID × 2 cm	AR-6592-10-20
10 mm ID × 3 cm	AR-6592-10-30
10 mm ID × 4 cm	AR-6592-10-40
10 mm ID × 5 cm	AR-6592-10-50
12 mm ID × 3 cm, 1/box	AR-6592-12-30
12 mm ID × 4 cm, 1/box	AR-6592-12-40
12 mm ID × 5 cm, 1/box	AR-6592-12-50
Divider, 12 mm, 1/box (a)	AR-6592-12D
PassPort Hemostat	AR-6592
PassPort Selection Guide	AR-6592M
8 mm PassPort Inserter	AR-6592-08PI
10 mm PassPort Inserter (b)	AR-6592-10PI
12 mm PassPort Inserter	AR-6592-12PI



Gemini™* SR8 Cannula



Gemini SR8 cannulas incorporate a deployable low-profile wing to help prevent cannula “fallout” while inserting and removing instruments during Bankart, SLAP, and rotator cuff repairs. Additionally, the cannula’s inner sleeve can telescope relative to the outer portion of the cannula, providing surgeons a clear, unobstructed channel to the pathology. A disposable cannulated obturator is included for easy insertion over a 2.6 mm switching stick.

Product Description	Item Number
Gemini Cannula, 8.25 mm × 9 cm	AR-6572
Switching Stick, 2.6 mm × 305 mm	AR-6572S

*Gemini is a trademark of Frantz Medical Development Ltd.

Expanula™ Cannula



Arthroscopic rotator cuff repair is easier with the Expanula cannula, which expands the workspace and view. Upon insertion, rotate its outer sheath to expand the distal end of the cannula beneath the deltoid. This design creates an extremely stable portal for inserting and removing instruments while preventing cannula loss. Each cannula is supplied with a “no-squirt” cap and disposable obturator.

Product Description	Item Number
Expanula Cannula, 8.25 mm ID × 7.5 cm	AR-6569


Shoehorn Cannula



The Shoehorn cannula has a longitudinal slot on top of the cannula for introducing oversized instruments. Each cannula is supplied with a disposable obturator.

Product Description	Item Number
Shoehorn Cannula, 6 mm ID × 9 cm	AR-6565

Triple-Dam Twist-In™ Cannulas



This cannula's unique, 3-dam design allows instrumentation and implants to be passed arthroscopically. Additionally, it maintains pressurization for a clear visual fluid environment while minimizing squirting. The side port offers a luer lock connection and a valve switch for manual adjustment of fluid outflow and inflow.


Product Description	Item Number
Triple-Dam Twist-In Cannula, 8.25 mm ID × 7 cm	AR-6530TD
Triple-Dam Twist-In Cannula, 7 mm ID × 7 cm (a)	AR-6570TD
Triple-Dam Twist-In Cannula, 7 mm ID × 7 cm, no valve (b)	AR-6570TDNV
Triple-Dam Instrument Cannula, 7 mm ID × 7 cm (c)	AR-6550TD
Triple-Dam Twist-In Cannula, 6 mm ID × 7 cm	AR-6535TD

Cannula Accessories

Product Description	Item Number
Wissinger Rod, 4 mm × 12.7 in	AR-3025
Extra Long Switching Stick, 4 mm × 12 in	AR-3026
Switching Stick Inserter	AR-3027
Portal Dilatation Set (4 mm, 6 mm, and 8 mm dilators)	AR-6520
1.5 mm Pins for AR-6520, qty. 5	AR-6521
Quick Connect Handle, non-ratcheting, small	AR-2005NR

Product Description	Standard	Quick Connect
Reusable Obturator for AR-6530, AR-6530N, AR-6566 Twist-In Cannula	AR-6531	AR-6531Q
Reusable Obturator for AR-6535 Twist-In Cannula	AR-6536	AR-6536Q
Reusable Obturator for AR-6540, AR-6575-09 Twist-In Cannula	AR-6541	AR-6541Q
Reusable Obturator for AR-6545 Twist-In Cannula	AR-6546	AR-6546Q
Reusable Obturator for AR-6550, AR-6567, AR-6570, AR-6570F Twist-In Cannula	AR-6549	AR-6549Q
Reusable Obturator for AR-6560, AR-6562, AR-6564, AR-6560F Instrument Cannula	AR-6563	AR-6563Q

Twist-In™ Cannulas



The translucent Twist-In cannula allows direct arthroscopic visualization of instruments and suture passing through the cannula. A flexible option easily conforms to large or curved instruments. Each cannula is supplied with a “no-squirt” cap and disposable obturator.

Product Description	Item Number
Twist-In Cannula, 8.25 mm ID × 7 cm, notched	AR-6530N
Twist-In Cannula, 8.25 mm ID × 7 cm	AR-6530
Twist-In Cannula, 8.25 mm ID × 9 cm	AR-6540
Twist-In Cannula, 7 mm ID × 7 cm (a)	AR-6570
Twist-In Cannula, 7 mm ID × 7 cm, flexible	AR-6570F
Twist-In Cannula, 6 mm ID × 7 cm	AR-6535
Twist-In Cannula, 6 mm ID × 9 cm	AR-6545
Partially Threaded Cannula, 8.25 mm ID × 7 cm	AR-6566
Partially Threaded Cannula, 8.25 mm ID × 9 cm	AR-6575-09
Partially Threaded Cannula, 8.25 mm ID × 11 cm	AR-6575-11
Partially Threaded Cannula, 7 mm ID × 7 cm (b)	AR-6567
Instrument Cannula, 5.5 mm ID × 9 cm	AR-6532
Instrument Cannula, 7 mm ID × 7 cm (c)	AR-6550



Crystal Cannula® Cannulas



The Crystal Cannula® cannula offers a unique barrel-shaped retention bowl that pools fluid to eliminate cannula “squirt.” Three atraumatic distal retaining options help to prevent cannula “fallout.” A flexible option easily conforms to large or curved instruments. Each cannula is supplied with a disposable orange obturator.

Product Description	Item Number
5.75 mm ID × 7 cm, distal ring	AR-6560
5.75 mm ID × 7 cm, smooth	AR-6562
5.75 mm ID × 7 cm, partially threaded	AR-6564
5.75 mm ID × 7 cm, flexible, distal ring	AR-6560F
5.75 mm ID × 5 cm, partially threaded, short	AR-6568


Low-Profile 5 mm Cannula



This 5 mm cannula is the lowest profile cannula on the market that allows for direct visualization of instruments and suture passing. The proximal portion of the cannula is similar to the Crystal Cannula cannula. Each cannula is supplied with a disposable blue obturator.

Product Description	Item Number
5 mm ID × 7 cm	AR-6548

Low-Profile Percutaneous Cannula Kit




This 5 mm low-profile percutaneous cannula instrument kit provides a solution for precise anchor and suture placement and is universally compatible with multiple anchors, including the Knotless FiberTak® anchor.

Product Description	Item Number
Low-Profile Percutaneous Instrument Kit Includes 5 mm x 7 cm low profile cannula, obturator, needle w/ stylet, and K-wire	AR-6548PK

Instruments

SutureLasso™ Suture Passers	42
Rotation Lasso Suture Passer	43
Transosseous Suture Passer	43
SutureLasso Suture Passer	43
Micro SutureLasso Suture Passer	43
FastPass Scorpion™ SL Suture Passer	44
SCORPION-multifire Suture Passer	44
Scorpion Suture Passer	45
Labral Scorpion Suture Passer	45
BirdBeak® Suture Passer	46
Rhino Suture Passers	46
Penetrator™ Suture Retriever	47
Suture Retriever	47
FiberTape® Retriever	48
KingFisher® Retriever/Grasper	48
Rotator Cuff Grasper	49
Back Grasper	49
FiberWire® Grasper	49
SCR Guide	50
Suture Cutters	50
CrabClaw™ Knot Pusher/Suture Retriever	51
Knot Pusher	51
6th Finger™ Knot Pusher	51
Crochet Hook	52
Curettes	52
Probes	52
Tissue Elevator and Bone Preparation System	53
Shoulder Arthroscopy Instrument System	54
Shoulder Express Instrument Systems	55


SutureLasso™ Suture Passers




The SutureLasso SD suture passers feature a thumb pad and the QuickPass™ SutureLasso suture passers use thumbwheels to quickly and easily advance the supplied nitinol wire loop, a #2 FiberStick™ suture, or a monofilament (PDS) suture. The small, 1.8 mm-diameter tip provides atraumatic suture passage and is sterile for single-patient use, which ensures a sharp instrument every time.

ReelPass SutureLasso™ suture passers contain 10 yd of blue #1 nonabsorbable monofilament to provide a continuous supply without the need to reload. The ergonomic handle and thumbwheels were designed to easily advance the monofilament. This single-use device has a 1.5 mm-diameter tip that is the smallest on the market and is always sharp.

		QuickPass™	SD	ReelPass
Product Description		Item Number	Item Number	Item Number
	30° Straight	AR-6068-30	AR-4068-30	N/A
	90° Straight	AR-6068-90	AR-4068-90	AR-6069-90 (e)
	25° Tight Curve, left	AR-6068-25TL (a)	AR-4068-25TL	N/A
	25° Tight Curve, right	AR-6068-25TR (b)	AR-4068-25TR	N/A
	45° Curve, left	AR-6068-45L	AR-4068-45L (c)	AR-6069-45L
	45° Curve, right	AR-6068-45R	AR-4068-45R (d)	AR-6069-45R
	90° Curve, left	AR-6068-90L	AR-4068-90L	N/A
	90° Curve, right	AR-6068-90R	AR-4068-90R	N/A
	90° Tight Curve	AR-6068-90T	N/A	N/A
	Crescent	AR-6068C	AR-4068C	N/A
	25° Tight Curve, w/ FiberStick suture, left	N/A	AR-4068-25TLF	N/A
	25° Tight Curve, w/ FiberStick suture, right	N/A	AR-4068-25TRF	N/A
	Wire Loop, 5/box	AR-6068-01	AR-4068-05SD	N/A




Rotation Lasso Suture Passer



The Rotation Lasso suture passer features a rotating and steerable loop for arthroscopic and open soft-tissue repair. The simplicity of the design helps reduce procedural steps, prevents twists and tangles inside cannulas, and eliminates the need for outside suture shuttling.

Product Description	Item Number
Rotation Lasso Suture Passer, 90°	AR-6065-90

Transosseous Suture Passer




The transosseous suture passer comes sterile and features a 6-inch malleable shaft, a durable nitinol capture loop, and a comfortable handle. Applications for this instrument include any passing of grafts or suture through bony sockets and tunnels.

Product Description	Item Number
Transosseous Suture Passer	AR-1000

Literature

Product Description	Item Number
Lasso Suture Passers Brochure	LB1-000382-en-US


SutureLasso™ Suture Passer



These passers have various curved-tip configurations for arthroscopic Bankart, SLAP, and rotator cuff repairs. Each suture passer has a 2.3 mm outer diameter tip and comes preloaded with a nitinol loop to simplify shuttling and passing suture through tissue. The Banana SutureLasso™ suture passer is designed for passing sutures through the rotator cuff via a superior, percutaneous approach (modified Neviaser portal) or along the acromial border.

Product Description	Item Number
Banana SutureLasso™ Suture Passer (a)	AR-4065B
45° w/ wire loop	AR-4065W
90° w/ wire loop	AR-4065-90W
45° curve, left (b)	AR-4065-45L
45° curve, right (c)	AR-4065-45R
45° w/ #2 monofilament	AR-4065S
90° w/ #2 monofilament	AR-4065-90S


Micro SutureLasso™ Suture Passer



These 1.25 mm-diameter instruments work well for percutaneous suture passing during rotator cuff and glenoid labrum repairs.

Product Description	Item Number
Minor bend (b)	AR-8701
Major bend (a)	AR-8702
Straight (c)	AR-8703
Micro SutureLasso Retriever (d)	AR-8701SR

FastPass Scorpion™ SL Suture Passer




The FastPass Scorpion SL suture passer simplifies antegrade suture passing by loading the suture from the side, preventing suture and soft tissue interference during suture passing. The low-profile design fits down a 6 mm cannula and the top jaw mechanism helps pass and retrieve FiberWire® sutures in one step.

Product Description	Item Number
FastPass Scorpion SL Suture Passer	AR-13999MF
HD Scorpion Needle w/ MegaLoader, 5/box	AR-13999HDN




SCORPION-multifire Suture Passer



With the FastPass SCORPION-multifire suture passer, surgeons can load 2 sutures outside and independently pass 2 sutures inside. The low-profile design of the straight suture passer fits down 6 mm cannulas. The FastPass upper jaw helps pass and retrieve FiberWire® suture in one step.

Product Description	Item Number
FastPass SCORPION-multifire Suture Passer (a)	AR-13997MF
SCORPION-multifire Suture Passer, straight (b)	AR-13996
HD Scorpion Needle w/ MegaLoader, 5/box	AR-13999HDN



Heavy Duty Scorpion Needle w/ MegaLoader


Scorpion™ Suture Passer



The FastPass Scorpion suture passer grasps cuff tissue then directly passes and retrieves FiberWire® suture. Scorpion suture passers grasp 16 mm of tissue and fit through a 6 mm cannula.


Product Description	Item Number
Scorpion Suture Passer	AR-13990
SureFire® Scorpion Needle, 5/box	AR-13991N

Labral Scorpion Suture Passer




Pass and retrieve #2 FiberWire® suture in the labrum using the Labral Scorpion QL suture passer. Its quickloading bottom jaw and suture-retrieving top jaw are designed for simplicity and easy use, supporting more successful suture passes. The Labral Scorpion suture passer's streamlined, low-profile jaws are designed to maneuver around the glenohumeral space and fit through a 6 mm cannula. Use the Labral Scorpion suture passer through a single portal to self-retrieve sutures with the upper jaw's capturing mechanism.

Product Description	Item Number
Labral Scorpion QL Suture Passer	AR-13998QL
HD Scorpion Needle w/ MegaLoader, 5/box	AR-13999HDN



BirdBeak® Suture Passer




BirdBeak suture passers have an extremely sharp tip to penetrate soft tissue easily and a stiff shaft that resists bending during tissue-shifting procedures. BirdBeak suture retriever devices are an essential tool for arthroscopic labral, SLAP, and rotator cuff repairs. The BirdBeak Evolution suture passer has a uniquely designed handle that allows for easy operation from virtually any hand position.

Product Description	Item Number
BirdBeak Suture Passer, 45° up tip (d)	AR-11800
BirdBeak Suture Passer, 22° up tip (c)	AR-11890
BirdBeak Suture Passer, straight (b)	AR-11880
BirdBeak Suture Passer, straight, w/ 45° handle, left	AR-11887
BirdBeak Suture Passer, straight, w/ 45° handle, right	AR-11886
BirdBeak Evolution Suture Passer, 45° up tip	AR-11800E
BirdBeak Evolution Suture Passer, 22° up tip	AR-11890E
BirdBeak Evolution Suture Passer, straight	AR-11880E
BirdBeak Evolution Suture Passer, 15° up curve (a)	AR-11881E
Grasping BirdBeak Suture Passer, 22° up curve w/ SR handle (f)	AR-11810SR

Inverted BirdBeaks

Product Description	Item Number
BirdBeak Suture Passer, 45° up tip, inverted jaw w/ WishBone™ handle (e)	AR-11805W
BirdBeak Suture Passer, 22° up tip, inverted jaw w/ WishBone™ handle	AR-11895W


Rhino Suture Passers



Rhino suture passers have an extremely sharp, small-diameter upturned tip that will easily penetrate soft tissue. The ridged 3.4 mm shaft resists bending and flexing during tissue-shifting procedures. The novel bottom opening jaw design is conveniently positioned to capture the suture from an anchor, eliminating the need to rotate the instrument.

Product Description	Item Number
Straight (a)	AR-11850SR
Left curve (c)	AR-11852SR
Right curve (b)	AR-11851SR
Straight, w/ WishBone handle	AR-11850W


Penetrator™ Suture Retriever



This unique instrument combines a small penetrating tip with a suture grasper to allow suture delivery or extraction in one step. The 2.7 mm-diameter tip slides easily through the tissue with the suture either sliding or grasped within the self-ratcheting mechanism. This device is ideal for instability and rotator cuff repairs.

Product Description	Item Number
Penetrator Suture Retriever, 15° up curved	AR-2167-2
Penetrator Suture Retriever, straight (a)	AR-2167ST-2
Penetrator FiberTape® Suture Retriever, 15° up curved	AR-2167-3
Penetrator FiberTape Suture Retriever, straight (b)	AR-2167ST-3
Penetrator FiberTape Suture Retriever, 15° up curved w/ WishBone™ handle	AR-2167W-3
Penetrator FiberTape Suture Retriever, straight w/ WishBone handle	AR-2167STW-3


Suture Retriever



The suture retriever is designed for atraumatic suture retrieval and manipulation during arthroscopic procedures. The jaws create a closed loop that allows suture to slide freely during suture extraction. The uniquely designed tip is made to spread parallel strands of suture and facilitate retrieval. Its small diameter and low profile allow access into the tightest joint spaces. The tip angle on the 45° suture retriever allows access to suture strands in deep, hard-to-reach spaces.

Product Description	Item Number
Suture Retriever, 3.4 mm, straight (a)	AR-12540
Suture Retriever, 3.4 mm, 15° up	AR-12550
Suture Retriever, 3.4 mm, 45° right (b)	AR-12580
Suture Retriever, 3.4 mm, 45° left	AR-12590
Suture Retriever w/ WishBone™ handle, 3.4 mm, straight	AR-12540W
Suture Retriever w/ WishBone handle, 3.4 mm, 15° up	AR-12550W

FiberTape® Retriever




The FiberTape retriever has a wide jaw specifically designed for easy FiberTape suture management during a SpeedBridge™ or SpeedFix™ repairs. The serrated tip also allows the grasper to be used to securely grab sutures.

Product Description	Item Number
FiberTape Retriever w/ SR handle (a)	AR-13974SR
FiberTape Retriever w/ NR handle	AR-13974NR
FiberTape Retriever w/ WishBone™ handle	AR-13974W

Product Description	Item Number
Mini FiberTape Retriever w/ SR handle (b)	AR-12974SR
Mini FiberTape Retriever w/ NR handle	AR-12974NR
Mini FiberTape Retriever w/ WishBone handle	AR-12974W

The Mini FiberTape retriever is designed for use with LabralTape™ suture.

KingFisher® Retriever/Grasper



Perform multiple tasks and improve procedural speed and efficiency with KingFisher suture retriever and tissue grasper. The tip design is optimal for arthroscopic tissue grasping/reduction, foreign body removal, and suture retrieval/management.

Product Description	Item Number
KingFisher Retriever/Grasper, 4.2 mm w/ SR handle (a)	AR-13970SR
KingFisher Retriever/Grasper, 4.2 mm w/ NR handle	AR-13970NR
KingFisher Retriever/Grasper, 4.2 mm w/ WishBone™ handle	AR-13970W
KingFisher FiberTape Retriever/Grasper, 4.2 mm w/ SR handle (b)	AR-13971SR
KingFisher FiberTape Retriever/Grasper, 4.2 mm w/ NR handle	AR-13971NR
KingFisher FiberTape Retriever/Grasper, 4.2 mm w/ WishBone handle	AR-13971W
Mini KingFisher Retriever/Grasper, 3.4 mm w/ SR handle (c)	AR-13972SR
Mini KingFisher Retriever/Grasper, 3.4 mm w/ NR handle	AR-13972NR
Mini KingFisher Retriever/Grasper, 3.4 mm w/ WishBone handle	AR-13972W



Self-releasing/
ratcheting handle



Non-ratcheting
handle



WishBone
handle

Rotator Cuff Grasper



The low-profile 4 mm tip was specifically designed to securely grasp the tendon and pull it into the proper anatomic position. This instrument also features an ergonomic ring handle with a tip-locking mechanism.

Product Description	Item Number
Rotator Cuff Grasper w/ SR handle	AR-13960SR
Rotator Cuff Grasper w/ NR handle	AR-13960NR
Rotator Cuff Grasper w/ WishBone™ handle	AR-13960W

Back Grasper



The back grasper is a simple but effective tool to aide introducing grafts, implants, and instruments into the body.

Product Description	Item Number
Back Grasper w/ NR handle	AR-12531NR
Back Grasper w/ SR handle	AR-12531SR
Back Grasper w/ WishBone™ handle	AR-12531W

FiberWire® Grasper



The FiberWire grasper is designed to grasp and retrieve multiple strands of FiberWire suture. This fully toothed design and long jaw assembly provides a substantial working range for easily grabbing suture within the glenohumeral or subacromial joint space.

Product Description	Item Number
FiberWire Grasper w/ SR handle	AR-13975SR
FiberWire Grasper w/ NR handle	AR-13975NR
FiberWire Grasper w/ WishBone™ handle	AR-13975W

SCR Guide



The SCR Guide helps simply and reproducibly estimate the distances between suture anchors during superior capsular reconstruction (SCR). This reusable hand instrument's suture-retrieving jaw is designed for both suture and FiberTape® suture. The highlighted black sleeve clearly shows graduations for distances referenced off of sutures.

Product Description	Item Number
SCR Guide	AR-16950SR

Suture Cutters



The suture cutter was designed to facilitate arthroscopic cutting of FiberWire® and braided sutures. The uniquely designed cutting jaws remain sharp throughout repeated use. The suture cutter is available in a closed and open end, left notch version. The closed-end suture cutter allows the surgeon to leave a 3 mm suture tail without the possibility of cutting the knot. The open end, left notch version facilitates suture cutting inside the joint without having to top load the cutter. The FiberTape cutter was designed to allow for flush cutting of 2 FiberTape sutures.

Product Description	Item Number
Suture Cutter, 4.2 mm, straight, closed-end (used w/ #2 and #5 suture) (a)	AR-12250
Suture Cutter, 4.2 mm, straight, closed-end w/ WishBone™ handle	AR-12250W
Suture Cutter, 4.2 mm, open end, left notch (used w/ all suture) (b)	AR-11794L
Suture Cutter, 4.2 mm, w/ WishBone handle, open end, left notch (used w/ all suture)	AR-11794LW
Mini Suture Cutter, 3.4 mm, straight, closed-end (used w/ #2 FiberWire and LabralTape sutures)	AR-13255
Mini Suture Cutter, 3.4 mm, straight, closed-end w/ WishBone handle (used w/ #2 FiberWire and LabralTape™ sutures)	AR-13255W
FiberTape Cutter, 4.2 mm, straight, closed-end	AR-13250
FiberTape Cutter, 4.2 mm, straight, closed-end w/ WishBone handle	AR-13250W
FiberWire Scissor (for open procedures) (c)	AR-11796

CrabClaw™ Knot Pusher/Suture Retriever



The opening rotary jaw of the CrabClaw instrument makes it an effective suture retriever and knot pusher. It allows for intra-articular capture of suture, which dramatically simplifies arthroscopic knot tying. Suture posts are alternated simply by opening the jaws and engaging the parallel suture limb outside the cannula or in the joint. When used as a knot pusher, the ratcheting handle secures the jaw in a closed position to allow for smooth knot advancement.

Product Description	Item Number
CrabClaw Knot Pusher/Suture Retriever	AR-12960

6th Finger™ Knot Pusher



The unique double-tube design allows surgeons to apply and maintain tension to the first throw, while advancing subsequent throws with the sliding plastic outer tube. The inner tube allows subsequent “past pointing” to apply opposite suture tension to the knot similar to open knot tying techniques. A wire loop is incorporated inside the disposable, sterile 6th Finger knot pusher for easy suture loading.

Product Description	Item Number
6th Finger Knot Pusher w/ suture passer, sterile, 5/box	AR-1930S


Knot Pusher



The single-hole knot pusher provides a simple method for advancing sliding knots and half-hitches. The closed-end knot pusher has a modified handle that provides an ergonomic feel. The distal tip has also been modified for easier advancement of slip knots and half-hitches.

Product Description	Item Number
Knot Pusher, closed end (a)	AR-1305
Single-Hole Knot Pusher	AR-1299


Crochet Hook



The crochet hook is a simple tool that performs well in tight spaces to retrieve suture loops during arthroscopic Bankart, SLAP, and rotator cuff repairs, or any suturing procedure. The smooth tip prevents abrasion of suture strands and the ergonomic, knurled handle facilitates instrument manipulation in the wet arthroscopic environment. The push/pull crochet hook was designed to push suture knots and/or retrieve suture with the same instrument.

Product Description	Item Number
Crochet Hook (a)	AR-5008H
Push/Pull Crochet Hook (b)	AR-5009H
Suture Hook (c)	AR-5007H

Curettes



Ring curettes are available with one or both sides cut. D-curettes feature a D-shaped, dual-sided cutting ring designed to prepare the articular margin prior to anchor insertion during rotator cuff repair. D-curettes are also an excellent choice for cartilage removal, particularly prior to glenoid placement in total shoulder arthroplasty procedures.

Product Description	Item Number
D-curette, 5.4 mm × 150 mm, both sides cut (a)	AR-22020
D-curette, 3.4 mm × 150 mm, both sides cut	AR-23020
Curette, 5.4 mm ring, one side cut (b)	AR-20010
Curette, 5.4 mm ring, both sides cut	AR-20020

Probes




With the articulating probe, surgeons can easily insert an instrument through a small stab incision in soft tissue with the tip in a straight position. The tip of the probe can then be changed to a locked 90° position to test the strength of a soft tissue repair or to probe defects.

Product Description	Item Number
Articulating Probe (a)	AR-10100
Hook Probe, 5.4 mm w/ 5 mm markings (b)	AR-10000
Hook Probe, 3.4 mm (c)	AR-10010


Tissue Elevator and Bone Preparation System

The Tissue Elevator and Bone Preparation System combines new tip designs and a modular handle to give you the ultimate selection of tools for work in the glenohumeral joint and subacromial space. A variety of rasps, curettes, tissue elevators, and a chondro pick are provided in one convenient system. The shaft length of each tip is long enough to effectively span the distance from an anterior working portal to the posterior rim of the glenoid, allowing the best possible angle of approach to work with ease and efficiency in this area.

Product Description	Item Number
Tissue Elevator Instrument Set	AR-1344S
Quick Connect Handle, ¼ sq. nonratcheting	AR-2003NR
Tissue Elevator Instrument Case	AR-1344C
Attachment Options	
Tissue Elevator, 15° down (a)	AR-1344-15
Tissue Elevator, 15° down, 30° left (b)	AR-1344-1530L
Tissue Elevator, 15° down, 30° right (c)	AR-1344-1530R
Tissue Elevator, 30° down (d)	AR-1344-30
Debridement Rasp (e)	AR-1344-DR
D-curette, 5.4 mm (f)	AR-1344-LC
Ball-Tip Rasp, 40° (g)	AR-1344-MR
Chondro Pick, angled (h)	AR-1344-P
Elevator, standard design, 15° (i)	AR-1344SE-15
Elevator, standard design, 30° (j)	AR-1344SE-30
Curette, round, double-sided, 5.4 mm (k)	AR-1344-RC



Shoulder Arthroscopy Instrument System



The Shoulder Arthroscopy Instrument System is a modular and customizable instrument system that can accommodate a comprehensive selection of instruments to facilitate arthroscopic shoulder repairs. The configured set contains a selection of the most popular and latest instruments, including a quick-connect system for portal dilation, cannula insertion, and tissue preparation, as determined by leading upper extremity surgeons.

Product Description	Item Number
Shoulder Arthroscopy Instrument System includes	AR-8405S
Hook Probe, 5.4 mm tip w/ 5 mm markings	AR-10000
Suture Cutter, open-ended, left notch	AR-11794L
FiberWire® Scissor	AR-11796
Suture Cutter, 4.2 mm, straight, closed-end	AR-12250
Suture Retriever, 3.4 mm, straight	AR-12540
Knot Pusher, closed end	AR-1305
FiberTape® Cutter, 4.2 mm, straight, closed-end	AR-13250
Debridement Rasp	AR-1345-DR
Ring Curette	AR-1345-RC
Elevator, 15°	AR-1345SE-15
Elevator, 30°	AR-1345SE-30
KingFisher® FiberTape® Retriever/Grasper w/ SR handle	AR-13971SR
FiberTape Retriever w/ SR handle	AR-13974SR
FiberWire® Grasper w/ SR handle	AR-13975SR
FastPass Scorpion™ SL Suture Passer	AR-13999MF
Punch for 5.5 Corkscrew® FT and 4.75 and 5.5 SwiveLock anchor	AR-1927PB
Quick-Connect Handle, small, non-ratcheting, qty. 2	AR-2005NR
Kirk Mallet	AR-2966
Extra-Long Switching Stick, qty. 2	AR-3026

Product Description	Item Number
Quick-Connect Reusable Obturator for AR-6530, AR-6530N, and AR-6566 Twist-In™ Cannula	AR-6531Q
Quick-Connect Reusable Obturator for AR-6535 Twist-In Cannula	AR-6536Q
Quick-Connect Reusable Obturator for AR-6540 and AR-6575-09 Twist-In Cannula	AR-6541Q
Quick-Connect Reusable Obturator for AR-6550, AR-6567, AR-6570, and AR-6570F Twist-In Cannula	AR-6549Q
Quick-Connect Reusable Obturator for AR-6560, AR-6562, AR-6564, and AR-6560F Instrument Cannula	AR-6563Q
PassPort Cannula Inserter, 8 mm	AR-6592-08PI
PassPort Cannula Inserter, 10 mm	AR-6592-10PI
Shoulder Repair Set Instrumentation Case	AR-8405C

Shoulder Express Instrument Systems



The Rotator Cuff and Instability Express Instrument System cases feature single-level, low-profile, and customizable designs that house the latest instrumentation for Arthrex next-generation implants.

Rotator Cuff Express Instrument System

Product Description	Item Number
Rotator Cuff Repair Express Case (a)	AR-8405CRC

Recommended Rotator Cuff Repair Express Instruments (must be purchased separately)


Product Description	Item Number
FastPass Scorpio™ SL Suture Passer	AR-13999MF
FiberTape® Retriever w/ SR handle	AR-13974SR
KingFisher® FiberTape Retriever/Grasper w/ SR handle	AR-13971SR
FiberTape Cutter, 4.2 mm, straight, closed end	AR-13250
Suture Retriever, ø3.4 mm, straight	AR-12540
FiberWire® Scissor	AR-11796
Punch for 2.6 FiberTak anchors	AR-3656
Angled Circumferential Teeth Spear for 2.6 FiberTak anchors	AR-3655
Sharp Obturator for 2.6 FiberTak Spear (AR-3655)	AR-3658T
Punch for 4.5 mm PushLock® and 4.5 mm Corkscrew® FT anchors	AR-1922P
Punch for 5.5 mm Corkscrew FT, 4.75 mm and 5.5 mm SwiveLock® anchors	AR-1927PB
Circumferential Teeth Spear w/ sharp obturator	AR-1941CT
Fishmouth Spear w/ sharp obturator	AR-1941DGF

Instability Express Instrument System

Product Description	Item Number
Instability Repair Express Case (b)	AR-8405CIN

Recommended Instability Repair Express Instruments (must be purchased separately)

Product Description	Item Number
Extra Long Switching Stick, 4 mm x 12 in (2 required)	AR-3026
Kirk Mallet	AR-2966
Quick-Connect Handle, small, non-ratcheting (2 required)	AR-2005NR
Elevator, standard design, 15°	AR-1345SE-15
Elevator, standard design, 30°	AR-1345SE-30
Debridement Rasp	AR-1345-DR
Curette, Round, double-sided, 5.4 mm	AR-1345-RC
Quick-Connect Reusable Obturator (gold) for AR-6530, AR-6530N, and AR-6566 Twist-In Cannula	AR-6531Q
Quick-Connect Reusable Obturator (purple) for AR-6550, AR-6567, AR-6570, and AR-6570F Twist-In Cannula	AR-6549Q
Suture Retriever, ø3.4 mm, straight	AR-12540
Mini Suture Cutter, 3.4 mm, straight, closed end	AR-13255
Suture Cutter, 4.2 mm, open-end, left notch (used w/ all suture)	AR-11794L
Mini KingFisher Retriever/Grasper, 3.4 mm w/ NR handle	AR-13972NR
Circumferential Teeth Curved Spear w/ flexible sharp obturator for FiberTak and Knotless FiberTak anchors	AR-3610CTC
Slotted Spear for FiberTak and Knotless FiberTak anchors	AR-3610ST
FiberWire Scissor	AR-11796



Biceps and Pectoralis Repair

Loop 'N' Tack™ Biceps Tenodesis	58
2.6 FiberTak® Button	59
FiberTak Biceps Implant System	59
Knotless FiberTak Biceps Implant System.....	60
Knotless TensionTight™ Implant System	60
Distal and Subpectoral Biceps Repair	61
Proximal Tenodesis Button.....	62
Two-Incision Distal Biceps Implant System	62
Pec Button.....	63
Birmingham Biceps™ Implant System	63
SwiveLock® Tenodesis Implants.....	64
Tenodesis Screw System	65

Loop 'N' Tack™ Biceps Tenodesis



The Loop 'N' Tack technique is a simple approach to arthroscopic biceps tenodesis. The technique allows for suture to be passed and retrieved under direct visualization through the biceps tendon while still attached to the superior labrum. The biceps can be fixated to the tuberosity or incorporated into a rotator cuff repair.

Product Description	Item Number
4.75 Knotless BC SwiveLock w/ HD Scorpion Needle Loop 'N' Tack Tenodesis Implant System Includes 4.75 mm × 19.1 mm Knotless BioComposite SwiveLock anchor, punch, FiberLink SutureTape, and HD Scorpion needle w/MegaLoader	AR-1665KBCSLSN
4.75 Knotless BC Loop 'N' Tack Tenodesis Implant System Includes 4.75 mm × 19.1 mm Knotless BioComposite SwiveLock® anchor, punch, Loop 'N' Tack SwiftStitch™ passer, and FiberLink™ SutureTape	AR-1665KBCSL
4.75 BC SwiveLock w/ HD Scorpion Needle Loop 'N' Tack Tenodesis Implant System Includes 4.75 mm × 19.1 mm BioComposite SwiveLock anchor, punch, FiberLink SutureTape, and HD Scorpion needle w/MegaLoader	AR-1665BCSLSN
4.75 BC Loop 'N' Tack Tenodesis Implant System Includes 4.75 mm × 19.1 mm BioComposite SwiveLock anchor, punch, Loop 'N' Tack SwiftStitch passer, and FiberLink SutureTape	AR-1665BCSL
3.9 BC SwiveLock w/ HD Scorpion Needle Loop 'N' Tack Tenodesis Implant System Includes 3.9 mm × 17.9 mm BioComposite SwiveLock anchor, punch, FiberLink SutureTape, and HD Scorpion needle w/MegaLoader	AR-1665BCSLSN-39
3.9 BC Loop 'N' Tack Tenodesis Implant System Includes 3.9 mm × 17.9 mm BioComposite SwiveLock anchor, punch, Loop 'N' Tack SwiftStitch passer, and FiberLink SutureTape	AR-1665BCSL-39
2.9 BC PushLock w/HD Scorpion Needle Loop 'N' Tack Tenodesis Implant System Includes 2.9 mm × 12.5 mm BioComposite short PushLock anchor, spear, trocar, drill, FiberLink SutureTape, and HD Scorpion needle w/ MegaLoader	AR-1665BCPLSN
2.9 BC Loop 'N' Tack Tenodesis Implant System Includes 2.9 mm × 12.5 mm BioComposite short PushLock® anchor, spear, trocar, drill, Loop 'N' Tack SwiftStitch passer, and FiberLink SutureTape	AR-1665BCPL


Product Description	Item Number
4.75 Knotless PEEK Loop 'N' Tack Tenodesis Implant System Includes 4.75 mm × 19.1 mm Knotless PEEK SwiveLock anchor, punch, Loop 'N' Tack SwiftStitch passer, and FiberLink SutureTape	AR-1665KPSL
4.75 PEEK Loop 'N' Tack Tenodesis Implant System Includes 4.75 mm × 19.1 mm PEEK SwiveLock anchor, punch, Loop 'N' Tack SwiftStitch passer, and FiberLink SutureTape	AR-1665PSL
3.9 PEEK Loop 'N' Tack Tenodesis Implant System Includes 3.9 mm × 17.9 mm PEEK SwiveLock anchor, punch, Loop 'N' Tack SwiftStitch passer, and FiberLink SutureTape	AR-1665PSL-39
2.9 PEEK Loop 'N' Tack Tenodesis Implant System Includes 2.9 mm × 12.5 mm PEEK short PushLock anchor, spear, trocar, drill, Loop 'N' Tack SwiftStitch passer, and FiberLink SutureTape	AR-1665PPL
Loop 'N' Tack SwiftStitch passer	AR-4068LNT

Product Description	Item Number
Instruments	
Suture Cutter, 4.2 mm, open end, left notch	AR-11794L

Literature

Product Description	Item Number
Loop 'N' Tack™ Biceps Tenodesis Surgical Technique	LT1-00136-EN

2.6 FiberTak® Button




The FiberTak button is an all-suture button that can be used for both bicortical and unicortical tension slide tenodesis techniques. The button is a reinforced 2.6 mm diameter sheath preloaded with two 2-0 FiberLink™ sutures that are used to shuttle the whipstitch limbs through the sheath.

Product Description	Item Number
2.6 FiberTak Button Implant System Includes FiberTak button w/ inserter and double-loaded FiberLink, 2.6 mm spade tip drill, drill guide, #2 FiberLoop® suture, and free needle	AR-3680
2.6 FiberTak Button	AR-3681
2.6 mm Spade-Tip Drill, sterile	AR-3682
Drill Guide for FiberTak Button, circumferential teeth, reusable	AR-3683
SutureTag w/ double 14 in #2 tails and 36.6 mm ½ circle tapered needle (white/black)	AR-7274T
SutureTag w/ double 14 in #2 tails and 36.6 mm ½ circle tapered needle (white/blue/black)	AR-7274TT

Literature

Product Description	Item Number
Subpectoral Biceps Tenodesis Using the FiberTak® Button	LT1-000166-en-US

FiberTak Biceps Implant System




The FiberTak Biceps Implant System delivers an all-suture anchor designed for use in open, onlay-tissue fixation procedures, particularly proximal and distal biceps repairs.

Product Description	Item Number
FiberTak Biceps Implant System Includes FiberTak biceps implant w/ 4 needles, slotted drill guide, 1.9 mm drill, and free curved needle	AR-3670
FiberTak Biceps Implant, double-loaded w/ white/blue and white/black SutureTape and 4 needles	AR-3671
FiberTak Biceps Disposables Kit Includes slotted drill guide and 1.9 mm drill	AR-3671DS
1.9 mm Drill for FiberTak Biceps Implant	AR-3672
Drill Guide for FiberTak Biceps Implant, circumferential teeth, slotted, reusable	AR-3673

Literature

Product Description	Item Number
FiberTak® Biceps Implant System Biceps Tenodesis Surgical Technique	LT1-000046-en-US

Knotless FiberTak® Biceps Implant System



The Knotless FiberTak Biceps Implant System offers a knotless, all-suture anchor solution to perform an onlay proximal biceps tenodesis at any location, depending on surgeon preference. The implant uniquely integrates two individual tensionable knotless mechanisms within a single FiberTak sheath.

Product Description	Item Number
Knotless FiberTak Biceps Implant System Includes: Knotless FiberTak biceps implant, 3 mm drill, drill guide, Loop 'N' Tack™ SwiftStitch™ suture passer	AR-3688
Knotless FiberTak Biceps Implant	AR-3687

Literature

Product Description	Item Number
Suprapectoral Biceps Tenodesis with the Knotless FiberTak Biceps Anchor Surgical Technique	LT1-000251-en-US



Knotless TensionTight™ Implant System





The TensionTight button is used to perform a knotless onlay biceps tenodesis, using either an arthroscopic suprapectoral or open subpectoral approach. A locking jaw in the button allows for a knotless biceps tenodesis. Securing the biceps tendon repair with a Loop 'N' Tack stitch enables an all-arthroscopic suprapect onlay tenodesis.

Product Description	Item Number
Knotless TensionTight Button Implant System Includes Knotless TensionTight Button with inserter, two #5 FiberLink™ sutures, 3.7 mm drill pin, curved needle with nitinol loop, and shoehorn cannula	AR-2350
Loop 'N' Tack SwiftStitch Suture Passer	AR-4068LNT

Literature

Product Description	Item Number
Suprapectoral Biceps Tenodesis With the TensionTight™ Locking Button Surgical Technique	LT1-000232-en-US

Distal and Subpectoral Biceps Repair



The tension-slide technique with the BicepsButton™ implant system provides surgeons a simple, reproducible, and biomechanically stable distal and proximal biceps repair. This “tensioning” technique draws the tendon against the distal cortex of the bone socket. A tenodesis screw can be added to augment the repair and place the tendon in a more anatomic position.

Literature

Product Description	Item Number
Distal Biceps Repair Implant Kit Includes BicepsButton implant, button inserter, FiberLoop® suture, 7 mm × 10 mm PEEK tenodesis screw, and BicepsButton drill pin	AR-2260
BioComposite Distal Biceps Implant Kit Includes BicepsButton implant, button inserter, FiberLoop suture, 7 mm × 10 mm biocomposite tenodesis screw, disposable tenodesis driver, BicepsButton drill pin, suture passing wire, and free needle	AR-2260BC
BicepsButton Implant, 2.6 mm × 12 mm	AR-2261
Knotless Distal Biceps Implant Kit Includes Knotless BicepsButton implant, button inserter, BicepsButton drill pin, FiberTape® loop, FiberLoop suture, and curved free needle	AR-2280
Knotless BicepsButton Implant, 2.6 mm × 12 mm	AR-2281
FiberTape Loop, 2 mm, (blue)	AR-7275
#2 FiberLoop® Suture w/ straight needle	AR-7234
PEEK Tenodesis Screw, 7 mm × 10 mm	AR-1670PS
Required Instruments	
Tenodesis Screw System includes	AR-1675S
Driver for tenodesis screws	AR-1670DB
Tear-Drop Handle w/ suture cleat	AR-2001BT
Headed Reamer, 8 mm, cannulated	AR-1408
Required Disposables	
Button Inserter	AR-2262
BicepsButton Drill Pin, 3.2 mm, sterile	AR-2263
Tapered Needle w/ Nitinol loop, 26.5 mm, ½ circle	AR-7281
Nitinol Suture-Passing Wire	AR-1255-18
Optional Disposables	
Low-Profile Reamer, 8 mm	AR-1408LP

Product Description	Item Number
Distal Biceps Repair Surgical Technique	LT1-0592-EN
Subpectoral Biceps Tenodesis Surgical Technique	LT1-0591-EN


Proximal Tenodesis Button



The proximal tenodesis button has an angled face to promote a toggle effect when the button contacts the far cortex, making it ideal for unicortical biceps tenodesis.

Product Description	Item Number
Proximal Tenodesis Implant Kit Includes Proximal BicepsButton™ implant, button inserter, FiberLoop® suture, free needles, drill pin, and a shoehorn cannula	AR-2290
Proximal Tenodesis Button, 2.6 mm × 8.5 mm, w/ inserter	AR-2291
#2 FiberLoop® Suture (blue) w/ straight needle	AR-7234
Tapered Needle w/ Nitinol Loop, 26.5 mm, ½ circle	AR-7281
Drill Pin, 3.2 mm, sterile	AR-2263
Shoehorn Cannula	AR-6565

Two-Incision Distal Biceps Implant System



The Two-Incision Distal Biceps Implant System provides everything required to perform a distal biceps repair from a posterior approach using two 7 mm unicortical buttons.

Product Description	Item Number
Two-Incision Distal Biceps implant System Includes two 2.6 mm × 7 mm buttons w/ inserters, two #2 FiberWire® w/ needles (1 blue/white, 1 black/white) 3.2 mm drill pin, flexible obturator, and curved obturator w/ nitinol wire	AR-2295
Two-Incision Tenodesis Button, 2.6 mm × 7 mm	AR-2296
Curved Obturator, cannulated, 3 mm	AR-2294


Pec Button



The pec button and large pec button have an angled face to promote a toggle effect when the button contacts the opposite cortex, making these buttons ideally suited for unicortical repairs, particularly pectoralis major injuries. The large pec button features larger eyelets to accommodate FiberTape® and #5 FiberWire® sutures.

Product Description	Item Number
Pec Button Implant Repair Kit Includes 3 pec buttons w/ FiberWire suture and needles, button inserter, and 3.2 mm drill pin	AR-2268
Pec Button, 2.6 mm × 10.9 mm (a)	AR-2266
Button Inserter	AR-2262
3.2 mm Drill, sterile	AR-2263
Spear w/ Circumferential Teeth w/ Obturator	AR-1906
#2 FiberWire Suture, 38 in (1 blue, 1 black/white)	AR-7201
Pec Button, Larger Eyelet Repair Kit Includes 3 large pec buttons w/ attached inserter, three 17 in FiberTape sutures w/ needles, three 38 in #5 FiberWire suture w/ needles, and 3.7 mm drill pin	AR-2269
Large Pec Button w/ Inserter, 3.2 mm × 10.9 mm, larger eyelets (b)	AR-2267
3.7 mm Drill	AR-2272
FiberTape Suture, 2 mm, 17 in (blue), each end tapered to #2 FiberWire suture w/ tapered needle, 26.5 mm, ½ circle	AR-7237-17N


Birmingham Biceps™ Implant System



The Birmingham Biceps Implant System for subpectoral tenodesis was designed to minimize the potential for postoperative humeral fractures by placing a small 5.5 mm PEEK SwiveLock® tenodesis implant near the top of the pectoralis at a newly defined anatomical landmark, the pallium pectoralis. The implant system includes everything needed for the procedure, making it very convenient for surgeon and staff.

Product Description	Item Number
Birmingham Biceps Implant Kit Includes 5.5 mm PEEK SwiveLock anchor, FiberLink Plus™ suture, free needle, 5.5 mm pilot-headed reamer, and 5.8 mm tap	AR-1662S-1

SwiveLock® Tenodesis Implants



The family of SwiveLock tenodesis implants has been designed to support a multitude of proximal biceps tenodesis techniques and maximize OR efficiency. Implants are available with closed eyelets for surgeons who prefer a whipstitching technique and forked eyelet versions are offered for an all-arthroscopic approach.

Product Description	Item Number
Forked Eyelet Implant Options	
Biocomposite, 7 mm × 19.5 mm	AR-1662BC-7
Biocomposite, 8 mm × 19.5 mm (a)	AR-1662BC-8
PEEK, 4.75 mm × 14 mm	AR-7324PSL
PEEK, 7 mm × 19.5 mm	AR-1662PSL-7
PEEK, 8 mm × 19.5 mm	AR-1662PSL-8
Closed Eyelet Implant Options	
Biocomposite, 6.25 mm × 19.1 mm	AR-1662BC
Biocomposite, 7 mm × 19.1 mm	AR-1662BCC-7
Biocomposite, 8 mm × 19.1 mm (b)	AR-1662BCC-8

Pilot-Headed Reamers (Supplied Sterile)

Product Description	Item Number
6 mm	AR-1450
6.5 mm	AR-1451
7 mm	AR-1452
7.5 mm	AR-1453
8 mm (c)	AR-1454
8.5 mm	AR-1455
9 mm	AR-1456


SwiveLock Tenodesis Sizers

Product Description	Item Number
SwiveLock Tenodesis Sizer, 7 m	AR-1662T-7
SwiveLock Tenodesis Sizer, 8 mm (d)	AR-1662T-8

Literature

Product Description	Item Number
Biceps Tenodesis SwiveLock® System Surgical Technique	LT1-0595-EN

Tenodesis Screw System



The Tenodesis Screw System was designed specifically for reattaching both ligament and tendon to bone, making it ideal for the reattachment of soft tissue to bone in upper extremity procedures including rotator cuff repairs, proximal/distal biceps tenodesis, and acromioclavicular joint reconstruction.


Product Description	Item Number
Tenodesis Master Set includes	AR-1675S
Cannulated Drill, 4 mm	AR-1204L
Cannulated Drill, 4.5 mm	AR-1204.5L
Cannulated Headed Reamers Sizes: 5, 5.5, 6, 6.5, 7, 7.5, 8, 8.5, 9, and 10 mm	AR-1405 - AR-1410
Tear Drop Handle w/ suture cleat	AR-2001BT
Driver for 10 mm tenodesis screws (AR-1540BC)	AR-1540DB
Driver for 10 mm and 12 mm tenodesis screws (AR-1670BC, AR-1680BC)	AR-1670DB
Driver for 15 mm tenodesis screws (AR-1547BC, AR-1555BC, AR-1562BC, AR-1655PS, AR-1655PS-10, and AR-1655PS-12)	AR-1350D
Driver for 23 mm tenodesis screws (AR-1570BC, AR-1580BC, and AR-1590BC)	AR-1570DB
Tenodesis Screw Instrumentation Case	AR-1675C
Optional Accessories	
Tenodesis Tap, 4 mm × 10 mm	AR-1540T
Tenodesis Tap, 4.75 mm × 15 mm	AR-1547T
Tenodesis Tap, 5.5 mm × 15 mm	AR-1555T
Tenodesis Tap, 6.25 mm × 15 mm	AR-1562T
Tenodesis Tap, 7 mm × 10 mm	AR-1670T
Tenodesis Tap, 7 mm × 23 mm	AR-1570T
Tenodesis Tap, 8 mm × 12 mm	AR-1680T
Driver Shaft for 10 mm tenodesis screw	AR-1540DB-01
Driver Shaft for 12 mm tenodesis screw	AR-1670DB-01
Driver Shaft for 15 mm tenodesis screw	AR-1350D-01
Driver Shaft for 23 mm tenodesis screw	AR-1570DB-01
Tear-Drop Handle (required for starter awl)	AR-2001

Product Description	Item Number
Implants	
BioComposite Tenodesis Screw w/ handled inserter, 3 mm × 8 mm	AR-1530BC
BioComposite Tenodesis Screw, 4 mm × 10 mm	AR-1540BC
Disposable Tenodesis Driver w/ 4 mm × 10 mm BioComposite screw and #2 FiberWire® suture	AR-1540CDS
BioComposite Tenodesis Screw, 4.75 mm × 15 mm	AR-1547BC
Disposable Tenodesis Driver w/ 4.75 mm × 15 mm BioComposite screw and #2 FiberWire suture	AR-1547CDS
BioComposite Tenodesis Screw, 5.5 mm × 15 mm	AR-1555BC
Disposable Tenodesis Driver w/ 5.5 mm × 15 mm BioComposite screw and #2 FiberWire suture	AR-1555CDS
BioComposite Tenodesis Screw, 6.25 mm × 15 mm	AR-1562BC
Disposable Tenodesis Driver w/ 6.25 mm × 15 mm BioComposite screw and #2 FiberWire suture	AR-1562CDS
BioComposite Tenodesis Screw, 7 mm × 10 mm	AR-1670BC
BioComposite Tenodesis Screw, 7 mm × 23 mm	AR-1570BC
BioComposite Tenodesis Screw, 8 mm × 12 mm	AR-1680BC
BioComposite Tenodesis Screw, 8 mm × 23 mm	AR-1580BC
BioComposite Tenodesis Screw, 9 mm × 23 mm	AR-1590BC
PEEK Tenodesis Screw w/ handled inserter, 3 mm × 8 mm	AR-1530PS
PEEK Tenodesis Screw, 4 mm × 10 mm	AR-1540PS
PEEK Tenodesis Screw, 4.75 mm × 15 mm	AR-1547PS
PEEK Tenodesis Screw, 5.5 mm × 8 mm	AR-1655PS
PEEK Tenodesis Screw, 5.5 mm × 15 mm	AR-1555PS
PEEK Tenodesis Screw, 6.25 mm × 15 mm	AR-1562PS
PEEK Tenodesis Screw, 7 mm × 10 mm	AR-1670PS
PEEK Tenodesis Screw, 7 mm × 23 mm	AR-1570PS
PEEK Tenodesis Screw, 8 mm × 12 mm	AR-1680PS
PEEK Tenodesis Screw, 8 mm × 23 mm	AR-1580PS
PEEK Tenodesis Screw, 9 mm × 23 mm	AR-1590PS
Disposables	
Tenodesis Disposables Kit Includes 2.4 mm short guide pin, suture passing wire, #2 FiberWire suture w/ needle, 2-0 FiberWire suture w/ needle, #2 FiberWire suture, 2-0 FiberWire suture, #2 FiberLoop® suture w/ straight needle, skin marking pen, and ruler	AR-1676DS
Tenodesis Disposables Kit for 3.0 mm × 8 mm screw	AR-1530DS
#2 FiberSnare® Suture, #2 FiberWire suture, 26 in w/ closed loop, one end stiffened, 12 in	AR-7209SN
#2 FiberLoop Suture w/ straight needle	AR-7234
Short Guide Pin, 2.4 mm, qty. 6	AR-1250SB
Suture-Passing Wire	AR-1255-18
Guidewire, 0.86 mm w/ laser line, SS	AR-8737-21

Elbow Arthroscopy

Elbow Arthroscopy Set.....	68
Elbow Ulnar Collateral Ligament (UCL) Reconstruction Set.....	69

Elbow Arthroscopy Set



The Elbow Arthroscopy Instrument Set is the first and only comprehensive set of instruments available for elbow arthroscopy. Precision manufactured to the highest-quality standards and thoughtfully designed around the feedback of leading elbow arthroscopists, the set enables improved portal placement and includes unique instruments designed for tissue retraction/ elevation, microfracture, precise soft-tissue resection with arthroscopic scissors and punches, loose body graspers, osteotomes, and curettes. The previously unreachable and unfixable can now be addressed, opening up new horizons in patient care for elbow pathology.

Product Description	Item Number
Elbow Arthroscopy Set includes	AR-8650S
Osteotome, angled up (a)	AR-1770
Osteotome, straight	AR-1771
Extra Long Switching Stick, 4 mm	AR-3026
Graduated Black Probe, 2.5 mm Tip	AR-5007
Portal Dilating Set	AR-8620
Articulating Paddle Elevator	AR-8630
Articulating Hook Elevator	AR-8631
Howard Elevator	AR-8650-06
Cup Curette, 45°	AR-8650-07
Osteotome, 5.5 mm, angled up	AR-8650-08
Osteotome, 5.5 mm, straight	AR-8650-09
Ring Curette, reverse angled	AR-8655-04
Chondral Pick, 30°	AR-8655-05
Chondral Pick, 60°	AR-8655-06
Osteotome, curved	AR-8655-09
Ring Curette, angled (c)	AR-8655-15
Probe (b)	AR-8655-16
Cup Curette, 100 mm long, curved shaft	AR-8660
Cup Curette, 100 mm long, straight shaft (d)	AR-8661
Chondro Pick, 40°, straight	AR-8670
Obturator	AR-8680
Hook Probe, 5.4 mm, tip w/ 5 mm markings	AR-10000

Product Description	Item Number
Articulating Probe	AR-10100
Grasper, mini straight tip, ø2.75 mm, straight shaft w/ SR handle	AR-11900SR
Scissor, serrated tooth, straight tip, ø3.4 mm, straight shaft	AR-12140
Scissor, ø3.4 mm, straight shaft, right curved tip (e)	AR-12150
Scissor, ø3.4 mm, straight shaft, left curved tip	AR-12160
Punch, 3.4 mm, capsular, inverted scoop tip	AR-12170
WideBiter™ Punch, 3.4 mm, capsular, inverted scoop tip	AR-12180
Punch, 3.4 mm, straight shaft, large straight tip	AR-12200
Punch, 3.4 mm, 15° up curved shaft, large straight tip	AR-12210
Suture Retriever, 3.4 mm, straight	AR-12540
Graft-Passing Forceps w/ SR handle	AR-13400SR
Grasper, ø4.2 mm, straight shaft w/ SR handle, blunt straight tip	AR-13500SR
Grasper, ø4.2 mm, straight shaft w/ SR handle, alligator hook tip	AR-13600SR
KingFisher® Suture Retriever (f)	AR-13970SR
Elbow Instrument Case	AR-8650C
Required Disposables	
1.5 mm Pins for AR-6520, qty. 5	AR-6521
Blunt Needle for Portal Dilatation Set, qty. 5	AR-6523

Elbow Ulnar Collateral Ligament (UCL) Reconstruction Set


This set includes all the instrumentation needed to perform elbow UCL reconstructions and accommodates the modified Jobe, docking, and DANE TJ techniques in addition to techniques using tenodesis screws or button fixation. This set's guided instruments, drills, and reamers allows precise placement of bone tunnels and sockets in the ulna and humeral epicondyle. The recommended disposable kits provide novel instruments to easily pass sutures through the prepared bone tunnels and sockets, making graft passage quick and easy.

Product Description	Item Number
Elbow UCL Reconstruction Set (a) includes	AR-7700S
Drill, 4.5 mm, cannulated	AR-1204.5L
Drill, 4 mm, cannulated	AR-1204L
Drill, 5 mm, cannulated	AR-1205L
UCL Humeral Tunnel Guide Sleeve, 5 mm × 15 mm (b)	AR-7710
UCL Humeral Tunnel Drill Bit, 4.5 mm	AR-7710-4.5
UCL Humeral Tunnel Drill Bit, 5.0 mm	AR-7710-5.0
UCL Humeral Socket Drill Guide, 4.5 mm × 15 mm	AR-7711
Adjustable Humeral Guide, 4.5 mm (c)	AR-7720
Humeral Drill, 2.0 mm	AR-7720-2.0
Humeral Drill, 3.5 mm	AR-7720-3.5
Humeral Drill 4.0 mm (h)	AR-7720-4.0
UCL Graft Sizing Block (f)	AR-7730
V-guide Drill, 3.5 mm	AR-7750-3.5
Ulna V-guide, 55° (d)	AR-7750-55
Ulna Partial V-guide Drill, 55°	AR-7750-55M
Intersecting V-guide Obturator	AR-7751
Drill Guide, 3.5 mm	AR-7752
Ulna Offset Guide, 7.0 mm (e)	AR-7755
#2 Curette, 7 in long, 45° tip	AR-7799-2
3-0 Curette, 7 in long, 45° tip	AR-7799-3-0
4-0 Curette, 7 in long, 45° tip (g)	AR-7799-4-0
Elbow UCL Reconstruction Case	AR-7700C

Product Description	Item Number
UCL suture Passing Disposables Kit (i) Includes #2 TigerLoop™ suture, #2 FiberLoop® suture, curved Micro SutureLasso™ suture passer w/ wire loop, 4.5 mm wire skid, chamfer tool, 2 suture-passing wires, suture passer, blunt curved needle w/ nitinol loop, and 6-in ruler	AR-7715-4.5
UCL InternalBrace Ligament Augmentation Kit (j) Includes two 3.5 mm PEEK SwiveLock® anchors, 8-in collagen-coated FiberTape® suture (blue), 0 FiberWire® suture (white), 2.7 mm drill, drill guide, 3.8 mm punch/tap (gray), 4.2 mm tap (red), and free needle (½ circle, tapered)	AR-7715

Literature

Product Description	Item Number
Elbow UCL Repair With the UCL <i>InternalBrace</i> ™ System	LT1-000230-en-US



AC Reconstruction Systems

Knotless AC TightRope® Implant.....	72
Low-Profile AC Implant.....	72
AC FiberTape® Cerclage Implant System.....	72
Dog Bone™ Button	73
AC TightRope Implant.....	73
Twin Tail TightRope Implant.....	73
AC <i>Internal</i> Brace™ Ligament Augmentation Repair	74
Adjustable AC Guides.....	74
Coracoid Passing Instruments.....	74
AC Joint Master Set.....	75
Open Chronic AC Joint Reconstruction.....	75

Knotless AC TightRope® Implant



The Knotless AC TightRope implant with Dog Bone™ button combine to form a knotless tensionable device designed for the reduction and stabilization of acute and chronic acromioclavicular joint separations.

Product Description	Item Number
Knotless AC TightRope Implant w/ Dog Bone Button	AR-2370BL
Knotless AC TightRope System Includes implant and 3 mm drill	AR-2371BL
Knotless AC TightRope Open Repair Implant w/ Pec Button (a)	AR-2372BLO

Literature

Product Description	Item Number
Knotless AC Repair System	LT1-000162-en-US

Low-Profile AC Implant



The Low-Profile AC implant and Dog Bone button combine to form a tensionable device optimized for the stabilization of acute and chronic AC joint separations.

Product Description	Item Number
Low-Profile AC Implant Kit Includes Knotless AC implant, Dog Bone button, 3 mm drill, SutureLasso™ SD wire loop, and 5.1 mm cannulated reamer	AR-2371
Knotless AC Implant	AR-2370
5.1 mm Cannulated Reamer w/ stop	AR-2374
Required Instruments and Disposables	
Acromioclavicular Joint Master Set	AR-2255MS
SutureLasso SD Wire Loop, 5/box	AR-4068-05SD
Cannulated Drill for AC repair, 3 mm or	AR-2257D-30
Cannulated Drill for AC repair, 2.4 mm	AR-2257D-24

Literature

Product Description	Item Number
Low-Profile AC Repair System	LT1-00059-EN

AC FiberTape® Cerclage Implant System




The FiberTape cerclage suture can eliminate the need to drill clavicle and/or coracoid tunnels for AC repair. The system offers a convenient kit that includes a complete set of instruments and implants for performing various AC cerclage techniques.

Product Description	Item Number
AC FiberTape Cerclage Implant System (b) Includes FiberTape cerclage, TigerTape™ sternal closure w/ blunt needle, TigerLink™ suture, FiberLink™ SutureTape, tensioner, MegaLoader, small cerclage passer, transosseous suture passer, flexible obturator, 2.4 mm drill bit, AC wire passer, SutureLasso SD wire loop, and Dog Bone button	AR-7820AC

Literature

Product Description	Item Number
AC FiberTape Cerclage Surgical Technique	LT1-000298-en-US



Dog Bone™ Button




The precontoured, titanium Dog Bone button allows the use of knotless TightRope® technology or multiple FiberTape® sutures for AC joint reduction. The buttons are attached to the Knotless TightRope® implant or FiberTape sutures independently, only suture material will be passed through the clavicle and coracoid tunnels, allowing the repair to be completed through smaller tunnels.

Product Description	Item Number
Acute AC Repair Kit Includes Dog Bone button preloaded onto FiberTape and TigerTape™ loops, free Dog Bone button, 3 mm cannulated drill, and nitinol suture-passing wire	AR-2271
Dog Bone Button	AR-2270
FiberTape Loop	AR-7275
TigerTape Loop	AR-7275T
Required Instruments and Disposables	
Acromioclavicular Joint Master Set	AR-2255MS
SutureLasso™ SD Wire Loop, 5/box	AR-4068-05SD
Cannulated Drill for AC repair, 3 mm or	AR-2257D-30
Cannulated Drill for AC repair, 2.4 mm	AR-2257D-24
Optional	
Clavicle Button	AR-2257-01
Coracoid Button	AR-2257-02
Clavicle Washer	AR-2258-03

Literature

Product Description	Item Number
AC Repair - Dog Bone Button Surgical Technique	LT1-0531-EN

AC TightRope® Implant



This implant allows surgeons to easily reconstruct acute AC joint separations in a minimally invasive manner, either open or arthroscopic, using a 4-strand continuous loop of #5 FiberWire® suture interlaced between 2 titanium buttons.

Product Description	Item Number
AC TightRope Repair Kit Includes AC TightRope implant, suture-passing wire, and knot pusher	AR-2257
Required Instruments	
AC Joint Reconstruction Master Set	AR-2255MS
Required Disposables	
2.4 mm Drill-Tip Guide Pin	AR-1250L
Optional Disposables	
SutureLasso™ SD Wire Loop, 5/box	AR-4068-05SD
Button Inserter	AR-2262


Twin Tail TightRope® Implant



This implant features 2 independent clavicle button tails and is designed to help reduce and stabilize the AC joint for open, acute AC injuries.

Product Description	Item Number
Twin Tail TightRope	AR-2264
Required Instruments	
Drill-Tip Guide Pin, 2.4 mm	AR-1250L
Cannulated Drill, 4.5 mm	AR-1204.5L
Cannulated Long Drill, 4 mm	AR-1204LX
Button Inserter	AR-2262
SutureLasso SD Suture Passer, crescent or	AR-4068C
Micro SutureLasso™ Suture Passer, minor bend	AR-8701
Optional Instruments	
Guide Pin Sleeve	AR-2255CG-02

AC *InternalBrace*™ Ligament Augmentation Repair



InternalBrace ligament augmentation repair combined with CC ligament fixation may help limit posterior displacement of the clavicle. A 5 mm offset guide can help drill the acromion and clavicle anteriorly. The guide accepts both the existing 3.5 mm and 4.75 mm SwiveLock® drills.


Product Description	Item Number
5 mm Offset Guide for AC <i>InternalBrace</i> Ligament Augmentation repair (a)	AR-2276-05
Twist Drill for 4.75 mm SwiveLock anchor (b)	AR-2276D
Drill for 3.5 mm SwiveLock anchor, sterile	AR-2325D-ST
BioComposite SwiveLock C Anchor, 4.75 mm × 19.1 mm	AR-2324BCC
BioComposite SwiveLock Anchor, 3.5 mm × 15.8 mm	AR-2325BCC
FiberTape® Suture, 2 mm, 7 in (blue), each end tapered to #2 FiberWire® suture	AR-7237-7

The *InternalBrace* surgical technique is intended only to augment the primary repair/reconstruction by expanding the area of tissue approximation during the healing period and is not intended as a replacement for the native ligament. The *InternalBrace* technique is for use during soft tissue-to-bone fixation procedures and is not cleared for bone-to-bone fixation.

Adjustable AC Guides



Coracoid Passing Instruments




The AC wire passer and flexible obturator allow surgeons to arthroscopically pass a graft around the coracoid to address chronic AC injuries. This system is intended to be used in conjunction with mechanical reduction and fixation of the CC ligaments, particularly with Dog Bone™ buttons and FiberTape® sutures. The coracoid passer is designed to fit through a small incision on the clavicle and curve around the coracoid, allowing the surgeon to pass a FiberStick™ suture through the cannulation and retrieve it on the lateral side of the coracoid. The flexible obturator is designed to shuttle around the coracoid and dilate the soft tissue prior to passing an allograft or autograft.

Product Description	Item Number
AC Wire Passer	AR-2252
Flexible Obturator	AR-2275

The adjustable AC guides give surgeons the ability to arthroscopically hook or rest the tip of the guide on the medial side of the coracoid, as well as independently adjust the guide medially and laterally on the clavicle.

Product Description	Item Number
AC Target Guide, left	AR-2253L
AC Target Guide, right	AR-2253R
AC Drill Guide Assembly	AR-2273
ACL Guide Frame Handle Assembly	AR-1510H
Drill Sleeve, 3 mm	AR-2255CG-05

AC Joint Master Set



The Acromioclavicular Joint Reconstruction System contains the instrumentation necessary to complete the acute and chronic AC joint repair techniques.

Product Description	Item Number
AC Joint Reconstruction System includes	AR-2255MS
Cannulated Drill, 4 mm	AR-1204L
Cannulated Drill, 4.5 mm	AR-1204.5L
Cannulated Headed Reamer, 5 mm	AR-1405
Cannulated Headed Reamer, 5.5 mm	AR-1405.5
Cannulated Headed Reamer, 6 mm	AR-1406
Cannulated Headed Reamer, 6.5 mm	AR-1406.5
ACL Guide Frame handle Assembly	AR-1510H
AC Guide, left	AR-2254L
AC Guide, right	AR-2254R
Guide Pin Sleeve, 2.4 mm	AR-2255CG-02
Clavicle Drill Positioner	AR-2255CG-03
Drill Stop	AR-2255CG-04
Drill Sleeve, 3 mm	AR-2255CG-05
AC Tenodesis Screw Driver	AR-2255D
Coracoid Graft Passer, left	AR-2256L
Coracoid Graft Passer, right	AR-2256R
AC GraftRope Graft Sizer	AR-2265
Forked Probe	AR-6002
AC Joint Instrumentation Case	AR-2255MC
Optional	
Guide Arm, left	AR-2253L
Guide Arm, right	AR-2253R
AC Drill Guide Assembly	AR-2273

Open Chronic AC Joint Reconstruction

This system enables surgeons to reconstruct chronic AC joint separations anatomically using a tendon graft with backup mechanical fixation. The technique uses 5.5 mm × 8 mm PEEK tenodesis screws to achieve strong graft fixation in the clavicle and allows surgeons to recreate the conoid and trapezoid coracoclavicular ligament bundles. The graft is protected during the healing process by using a strand of FiberWire® suture as backup mechanical fixation.

Product Description	Item Number
PEEK Tenodesis Screw, 5.5 mm × 8 mm	AR-1655PS
Required	
AC Joint Reconstruction Master Set	AR-2255MS
Disposables	
SutureLasso™ SD Wire Loop, 5/box	AR-4068-05SD
#2 FiberWire Suture, 38 in w/ tapered needle	AR-7200
#2 FiberWire Suture, 38 in	AR-7233
Optional	
Bio-Tenodesis Tap, 5.5 mm × 15 mm	AR-1555T
Tear-Drop Handle (required w/ tap)	AR-2001
Nitinol Guide Pin, 1.1 mm, 6/box	AR-1249
Bio-Tenodesis Disposable Kit, sterile, 5/box	AR-1676DS

Glenoid Bone Loss

Glenoid Bone Loss Set.....	78
FiberTape® Cerclage for Glenoid Bone Loss.....	79
Distal Tibia Allograft Workstation.....	79

Glenoid Bone Loss Set



The Glenoid Bone Loss Set helps surgeons address the complexity of shoulder instability caused by bony pathology such as anterior glenoid bone loss, bony Bankart lesions, glenoid fractures, and engaging Hill-Sachs lesions. Unique open Latarjet instruments help make the Latarjet technique more consistent and reproducible. The included bony Bankart/glenoid fracture instruments have been designed to allow placement of the partially threaded 3.75 mm cannulated screws, either percutaneously or through a standard arthroscopic cannula.


Product Description	Item Number
Glenoid Bone Loss Set includes	AR-7000S
Osteotome Blade Shield	AR-7000-02
Parallel Drill Guide, 4 mm offset	AR-7000-03
Parallel Drill Guide, 6 mm offset	AR-7000-04
Parallel Drill Guide, 8 mm offset	AR-7000-05
Screw Length Sizer	AR-7000-06
Coracoid Drill Guide	AR-7000-07
Fukuda Retractor, small	AR-7000-08
Kolbel Glenoid Retractor	AR-7000-09
Nesting Guide Sleeves	AR-7000-12
Cannulated Hex Driver, 2.5 mm	AR-7000-13
Cannulated Drill, 2.75 mm	AR-7000-14
“No Offset” Parallel Drill Guide	AR-7000-16
Drill, 4 mm	AR-1204D
Osteotome Handle	AR-2961
Drill Guide Handle	AR-9215-1-01
Cannulated Driver Handle w/ AO Connection	AR-13221AOC
Gelpi Subscapularis Retractor	AR-8104
Glenoid Bone Loss Instrument Case	AR-7000C
Screw Caddy, 3.75 mm, fully threaded screw	AR-7000SC-1
Screw Caddy, 3.75 mm, partially threaded screw	AR-7000SC-2
Optional Instruments	
Angled Saw Blade	AR-300-450S
Coracoid Drill Guide, shallow jaw	AR-7000-17
Bony Bankart Drill Guide	AR-7050
Mini-Open Shoulder Retractor	AR-8100-21
Bone Block Grasper w/ drill stop	AR-8103
Arthroscopic Chisel, 15 mm, large wedge, J-Span	AR-7100-03NL

Product Description	Item Number
Modular Soft-Tissue Retractor Body	AR-8170
Modular Soft-Tissue Retractor Set (includes body)	AR-8170S
Modular Soft-Tissue Retractor Atraumatic Set	AR-8170DS
Modular Soft-Tissue Retractor Atraumatic Paddle Set, 75 mm	AR-8170-75DS
Modular Retractor Replacement Paddles - Atraumatic	
50 mm, left	AR-8170-50DL
50 mm, right	AR-8170-50DR
75 mm, left	AR-8170-75DL
75 mm, right	AR-8170-75DR
Implants and Disposables	
LatarJet Implant System Includes: osteotome blade, a cannulated and non-cannulated drill, three guidewires (6, 7, and 12 in), and wedged profile plate	AR-7000KL
LatarJet Disposable Kit Includes: osteotome blade, a cannulated and non-cannulated drill, and three guidewires (6, 7, and 12 in)	AR-7000KL-1
3.75 mm Partially Threaded Cannulated screw, 26 mm – 42 mm (2 mm increments)	AR-7000-26 to -42
3.75 mm Fully Threaded Cannulated screws, 26 mm – 42 mm (2 mm increments)	AR-7000-26FT to 42FT
Titanium Washer for 3.75 mm screw	AR-7000-15
Wedged Profile Plate	AR-8111
Titanium Suture Washer, w/ #2 FiberWire® suture and curved needle	AR-7000-18T
Osteotome Blade Latarjet	AR-7000-01
Guidewires	
0.062 in (1.6 mm), 6 in long	AR-8941-6
0.062 in (1.6 mm), 7 in long	AR-8941-7
0.062 in (1.6 mm), 12 in long	AR-8941-12

Literature

Product Description	Item Number
Congruent-Arc Latarjet Surgical Technique	LT1-0556-EN

FiberTape® Cerclage for Glenoid Bone Loss



Bone block cerclage is an innovative arthroscopic surgical technique that preserves the integrity of the subscapularis and provides a metal-free fixation of the bone graft.¹ Two interconnected FiberTape and TigerTape™ cerclage sutures provide firm compression of the graft and strong fixation that minimizes construct displacement at high loads.²

The FiberTape Cerclage for Glenoid Bone Loss Implant system contains all the necessary implants and disposables for the procedure (except capsulolabral repair anchors). To help provide an easy and reproducible surgical technique, the kit includes an arthroscopic measuring probe, two 3-mm cannulated drills of different lengths to facilitate unobstructed drilling, and the MegaLoader suture-passing loop to aid in the passage of sutures into instruments and bone grafts.

Product Description	Item Number
FiberTape Cerclage for Glenoid Bone Loss Implant System	AR-7268GBL
FiberTape Cerclage w/ TigerLink™ shuttle suture	
TigerTape Cerclage w/ FiberLink™ shuttle suture	
3 mm Cannulated Drill	
3 mm short Cannulated Drill	
Measurement Probe, 60°, 220 mm	
MegaLoader suture passing loop	
Instruments	
TransGlenoid Drill Guide Hook for anterior instability	AR-7800TGH
Parallel Drill Guide Hook for posterior instability	AR-7800PGH
Drill Guide Sleeve	AR-7800DGS
Tensioner, reusable	AR-7800
Tensioner Handle, ratcheting, reusable	AR-7801

Literature

Product Description	Item Number
Bone Block Cerclage Surgical Technique	LT1-000202-en-US

References

- Hachem AI, Del Carmen M, Verdalet I, Rius J. Arthroscopic bone block cerclage: a fixation method for glenoid bone loss reconstruction without metal implants. *Arthrosc Tech.* 2019;8(12):e1591-e1597. doi:10.1016/j.eats.2019.08.014
- Arthrex GmbH. Data on file (LA2-000140-en-US_A). Munich, Germany; 2021.

Distal Tibia Allograft Workstation



The workstation facilitates preparation of fresh osteochondral distal tibia allograft to help address the complex issue of shoulder instability caused by bony pathology in conjunction with the Glenoid Bone Loss Set.

Product Description	Item Number
Distal Tibia Allograft Workstation includes	AR-7001S
Cannulated Drill, 4 mm	AR-7000D
Sizing Template, 7 mm × 5°	AR-7000-7-5
Sizing Template, 7 mm × 10°	AR-7000-7-10
Sizing Template, 7 mm × 15°	AR-7000-7-15
Sizing Template, 10 mm × 5°	AR-7000-10-5
Sizing Template, 10 mm × 10°	AR-7000-10-10
Sizing Template, 10 mm × 15°	AR-7000-10-15
Allograft Workstation	AR-7001
Cutting Block Assembly, 7 mm × 5°	AR-7001-7-5
Cutting Block Assembly, 7 mm × 10°	AR-7001-7-10
Cutting Block Assembly, 7 mm × 15°	AR-7001-7-15
Cutting Block Assembly, 10 mm × 5°	AR-7001-10-5
Cutting Block Assembly, 10 mm × 10°	AR-7001-10-10
Cutting Block Assembly, 10 mm × 15°	AR-7001-10-15
Holding Post Assembly	AR-7002
Cutting Block Base Assembly	AR-7003
Cutting Block Assembly, graft thickness	AR-7004
Graft Post Assembly	AR-7006
Parallel Drill Guide Assembly, w/ extensions	AR-7007
Parallel Drill Guide Assembly, w/o extensions	AR-7008
Cutter Stop	AR-7009
Cutting Block Assembly, length	AR-7010
Guidewire, 0.090 in × 4 in	AR-7011
Glenoid Drill Guide Handle, long	AR-9215-1-02
Distal Tibia Allograft Workstation Case	AR-7001C
Implants and Disposables	
DTA Implant System w/ Suture Washer	AR-7000KD
Includes washer, 8 guidewires, 2 drills	
DTA Disposable Kit (No washer)	AR-7000KD-1
Includes 8 guidewires, 2 drills	

Literature

Product Description	Item Number
Distal Tibia Allograft Workstation for Glenoid Bone Loss Surgical Technique	LT1-00062-EN

Fracture Management

FiberTape® Cerclage System	82
QuickPass™ FiberTape Cerclage	83
Proximal Humeral Plating System.....	84
Humeral Nail System.....	85
Clavicle Plate and Screw System.....	86
Elbow Fracture Plating System.....	88

FiberTape® Cerclage System



The FiberTape cerclage system consists of specially designed FiberTape sutures, passing needles, passing hooks, and a tensioner. The system is designed to replace standard cable and wire systems commonly used in upper and lower extremity trauma, arthroplasty, and sternal closure indications.

- 100% radiolucent – no radiographic interference from cerclage
- No concerns of metal mismatch with adjacent implants
- No sharp ends or broken wires to potentially harm the patient or medical staff
- No metal debris migration or concerns of metallosis

Product Description	Item Number
FiberTape Cerclage (blue), w/ TigerLink™ shuttle suture and loader	AR-7268
TigerTape™ Cerclage (white/black), w/ FiberLink™ shuttle suture and loader	AR-7268T
Radiopaque FiberTape Cerclage w / TigerLink shuttle suture (no needles)	AR-7268R
FiberTape Cerclage System with Small Passer Includes: two FiberTape cerclage (2 mm, 48 in, blue), TigerLink suture, FiberTape cerclage tensioner, MegaLoader, and small FiberTape cerclage passer	AR-7820SB
FiberTape Cerclage Tension Band System Includes FiberTape Cerclage (2 mm, 48", blue), TigerLink suture, FiberTape cerclage tensioner, MegaLoader, 2.5 mm drill, transosseous suture passer, suture passing wire and two 1.6 mm K-wires	AR-7820TB
FiberTape Cerclage Modular Instrument Set includes	AR-7830S
Tensioner, reusable	AR-7800
Tensioner Handle, ratcheting, reusable	AR-7801
Quick-Connect Handle, 1/4 in square, Nonratcheting handle, reusable	AR-2003NR
Straight Medium Passing Hook, reusable	AR-7831
Straight Large Passing Hook, reusable	AR-7832
Trochanteric Passing Hook, reusable	AR-7833
Left Offset Passing Hook, reusable	AR-7834L
Right Offset Passing Hook, reusable	AR-7834R


Curved Passing Hook, reusable	AR-7835
Blunt Medium Passing Hook, reusable	AR-7836
Blunt Large Passing Hook, reusable	AR-7837
Cerclage Modular Instrument Case	AR-7830C
Disposables	
Cerclage Passing Needle, medium, disposable, qty 5	AR-7816
FiberTape Cerclage Disposable Tensioner (a)	AR-7820
FiberTape Cerclage Straight Passer, medium, disposable (white) (d)	AR-7821
FiberTape Cerclage Straight Passer, large, disposable (black) (c)	AR-7822
FiberTape Cerclage Trochanter Passer, disposable (yellow) (b)	AR-7823
FiberTape Cerclage Offset Passer, right, disposable (red) (f)	AR-7824R
FiberTape Cerclage Offset Passer, left, disposable (green) (g)	AR-7824L
FiberTape Cerclage Straight Passer, small, disposable (blue) (e)	AR-7825

Literature

Product Description	Item Number
FiberTape Cerclage System Surgical Technique	LT1-00130-EN



QuickPass™ FiberTape® Cerclage





The QuickPass FiberTape Cerclage System includes the FiberTape Cerclage with QuickPass tube and cannulated, slotted FiberTape Cerclage QuickPass suture passers. The QuickPass tube used in conjunction with the cannulated, slotted hooks enable the suture to be passed around the bone twice without the need for a separate shuttle suture.

Product Description	Item Number
FiberTape Cerclage, 2 mm, 48 in w/ QuickPass tube, qty. 5 (a)	AR-7269
QuickPass FiberTape Cerclage System – For Humerus (b) includes a medium QuickPass passer, FiberTape cerclage tensioner, MegaLoader, and two FiberTape cerclage w/ QuickPass tube	AR-7820HF
QuickPass FiberTape Cerclage System – For Femur includes a medium QuickPass passer, FiberTape cerclage tensioner, MegaLoader, and two FiberTape cerclage w/ QuickPass tube	AR-7820FF
Disposables	
QuickPass FiberTape Cerclage Passer, medium, white, disposable (c)	AR-7841
QuickPass FiberTape Cerclage Passer, large, black, disposable (d)	AR-7842
QuickPass FiberTape Cerclage Passer, curved, gray, disposable (e)	AR-7843
QuickPass FiberTape Cerclage Passer, right offset, red, disposable (f)	AR-7844R
QuickPass FiberTape Cerclage Passer, left offset, green, disposable (g)	AR-7844L

Literature

Product Description	Item Number
QuickPass FiberTape Cerclage System Surgical Technique	LT1-000326-en-US



Proximal Humeral Plating System



The ALPHA plate is a unique anatomical side-specific plate with a proximal contour that allows for sparing of the deltoid insertion and avoidance of dissection/soft-tissue stripping at the subdeltoid space.

Proximal Humeral Plate Main Tray (9912-100S)

Product Description	Item Number
Instruments	
Trocar Tip Guidewire, 2.0 mm × 150 mm	0102-150
Cortical Tap, 3.5 mm, AO style	0226-100
Calibrated Drill, AO style, 2.8 mm × 95 mm	0237-200
Calibrated Drill, long, 2.5 mm × 130 mm	0250-000
Calibrated Step Drill, AO style, 3.5/2.5 mm × 130 mm	0254-000
Cancellous Tap, 3.5 mm, AO style	0255-200
Screw Sheath	0316-100
Drill Guide, nonlocking/locking, 2.8 mm	0318-200
Drill Guide, 2.8 mm	0324-500
Drill Guide, 2.5 mm	0327-000
Drill Guide, 3.5 mm	0331-000
T15 Driver, AO style	0432-300
T15 Driver, AO style, retaining	0432-500
Screw Gripper, 3.5 mm/4.0 mm screws	0433-000
Hex Driver, 2.5 mm	0448-000
T15 Solid Hex Driver	0448-200
T15 Solid Screwdriver, retaining	0448-400
T15 Screwdriver, large Hudson	0452-100
T15 Screwdriver, large Hudson, retaining	0452-400
Hex Driver, 5.0 mm	5001-000
Ratcheting Small Axial Handle, AO quick connect, cannulated	5026-100
Proximal Screw Depth Gauge	0522-100
Hook Tip Depth Gauge, locking/nonlocking	0523-100
Alignment Guide, proximal humeral plate	1250-200
Bolt, T15, alignment guide, proximal humeral plate	1251-200
Proximal Humeral Plate Case Assembly	9912-100

Product Description	Item Number
Plates and Screws	
Proximal Humeral Plate, 3H	3012-003
Proximal Humeral Plate, 5H (a)	3012-005
Proximal Humeral Plate, 8H	3012-008
Proximal Humeral Plate, 11H	3012-011
Proximal Humeral Plate, 14H	3012-014
Cortical Screw, 3.5 mm × 20 mm-32 mm (2 mm increments), 35 mm, 40 mm, 45 mm	8010-200-450
Cortical Bone Screw, 3.5 mm × 20 mm-32 mm (2 mm increments), 35 mm, 40 mm, 45 mm	8110-020-045
Cortical Screw, double lead lock, 3.5 mm × 20 mm-36 mm (2 mm increments)	8114-020-036
Cancellous Locking Screw, fully threaded, 4.0 mm × 22 mm-60 mm (2 mm increments), 65 mm, 70 mm	8124-022-070
Cancellous Locking Screw, partially threaded, 4.0 mm × 22 mm-60 mm (2 mm increments), 65 mm, 70 mm	8128-022-070

Proximal Humeral Plate Expansion Kit (9941-0001)

Product Description	Item Number
Instruments	
Alignment Guide Locking Bolt, proximal humeral plate	1251-200
Alignment Guide, 95° PHP	1288-000
Alignment Guide, ALPHA plate, left	1296-000
Alignment Guide, ALPHA plate, right	1297-000
AOS Proximal Humeral Plate Expansion Kit Tray	9941-000
Guide Pin, 1.5 mm	0109-150
Implants	
95° Proximal Humeral Plate, 3H	3022-003
95° Proximal Humeral Plate, 5H (b)	3022-005
ALPHA Plate, left, 4H	3024-004
ALPHA Plate, left, 5H	3024-005
ALPHA Plate, left, 7H	3024-007
ALPHA Plate, left, 11H	3024-011
ALPHA Plate, left, 14H	3024-014
ALPHA Plate, left, 16H	3024-016
ALPHA Plate, left, 18H	3024-018
ALPHA Plate, right, 4H	3025-004
ALPHA Plate, right, 5H	3025-005
ALPHA Plate, right, 7H (b)	3025-007
ALPHA Plate, right, 11H	3025-011
ALPHA Plate, right, 14H	3025-014
ALPHA Plate, right, 16H	3025-016
ALPHA Plate, right, 18H	3025-018

Literature

Product Description	Item Number
Proximal Humerus Plate System Surgical Technique	LT1-000234-en-US

Humeral Nail System



The Arthrex Humeral Nail System consists of both short and long intramedullary nails. The 15 cm short nail is best suited for simple two-part fractures of the proximal humerus involving the surgical neck of the humerus and for proximal third humeral fractures without comminution. The long humeral nail is primarily intended for humeral shaft fractures, which are inherently not prone to varus malposition.

Humeral Nail System (9903-200S)


Product Description	Item Number
Instruments	
Straight Awl, cannulated	0218-200
Entry Reamer, cannulated, 10.5 mm	0221-100
Flexible Shaft	0222-200
Reamer Head, 6.0 mm	0223-060
Reamer Head, 6.5 mm	0223-065
Reamer Head, 7.0 mm	0223-070
Reamer Head, 7.5 mm	0223-075
Reamer Head, 8.0 mm	0223-080
Reamer Head, 8.5 mm	0223-085
Reamer Head, 9.0 mm	0223-090
Reamer Head, 9.5 mm	0223-095
Reamer Head, 10.0 mm	0223-100
Drill Guide, green, 2.5 mm	0307-000
Drill Guide, gold, 3.5 mm	0309-000
Driver, hex	0407-100
T-handle	0411-100
Screw Gripper	0413-000
Screw Driver, hex, 2.5 mm	0416-200
Screw Driver, hex, 3.5 mm	0417-100
Guidewire Gripper, 2.0 mm / 3.0 mm	0419-200
Depth Gauge, hook tip	0513-000
Guidewire Depth Gauge	0515-000
Screw Sheath	0601-100

Product Description	Item Number
Distal Sheath	0602-000
Drill Guide, 2.9 mm	0604-000
Screw Trocar	0607-000
Extractor	0803-000
Reduction Tool	0807-000
Impactor Pad	0808-200
Combination Wrench	0813-000
Targeting Module	1228-200
Insertion Guide	1230-100
Locking Bolt	1231-100
Humeral Nail Case	9903-200
Implants	
Proximal Humeral Nail, 8 mm × 15 cm	1170-150
Humeral Nail, 7 mm × 20 cm - 30 cm (2.5 increments)	1171-200-300
Humeral Nail, 8 mm × 20 cm - 30 cm (2.5 increments)	1172-200-300
Humeral Nail, 9 mm × 20 cm - 30 cm (2.5 increments)	1173-200-300
Cancellous Screw, 5.0 mm × 25 mm-60 mm (5 mm increments)	8022-025-060
Cortical Screw, 3.5 mm × 17.5 mm-40 mm (2.5 mm increments)	8010-175-400
End Cap, 2 mm	1074-020
End Cap, 5 mm	1074-050
Fracture Reduction Screw	1175-000
Disposables	
Guide Pin, 3.2 mm	0100-000
Guidewire w/ trocar tip, 2.0 mm	0102-500
Guidewire, ball nose, 2.0 mm	0103-800S
Drill Bit, calibrated, green, 2.5 mm	0217-100
Drill Bit, calibrated, gold, 3.5 mm	0220-100
Drill Bit, short, 2.9 mm	0224-100
Tap, cortical, 3.5 mm	0226-100

Literature

Product Description	Item Number
Humeral Nail System Surgical Technique	LT1-000241-en-US

Clavicle Plate and Screw System



This system is a comprehensive set of plates, screws, and instrumentation designed to treat central third and distal clavicle fractures. The low-profile, precontoured, stainless steel plates combine locking and nonlocking options.

Product Description	Item Number
Clavicle Plate Set includes	AR-2650S
Cannulated Driver Handle	AR-13221AOC
BB-Tak, disposable	AR-13226
Distal Clavicle Plate Positioning Handle	AR-2659
Plate Bending Pliers	AR-2660
Distal Clavicle Screw and Drill Guide, left	AR-2664L
Distal Clavicle Screw and Drill Guide, right	AR-2664R
Verbrugge Forceps w/ pivoting jaw	AR-2662
K-wire, 0.062 in, 3 in, disposable	AR-2663
T10 Hexalobe Driver	AR-2665-T10
T15 Hexalobe Drive Shaft, long	AR-2666-T15
Drill/screw Guide Sleeve	AR-2668
Drill 2 mm	AR-2669
Screw Depth Indicator, 2.7 mm	AR-2670
Step Drill for 2.5 mm fragment screw	AR-2752
T8 Hexalobe Driver	AR-8916-22
T15 Hexalobe Driver	AR-8941DH
Drill Guide, 2.6 mm/1.35 mm	AR-8943-03
Bone Reduction Forceps	AR-8943-07
Screwdriver, T10 hexalobe	AR-8943-08
Screwdriver, T15 hexalobe	AR-8943-10
Screw Holding Sleeve	AR-8943-11
Drill Bit, 2.5 mm, disposable	AR-8943-13
Drill Guide, 3.5 mm/2.5 mm	AR-8943-14
Depth Device	AR-8943-15
Drill Bit, 2 mm, disposable	AR-8943-16
Drill Guide, 2.7 mm, threaded	AR-8943-17
Plate Bending Iron	AR-8943-18
Freer Elevator	AR-8943-19
Periosteal Elevator, 6 mm curved blade	AR-8943-20
Sharp Hook	AR-8943-21
Hohmann Retractor, 15 mm	AR-8943-22

Lobster Claw Forceps	AR-8943-23
Drill Guide, 3.5 mm, threaded	AR-8943-26
Locking Bending Guide, 3.5 mm	AR-8954-07
Clavicle Plating System Instrument Case	AR-2650C
Additional Instruments	
Cannulated Drill, 3 mm	AR-2257D-30
Cannulated Drill, 2.4 mm	AR-2257D-24
Drill Pin, 3.7 mm for AR-2658TR	AR-2272
AC Guide, left	AR-2254L
AC Guide, right	AR-2254R
ACL Guide Frame Handle Assembly	AR-1510H
Drill Sleeve, 3 mm	AR-2255CG-05
T8 Hexalobe Driver	AR-2665-T8

Clavicle Plate Sizing Template for

Product Description	Item Number
AR-2685DL/DR-10	AR-2685DT-10
AR-2622DL/DR	AR-2622DT
AR-2680ST	AR-2680STT

Low-Profile Screws (2 mm increments)


Product Description	Item Number
3.5 mm Nonlocking Screw, cortical, 10 mm - 30 mm	AR-8835-10 to -30
3.5 mm Locking Screw, 8 mm - 24 mm	AR-8835L-08 to -24
2.7 mm Nonlocking Screw, cortical, 10 mm - 24 mm	AR-8827-10 to -24
2.7 mm Locking Screw, 10 mm - 24 mm	AR-8827L-10 to -24
2.7 mm Screw Caddy	AR-8827C-10
4 mm Nonlocking Screw, cancellous, 10 mm - 24 mm	AR-8840-10 to -24
3 mm Nonlocking Screw, cancellous, fully threaded, 10 mm - 24 mm	AR-8830-10 to -24

Fragment Screws (2 mm increments)

Product Description	Item Number
2.5 mm Fragment Screw, 10 mm - 24 mm	AR-2665-10H to -24H
Required Instruments for Fragment Screw	
Step Drill	AR-2752
T8 Hexalobe Driver	AR-2665-T8

Clavicle plates and screws are stainless steel

Clavicle Plate and Screw System (Cont.)



Clavicle Fracture Plates, SS, sterile

Product Description	Item Number
Central Third, left (a)	AR-2650CL
Central Third, right (g)	AR-2650CR
Central Third, left (b)	AR-2651CL
Central Third, right (h)	AR-2651CR
Central Third, left (c)	AR-2652CL
Central Third, right (i)	AR-2652CR
Central Third, left (d)	AR-2653CL
Central Third, right (j)	AR-2653CR
Central Third, left (e)	AR-2654CL
Central Third, right (k)	AR-2654CR
Central Third, left (f)	AR-2655CL
Central Third, right (l)	AR-2655CR
Central Third, 7 hole, straight* (m)	AR-2680ST

Distal Clavicle Fracture Plates, SS, sterile

Product Description	Item Number
Short, left (p)	AR-2656DL
Short, right (q)	AR-2656DR
Long, left (r)	AR-2657DL
Long, right (s)	AR-2657DR
10 Hole, left* (t)	AR-2685DL-10
10 Hole, right* (u)	AR-2685DR-10
22°, left* (n)	AR-2622DL
22°, right* (o)	AR-2622DR

Implants


Product Description	Item Number
Knotless Distal Clavicle Plate Button TightRope® Implant (Knotless distal clavicle plate button attached to large pec button preloaded on inserter)	AR-2658TR
Knotless Distal Clavicle Plate Button (Knotless TightRope suture preloaded on clavicle button for use w/ Dog Bone™ button)	AR-2658T
Distal Clavicle Plate Button	AR-2658
Dog Bone Button	AR-2270
FiberTape® Loop	AR-7275
TigerTape™ Loop	AR-7275T

Literature

Product Description	Item Number
Clavicle Plate and Screw System Surgical Technique	LT1-0255-EN
Inventory Control Form	LF1-0255-EN

*Plates are sterile packaged.

Elbow Fracture Plating System



The Arthrex Elbow Plating System is a comprehensive distal humerus and olecranon plating system. For distal humerus, both 90° and 180° plating options are available, as well as extra-articular plates that go up to 294 mm. Dorsal olecranon and olecranon osteotomy plates complete the portfolio.

Elbow Plating System (AR-7650S)


Product Description	Item Number
2.7 mm Screw Instrumentation	
Drill Guide, 2.0 mm/2.7 mm	AR-18800-22
Drill/Depth Guide, 2.7 mm, long	AR-18800-29
VAL Guide, cone style, 2.7 mm	AR-18800-31
Locking Tower, 2.7 mm	AR-18800-30
Countersink, 2.7 mm	AR-18800-23
Profile Drill, 2.7 mm KreuLock™	AR-18900-03
Bone Tap, 2.7 mm, AO	AR-18800-34
3.5 mm Shaft Hybrid Screw Instrumentation	
3.5 mm/2.5 mm Straight Drill Guide	AR-7650-33
Eccentric Drill Guide, T10	AR-7650-04
Drill Guide, VAL cone/straight	AR-7650-02
Locking Drill Guide, 3.5 mm shaft-hybrid	AR-8933HGL
Profile Drill, 3.5 mm shaft hybrid KreuLock	AR-18900-04
3.5 mm Screw Instrumentation	
Driver Shaft, T15, retaining, AO	AR-7650-09
Driver Shaft, T15, nonretaining, AO	AR-7650-10
Holding Sleeve, 3.5 mm	AR-7650-03
Eccentric Drill Guide	AR-7650-12
Driver Shaft, short, T15	AR-7650-14
PEEK Aiming Guides	
PEEK Aiming Guide - olecranon fracture, left	AR-7650-15-L
PEEK Aiming Guide - olecranon fracture, right	AR-7650-15-R
PEEK Aiming Guide - distal humerus, medial, left	AR-7650-16-L

PEEK Aiming Guide - distal humerus, medial, right	AR-7650-16-R
PEEK Aiming Guide - distal humerus, posterolateral, left	AR-7650-17-L
PEEK Aiming Guide - distal humerus, posterolateral, right	AR-7650-17-R
PEEK Aiming Guide - distal humerus, lateral biplanar, left	AR-7650-18-L
PEEK Aiming Guide - distal humerus, lateral biplanar, right	AR-7650-18-R
PEEK Aiming Guide - extra-articular, left	AR-7650-19-L
PEEK Aiming Guide - extra-articular, right	AR-7650-19-R
PEEK Aiming Guide - true lateral, left	AR-7650-20-L
PEEK Aiming Guide - true lateral, right	AR-7650-20-R
PEEK Aiming Guide - olecranon osteotomy, left	AR-7650-21-L
PEEK Aiming Guide - olecranon osteotomy, right	AR-7650-21-R
General Instrumentation	
Driver Shaft, T10, long	AR-18800-26
Driver Shaft, T10	AR-18800-24
Driver Shaft, T10, nonretaining	AR-18800-25
Short Driver, T10	AR-7650-13
Solid Screwdriver, T10	AR-18800-27
Ratcheting Handle, AO QC	AR-8700RH
Depth Device, 2.7/3.5 mm, short	AR-18800-39
Screw Holding Sleeve, T10	AR-18800-32
Screw Holding Forceps	AR-8941F
Depth Device, 2.7/3.0/3.5	AR-7650-01
Reduction, Bending, and Cutting Tools	
Toothed Kocher Forceps, large	AR-7650-28
Langenbeck/Freer Elevator	AR-7650-29
Cup Curette/Dental Pick	AR-7650-26
Lobster Claw, large	AR-7650-36
Weber Clamp	AR-8943-24
Bending Iron	AR-7650-30
Torque-Limiting Attachment, 1.5 Nm	AR-8963TL-01
French Bending Iron	AR-7650-27
Verbrugge	AR-8943-39
Plate Cutter, straight	AR-8957-06
Compatible Sutures With 22.2 or 17.9 mm Needles	
#0 FiberWire®, 97 cm, with tapered needle	AR-7250
#0 FiberWire, with diamond point needle, 22.2 mm, 1/2 circle, qty. 12	AR-7251
#2-0 FiberWire, 46 cm, with tapered needle, 17.9 mm, 3/8 circle, qty. 12	AR-7220
SutureTape, 0.9 mm, with two 22.2 mm curved needles, (white/blue)	AR-7570-02
Elbow Fracture System Cases	
Elbow Fracture System Instrument Case (a)	AR-7650C-02
Elbow Fracture System Implant Case (b)	AR-7650C-01
Plates	
DH Plate, medial, Ti, left, short	AR-7654ML-S
DH Plate, medial, Ti, left, medium	AR-7654ML-M

Elbow Fracture Plating System (cont.)

DH Plate, medial, Ti, left, long	AR-7654ML-L
DH Plate, medial, Ti, right, short	AR-7654MR-S
DH Plate, medial, Ti, right, medium	AR-7654MR-M
DH Plate, medial, Ti, right, long	AR-7654MR-L
Distal Humerus Biplanar Plates	
DH Plate lateral biplanar, Ti, left, short	AR-7654BLL-S
DH Plate, lateral biplanar, Ti, left, medium	AR-7654BLL-M
DH Plate, lateral biplanar, Ti, left, long	AR-7654BLL-L
DH Plate, lateral biplanar, Ti, right, short	AR-7654BLR-S
DH Plate, lateral biplanar, Ti, right, medium	AR-7654BLR-M
DH Plate, lateral biplanar, Ti, right, long	AR-7654BLR-L
Distal Humerus Posterolateral Plates	
DH Plate, posterolateral, Ti, left, short	AR-7654PL-S
DH Plate, posterolateral, Ti, left, medium	AR-7654PL-M
DH Plate, posterolateral, Ti, left, long	AR-7654PL-L
DH Plate, posterolateral, Ti, right, short	AR-7654PR-S
DH Plate, posterolateral, Ti, right, medium	AR-7654PR-M
DH Plate, posterolateral, Ti, right, long	AR-7654PR-L
Distal Humerus Lateral Plates	
DH Plate, lateral, Ti, left, short	AR-7654LL-S
DH Plate, lateral, Ti, left, medium	AR-7654LL-M
DH Plate, lateral, Ti, left, long	AR-7654LL-L
DH Plate, lateral, Ti, right, short	AR-7654LR-S
DH Plate, lateral, Ti, right, medium	AR-7654LR-M
DH Plate, lateral, Ti, right, long	AR-7654LR-L
Distal Humerus Extraarticular Plates	
DH Plate, EA, Ti, left, 7H, sterile	AR-7654EL-07S
DH Plate, EA, Ti, left, 9H, sterile	AR-7654EL-09S
DH Plate, EA, Ti, left, 11H, sterile	AR-7654EL-11S
DH Plate, EA, Ti, left, 13H, sterile	AR-7654EL-13S
DH Plate, EA, Ti, left, 15H, sterile	AR-7654EL-15S
DH Plate, EA, Ti, left, 17H, sterile	AR-7654EL-17S
DH Plate, EA, Ti, right, 7H, sterile	AR-7654ER-07S
DH Plate, EA, Ti, right, 9H, sterile	AR-7654ER-09S
DH Plate, EA, Ti, right, 11H, sterile	AR-7654ER-11S
DH Plate, EA, Ti, right, 13H, sterile	AR-7654ER-13S
DH Plate, EA, Ti, right, 15H, sterile	AR-7654ER-15S
DH Plate, EA, Ti, right, 17H, sterile	AR-7654ER-17S
Dorsal Olecranon Plates	
Olecranon Plate, Ti, left, 2H	AR-7651L-02
Olecranon Plate, Ti, left, 6H	AR-7651L-06
Olecranon Plate, Ti, left, 10H	AR-7651L-10
Olecranon Plate, Ti, left, 14H	AR-7651L-14
Olecranon Plate, Ti, right, 2H	AR-7651R-02
Olecranon Plate, Ti, right, 6H	AR-7651R-06
Olecranon Plate, Ti, right, 10H	AR-7651R-10
Olecranon Plate, Ti, right, 14H	AR-7651R-14

Olecranon Osteotomy Plates	
Olecranon Osteotomy Plate, Ti, left	AR-7651AL
Olecranon Osteotomy Plate, Ti, right	AR-7651AR
Screws	
Cortical Screw, 2.7 mm × 10 mm-60 mm (2 mm increments), 65 mm, 70 mm, 75 mm, 80 mm	AR-18827-10-80
KreuLock™ Screw, VAL, Ti, 2.7 mm × 10 mm-60 mm (2 mm increments) 65 mm, 70 mm, 75 mm, 80 mm	AR-8927VCL-10-80
Low Profile Cortical Screw, LPS, hybrid, 3.5 mm shaft × 10 mm - 60 mm (2 mm increments)	AR-8933HY-10-60
VAL KreuLock Screw, hybrid, 3.5 mm × 10 mm-60 mm (2 mm increments)	AR-8933HVCL-10-60
Disposables	
Guidewire w/ trocar tip, .062 in	AR-8941K
BB-Tak, MTP	AR-13227
2.7 mm Screw Disposables	
Drill Bit, 2.0 mm, short, AO	AR-18800-18
Drill Bit, 2.0 mm, long, AO	AR-18800-19
Drill Bit, 2.7 mm, AO	AR-18800-20
3.5 mm Shaft Hybrid Screw Disposables	
Drill Bit, 2.5 mm, long, AO	AR-8933HD
Drill Bit, calibrated, 2.5 mm, AO	AR-7650-05
Drill Bit, 3.5 mm, AO	AR-7650-34
Drill Bit, short, 2.5 mm, AO	AR-7650-35



Cardiothoracic

FiberTape® Sternal Closure	92
DynaNite® CT Nitinol Staples Sternal Closure.....	93
Cardiothoracic Biologics.....	93
NanoNeedle™ Thoracic Access Set.....	93
NanoNeedle Scope CT	93

FiberTape® Sternal Closure



The FiberTape Sternal Closure System is a 100% nonmetallic, suture-based alternative to metal wires traditionally used to close the sternum after cardiac surgery. The FiberTape sternal closure implant, which is implanted similarly to metal wires, provides strong, reproducible compression¹ that helps support sternal approximation and stabilization.

Product Description	Item Number
FiberTape Sternal Closure (white/blue), w/ tapered cutting needle	AR-7288
FiberTape Sternal Closure (white/black), w/ conventional cutting needle	AR-7288T
Radiopaque FiberTape Sternal Closure, w/ tapered cutting needle	AR-7288R
FiberTape Sternal Closure (white/blue), w/ blunt needle	AR-7289
FiberTape Sternal Closure (white/black), w/ blunt needle	AR-7289T
Radiopaque FiberTape Sternal Closure, w/ blunt tapered needle	AR-7289R

FiberTape Sternal Closure Implant System Options

Implant systems include FiberTape Sternal Closure sutures, 2 mm, 36", blue with attached conventional cutting needle, a FiberTape cerclage disposable tensioner, 60 mm free cutting needle and MegaLoader


FiberTape Sternal Closure X4	AR-7820SC-4
FiberTape Sternal Closure X5	AR-7820SC-5
FiberTape Sternal Closure X6	AR-7820SC-6
Disposables	
Cerclage Passing Needle, medium, disposable, qty 5	AR-7816
FiberTape Cerclage Disposable Tensioner	AR-7820
Reverse Cutting Spring-Eye Needle, 60 mm	AR-7283
Instrumentation	
FiberTape Cerclage Tensioner, reusable	AR-7800
FiberTape Cerclage Tensioner Handle, reusable	AR-7801
Sternal Needle Driver	AR-5050-04
FiberWire® Scissor	AR-11796
FiberTape Cerclage Tray	AR-7800C

Literature


Product Description	Item Number
FiberTape® Sternal Closure System Technique	LT1-000193-EN-US

Reference

1. Arthrex, Inc. Data on file (APT 5056). Naples, FL; 2021.



Stronger Than Conventional Methods¹
Load-to-Failure Comparison



DynaNite® CT Nitinol Staples Sternal Closure



DynaNite CT nitinol staples used in conjunction with FiberTape® sternal closure to augment sternotomy fixation. These shape-changing, dynamic implants increase sternal stability by producing constant compression across the sternotomy.

Product Description	Item Number
DynaNite CT Nitinol Staple Kit, 13 mm (W) × 10 mm (L), w/ Instrumentation Includes preloaded staple, drill, guide, and tamp	AR-8718CTDS-1310
DynaNite CT Nitinol Staple Kit, 15 mm (W) × 12 mm (L), w/ Instrumentation Includes preloaded staple, drill, guide, and tamp	AR-8718CTDS-1512
DynaNite CT Nitinol Staple Kit, 13 mm × 10 mm	AR-8718CT-1310
DynaNite CT Nitinol Staple Kit, 15 mm × 12 mm	AR-8718CT-1512

Literature

Product Description	Item Number
DynaNite CT Nitinol Staples Sternotomy Closure Surgical Technique	LT1-000340-EN-US

Cardiothoracic Biologics



CT biologic sternal sets combine the benefits of JumpStart® microcurrent and antibacterial technology and growth factor support from Angel® PRP in one convenient kit.

Product Description	Item Number
JumpStart Angel PRP Sternal Set	ABS-4006L-61T
JumpStart Sternal Set	ABS-4006L-T
FlexEFit® Angel PRP Sternal Set, 1.5 × 4.5 in	ABS-4060L-61T
FlexEFit Sternal Set, 1.5 × 4.5 in	ABS-4060L-T
Viscous Spray Applicator w/ 3 cm mixing tip, 5/box	ABS-10051
Blending Connector w/ single spray, 10/box	SA-3674
Viscous Spray II, 5/box	ABS-10052
Dual Spray Tip, 10/box	SA-3660
Dual Spray Procedure Kit, 11:1 ratio, 5/box	SA-4460

NanoNeedle™ Thoracic Access Set



Featuring NanoNeedle Scope visualization, the Thoracic Access Sets combine the PassPort Button™ cannula with auxiliary lumen or the NanoNeedle Scope High-Flow Operative Sheath Kit.

Product Description	Item Number
Thoracic Access Set w/ 125 mm NanoNeedle scope 2.0 and 3 PassPort Button cannulas	AR-0125-P-T
Thoracic Access Set w/ 180 mm NanoNeedle scope 2.0 and 3 PassPort Button cannula	AR-0180-P-T
Thoracic Access Set w/ 250 mm NanoNeedle scope 2.0 and 3 PassPort Button cannulas	AR-0250-P-T

NanoNeedle Scope CT



Product Description	Item Number
NanoNeedle Scope 2.0, 125 mm	AR-3210-0070-T
NanoNeedle Scope 2.0, 180 mm	AR-3210-0071-T
NanoNeedle Scope 2.0, 250 mm	AR-3210-0072-T

Limb Positioning Devices

TRIMANO FORTIS® Support Arm	96
Beach Chair Positioner	97
Shoulder Suspension System (S3)	98
ElbowLOC® Arm Positioning System	98
3-Point Shoulder Distraction System	99
Arm Sleeves.....	99
Lateral Decubitus Shoulder Traction Tower	99
Accessory Items for Limb Positioners	99

TRIMANO FORTIS® Support Arm

The TRIMANO FORTIS support arm acts as the surgical team's "third hand" by securely and safely holding the patient's arm in any desired position during arthroscopic or open shoulder surgery performed in the beach chair position. The compact and lightweight device is easily attached to any OR table Clark rail and is ready for immediate use. No additional power or air connections are required. The patient's arm can be moved in any direction by simply pressing the support arm handle. Releasing the handle locks it into the desired position. The convenient TRIMANO Beach Chair Kit includes a sterile drape for the support arm, an ergonomically designed foam arm holder and Coban™ self-adherent wrap, and allows for quick and easy patient preparation, while gently protecting the operative arm.

Product Description	Item Number
TRIMANO FORTIS Support Arm (a)	AR-1740
TRIMANO FORTIS Adapter (b)	AR-1741
TRIMANO Beach Chair Kit, 6/box (c)	AR-1644
TRIMANO Elbow Kit, 6/box (d)	AR-1646
TRIMANO Drape, 10/box	AR-1648
TRIMANO Wrist Positioner	AR-1647
TRIMANO FORTIS Leg Holder, reusable	AR-1747-B
TRIMANO FORTIS Knee Kit, 6/box	AR-1747

*TRIMANO FORTIS is a registered trademark of Maquet GmbH.
Coban is a registered trademark of 3M.



Beach Chair Positioner

The “Lift-Assist” Beach Chair Positioner System allows unobstructed posterior access and makes repositioning the patient quick and easy. The free-sliding head positioner feature protects the patient’s neck when raising and lowering the table and the fold-away shoulder wings remove completely for unobstructed shear-free access to the operative shoulder. The lift-assist design enables nearly effortless positioning as the piston supports most of the patient’s weight.

Product Description	Item Number
“Lift-Assist” Beach Chair Positioner Includes reusable pad set, intubation pad and plate, beach chair clamps (2), and reusable security strap (a)	AR-1627
Required	
LPS Arm Support, required (d)	AR-1627-01
Counter Traction Strap, required (e)	AR-1627-03
Easy Lock Socket, 2 required (d)	AR-1627-12
Universal Head Positioner w/ handles, required	AR-1627-05
Universal Head Positioner w/o handles (b)	AR-1627-20
Head Positioner Disposable Kit, 10/box, required (c)	AR-1627-06
Beach Chair Cart, recommended	AR-1627-13
Replacement	
Reusable Pad Set (replacement)	AR-1627-09
Intubation Pad and Plate (replacement)	AR-1627-10
Beach Chair Rail Clamp US replacement, qty. 1	AR-1627-11
Reusable Security Strap for AR-1627-01 (replacement)	AR-1627-15
UHP Velcro Strap (replacement), qty. 1	AR-1627-17

Literature

Product Description	Item Number
Lift-Assist Beach Chair Positioner System	FL1-000368-en-US



Shoulder Suspension System (S3)




The Arthrex Shoulder Suspension System for lateral decubitus position combines simple and intuitive functionality, providing the surgeon with optimum arthroscopic visualization for all shoulder pathology. The tower is lightweight and folds into a small vertical footprint for easy storage and setup by the clinical staff. Intraoperative control of abduction, internal/external rotation, and lateral distraction can be easily performed in a sterile fashion, allowing for the utmost in direct control.

Product Description	Item Number
Shoulder Suspension System (S3) includes	AR-1650S
Shoulder Suspension Tower (a)	AR-1650-01
Arm Sleeve Connector, S3 (b)	AR-1650-02
Lateral Traction Sling Connector, S3 (c)	AR-1650-03
Weight Hanger, 10 lb (4.5 kg), S3	AR-1650-04
Weight, 2.5 lb (1.1 kg)	AR-1650-05
Storage Cart, S3, recommended	AR-1650-06
US Clamp, S3	AR-1650-07
Arm Sleeve Disposable, S3, qty. 6 (d)	AR-1651
Lateral Traction Sling Disposable, S3, qty. 6	AR-1652

Literature

Product Description	Item Number
Shoulder Suspension System Brochure	FL1-000322-en-US
Shoulder Suspension System Quick Setup Guide	LR1-00037-EN

ElbowLOC® Arm Positioning System



The Arthrex ElbowLOC Arm Positioning System is an upper-extremity positioning device designed specifically for surgery from the mid-humerus to the fingers in one self-contained sterile system. It is easy to assemble and intuitive to use. The system enables traction across the elbow or wrist and unhindered intraoperative elbow motion and forearm manipulation, depending on the fracture and reduction needed. The system offers the surgeon the choice of supine, lateral, or supine suspended positioning for elbow surgery as well as a wrist tower traction system for wrist or hand surgery.

Product Description	Item Number
ElbowLOC® Arm Positioning System includes	AR-1600S
Clamp Base	AR-1600-1
Supine Positioner	AR-1600-2
Lateral Positioner	AR-1600-3
Wrist Support/Capture	AR-1600-4
Humeral Support/Capture	AR-1600-5
Reducer/Distractor	AR-1600-6
Wrist Tower	AR-1600-7
Traction Bar	AR-1600-8
Traction Bar Capture	AR-1600-9
Traction Gauge	AR-1600-10
Distraction Nut	AR-1600-11
Instrument Case	AR-1600C
Disposables	
Double Finger Trap, small	AR-1600-FTS-S
Double Finger Trap, medium	AR-1600-FTS-M
Double Finger Trap, large	AR-1600-FTS-L
Double Finger Trap, x-large	AR-1600-FTS-X
ElbowLOC Sterile Field kit w/ Foam Tubes and Adherent Wrap	AR-1600-FKS

Literature

Product Description	Item Number
ElbowLOC Positioning System Surgical Technique Guide	LT1-000303-en-US

3-Point Shoulder Distraction System



This versatile system provides safe, effective, and easy positioning of the shoulder during all types of arthroscopic or open shoulder surgery performed in the lateral decubitus position.

Product Description	Item Number
3-Point Shoulder Distraction System	AR-1600M

Arm Sleeves



The lateral traction arm sleeve provides a simple, economical “roll-down” solution for traction in the lateral decubitus position.

The STaR Sleeve™ arm holder is a soft foam traction boot designed to gently cradle the arm, forearm, and wrist during distal distraction in any desired position of abduction.

Product Description	Item Number
Lateral Traction Arm Sleeve, sterile, qty. 5 (a)	AR-1635
STaR Sleeve, compact, Coban, sterile, qty. 6 (b)	AR-1606C
STaR Sleeve, Coban, sterile, qty. 6	AR-1606
STaR Sleeve, velcro, sterile, qty. 6	AR-1606V
STaR Sleeve, velcro, extra long, sterile, qty. 6	AR-1606LV
STaR Sleeve, velcro, small, sterile, qty. 6	AR-1606SV
Cushioned Lateral Traction Sling for Instability, qty. 6	AR-1603
Atraumatic Hand Holder Traction Attachment	AR-1602D


Lateral Decubitus Shoulder Traction Tower



This traction tower features a single boom arm that allows for abduction, forward flexion, and traction during open and arthroscopic procedures.

Product Description	Item Number
Lateral Decubitus Shoulder Traction Tower	AR-1630

Accessory Items for Limb Positioners



Product Description	Item Number
Optional Accessories	
Traction Scale Attachment (a)	AR-1604
Weight Hanger, 10 lb (4.5 kg) (b)	AR-1650-04
Weight, 2.5 lb (1.1 kg) (c)	AR-1650-05
Storage Hook and Wall Mount	AR-1605S
Storage Stand for Limb Positioners (d)	AR-1600SS
Replacement Parts	
White Plastic Pulley Wheel	AR-16000602
Pulley Post Assembly for AR-16000602	AR-16000510
Lateral Traction Hub and Handle for AR-1600M	AR-16001400
handle for AR-1600M Lateral Traction Hub	AR-16001408
Extension Arm Elevation Shaft Cap	AR-16001607
Small Black Knob for AR-1600M	AR-1600A-01
Base and Small Knobs Assembly for AR-1600M	AR-16001500
Base Assembly 5/16-18 screw for AR-16001500	AR-16110407
Large Black Knob	AR-1600A-02
Swivel Lock Base for Large Knob Assembly for AR-1600M	AR-16001200


Anchor Charts

Rotator Cuff Anchor Chart.....	102
Instability Anchor Chart.....	108
Knotless Tensionable Anchor Chart.....	110
Soft Anchors.....	112

Rotator Cuff Anchor Chart

	Dimensions	Anchor	Material	Suture	Item Number
	3.9 mm × 11.2 mm	Knotless Corkscrew®	β-TCP	w/ #5 FiberWire® CL	AR-1941BC
	3.9 mm × 11.2 mm	Knotless Corkscrew	PEEK	w/ #5 FiberWire CL	AR-1941PSV
	3.9 mm × 17.9 mm	SwiveLock C	β-TCP/PEEK	Closed eyelet	AR-2326BCC
	3.9 mm × 17.9 mm	SwiveLock C	PEEK	Closed eyelet	AR-2326PSLC
	4.5 mm × 14 mm	Corkscrew FT	β-TCP	w/ two #2 FiberWire	AR-1927BCF-45
	4.5 mm × 14 mm	Corkscrew FT	β-TCP	w/ two #2 TigerTail®	AR-1927BCFT-45
	4.5 mm × 14 mm	Corkscrew FT	β-TCP	w/ two #2 FiberWire and 4 needles	AR-1927BCNF-45
	4.5 mm × 14 mm	Corkscrew FT	PEEK	w/ two #2 FiberWire and 4 needles	AR-1927PNF-45
	4.5 mm × 14 mm	Corkscrew FT	PEEK	w/ two #2 FiberWire	AR-1927PSF-45
	4.5 mm × 14 mm	Corkscrew FT	Titanium	w/ two #2 FiberWire	AR-1928SF-45
	4.75 mm × 14 mm	Corkscrew FT	β-TCP	w/ three #2 TigerTail	AR-1927BCF-475
	4.75 mm × 14 mm	Corkscrew FT	β-TCP	w/ two 1.3 mm SutureTapes (wht/blk, wht/blu)	AR-1927BCT-475
	4.75 mm × 14 mm	Corkscrew FT	PEEK	w/ three #2 TigerTail	AR-1927PSF-475
	4.75 mm × 14 mm	Corkscrew FT	PEEK	w/ two 1.3 mm SutureTapes (wht/blk, wht/blu)	AR-1927PST-475
	4.75 mm × 14.7 mm	Knotless Corkscrew	β-TCP	w/ two #2 knotless repair sutures	AR-1943BC
	4.75 mm × 14.7 mm	Knotless Corkscrew	PEEK	w/ two #2 knotless repair sutures	AR-1943PS
	4.5 mm × 24 mm	PushLock®	β-TCP/PEEK	Closed eyelet	AR-1922BC
	4.5 mm × 24 mm	PushLock	PEEK	Closed eyelet	AR-1922PS
	4.5 mm × 28 mm	PushLock SP	β-TCP/Titanium	Self-punching	AR-1922BCM
	4.5 mm × 28 mm	PushLock SP	PEEK/Titanium	Self-punching	AR-1922PSM
	4.75 mm × 19.1 mm	Knotless SwiveLock C	β-TCP/PEEK	w/ #2 suture (blu) and SutureTape Link (wht/blk)	AR-2324KBCC
	4.75 mm × 19.1 mm	Knotless SwiveLock C	β-TCP/PEEK	w/ FiberTape® Loop (blu), #2 suture (wht/blk), and SutureTape Link (wht/blu)	AR-2324KBCCT
	4.75 mm × 19.1 mm	Knotless SwiveLock C	β-TCP/PEEK	w/ TigerTape™ Loop (wht/blk), #2 suture (blu), and SutureTape Link (wht/blk)	AR-2324KBCCTT
	4.75 mm × 19.1 mm	SwiveLock C	β-TCP/PEEK	Closed Eyelet	AR-2324BCC
	4.75 mm × 19.1 mm	SwiveLock C	β-TCP/PEEK	Closed Eyelet w/ FiberTape® Loop (blu)	AR-2324BCCT
	4.75 mm × 19.1 mm	SwiveLock C	β-TCP/PEEK	Closed Eyelet w/ TigerTape™ Loop (wht/blk)	AR-2324BCCTT
	4.75 mm × 22 mm	SwiveLock C	β-TCP/PEEK	Closed Eyelet w/ two #2 TigerWire® CL	AR-2324BCC-2
	4.75 mm × 22 mm	SwiveLock C	β-TCP/PEEK	Closed Eyelet w/ two #2 TigerTail (wht/blk, blu/blk)	AR-2324BCT-2
	4.75 mm × 19.1 mm	Knotless SwiveLock C	PEEK	w/ #2 suture (blu) and SutureTape Link (wht/blk)	AR-2324KPSLC
	4.75 mm × 19.1 mm	Knotless SwiveLock C	PEEK	w/ FiberTape Loop (blu), #2 suture (wht), SutureTape Link (wht/blu)	AR-2324KPCT
	4.75 mm × 19.1 mm	Knotless SwiveLock C	PEEK	w/ TigerTape Loop (wht/blk), #2 suture (blu), SutureTape Link (wht/blk)	AR-2324KPCTT


	Punch (required)	Punch/Tap (for hard bone)	Drill	Other
	AR-1941P (a)	AR-1941PTB	AR-1941D	AR-1941DG (b) AR-1941DS AR-1941PK
	AR-2326P (c)	AR-2326PTB	AR-2326D	N/A
	AR-1922P (d) or AR-1927PBS (disposable)	AR-1927PTB-45 (e) or AR-1927PTB-45-2	N/A	N/A
	Metal anchors do not require instrumentation			
		AR-1927PTB-475 (required)	N/A	N/A
	AR-1927KPB-475	N/A	N/A	N/A
	AR-1922P (f) or AR-1922PBS (disposable)	N/A	N/A	N/A
	Self-punching PushLock® anchors do not require instrumentation			
	AR-1927PB (g) or AR-1927PBS (disposable)	AR-2324PTB	AR-1927D	N/A



Rotator Cuff Anchor Chart (Cont.)

	Dimensions	Anchor	Material	Suture	Item Number
	4.75 mm × 19.1 mm	SwiveLock® C	PEEK	Closed eyelet	AR-2324PSLC
	4.75 mm × 19.1 mm	SwiveLock C	PEEK	Closed eyelet w/ FiberTape® Loop (blu)	AR-2324PCT
	4.75 mm × 19.1 mm	SwiveLock C	PEEK	Closed eyelet w/ TigerTape™ Loop (wht/blk)	AR-2324PCTT
	4.75 mm × 22 mm	SwiveLock C	PEEK	Closed eyelet w/ two #2 TigerTail® (wht/blk, blu/blk)	AR-2324PSLC-2
	4.75 mm × 19.1 mm	SwiveLock C	PLLA/PEEK	Closed eyelet	AR-2324BSLC
	4.75 mm × 24.5 mm	Knotless SwiveLock SP	β-TCP/ PEEK	w/ #2 suture (blu), self-punching PEEK eyelet	AR-2324KBCSP
	4.75 mm × 24.5 mm	SwiveLock SP	β-TCP/ PEEK	w/ 1.3 mm SutureTape (wht/blu), self-punching PEEK eyelet	AR-2324BCSP
	4.75 mm × 24.5 mm	SwiveLock SP	β-TCP/Titanium	Self-punching	AR-2324BCM
	4.75 mm × 24.5 mm	Knotless SwiveLock SP	PEEK	w/ #2 suture (blu), self-punching PEEK eyelet	AR-2324KPSP
	4.75 mm × 24.5 mm	SwiveLock SP	PEEK	w/ 1.3 mm SutureTape (wht/blu), self-punching PEEK eyelet	AR-2324PSP
	4.75 mm × 24.5 mm	SwiveLock SP	PEEK/Titanium	Self-punching	AR-2324PSLM
	4.75 mm × 24.5 mm	SwiveLock SP	Titanium	Self-punching	AR-2324SLM
	5.5 mm × 14.7 mm	Corkscrew® FT	β-TCP	w/ two 1.3 mm SutureTapes (wht/blk, wht/blu)	AR-1927BCT
	5.5 mm × 14.7 mm	Corkscrew FT	β-TCP	w/ three 1.3 mm SutureTapes (wht/blu, blk/wht, blu/blu)	AR-1927BCT-3
	5.5 mm × 14.7 mm	Corkscrew FT	β-TCP	w/ two #2 FiberWire®	AR-1927BCF
	5.5 mm × 14.7 mm	Corkscrew FT	β-TCP	w/ three #2 FiberWire	AR-1927BCF-3
	5.5 mm × 14.7 mm	Corkscrew FT	β-TCP	w/ two #2 TigerTail®	AR-1927BCFT
	5.5 mm × 14.7 mm	Corkscrew FT	β-TCP	w/ two #2 FiberWire and 4 needles	AR-1927BCNF
	5.5 mm × 14.7 mm	Corkscrew FT	PEEK	w/ two 1.3 mm SutureTapes	AR-1927PST
	5.5 mm × 14.7 mm	Corkscrew FT	PEEK	w/ three 1.3 mm SutureTapes	AR-1927PST-3
	5.5 mm × 14.7 mm	Corkscrew FT	PEEK	w/ two 1.3 mm SutureTapes and 4 needles	AR-1927PNT
	5.5 mm × 14.7 mm	Corkscrew FT	PEEK	w/ two #2 FiberWire and 4 needles	AR-1927PNF
	5.5 mm × 14.7 mm	Corkscrew FT	PEEK	w/ two #2 FiberWire	AR-1927PSF
	5.5 mm × 14.7 mm	Corkscrew FT	PEEK	w/ three #2 FiberWire	AR-1927PSF-3
	5.5 mm × 14.7 mm	Corkscrew FT	PLLA	w/ two #2 FiberWire	AR-1927BF
	5.5 mm × 16.3 mm	Corkscrew FT	Titanium	w/ two 1.3 mm SutureTapes	AR-1928ST
	5.5 mm × 16.3 mm	Corkscrew FT	Titanium	w/ three 1.3 mm SutureTapes	AR-1928ST-3
	5.5 mm × 16.3 mm	Corkscrew FT	Titanium	w/ two 1.3 mm SutureTapes and 4 needles	AR-1928STN
	5.5 mm × 16.3 mm	Corkscrew FT	Titanium	w/ two #2 FiberWire	AR-1928SF-2
	5.5 mm × 16.3 mm	Corkscrew FT	Titanium	w/ three #2 FiberWire	AR-1928SF-3
	5.5 mm × 16.3 mm	Corkscrew FT	Titanium	w/ two #2 TigerTail	AR-1928SFT-2
	5.5 mm × 16.3 mm	Corkscrew FT	Titanium	w/ two #2 FiberWire and 4 needles	AR-1928SNF-2

	Punch (required)	Punch/Tap (for hard bone)	Drill	Other
	AR-1927PB (a) or AR-1927PBS (disposable)	AR-2324PTB	AR-1927D	N/A
	Self-punching SwiveLock® anchors do not require instrumentation			
	AR-1927PB or AR-1927PBS (disposable)	AR-1927CTB (b)	AR-1927D	N/A
	Titanium Corkscrew® FT anchors do not require instrumentation			



Rotator Cuff Anchor Chart (Cont.)

	Dimensions	Anchor	Material	Suture	Item Number
	5.5 mm × 19.1 mm	Knotless SwiveLock®	β-TCP/PEEK	w/ #2 suture (blu) and SutureTape Link (wht/blk)	AR-2323KBCC
	5.5 mm × 19.1 mm	SwiveLock® C	β-TCP/PEEK	Closed eyelet	AR-2323BCC
	5.5 mm × 22 mm	SwiveLock C	β-TCP/PEEK	w/ two #2 TigerTail® (wht/blk, blu/blk)	AR-2323BCT-2
	5.5 mm × 19.1 mm	Knotless SwiveLock	PEEK	w/ #2 suture (blu) and SutureTape Link (wht/blk)	AR-2323KPSLC
	5.5 mm × 19.1 mm	SwiveLock C	PEEK	Closed eyelet	AR-2323PSLC
	5.5 mm × 22 mm	SwiveLock C	PEEK	w/ two #2 TigerTail (wht/blk, blu/blk)	AR-2323PSLC-2
	5.5 mm × 24.5 mm	Knotless SwiveLock SP	β-TCP/ PEEK	w/ #2 suture (blu) and self-punching PEEK eyelet	AR-2323KBCCSP
	5.5 mm × 24.5 mm	SwiveLock SP	β-TCP/ PEEK	w/ 1.3 mm SutureTape (wht/blu), self-punching PEEK eyelet	AR-2323BCSP
	5.5 mm × 24.5 mm	SwiveLock SP	β-TCP/Titanium	Self-punching	AR-2323BCM
	5.5 mm × 24.5 mm	Knotless SwiveLock SP	PEEK	w/ #2 suture (blu) and self-punching PEEK eyelet	AR-2323KPSP
	5.5 mm × 24.5 mm	SwiveLock SP	PEEK	w/ 1.3 mm SutureTape (wht/blu), self-punching PEEK eyelet	AR-2323PSP
	5.5 mm × 24.5 mm	SwiveLock SP	PEEK/Titanium	Self-punching	AR-2323PSLM
	6.5 mm × 14.7 mm	Corkscrew® FT	β-TCP	w/ two #2 FiberWire®	AR-1927BCF-65
	6.5 mm × 14.7 mm	Corkscrew FT	PEEK	w/ two #2 FiberWire	AR-1927PSF-65
	6.5 mm × 16.3 mm	Corkscrew FT	Titanium	w/ three #2 FiberWire	AR-1929SF-3


The 6.5 mm anchors are intended for backup; no instrumentation required.

	Punch (required)	Punch/Tap (for hard bone)	Drill	Other
	AR-1927PB or AR-1927PBS (disposable)	AR-1927CTB	AR-1927D	N/A
	Self-punching SwiveLock® anchors do not require instrumentation			
	The 6.5 anchors are intended for backup, no instrumentation required.			

Instability Anchor Chart

	Dimensions	Anchor	Material	Suture	Item Number
	2.4 mm × 12 mm	SutureTak	β-TCP	w/ two #2 FiberWire	AR-1934BCF-24-2
	2.4 mm × 11.3 mm	PushLock®	β-TCP	Use w/ #1 FiberWire, #2 FiberLink™ or 1.3 mm SutureTape	AR-2922BC
	2.4 mm × 11.3 mm	PushLock	PEEK	Use w/ #1 FiberWire, #2 FiberLink, or 1.3 mm SutureTape	AR-2922PS
	2.9 mm × 12.5 mm	PushLock (short)	β-TCP		AR-2923BC
	2.9 mm × 12.5 mm	PushLock (short)	PEEK		AR-2923PS
	3 mm × 12.7 mm	Knotless SutureTak	β-TCP	w/ #2 FiberWire	AR-1936BC
	3 mm × 12.7 mm	Knotless SutureTak	β-TCP	w/ #2 FiberWire CL suture	AR-1938BC
	3 mm × 12.7 mm	Knotless SutureTak	PEEK	w/ #2 FiberWire	AR-1936PS
	3 mm × 12.7 mm	Knotless SutureTak	PEEK	w/ #2 FiberWire	AR-1938PS
	3 mm × 14.5 mm	SutureTak	β-TCP	w/ 1.3 mm SutureTape (wht/blu)	AR-1934BCT
	3 mm x 14.5 mm	SutureTak	β-TCP	w/ two 1.3 mm SutureTape (blu, blk/wht)	AR-1934BCT-2
	3 mm × 14.5 mm	SutureTak	β-TCP	w/ two #2 FiberWire	AR-1934BCF-2
	3 mm x 14.5 mm	SutureTak	PEEK	w/ two 1.3 mm SutureTape (blu, blk/wht)	AR-1934PST-2
	3 mm × 12.2 mm	SutureTak	PEEK	w/ two #2 FiberWire	AR-1934PS-2
	3.5 mm × 15.8 mm	SwiveLock® C	β-TCP		AR-2325BCC
	3.5 mm × 15.8 mm	SwiveLock C	PEEK		AR-2325PSLC
	3.5 mm × 19.5 mm	PushLock	β-TCP		AR-1926BC
	3.5 mm × 19.5 mm	PushLock	PEEK		AR-1926PS
	3.5 mm × 19.5 mm	PushLock	PLLA		AR-1926B
	3.9 mm x 11.2 mm	Knotless Corkscrew®	β-TCP	w/ #5 suture	AR-1941BC
	3.9 mm x 11.2 mm	Knotless Corkscrew	PEEK	w/ #5 suture	AR-1941PSV






	Drill Options	Spear/Guide Options (see page 24)	Disposable Options	Perc Kits	Shaver Drill
	AR-1934D-24 AR-1934D-24-1 (hard bone) AR-1934D-24-2 (very hard bone)	AR-1948 (a) AR-1948R AR-1948CT	AR-1934-24DS AR-1945S AR-1934GS	AR-1934PI	N/A
	AR-2922D-24-2 (sterile) AR-2922D-24-3 (hard bone, sterile)	AR-1949 (b) AR-1934R AR-1946 AR-1317	AR-1949S	AR-2922PK	AR-2922DSR-24 AR-2922DTSR-24
	AR-2923D AR-2923D-ST (sterile) AR-2923DT (hard bone) AR-2923DT-ST (hard bone, sterile)	AR-1949 AR-1946 AR-1317 AR-1949	AR-2923DS AR-1949S AR-1934GS AR-1949S	AR-1923PK AR-1923PK	AR-2923DSR AR-2923DSR
	AR-1938D AR-1938DL (very hard bone, sterile)	AR-1949	AR-1938DS AR-1949S	AR-1938PK	AR-1938DSR AR-1938DSSK
	AR-1250LT AR-1257 (spade tip)	AR-1949 AR-1934R AR-1946 AR-1317	AR-1934DS-2 AR-1949S AR-1934GS	AR-1934PI-30	AR-1934DSSK AR-1250LTSR
	AR-2325D AR-2325D-ST (sterile)	AR-1907 (c) AR-1909R AR-1906	N/A	N/A	N/A
	AR-1912	AR-1907 AR-1909R AR-1906	AR-1926DS AR-1926DS-2	N/A	N/A
	AR-1941D	AR-1941DG (d) AR-1941DGF AR-1941CT	AR-1941DS	AR-1941PK	N/A



Knotless Tensionable Anchor Chart

Knotless FiberTak®	Dimensions	Anchor	Material	Suture	Item Number
	N/A	Knotless 1.8 FiberTak® Gen2	Suture	w/ #2 suture	AR-3636
	N/A	Knotless 2.6 FiberTak	Suture	w/ #5 suture	AR-3641
	N/A	Self-Punching Knotless 2.6 FiberTak	Suture	w/ #5 suture	AR-3641SP
	N/A	Self-Punching Knotless 2.6 FiberTak	Suture	w/ two #2 knotless repair suture (wht/blu, wht/blk)	AR-3642SP
	N/A	Self-Punching Knotless 2.6 FiberTak RC	Suture	w/ 1.7 mm TigerTape™ loop (wht/blk) and #2 suture (blu)	AR-3653SP
	N/A	Self-Punching Knotless 2.6 FiberTak RC	Suture	w/ 1.7 mm FiberTape® loop (blu) and #2 suture (wht/blk)	AR-3653TSP
	N/A	Self-Punching Knotless 2.6 FiberTak RC	Suture	w/ 1.7 mm FiberTape loop (blk/blu) and #2 suture (blu/blk)	AR-3653TTSP
Knotless SutureTak®	Dimensions	Anchor	Material	Suture	Item Number
	3 mm × 12.7 mm	Knotless SutureTak®	β-TCP	w/ #2 suture	AR-1936BC
	3 mm × 12.7 mm	Knotless SutureTak	PEEK	w/ #2 suture	AR-1936PS
	3 mm × 12.7 mm	Knotless SutureTak®	β-TCP	w/ #2 suture	AR-1938BC
	3 mm × 12.7 mm	Knotless SutureTak	PEEK	w/ #2 suture	AR-1938PS
Knotless Corkscrew®	Dimensions	Anchor	Material	Suture	Item Number
	3.9 mm × 11.2 mm	Knotless Corkscrew®	β-TCP	w/ #5 FiberWire® CL	AR-1941BC
	3.9 mm × 11.2 mm	Knotless Corkscrew	PEEK	w/ #5 FiberWire CL	AR-1941PSV
	4.75 mm × 14 mm	Knotless Corkscrew	β-TCP	w/ #5 suture	AR-1927KBC
	4.75 mm × 14 mm	Knotless Corkscrew	PEEK	w/ #5 suture	AR-1927KPS
Knotless SwiveLock®	Dimensions	Anchor	Material	Suture	Item Number
	4.75 mm × 24.5 mm	Knotless SwiveLock® SP	β-TCP/ PEEK	w/ #2 suture (blu), self-punching PEEK eyelet	AR-2324KBCSP
	4.75 mm × 19.1 mm	Knotless SwiveLock C	β-TCP/ PEEK	w/ #2 suture (blu) and SutureTape Link (wht/blu)	AR-2324KBCC
	4.75 mm × 19.1 mm	Knotless SwiveLock C	β-TCP/ PEEK	w/ FiberTape loop (blu), #2 suture (wht/blk), SutureTape Link (wht/blu)	AR-2324KBCCT
	4.75 mm × 19.1 mm	Knotless SwiveLock C	β-TCP/ PEEK	w/ TigerTape loop (wht/blk), #2 suture (blu), and SutureTape Link (wht/blk)	AR-2324KBCCTT
	4.75 mm × 24.5 mm	Knotless SwiveLock SP	PEEK	w/ #2 suture (blu), self-punching PEEK eyelet	AR-2324KPSP
	4.75 mm × 19.1 mm	Knotless SwiveLock C	PEEK	w/ #2 suture (blu) and SutureTape Link (wht/blk)	AR-2324KPSLC
	4.75 mm × 19.1 mm	Knotless SwiveLock C	PEEK	w/ FiberTape loop (blu), #2 suture (wht), SutureTape Link (wht/blu)	AR-2324KPCT
	4.75 mm × 19.1 mm	Knotless SwiveLock C	PEEK	w/ TigerTape loop (wht/blk), #2 suture (blu), and SutureTape Link (wht/blk)	AR-2324KPCTT
	5.5 mm × 24.5 mm	Knotless SwiveLock SP	β-TCP/ PEEK	w/ #2 suture (blu), self-punching PEEK eyelet	AR-2323KBCSP
	5.5 mm × 19.1 mm	Knotless SwiveLock C	β-TCP/ PEEK	w/ #2 suture (blu) and SutureTape Link (wht/blu)	AR-2323KBCC
	5.5 mm × 24.5 mm	Knotless SwiveLock SP	PEEK	w/ #2 suture (blu), self-punching PEEK eyelet	AR-2323KPSP
	5.5 mm × 19.1 mm	Knotless SwiveLock C	PEEK	w/ #2 suture (blu) and SutureTape Link (wht/blk)	AR-2323KPSLC



Punch/Tap	Drill Options	Spear/Guide Options (see page 24)	Disposable Options	Perc Kits	ShaverDrill
N/A	AR-3610ND-2 (1.8 mm, flexible w/ hub, sterile) AR-3600D-2 (1.8 mm, rigid, sterile) AR-3600D-4 (1.9 mm, rigid, sterile) AR-3610ND-4 (1.9 mm, flexible w/ hub, sterile)	AR-3610F AR-3610CT (a) AR-3610ST AR-3610CTC AR-3610CTC-2	AR-3638DC AR-3638DS AR-3610DC-3	AR-3610PK-3	AR-3610NSD-2
AR-3656	AR-3657	AR-3655 (b) AR-1941CT AR-1941DGF	AR-3650DS	N/A	AR-3657SD
Punch/Tap	Drill Options	Spear/Guide Options (see page 24)	Disposable Options	Perc Kits	ShaverDrill
N/A	AR-1938D (sterile) AR-1938DL (very hard bone, sterile)	AR-1949 (c)	AR-1938DS	AR-1938PK	AR-1938DSR AR-1938DSSK
Punch/Tap	Drill Options	Spear/Guide Options (see page 24)	Disposable Options	Perc Kits	ShaverDrill
AR-1941P AR-1941PTB	AR-1941D (sterile)	AR-1941DG AR-1941DGF (d) AR-1941CT	AR-1941DS	AR-1941PK	N/A
AR-1927KPB-475	AR-1927KPTB-475	N/A	N/A	N/A	N/A
Punch/Tap	Drill Options	Spear/Guide Options (see page 24)	Disposable Options	Perc Kits	ShaverDrill™
AR-1927PB (e) AR-2324PTB	AR-1927D	N/A	AR-1927PBS	N/A	N/A
AR-1927P AR-1927CTB	AR-1927D	N/A	AR-1927PBS	N/A	N/A



Soft Anchors

For Glenoid Repair	Anchor	Suture	Item Number
	FiberTak®	w/ #2 suture	AR-3600
	FiberTak	w/ two #2 sutures	AR-3600-2
	FiberTak	w/ 1.3 mm SutureTape	AR-3602
	FiberTak	w/ two 1.3 mm SutureTape	AR-3602-2
	Knotless 1.8 FiberTak	w/#2 suture	AR-3636
For Rotator Cuff Repair	Anchor	Suture	Item Number
	Self-Punching 2.6 FiberTak RC	w/ 1.7 mm TigerTape™ Loop (wht/blk) and 1.3 mm SutureTape (wht/blu)	AR-3652SP
	Self-Punching 2.6 FiberTak RC	w/ 1.7 mm FiberTape® Loop (blu) and 1.3 mm SutureTape (wht/blk)	AR-3652TSP
	Self-Punching 2.6 FiberTak RC	w/ 1.7 mm FiberTape Loop (blk/blu) and 1.3 mm SutureTape (blu)	AR-3652TTSP
	Self-Punching Knotless 2.6 FiberTak RC	1.7mm TigerTape Loop (wht/blk) and #2 suture (blu)	AR-3653SP
	Self-Punching Knotless 2.6 FiberTak RC	1.7 mm FiberTape Loop (blu) and #2 suture (wht/blk)	AR-3653TSP
	Self-Punching Knotless 2.6 FiberTak RC	1.7 mm FiberTape Loop (blu/blk) and #2 suture (blk/blu)	AR-3653TTSP
	Knotless 2.6 FiberTak	w/ #5 suture	AR-3641
	Self-Punching Knotless 2.6 FiberTak	w/ #5 suture	AR-3641SP
	Self-Punching Knotless 2.6 FiberTak	w/ two #2 knotless repair suture	AR-3642SP
	Self-Punching 2.6 FiberTak	w/ two 1.3 mm SutureTapes (wht/blu, wht/blk)	AR-3632SP
	Self-Punching 2.6 FiberTak	w/ three 1.3 mm SutureTapes (wht/blu, blk/wht, blu)	AR-3633SP

Punch/ Tap	Drill Options	Spear/Guide Options (see page 24)	Disposable Options	Perc Kits	ShaverDrill
N/A	AR-3600D-1 (1.6 mm, rigid, sterile) AR-3610ND-1 (1.6 mm, flexible w/hub, sterile) AR-3600D-3 (1.7 mm, rigid, sterile) AR-3610ND-3	AR-1948R AR-3610F (a) AR-3610CT AR-3610ST AR-3610CTC AR-3610CTC-2	AR-3600D AR-3600DC	AR-3610PK-1 AR-6548PK	AR-3600NDSR-1
N/A	AR-3610ND-2 (1.8 mm, flexible w/hub, sterile) AR-3600ND-2 (1.8 mm, flexible, sterile) AR-3600D-2 (1.8 mm, rigid, sterile) AR-3610ND-4 (1.9 mm, flexible w/hub, sterile) AR-3600ND-4 (1.9, flexible, sterile) AR-3600D-4 (1.9 mm, rigid, sterile)	AR-3610F AR-3610CT AR-3610ST AR-3610CTC AR-3610CTC-2	AR-3638DC AR-3638DS AR-3610DC-3	AR-3610PK-3 AR-6548PK	AR-3610NSD-2
Punch/ Tap	Drill Options	Spear/Guide Options (see page 24)	Disposable Options	Perc Kits	ShaverDrill™
AR-3656	AR-3657 (2.6 mm)	AR-3655 (b) AR-1941CT AR-1941DGF	AR-3650DS	N/A	AR-3657SD





This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience, and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level and/or outcomes.

View U.S. patent information at www.arthrex.com/corporate/virtual-patent-marking

arthrex.com

© 2026-02 Arthrex, Inc. All rights reserved. LB1-0220-EN_AE