English

A. DEVICE DESCRIPTION

Arthrex saw blades made of high-grade stainless steel are optimally suited for their specific surgical use. They are designed, constructed and produced with utmost care. These quality instruments assure best working results and long service life provided they are used in the proper manner. Therefore, the following instructions for use and safety recommendations have to be observed.

Improper use may damage tissue and may lead to premature wear, destruction of the instruments and injury to the operator, patient or other persons.

B. PROPER USE

- All saw blades are to be used in the appropriate handpieces/power systems only.
- Before use, make absolutely sure that the saw blades to be used and their packaging are in sterile and in technically sound condition. Discard any damaged saw blades and sterile packaging.
- Reciprocating saw blades should be inserted as deeply as possible. Saw blades which are not inserted to their full chucking length cause high centrifugal forces and vibrations.
- Check instruments for their perfect fit prior to use.
- The saw blades are not suitable for work on metal (e.g. implants).
- The saw blades should be actuated after insertion in the saw guide and before contact is made with the bone.
- When used in combination with templates, the cutting thickness of the saw blade should not exceed the thickness of the template. Ensure the blade moves freely in the template before it is used on the patient.
- Jamming, leverage and bending of the saw blade must be avoided while in the template. Otherwise the saw blade or template will be overheated which leads to jamming of the saw blade or thermal necrosis, or even the risk of fracture.
- Contact of the saw blade with the template, guiding device or other metallic objects has to be avoided. Otherwise, damage to the saw blade, template and guiding device can occur as well as blunting of the blades in addition to injury to the bone and surrounding tissue.

Improper use leads to unacceptable results and increased risk of injury.

C. INDICATIONS

Arthrex saw blades have been developed for working on bone in the field of orthopaedics. The operator is responsible for choosing the appropriate saw blade according to the intended use.

D. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

E. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.

F. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.

- An internal fixation device must never be re-used.
- Do not re-sterilize this device.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

G. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

- In CE accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE accepting Countries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

M. CONTACT PRESSURE

Excessive contact pressure has to be avoided because this might lead to

- damage to the blade surfaces
- chipping of the blades
- reduced service life of saw blades
- excessive heat generation
- thermal necrosis

Excessive contact pressure may lead to overheating which in turn results in thermal necrosis; nicks in the blades may cause undesirable rough surfaces. At the worst, the instrument may break.

N. COOLING

Sufficient cooling should be provided when using saw blades. In case of insufficient cooling, the saw blade teeth will get clogged with bone chips, leading to increased heat generation. At worst, this may cause irreversible damage to the bone (thermal necrosis) and reduce service life of the saw blades considerably.

O. SAFETY AND LIABILITY

- The user is responsible for checking the product prior to use to ensure that it is suitable for the intended purpose.
- The user is responsible for the use of the instruments.
- In case of contributory negligence by the user, Arthrex partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.
- Do not regrind saw blades.**
- Multiple use can lead to the fracture of the instrument and constitute an increased risk for the patient and the operator!**

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Arthrex, Inc.
Arthrex GmbH