English

A. DESCRIPTION

The Apollo^{RF} probes are designed specifically for use with the Synerg^{RF}™ (Radio Frequency) Console only and are not compatible with any other electro-surgical generator.

Non-Aspirating RF probes provide a bipolar electro-surgical effect to target tissue. Located on the top of the RF probe handle, the black button adjusts the default power setting for the specific RF probe connected to the Synerg^{RF} Console. The yellow button invokes the Ablate function. The blue button that is farthest from the end of the RF probe invokes the Coagulation function. A Synerg^{RF} Footswitch is connected to the front panel of the Synerg^{RF} Console to override the RF probe’s control buttons. The user has the option to turn off the footswitch override through the Synerg^{RF} Console touch screen options. **Aspirating RF probes** provide a bipolar electro-surgical effect with control buttons as described above and can aspirate materials away from the surgical field.

The aspirating Apollo^{RF} probes come equipped with a suction connector capable of integrating with a standard luer lock syringe. The connector allows a surgeon to use a syringe to loosen any clots that the probes could experience under heavy use.

B. INDICATIONS

The RF probe are accessories to the Synerg^{RF} Console and are intended for use as a complete system in the resection, ablation, and coagulation of soft tissue and hemostasis of blood vessels and tissue in arthroscopic and orthopedic procedures. Specifically, the RF probes, Synerg^{RF} Console and their accessories are used for arthroscopic surgery of the shoulder, wrist, hand, elbow, hip, knee, foot and ankle.

C. CONTRAINDICATIONS

- Use of an RF probe in any non-arthroscopic surgical procedure.
- Use of an RF probe during an arthroscopic procedure without a conductive irrigation solution.
- Use of an RF probe in a surgical procedure for a patient where an arthroscopic procedure is contraindicated for any reason.
- Use of an RF probe in patients with cardiac pacemakers or other electronic device implants.

D. ADVERSE EFFECTS

- As a consequence of electro-surgery, damage to surrounding tissue through iatrogenic injury could occur.
- Additional damage to tissue from excessive force during blunt dissection is possible.
- Premature tip wear may be caused by vigorous use against bony surfaces.
- Thermal damage to joint tissue or dermal burns around portal incisions are possible, if fluid flow is not continuous.

E. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate sized device based on specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Prior to initial use, ensure that all package inserts, including warnings, precautions, *User’s Guides* and Directions for Use are read and understood.
- The RF probes are supplied sterile and intended for SINGLE USE ONLY. DO NOT clean, re-sterilize, or reuse the RF probe.
- DO NOT use excessive force when inserting or removing the RF probe. DO NOT insert the RF probe into an obstructed passageway. Patient injury and/or product damage may occur.
- DO NOT use an RF probe as a lever to enlarge the surgical site or gain access to tissue, which may result in a bent or damaged RF probe.
- Always keep the tip of the active RF probe completely immersed in conductive irrigation fluid (saline, Ringers’ lactate, etc.), regardless of the type of arthroscopic surgery being performed. DO NOT use pre-warmed conductive irrigation fluids as this may result in tissue damage or thermal injury.
- A continuous flow of irrigation fluid is necessary during use of the Apollo probes. Fluid flow assists in removing tissue debris and reducing the intra-joint temperature between activations. RF probes utilized in stagnant fluid can result in heated irrigant in the joint, which can result in tissue damage inside the joint, skin burns near portals or result in damage to the probe.
- For RF probes, failure to attach the Suction Adapter to an external suction source may cause thermal injury to the patient or user. The rate and depth of tissue ablation can be affected by a number of contributors, such as, but not limited to the style of RF probe selected, the Synerg^{RF} Console power setting, amount of pressure exerted on the target tissue by the RF probe tip, and the speed that the RF probe is moved over the target tissue to be ablated.
- DO NOT place cables from other equipment between the probe and Synerg^{RF} Console.
- DO NOT allow patient contact with grounded metal objects.
- Prolonged ablation should be avoided. Intermittent passes in ablation are recommended, to allow warm fluid and tissue debris to be evacuated from the joint. Prolonged probe activation while using the suction feature may result in elevated shaft or suction tubing temperatures.

F. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- Please be sure to read and understand the Synerg^{RF} Console *User’s Guide* (*Arthrex DFU-0221-XX*) prior to using any RF probes as they are designed to be operated as a complete system.
- DO NOT use an RF probes in the presence of flammable materials. Do not place activated RF probes near or in contact with flammable materials (gauze, surgical drapes, etc.). Sparking and heating associated with electro-surgery may be an ignition source and fire could result.
- DO NOT use non-conductive irrigation fluids (sterile water, air, gas, glycine, Purisole, sorbitol-mannitol, etc.).
- DO NOT touch the tip of the RF probe while activated.
- DO NOT contact metal objects while the RF probe is activated, as this may cause damage to the RF probe.

- DO NOT insert or withdraw the tip of the RF probe from the surgical site while the RF probe is activated.
- The tip of the RF probe can remain hot and cause burns after the RF probe is deactivated.
- Inadvertent activation of the RF probe or moving the RF probe tip outside of the field-of-view can cause unintentional injury to tissue.
- DO NOT make contact with an arthroscope during activation of the RF probe as it may damage the RF probe and/or scope.
- DO NOT use metal cannulas as they may damage the RF probes’ insulation or create an alternate current path (capacitive coupling) resulting in an unintended burn. Use of all plastic cannula systems will help to avoid this problem.
- Excessive wear of the RF probe tip may result from vigorous use against bony tissue.
- Probe wear will occur with normal use, depending on a variety of factors, including high ablation settings, length of use, prolonged use in tissues and conductive fluid, and use with minimal suction or fluid management. Occasionally inspect the electrode for wear. Replace the probe if excessive wear is noted.
- As with other electro-surgical units, electrodes and cables can provide paths for high frequency current. Position the cables to avoid contact with the patient or other leads. Other electrical equipment may experience interference when positioned near the Synerg^{RF} system.
- DO NOT bend or reshape an RF probe as insulation damage can result.

- DO NOT use an RF probe that has visible insulation damage. Insulation damage can cause unintentional injury to tissue.
- Always keep activated RF probes away from the patient and user when not in use. RF probes should always be placed in a clean, dry, non-conductive, and highly visible area within the sterile field when not in use.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-company representative, and to the health authority where the incident occurred.

G. PACKAGING AND LABELING

- RF probes should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbols/lossary.

H. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

I. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

J. INFORMATION

- In EU Only:** PROCEDES carried out using these devices may be used on the general population.
- In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

K. DIRECTIONS FOR USE

Please be sure to read and understand the Synerg^{RF} Console *User’s Guide* (*Arthrex DFU-0221-XX*) prior to using any RF probes as they are designed to be operated as a complete system.

Users of the RF probe are encouraged to contact their Arthrex representatives if, in their professional judgment, they require a more comprehensive surgical technique.

I. SAFETY AND SETUP

- Examine all devices to be connected to the Synerg^{RF} Console (i.e.: Synerg^{RF} Footswitch) and used during surgery (irrigation system, trocar/cannula, aspiration system, etc.) for visible damage and cross-compatibility.
- Prepare the Synerg^{RF} Console for use as described in the *User’s Guide* (*Arthrex DFU-0221-XX*), position it for surgery and press the main power switch (I/O) to activate the Synerg^{RF} system.
- Carefully open the outer packaging and place the sterile RF probe within the sterile field. Examine the RF probe to verify the integrity of its insulation.
- Fully insert the RF probe connector into the socket on the front of the Synerg^{RF} Console. Be sure not to allow any fluid to contact the end of the connector prior to insertion into the Synerg^{RF} Console socket.

NOTE: If using an aspirating RF probe, connect the RF probe to a length of tubing that leads to standard hospital-suction canisters and ultimately to wall suction. Optimal vacuum range targets should be between 200 mmHg (8 inHg) to 400 mmHg (16 inHg).

II. BEFORE SURGERY

- Once an RF probe is connected to the front panel of a Synerg^{RF} Console, intelligent device recognition (IDR) capabilities of the Synerg^{RF} Console will display the RF probe description information and default power setting. The user can change the Synerg^{RF} Console power setting manually by either depressing the black Power Button on the top of the RF probe, the black Power Button on the footswitch or by using the touch display on the front of the Synerg^{RF} Console. Although providing with this option, it is preferred that the default power setting displayed for the specific RF probe be used to achieve the desired therapeutic effect.
- Caution should be used when overriding the default power settings. Use the lowest power setting and the minimum tissue contact time necessary to achieve the appropriate surgical effect.

CAUTION: It is recommended to use the default power setting to create the desired therapeutic effect.

III. DURING SURGERY

- Insert the tip of the RF probe through the arthroscopic portal, into the visual surgical field and ensure that the tip of the probe is fully submerged in conductive fluid. Activate the RF probe slightly off of the tissue to be treated. It is always good practice to keep the tip of the RF probe moving along, slightly off the tissue being treated.
- WARNING:** An Active RF probe tip must never be outside of the user’s field-of-view during surgery, as this can result in unintentional damage to tissue.

- NOTE:** It is normal to view bubbles around the shaft of the RF probe as well as an illumination of its tip when active.
- NOTE:** When a Synerg^{RF} Footswitch is connected to the front panel of a Synerg^{RF} Console, its controls override the buttons of the RF probe’s handle. The user has the option of overriding this feature using the Synerg^{RF} Console touch screen. Please refer to the Synerg^{RF} Console *User’s Guide* (*Arthrex DFU-0221-XX*) for complete information.
- To use the Ablate function:** Press the yellow Ablate button on the top of the RF probe’s handle to invoke the Ablate function, or if using a Synerg^{RF} Footswitch, depress the yellow pedal on the left side of the Synerg^{RF} Footswitch.
- To use the Coagulation function:** Press the blue Coag button on the top of the RF probe’s handle to invoke the Coagulation function, or if using an Synerg^{RF} Footswitch, depress the blue pedal on the right side of the Synerg^{RF} Footswitch.
- To use aspirating RF probes:** After connecting to a suction source, open the roller clamp of an aspirating RF probe fully prior to and during use. To shut off the aspiration feature, close the clamp on the suction tube only when the RF probe is outside of the patient.

- To remove eschar from the RF probe tip:** Remove the RF probe from the patient, rub the RF probe tip against an abrasive pad, and then against a wet gauze pad.
- To remove a clog from an aspirating Apollo^{RF} probe:** first, remove the RF probe from the patient. Place the RF probe into a beaker of saline solution, pump tip facing downward into the container. Connect a 10-20cc saline-filled syringe to the suction connector and inject saline into the probe suction tubing to loosen the clog.

IV. TERMINATING SURGERY

- Push the main power switch to the “OFF” position on the front panel of the Synerg^{RF} Console to deactivate the Synerg^{RF} System.
- Disconnect the RF probe’s connector from the front panel of the Synerg^{RF} Console.
- Disconnect the suction tubing.
- Dispose of the RF probe by following standard hospital protocol.

L. ALARMS AND TROUBLESHOOTING

The Synerg^{RF} Console issues an audible continuous tone whenever the RF probe is active. It also issues an audible alarm when it detects a significant change in resistance. Refer to the Synerg^{RF} Console *User’s Guide* (*Arthrex DFU-0221-XX*) or contact your Arthrex representative for more information.

- Contact Customer Service if the package has been opened or altered.

All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbols/lossary.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

Users of the RF probe are encouraged to contact their Arthrex representatives if, in their professional judgment, they require a more comprehensive surgical technique.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

Users of the RF probe are encouraged to contact their Arthrex representatives if, in their professional judgment, they require a more comprehensive surgical technique.

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