<h1>English</h1>
<p><b>Symbols glossary can be found at</b>  <b>www.arthrex.com/symbolsglossary.</b></p>

#### A. DEVICE DESCRIPTION

These devices consist of the Bio-FASTak<sup>®</sup>, FASTak<sup>™</sup>, SutureTak<sup>®</sup> and FiberTak<sup>®</sup> Suture Anchors.

These suture anchor families consist of rigid and threaded implants and all-suture implants (i.e. FiberTak). FiberTak Suture anchors (including FiberTak Button) are comprised of a polyester sheath with single or multiple sutures assembled through the sheath. The FiberTak Button suture anchor is assembled with FiberLink™ shuttling sutures only. The FiberLink shuttling sutures are not intended to be implanted but are used for loading the implantable repair sutures. Arthrex recommends the use of #2 FiberWire<sup>®</sup> or SutureTape sutures with the FiberTak Button Suture Anchor.

All implants are provided preloaded on an inserter and may include eyelets which are preloaded with various types of suture, with or without needles.

#### B. INDICATIONS

The Bio-FASTak, FASTak, SutureTak and FiberTak Suture Anchors are intended to be used for suture or tissue fixation to the foot, ankle, knee, hand, wrist, elbow, shoulder, and hip.

The Bio-FASTak, FASTak, SutureTak and FiberTak Suture Anchors are intended for the following indications:

##### I. In Non CE Accepting Countries:

- Shoulder:** Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Acromio-Clavicular Separation Repair, Deltoïd Repair, Capsular Shift or Capsulolabral Reconstruction
- Foot/Ankle:** Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Mid-foot Reconstruction, Hallux Valgus Reconstruction, Metatarsal Ligament Repair, and Digital Tendon Transfers (**SutureTak** and **FiberTak Suture Anchors** only)
- Knee:** Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Patellar Tendon, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis, and Joint Capsule Closure (**with Knotless SutureTak** only)
- Hand/Wrist:** Scapholunate Ligament Reconstruction, Carpal Ligament Reconstructions (**FiberTak** and **Knotless SutureTak** Only), Repair/Reconstruction of Collateral Ligaments Digital Tendon Transfers, and Carpometacarpal joint arthroplasty (basal thumb joint arthroplasty) (**FiberTak** and **Micro BioComposite SutureTak Suture Anchors** only)
- Elbow:** Biceps Tendon Reattachment, Ulnar/Radial Collateral Ligament Reconstruction, and Lateral Epicondylitis Repair (**Knotless SutureTak** only)
- Hip:** acetabular labral repair, Capsular repair (**excludes FiberTak Suture Anchors**). Exclusions to hip indications include the 2.4 mm FASTak Suture Anchors and FiberTak DR and DX Suture Anchors.

The 2.0 mm Mini Hip SutureTak Suture Anchor is indicated for acetabular labral repair **ONLY**.

The FiberTak Button Suture Anchor is also used for fixation of bone to bone or soft tissue to bone, and is intended as fixation posts, a distribution bridge, or for distributing suture tension over areas of ligament or tendon repair in the knee, shoulder, and elbow and may include the following indications: ACL/PCL repair, pectoralis major (minor/major), biceps tendon repair and reattachment (distal/proximal), acromioclavicular repair, and ulnar collateral ligament reconstruction.

##### I. In CE Accepting Countries:

- Bio-FASTak:**
  - Shoulder:** Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair
- FASTak:**
  - Shoulder:** Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Capsular Shift, and Capsulolabral Reconstruction
  - Foot/Ankle:** Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Mid-foot Reconstruction
  - Elbow:** Biceps Tendon Reattachment
  - Hip:** acetabular labral repair. Exclusions to hip indications include the 2.4 mm FASTak Suture Anchors.
- SutureTak:**
  - Shoulder:** Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Labral Repair, Capsulolabral Reconstruction
  - Hip:** Acetabular Labral Repair
  - Foot/Ankle:** Lateral Stabilization, Medial Stabilization
- FiberTak:**
  - Shoulder:** Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Capsular Shift, and Capsulolabral Reconstruction
  - Foot/Ankle:** Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Mid-foot Reconstruction,
  - Elbow:** Biceps Tendon Reattachment
  - Hip:** acetabular labral repair

#### C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign Body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

#### D. ADVERSE EFFECTS

- Infection, both deep and superficial.
- Foreign body reactions.
- Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
- FASTak Suture Anchors** only: Shoulder dislocation/subluxation.

#### E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be re-used.
- Do not re-sterilize this device.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

#### F. MRI SAFETY INFORMATION

##### 1. MR Conditional

*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the metal (titanium) FASTak Suture Anchors are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:*

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the FASTak Suture Anchors are expected to produce a maximum temperature rise of up to 1.8°C after 15-minutes of continuous scanning.

*In non-clinical testing, the image artifact caused by the FASTak Suture Anchors can extend up to approximately 17 mm from this implant when imaged using a gradient eco pulse sequence and a 3-Tesla MR system.*

##### 2. MR Safe

*The Bio-FASTak, SutureTak FiberTak Suture Anchors and FiberTak Button devices manufactured from only polyester, polyetheretherketone (PEEK) poly (L-Lactide, PLLA), Poly (L-Lactide-co-D acid, PLDLA), and/or Poly (L-Lactide acid, PLLA), tri-calcium phosphate (TCP and Polyester (sheath)) are MR safe. Arthrex Sutures manufactured from UHMWPE and polyester with or without silicone elastomer coating, cyanoacrylate, and nylon are MR safe.*

#### G. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate suture-anchor type and size based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- FASTak Suture Anchors** only: Stopping and restarting the drill may result in excess torque being applied to the implant, which may cause the device to fail.
- FASTak Suture Anchors** only: It is important to stop drilling as soon as the chuck contacts the guide or grasper. Failure to do so may result in damage to the suture and/or implant.
- Insert the anchor with the same orientation as that of the prepared drilled bone hole to avoid damage to the anchor or inserter.
- SutureTak and FiberTak Suture Anchors** only: Anchors Loaded on Flexible Drivers: The drill guide tip must remain in contact with the bone surface during the drilling and implant impactation steps of the procedure. Failure to do so may result in difficulty seating the implant to its intended depth.
- Self-Punching Suture Anchors** only: Insertion in very hard bone may require pre-punching a bone socket to avoid damage to the implant.
- Self-Punching Suture Anchors** only: Ensure that the angle of anchor insertion is perpendicular to the bone.
- FiberTak Button Suture Anchor** only: Surgeons may elect to use a tenodesis screw during biceptal repairs with FiberTak Button for additional tendon fixation.

#### H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.

#### I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions. Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

#### J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

**Anchor:** The device is made of titanium, PEEK (polyetheretherketone), Poly (L-Lactide) or Poly (L-Lactide-co-D, L-Lactide). FiberTak is made with a polyester sheath.

Biocomposite anchors are made of Poly (L-Lactide) or Poly (L-Lactide-co-D, L-Lactide) and tricalcium phosphate (TCP).

**Suture:** See package label for size and type of suture provided with device.

FiberWire, FiberWire CL, TigerWire<sup>®</sup>, TigerWire CL, FiberTape<sup>®</sup>, and TigerTape™ sutures are made of Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester. The FiberWire and TigerWire sutures include FiberStick™ , TigerStick<sup>®</sup>, FiberChain™, FiberLink, TigerLink, TigerTail, FiberLoop and TigerLoop sutures. The FiberLink, TigerLink™ , FiberLoop<sup>®</sup> and TigerLoop™ sutures may also be braided from SutureTape sutures.

Additional materials may include silicone elastomer coating (except suture with the suffix -Tape), cyanoacrylate, and may include nylon. LabralTape™ suture is made from Ultra High Molecular Weight Polyethylene (UHMWPE) and may include polyester or nylon.

SutureTape suture is made from Ultra High Molecular Weight Polyethylene (UHMWPE), and polyester. Additional materials may include nylon and/or cyanoacrylate.

The #2 Suture and #5 Suture are a polyblend made from Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester.

### or

The suture is made from Polyester, polytetrafluoroethylene (PTFE) coating, cyanoacrylate.

The sutures supplied with the suture anchor meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements). The suture dyes may include: D&C Blue No. 6, D&C Green No. 6, Logwood Black and FD&C Blue No.2, Chromium-cobalt-aluminum oxide (Co-Cr-Al oxide) and D&C Black No. Suture strands that are dyed with Logwood black are made of nylon.

#### K. STORAGE CONDITIONS

Bioabsorbable devices must be stored in the original unopened packaging in a dry place at maximum 32° C / 90° F, and should not be used after the expiration date.

Non-Bioabsorbable devices must be stored in the original un-opened packaging, away from moisture and should not be used after the expiration date

#### L. INFORMATION

- In CE Accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE Accepting Countries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE Accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- For Russia:** registered names of the medical devices: “Screw”; “Anchor”

Under the scan conditions defined, the FASTak Suture Anchors are expected to produce a maximum temperature rise of up to 1.8°C after 15-minutes of continuous scanning.

*In non-clinical testing, the image artifact caused by the FASTak Suture Anchors can extend up to approximately 17 mm from this implant when imaged using a gradient eco pulse sequence and a 3-Tesla MR system.*

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