# English

## A. DEVICE DESCRIPTION

The Univers™ Shoulder Fracture System consists of: a stem with a height adjustment mechanism in the neck section for attachment to the humerus, a spherical head for replacing the humeral head, and a metal head support to connect the head to the stem.

## B. INDICATIONS

The Univers Shoulder Fracture System is indicated for severe pain or significant disability resulting from degenerative, rheumatoid, or traumatic disease or injury of the glenohumeral joint. This includes traumatic or pathological conditions of the shoulder resulting in fracture of the glenohumeral joint, including impression fractures, comminuted fracture, humeral head fracture, displaced 3- or 4-fragment proximal head fractures, avascular necrosis of the humeral head, and fractures of the anatomical neck.

The Univers Shoulder Fracture System is intended for cemented or non-cemented use. This device may be used for hemi or total shoulder repair, utilizing the appropriate Arthrex glenoid component.

**In non CE accepting Countries:** The glenoid components are intended for cemented fixation in the joint and must only be used with appropriate bone cement.

## C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections that may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection, including severe neuro-arthropathy.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated above.

## D. ADVERSE EFFECTS

- Infection, both deep and superficial.
- Allergies and other reactions to device materials.
- Loosening of the implant.
- Dislocation, subluxation, inadequate scope of movement, or unwanted shortening or lengthening of the limb concerned as a result of failure to achieve optimum positioning of the implant.
- Bone fractures as a result of one-sided overload or weakened bone substance.
- Temporary or permanent nerve damage as a result of pressure or hematoma.
- Cardiovascular diseases including venous thrombosis, pulmonary embolism, and cardiac arrest.
- Wound hematoma and delayed wound healing.
- Tissue reactions caused by allergic reactions to the implanted material, particularly metal, or caused by accumulations of wear particles or cement particles.

## E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device.
- The following operative situations may cause premature loosening and complications:
- Extreme weakening of the bone structure in preparing the bone bed;
- Unsuitable selection of the implant size;
- Detailed instructions on the use and limitations of this device, the patient leaflet ([www.arthrex.com/patientleaflets](http://www.arthrex.com/patientleaflets)) and the patient implant card should be given to the patient. Your surgeon will guide you in deciding what particular treatment is best for you and explain the benefits, risks, and contraindications associated with the treatment.
- Inadequate cleaning of the bone bed prior to implantation; and,
- Excessive use of force in placing or fastening the implant, provoking splintering fractures or causing the bone to tear.
- An internal fixation device must never be re-used.

- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Do not re-sterilize this device.
- The appropriate Arthrex delivery system is required for proper insertion of the implant.
- Only Arthrex measuring systems, instruments and trial prostheses should be used for the implantation procedure.
- The surfaces of the prosthesis may not be written on, nor come into contact with metallic or other hard objects (instruments, table top, etc.).
- Endoprostheses may not be altered mechanically, or changed in any other way.
- Do not implant any parts that that have been altered from their original state, scratched, or damaged.
- An artificial joint is subject to wear and can loosen over a period of time. Wear and loosening may make it necessary to re-operate on the artificial joint.
- An infection in an artificial joint may lead to implant removal.
- The patient is to be provided with a prosthesis card that outlines the use and limitations of the device.
- This device should only be used in conjunction with other implants designed specifically for use with this system.
- Follow your institutions policy for safe disposal of all needles and other sharps or medical waste.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

## F. FACTORS AND RISKS IMPACTING THE SAFETY AND SERVICE LIFE OF THE IMPLANT

- Patient weight. An overweight patient may present additional risk.
- Extreme stress or strain resulting from work or sport related activity.
- Patients with increased risk of fractures due to repeated strain or trauma, or medical conditions that increase the patients risk for trauma.
- Osteoporosis or osteomalacia.
- Exposure to infectious diseases with possible manifestation in the joints.
- Deformity of the operative site, which can prevent or impede anchoring of the implant.
- Tumors that weaken the supportive structure.
- Allergic reactions to implant materials.
- Thromboses and pulmonary infarction caused during the implant site preparation and subsequent treatment.

## G. MRI SAFETY INFORMATION

### 1. MR Conditional

*Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Univers Shoulder Fracture System is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:*

- Static magnetic field of 1.5-Tesla and 2.5-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the Univers Shoulder Fracture System is expected to produce a maximum temperature rise of 3 °C after 15-minutes of continuous scanning.

*In non-clinical testing, the image artifact caused by the Univers Shoulder Fracture System extends approximately 60 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*

## H. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

## I. PACKAGING AND LABELING

- Arthrex implants should be accepted only if the factory packaging and labeling are intact.
- All prosthesis parts and components should be stored unopened in the original packaging and should be kept in protective packaging until immediately prior to use.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbolsglossary](http://www.arthrex.com/symbolsglossary).

## J. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately

cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

## K. MATERIAL SPECIFICATIONS

The device consists of a stem body made of titanium, alloy and a head made of cobalt chromium (Co-Cr) alloy.  
**In non CE accepting Contries:** The head is also available as TiNbN coated cobalt-chromium (Co-Cr) alloy.

## L. STORAGE CONDITIONS

Joint prostheses must be treated with great care prior to implantation. Scratches or nicks on the joint surface can result in excessive warp and complications.

- Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.
- Do not remove the device from the protective packaging until immediately prior to use.

## M. INFORMATION

- In CE accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE accepting Countries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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