English

A. DESCRIPTION

Arthrex Footswitches are designed specifically to work with Arthrex consoles, be sure to match the correct footswitch to the correct console. AR-6483 is designed to control Arthrex DualWave™ Arthroscopy Pump (AR-6480). AR-8310, AR-8315C and AR-8315W are designed to control Synergy^{Resection}™ Console (AR-8305). AR-8310SP is designed to control both the Arthrex DualWave Arthroscopy Pump (AR-6480) and the Synergy^{Resection}™ Console (AR-8305).

The footswitch lets you adjust the settings at a distance from the control unit of the console.

For more details on the effect of each function refer to console's manual that the footswitch is being connected to.

Note: Each footswitch is designed to control specific functions for each console.

B. CAUTIONS

1. Caution: Federal law restricts this device to sale by or on the order of a physician.
2. This device is intended to be used by a trained medical professional.
3. Do NOT sterilize this device.
4. Do NOT clean and/or disinfect in a washer/disinfectant.
5. NEVER immerse this device in water, cleaning, or disinfectant solution.
6. Do NOT clean the device with abrasive cleaning or disinfectant compounds, solvents, or other materials that could scratch or damage the device.
7. Do NOT lift the device by the cord.
8. Keep this device away from magnetic and electromagnetic fields.
9. This unit is NOT suitable for use in the presence of flammable mixtures.
10. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
11. Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.
12. Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
13. Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical techniques information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

C. CLEANING AND DISINFECTATION

1. The footswitch is supplied non-sterile and is not to be sterilized. The footswitch must be adequately cleaned and disinfected prior to use or re-use.
2. Use a commercially available surface disinfectant. Arthrex recommends the use of a quaternary ammonium based product.
3. Do NOT clean the device with abrasive cleaning or disinfectant compounds, solvents, or other materials that could scratch or damage the device.
4. Always comply with the instructions issued by the manufacturer of the disinfectant regarding concentration, exposure times, temperature, and material compatibility.
5. To clean and disinfect the footswitch, unplug the cord from the console and add the protective cap to the shaver connector.

6. Using a disinfecting towelette or a clean, low-linting cloth dipped in disinfectant solution, gently wipe down all surfaces of the footswitch.
7. Ensure all surfaces remain visibly wet for the contact time recommended by the disinfectant manufacturer.
8. If required by the disinfectant manufacturer, rinse per instructions; otherwise, allow to air dry.
9. Check all surfaces for visible soil. Repeat cleaning/disinfecting steps if soil is visible and re-inspect.

CAUTION: Do NOT allow fluids to enter the connector or allow the receptacles to have any contact with liquids. Remove dust or moisture, if necessary, with dry compressed air. ONLY dry connectors should be plugged into the console.

D. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact customer service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

E. INSPECTION

Before use, the device should be inspected for signs of wear or damage. In particular look for any of the following problems:

- Damage to the cord
- Damage to the cord receptacle
- Damage to the connector
- Damage to the buttons
- Obvious damage to the unit
- Traces of debris or stains

For continued safety and performance of the footswitch, it is recommended that the functionality and electrical safety of the switch to be verified on an annual basis by a qualified person.

F. STORAGE CONDITIONS

The device should be stored at room temperature in a dry environment.

G. INFORMATION

1. **In EU Only:** Procedures carried out using these devices may be used on the general population.
2. **In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
3. **In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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