**A. DEVICE DESCRIPTION**

The Arthrex Mixing and Delivery System consists of a piston syringe with a movable plunger and cap to facilitate mixing and delivery; and may include straight and curved (tuohy) delivery needles; a mating obturator for delivery needles; luer connectors; guide pins; reamers; delivery syringes; and a funnel to facilitate filling of the syringe barrel. The system will be offered with either a 3 mL or 14 mL syringe barrel and may be provided either empty or pre-filled with allograft, autograft, or synthetic bone graft materials.

**B. INDICATIONS**

The Arthrex Mixing and Delivery System is indicated for the delivery of allograft, autograft, or synthetic bone graft materials to all orthopedic surgical sites. In addition, it is designed to facilitate pre-mixing of allograft, autograft, or synthetic bone graft materials with I.V. fluids, blood, plasma, platelet rich plasma, bone marrow or other specific blood component(s) as deemed necessary by the clinical use requirements.

**C. CONTRAINDICATIONS**

1. Insufficient quantity or quality of bone.
2. Blood supply limitations and previous infections, which may retard healing.
3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
4. Any active infection or blood supply limitations.
5. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
6. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
7. Do not use for surgeries other than those indicated.

**D. ADVERSE EFFECTS**

1. Allergies and other reactions to device materials.
2. Foreign body reactions.

**E. WARNINGS**

1. Caution: Federal law restricts this device to sale by or on the order of a physician.
2. Procedures carried out using these devices may be used on the general population.
3. The clinical benefits associated with the use of these devices outweigh the known clinical risks.
4. There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

5. This device is intended to be used by a trained medical professional.
6. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device.
7. Do not re-sterilize this device.
8. Do not re-use this device.
9. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.
10. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
11. Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

9. A mating obturator for delivery needles; luer connectors; guide pins; reamers; delivery syringes; and a funnel to facilitate filling of the syringe barrel. The system will be offered with either a 3 mL or 14 mL syringe barrel and may be provided either empty or pre-filled with allograft, autograft, or synthetic bone graft materials.
10. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
11. Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

**F. PRECAUTIONS**

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

**G. PACKAGING AND LABELING**

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at [www.arthrex.com/symbolsglossary](http://www.arthrex.com/symbolsglossary).

**H. STERILIZATION**

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or reuse. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

**I. MATERIAL SPECIFICATIONS**

Refer to the package label for the materials.  
**Needle:** Surgical-grade stainless steel, ABS Luer Lock  
**3 mL or 14 mL Syringe:** Polyamide, PBT, EDPM, PP  
**Cap:** PBT  
**Luer Connectors:** ABS

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**J. STORAGE CONDITIONS**

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

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