English

A. DEVICE DESCRIPTION

The iBalance® UKA System is a unicondylar knee arthroplasty system that consists of femoral components, tibial tray components and tibial bearing inserts. All components are available in a range of sizes to fit varying anatomical requirements. Femoral components and tibial tray components are available in left-medial/right-lateral and right-medial/left-lateral varieties. Tibial bearing inserts are symmetric for left or right use.

B. INDICATIONS

Indicated for use in unicompartmental knee arthroplasty as a result of:

- Moderately disabling joint disease of the knee resulting from painful osteoarthritis or posttraumatic arthritis
- Correction of functional deformity
- Revision of previous unsuccessful unicompartmental knee replacement or other procedure
- As an alternative to tibial osteotomy in patients with unicompartmental osteoarthritis

These components are single use only and are intended for implanta-tion with bone cement.

When used concurrently, the Arthrex iBalance UKA and PFJ systems create the Arthrex iBalance BiCompartmental Arthroplasty System. The Arthrex iBalance BiCompartmental Arthroplasty System is intended to be used as a multi-compartmental knee arthroplasty in patients with:

- Moderately disabling joint disease of the knee resulting from painful osteoarthritis or post traumatic arthritis;
- Correction of functional deformities;
- Revision of previous unsuccessful partial knee replacement or other procedure.

The BiCompartmental Arthroplasty System is not intended to be used as a dual-condyle or tri-compartmental knee.

These components are single use only and are intended for implan-tation with bone cement.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- The use of this device may not be suitable for patients with insufficient or immature bone.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- Do not use for surgeries other than those indicated.
- An overweight or obese patient can produce loads on the prosthesis which can lead to failure of fixation or failure of the device itself.
- Severe deformity and/or recurrent subluxation of the knee joint.
- Inflammatory arthritis conditions, particularly rheumatoid arthritis.
- Ligamentous instability or soft-tissue laxity such that the postoperative stability afforded by the prosthesis may be compromised or dislocate.
- Varus or valgus deformity greater than 15 degrees that is not passively correctable.

D. ADVERSE EFFECTS

- Loosening, migration, or loss of fixation of the implant.
- Foreign body reactions.
- Infections, both deep and superficial.
- Particulate wear debris and mild tissue discoloration from metallic components have been noted with other prosthetic devices constructed of similar materials. Some types of wear debris have been associated with osteolysis and implant loosening.
- Wear and damage to the articulating surfaces or soft tissue support structures.
- Postoperative symptoms include, but are not limited to: pain, persistent swelling, stiffness, limited range of motion, or incomplete resolution of preoperative symptoms.
- Periarticular calcification or ossification with or without impediment of joint mobility.
- Incomplete range of motion due to improper selection or positioning of components.
- Fatigue fracture of the implants.
- Intraoperative or postoperative bone fracture.
- Allergic reactions to materials.
- Other potential adverse effects of knee surgery may include: locking, neurovascular damage, dislocation, bending of the components, leg length discrepancies, and other less common adverse effects.
- Serious complications may be associated with any total joint replacement surgery. These complications include, but are not limited to genitourinary disorders; gastrointestinal disorders; vascular disorders, including thrombus; bronchopulmonary disorders, including emboli; myocardial infarction or death.
- Particles from metal and polyethylene components can be shed from the components during normal use and over time. These particles may be incorporated into the joint through natural biological processes or migrate to other parts of the body. The long-term effects of these particles may include, but are not limited to, cancer, lymphadenopathy and accumulation in other tissues and organs, and systemic disease.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- Care should be utilized in the handling of the components to minimize contamination and damage of the component surfaces.
- When using cement for fixation, the surgeon should use care to ensure complete cement contact and coverage on all parts of the prosthesis embedded in bone cement. Remove any excess cement to ensure that no cement remains on the articulating surfaces of the prosthesis.

- An internal fixation device must never be re-used.
- Do not re-sterilize this device.
- The surgeon must not allow damage to polished bearing surfaces because this may accelerate wear of the components. Any alteration or damage to a component may result in failure under load. Any prostheses so damaged must not be used.
- Components of the iBalance UKA System should not be used with those of another manufacturer since articular and dimensional compatibility cannot be assured.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Postoperatively and until healing is complete, the fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- This product’s useful life will vary between patients, and is dependent upon accumulated weight-bearing activities and lifestyle.
- Physicians should carefully assess the activity level of the patient prior to performing arthroplasty of the knee. Increased activity as well as increased weight can lead to accelerated wear of UHMWPE components.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex iBalance UKA is MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 2000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system

Under the scan conditions defined, the Arthrex iBalance UKA is expected to produce a maximum temperature rise of 2 °C after 15-minutes of continuous scanning.

Patient knee must be in direct contact with the patient table. Pillows should not be placed under the knee during scanning.

In non-clinical testing, the image artifact caused by the Arthrex iBalance UKA extends approximately 20 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

G. PRECAUTIONS

- Proper selection of the Arthrex iBalance UKA System depends on the judgment of the surgeon.
- For cemented use only, unless otherwise indicated.
- Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

NOTE: Ultra High Molecular Weight Polyethylene (UHMWPE) components may be damaged if resterilized. Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

Femoral Implant: Cobalt Chrome
Tibial Tray: Cobalt Chrome
Tibial Bearing: Ultra High Molecular Weight Polyethylene (UHMWPE) or Vitamin E blended UHMWPE (not available in all countries).

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

- In CE Accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE Accepting Countries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE Accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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