English

A. DESCRIPTION

The **FiberTape® Cerclage** and **TigerTape™ Cerclage** sutures are assembled from **FiberTape** and **TigerTape** sutures, in a loop configuration. The looped tip of the suture is tied as a hitch over a suture sheath securing a shuttling polyblend suture.

FiberTape and TigerTape are non-absorbable, flat braided sutures composed of UHMWPE and polyester yarns over a core of FiberWire® or TigerWire® sutures (each made of a polyblend of UHMWPE and polyester). Additional materials may include cyanoacrylate and/ or nylon.

These sutures meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements).

These sutures are supplied sterile, in pre-cut lengths, in various single and multi-loop configurations, and in some cases, with various swaged needles, with stiffened ends and with a pre-tied hitch. The sutures are available non-dyed, dyed, and fully or partially striped. Dyes may include: D&C Blue No. 6, D&C Green No. 6, and Logwood Black. Suture strands that are dyed black are made of nylon.

F. MRI SAFETY INFORMATION

1. MR Safe

FiberTape and TigerTape sutures are MR Safe.

G. PRECAUTIONS

1. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

2. Assure that all knots have been secured using accepted surgical knot tying techniques. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. Care should be taken to prevent damage to surrounding tissue or user puncture due to improper handling of the needlepoint.

3. Do not grasp the needle at the point or swage to avoid damage to these areas. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in “sharps” containers.

4. Avoid wrapping the suture over sharp metal or bone graft surfaces.

5. Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

6. Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.

7. Surgeons should be aware that the use of sutures may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in “sharps” containers.

8. Do not grasp the needle at the point or swage to avoid damage to these areas. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in “sharps” containers.

D. CONTRAINDICATIONS

1. These sutures are not for use in cardiac indications.

2. Spinal fixation applications.

3. Osteoporosis or insufficient bone quality.

4. Patients that are skeletally immature.

5. Cases where there is poor quality or osteoporotic bone at the repair.

6. Compromised vascularity that would inhibit adequate blood supply to the fracture or operative site.

E. WARNINGS

1. Caution: Federal law restricts this device to sale by or on the order of a physician

2. Procedures carried out using these devices may be used on the general population.

3. The clinical benefits associated with the use of these devices outweigh the known clinical risks.

4. There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

5. This device is intended to be used by a trained medical professional.

6. Do not re-sterilize.

7. Once open, discard unused suture.

8. Do not expose suture to heat.

9. Users should be familiar with surgical procedures and techniques involving non-absorbable sutures

before employing Arthrex sutures for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

10. As with any foreign body, prolonged contact of this or any other suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formations. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

11. This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.

12. The safety and effectiveness of these devices for use as an artificial ligament or tendon has not been established; hence, the device should not be used for those applications.

13. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.

14. Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Safe

FiberTape and TigerTape sutures are MR Safe.

G. PRECAUTIONS

1. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

2. Assure that all knots have been secured using accepted surgical knot tying techniques. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. Care should be taken to prevent damage to surrounding tissue or user puncture due to improper handling of the needlepoint.

3. Do not grasp the needle at the point or swage to avoid damage to these areas. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in “sharps” containers.

4. Avoid wrapping the suture over sharp metal or bone graft surfaces.

5. Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

6. Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.

7. Surgeons should be aware that the use of sutures may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in “sharps” containers.

8. Do not grasp the needle at the point or swage to avoid damage to these areas. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in “sharps” containers.

9. Do not grasp the needle at the point or swage to avoid damage to these areas. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in “sharps” containers.

H. ADVERSE REACTIONS

1. Adverse reactions have not been noted with the Arthrex FiberTape product in animal testing.

2. Common non-absorbable suture reactions may include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, enhanced bacterial infectivity, minimal acute inflammatory tissue reaction, pain, edema, and erythema at the wound site.

3. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

4. Loosening.

5. Need for revision surgery.

I. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.

2. Contact Customer Service if the package has been opened or altered.

3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at **www.arthrex.com/symbolsglossary**.

J. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

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