English

A. DESCRIPTION

The Arthrex Suture Family is comprised of non-absorbable sutures made of Ultra High Molecular Weight Polyethylene (UHMWPE), or a polyblend of UHMWPE and polyester.

FiberWire® and TigerWire™ is a braided suture made of UHMWPE and polyester over a core of UHMWPE yarns. Some braided polyblend sutures or UHMWPE sutures may be coreless. Additional materials to the Arthrex sutures include silicone elastomer coating, and may include cyanoacrylate and/or nylon. The coating acts as a lubricant for suture sliding, knot tying, and ease of passing suture through tissue. FiberWire CL and other coreless polyblend suture may be uncoated.

FiberTape®, TigerTape™, LabralTape™ and SutureTape™ are flat braided constructs. FiberTape and TigerTape suture is composed of UHMWPE and polyester yarns over a core of FiberWire or TigerWire. LabralTape and SutureTape are made of UHMWPE but may include polyester yarns and/or nylon. SailTape™ is a woven construct made of UHMWPE and polyester with tapered suture tails made of UHMWPE. Additional materials to the Arthrex tape constructs may include cyanoacrylate and/or nylon.

The sutures meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements). The diameters of the sutures typically span the diameter requirements of two USP sizes. For that reason, Arthrex tests the sutures to conform to the greater USP size specification requirement. SutureTape sutures meet or exceed USP standards for size 2 or size 2-0 non-absorbable surgical sutures, except for an oversize in diameter for the round tails and width of flat tape.

The Arthrex Suture Family is supplied sterile, in pre-cut lengths, in various single and multi-loop configurations, and in some cases, with various swaged needles and with stiffened ends. The Arthrex Suture family is available non-dyed, dyed, and fully or partially striped. Dyes may include: D&C Blue No. 6, D&C Green No. 6, and Logwood Black. Suture strands that are dyed black are made of nylon.

B. INDICATIONS

The Arthrex Suture Family is intended for use in soft tissue approximation and or ligation. These sutures may be incorporated, as components, into surgeries where constructs including those with allograft or autograft tissues are used for repair.

When used in conjunction with the Arthrex Blunt Tip Screws, FiberTape can be used to treat patella fractures.

C. ACTIONS

Arthrex sutures, when tested per ISO 10993, *Biological Evaluation of Medical Devices-Part 10: Tests for Irritations and Sensitization*, had no reactions of allergic or sensitive nature. The dyed suture and coating are pharmacologically inactive.

Arthrex FiberWire suture is not absorbed, but may become encapsulated in the surrounding connective tissues. The Arthrex FiberWire is not known to have significant change in tensile strength *in vivo*.

D. CONTRAINDICATIONS

Arthrex suture is not for use in cardiac indications.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- Procedures carried out using these devices may be used on the general population.
- The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.
- This device is intended to be used by a trained medical professional.
- Do not re-sterilize this device.
- Once open, discard unused suture.
- Do not expose suture to heat.
- Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing Arthrex FiberWire or Arthrex UHMWPE suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- As with any foreign body, prolonged contact of this or any other suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formations. Acceptable surgical practice

must be followed with respect to drainage and closure of infected or contaminated wounds.

- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- The safety and effectiveness of these devices for use as an artificial ligament or tendon has not been established.
- Detailed instructions on the use and limitations of this device, the patient leaflet and the patient implant card should be given to the patient.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Safe

Arthrex Suture manufactured from UHMWPE and polyester with or without silicone elastomer coating, cyanoacrylate, and nylon are MR safe.

G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
- Assure that all knots have been secured using accepted surgical knot tying techniques. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. Care should be taken to prevent damage to surrounding tissue or user puncture due to improper handling of the needlepoint.
- Do not grasp the needle at the point or swage to avoid damage to these areas. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in “sharps” containers.

H. ADVERSE REACTIONS

- Adverse reactions have not been noted with the Arthrex FiberWire product in animal testing.
- Common non-absorbable suture reactions may include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, enhanced bacterial infectivity, minimal acute inflammatory tissue reaction, pain, edema, and erythema at the wound site.
- Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.
- Silicone sensitivity, though very rare, has been reported.

I. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at **www.arthrex.com/symbolsglossary**.

J. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

Surgeons are advised to review the product specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.

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