<b>English</b>
<p><b>Symbols glossary can be found at <a href="http://www.arthrex.com/symbolsglossary">www.arthrex.com/symbolsglossary</a>.</b></p>

##### A. DEVICE DESCRIPTION

The Arthrex LoopLoc™ knotless suture is a knotless suture assembled with coreless FiberWire® suture looped through a FiberLink™ suture.

##### B. INDICATIONS

**In Non CE Accepting Countries:** The Arthrex LoopLoc knotless suture is intended for use in soft-tissue approximation for hip capsular closure.

**In CE Accepting Countries:** The LoopLoc knotless suture is intended for use in soft-tissue approximation for hip capsular closure.

##### C. CONTRAINDICATIONS

- Blood supply limitations and previous infections, which may retard healing.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- Do not use for surgeries other than those indicated.
- Arthrex LoopLoc knotless suture is not for use in cardiac indications.

##### D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.
- Silicone sensitivity, though very rare, has been reported.
- Adverse reactions have not been noted with the Arthrex FiberWire product in animal testing.
- Common non-absorbable suture reactions may include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occur, enhanced bacterial infectivity, minimal acute inflammatory tissue reaction, pain, edema, and erythema at the wound site.

##### E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- An internal fixation device must never be re-used.
- Do not re-sterilize this device.
- Once open, discard unused suture.
- Do not expose suture to heat.
- Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing Arthrex FiberWire suture for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/ or user.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.
- Arthrex FiberWire suture is not absorbed but may become encapsulated in the surrounding connective tissues. The Arthrex FiberWire suture is not known to have significant change in tensile strength *in vivo*.

##### F. MRI SAFETY INFORMATION

###### 1. MR Safe

The Arthrex LoopLoc knotless suture device is manufactured from UHMWPE and polyester with or without silicone elastomer coating, which are MR safe.

##### G. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

##### H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.

##### I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

##### J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

FiberWire suture is a braided suture made of Ultra High Molecular Weight Polyethylene (UHMWPE) and polyester, over a core of UHMWPE yarns. Some braided polyblend sutures or UHMWPE sutures may be coreless. Additional materials to the LoopLoc knotless suture include silicone elastomer coating and may include cyanoacrylate. The coating acts as a lubricant for suture sliding, knot tying, and ease of passing the suture through tissue.

The FiberWire suture includes FiberLink suture. The FiberLink suture is assembled from FiberWire sutures.

The sutures meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements). The diameters of the sutures typically span the diameter requirements of two USP sizes. For that reason, Arthrex tests the sutures to conform to the greater USP size specification requirement.

The Arthrex LoopLoc knotless suture is supplied sterile in a multi-loop configuration. The LoopLoc knotless suture is available dyed. The dye includes: D&C Blue No. 6.

##### K. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

##### L. INFORMATION

- In CE Accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE Accepting Countries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE Accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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