English

- technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Surgeons must apply their professional judgment when determining the appropriate size device based on the specific indication, preferred surgical technique, and patient history.
- Do not bend the plate near the locking hole. Bending the plate near the locking hole can distort the holes threading, which prohibits insertion of the screw.
- The Distal Clavicle Plate Button is designed to fit securely in the holes of the fracture plates. The Plate Button is available in one size only. The Distal Clavicle Plate Button is available as a single board or assembled in a construct with Arthrex suture. Washers are available in diameters of 6 – 13 mm.

A. DEVICE DESCRIPTION
Plates are available in different shapes, sizes and orientations (e.g. left and right types). The plates have specified sizes for screws to provide fixation. Some plates with multiple holes (i.e. Arthrex Mesh Plates) may cut to size.

The Distal Clavicle Plate Button is designed to fit securely in the holes of the fracture plates. The Plate Button is available in one size only. The Distal Clavicle Plate Button is available as a single board or assembled in a construct with Arthrex suture. Washers are available in diameters of 6 – 13 mm.

- B. INDICATIONS**
The Arthrex Mini CFS (Comprehensive Fixation System) Plates (1.4–2.4 mm) are intended for use in selective trauma, reconstructive procedures, and general surgery of the hand, wrist, and other small bones. The Arthrex Comprehensive Fixation System Plates are to be used with the Arthrex Comprehensive Fixation System Screws (1.4–2.4 mm solid).

The Arthrex Mini CFS Plates (2.0–2.4 mm) are intended for fixation of fractures, osteotomies, nonunions, replantations, and fusions of small bones and small bone fragments, particularly in osteopenic bone in the hand, wrist, foot, and ankle. The Arthrex Comprehensive Fixation System Plates are to be used with the Arthrex Comprehensive Fixation System Screws (2.0–2.4 mm solid).

The Arthrex Fracture Plates and Ankle Fusion Plates are intended to be used for internal bone fixation for bone fractures, fusions, osteotomies and non-unions in the ankle, foot, hand, wrist, clavicle, scapula, olecranon, humerus, radius, ulna, tibia, calcaneus and fibula. The Arthrex Fracture Plates are to be used with the Arthrex Lock Profile Plates.

The Distal Clavicle Plate Button is intended for use with the clavicle plates for clavicle indications such as for the treatment of syndesmotic trauma, such as fixation of acromioclavicular separations due to coracoclavicular ligament disruption, and this button may not be used alone. The button is intended to be used with #5 Fibres®.

The Arthrex Low Profile Plates, Arthrex Compression Plates and Calcaneal Fracture Plates are intended to be used for internal bone fixation for bone fractures, fusions, or osteotomies in the ankle, foot, hand, and wrist, such as opening wedge osteotomies of Hallux Valgus.

The Arthrex Extremity Extralyte Plate System and Arthrex Mesh Plates are intended for use in stabilization of fresh fractures, revision procedures, osteotomies, joint fusion and reconstruction of small bones and bone fragments of the hand/ wrist, foot/ankle, and osteopenic bone.

The Arthrex Distal Radius System is intended for fixation of fractures, fusions, osteotomies and non-unions of the radius, ulna, olecranon, metacarpal, metatarsal and metatarsals.

The Arthrex washers and certain plates are intended as adjunct devices.

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity rated out prior to implantation.
- Any active infection or blood supply limitations.
- Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with infection or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not delay, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
- Foreign body reactions.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or to the order of a physician
- This device is intended to be used by a trained medical professional and is not intended for home use.
- An internal fixation device must never be re-used.
- All metallic implant devices used for this surgical procedure should have the same metallic composition.
- Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
- Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Removal of supplemental fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Migration of implant position resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to the presence of the device; (6) Possible increased risk of infection; and (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
- Bioabsorbable waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex, Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

Non-clinical testing and in-vivo electromagnetic simulations demonstrated that the Arthrex Plates are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning in the Normal Operating Mode for operation for the MR system

- Under the scan conditions defined, the Arthrex Plates are expected to produce a maximum temperature rise of up to 6°C after 15-minutes of continuous scanning.
- In non-clinical testing, the image artifact caused by the Arthrex Plates can extend up to approximately 60 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.*

G. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical

- After preliminary cleaning is complete, load the devices in the washer-disinfector such that all device features of the device are accessible to cleaning and such that design features that might retain liquid can drain (for example, hinges should be open and cannulations/holes positioned to drain).
- If using alkaline cleaning agents, a neutralization step should be utilized as appropriate.
- Run an automated wash cycle with fundamentally approved efficiency of the washer-disinfector (for example, CE marking according to EN ISO 15883 or FDA approval/clearance/registration). The following minimum recommended automated wash cycle parameters were utilized by Arthrex during the validation of these instructions.

RECOMMENDED WASHING CYCLE PARAMETERS				
Phase	Recirculation Time	Temperature	Detergent	
Pre-Wash	3 Minutes	Cold Water	N/A	
Cleaning Wash	10 Minutes	Follow detergent manufacturer’s recommendation	Enzymatic or alkaline detergent	
Neutralization Rinse (optional)	2 Minutes	Follow detergent manufacturer’s recommendation	Neutralizing agent (if needed)	
Rinse	3 Minutes	Cold Water	N/A	
Thermal Disinfection Rinse	5 Minutes	90°C (194°F)	N/A	
Drying	Minimum 6 Minutes or until fully dry	Minimum 100°C (212°F)	N/A	

- Remove the devices from the washer-disinfector following the completion of the program and check devices for visible soil. Repeat cleaning if soil is visible and re-inspect; otherwise, proceed to Sterilization step.

IV. MANUAL CLEANING AND DISINFECTION

Following preliminary cleaning, the instructions for Manual Cleaning and Disinfection may be followed as an alternative cleaning method to Machine (Automated) Cleaning and Thermal Disinfection if an automated procedure is not available.

- After preliminary cleaning is complete, repeat steps 1-5 provided in the Preliminary Cleaning section within this DFU, including rinsing, immersion and sonication, and post-rinsing. Final rinsing should be completed with purified (critical, e.g. RO or DI) water.
- Check devices for visible soil. Repeat cleaning if soil is visible and re-inspect.
- Soak the devices for the given soaking time (provided by the disinfectant manufacturer) in disinfectant solution so that the devices are sufficiently covered. Make sure that there is no contact between the devices. Ensure that the device is in the open position during soaking. The device is to be soaked at least five times during disinfection, as applicable.
- Remove the devices from the disinfectant solution and rinse per disinfectant manufacturer’s instructions.
- Dry devices thoroughly utilizing filtered medical grade air or a soft, low-lining cloth. Proceed to sterilization stage.

K. STERILIZATION

This device may be provided either sterile or non-sterile. Check the package labeling for more information. For devices that are not provided in a terminally sterilized configuration, sterilization is to be performed following cleaning, disinfection, and sterile packaging prior to use, and may be re-sterilized (if unused) following cleaning, disinfection, and sterile packaging prior to use.

Devices that are provided in a terminally sterilized configuration should never be re-sterilized under any conditions.

Arthrex Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSIAAMI ST79 for specific information.

I. STERILE PACKAGING

Singly: Single devices should be packed as to ensure that the pack is large enough to contain the device without stressing the seals. Packaging should be completed utilizing a pop-top or wrap which conforms to the recommended specifications for steam sterilization as outlined below. If a wrap is utilized, it should be completed following AAMI double-wrap or equivalent guidelines with an appropriate wrap. An appropriate wrap is one that, for example, is cleared by the FDA or the local governing body at the point of use. Devices may also be placed into an approved reusable rigid sterilization container. Aesculap SteriContainer™ rigid containers with perforated bottoms and lids are approved for use with Arthrex, Inc. devices.

Sets: Where appropriate, cleaned, disinfected and inspected devices should be placed into trays/cases as provided or in general-purpose sterilization trays. The total weight of trays/cases should not exceed 11.4kg/25 lbs. (other local limits below 11.4kg/25 lbs. may apply). Trays/cases should be double wrapped following AAMI or equivalent guidelines with an appropriate wrap. An appropriate wrap is one that, for example, is cleared by the FDA or the local governing body at the point of use. Sets may also be placed into an approved rigid container with perforated bottoms and lids as approved for use with Arthrex, Inc. devices.

Areas designated for specific devices shall contain only devices intended for those areas. Devices should not be stacked or placed in close contact. Only Arthrex devices should be included in the trays or cases. These validated processing instructions are not applicable to trays or cases that include devices not intended to be used with Arthrex trays or cases.

II. STERILIZATION

Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in the table below. Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer’s instructions.

RECOMMENDED STEAM STERILIZATION PARAMETERS					
Cycle Type	Minimum Exposure Temperature	Minimum Exposure Time	Minimum Dwell Time ¹	Minimum Cooling Time ²	
US Preva-cure Cycle	132°C (270°F)	4 Minutes	30 Minutes	30 Minutes	
UK Pre-cure Cycle ³	134°C (273°F)	3 Minutes	30 Minutes	30 Minutes	
Preva-cure Cycle ⁴ (Prion Cycle)	134°C (273°F)	18 Minutes	30 Minutes	30 Minutes	

¹Drying times vary according to load size and should be increased for large loads.

²Cooling times vary according to the sterilizer used, device design, temperature and humidity of ambient environment, and type of packaging used. Cooling process should comply with ANSIAAMI ST79.

³Reprocessing parameters recommended by the World Health Organization (WHO) where there is concern regarding TSE/CJD contamination.

L. MATERIAL SPECIFICATIONS

Refer to the package label for the materials. This device is made of titanium or stainless steel.

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- N. INFORMATION**
- IN EU Only:** Procedures carried out using these devices may be used on the general population.
- IN EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- IN EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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