English

A. DEVICE DESCRIPTION

The Post, Bio-Post™, Bi-Cortical™, Bio-Post™, or Screws are threaded implants in a variety of designs and sizes. The Washers are available in smooth and spiked designs, and may be provided with a Post or Screw.

B. INDICATIONS

The Arthrex metal Post or Screw with or without a Washer, and the Arthrex bioabsorbable Bio-Post and Washers and the Bi-Cortical Bio-Post and Washers are used as anchoring devices for suture, or to secure soft tissue directly to bone in the knee in the following procedures:

Knee: Cruciate Ligament Repair, Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Patellar Tendon Repair, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis

C. CONTRAINDICATIONS

- Insufficient quantity or quality of bone.
- Blood supply limitations and previous infections, which may retard healing.
- Foreign-body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Foreign body Reactions. See Adverse Effects-Allergic Type Reactions.
- Any active infection or blood supply limitations.
- Conditions that tend to limit patient's ability or willingness to restrict activities or follow directions during the healing period.
- The use of this device may not be suitable for patients with sensitive or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants should not bridge, disturb or disrupt the growth plate.
- Do not use for surgeries other than those indicated.

D. ADVERSE EFFECTS

- Infections, both deep and superficial.
 - Foreign body reactions.
 - Allergic-like reactions to PLA materials (PLLA, PLDLA) have been reported. These reactions have sometimes necessitated the removal of the implant. Patient sensitivity to device materials must be considered prior to implantation.
- E. WARNINGS**
- Caution: Federal law restricts this device to sale by or on the order of a physician.
 - This device is intended to be used by a trained medical professional.
 - An internal fixation device must never be re-used.
 - All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
 - Patient sensitivity to the device materials should be considered prior to implantation. See Adverse Effects.
 - Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
 - Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device. The appropriate Arthrex delivery system is required for proper implantation of the device.

- Caution: Do not remove the device should take into consideration the potential risks to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
- Devices that have been implanted for a long period of time may require the use of screw removal instrumentation.
- Detailed instructions on the use and limitations of the device should be given to the patient.
- This is a single use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
- Removal of supplemental metal fixation after healing. If the supplemental fixation is not removed following the completion of its intended use, any of the following complications may occur:
 - (1) Corrosion, with localized tissue reaction or pain;
 - (2) Migration of any non-stainless metal resulting in injury;
 - (3) Risk of additional injury from post-traumatic trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult;
 - (5) Pain, discomfort, or abnormal sensations due to the presence of the device;
 - (6) Possible increased risk of infection; and
 - (7) Bone loss due to stress shielding. The surgeon should carefully weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid re-fracture.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex, Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

1. MR Conditional

Non-clinical testing and in vivo electromagnetic simulations demonstrated that the Post, Screw, and Washer are MR Conditional. A patient with this device can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3000 Gauss/cm or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 1 W/kg for 15 minutes of scanning in the Normal Operating Mode of operation for the MR system.

- Under the scan conditions specified, the Post, Screw, and Washers are expected to produce a maximum temperature rise of up to 6°C after 15-minutes of continuous scanning.
- In non-clinical testing, the image artifact caused by the Post, Screw, and Washers can extend up to approximately 120 mm from the implant, depending on the imaging sequence and pulse sequence and a 3-Tesla MR system.*

2. MR Safe

- The Post, Screw, and Washer devices manufactured from only Poly (L-Lactide) are MR Safe.*
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
 - Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
 - The Arthrex Washers are designed for use with specific Arthrex Posts and Screws. The use of other posts and screws is not recommended.
 - The Arthrex GablesSoft Scrap Tap is used to tap the proximal cortex prior to screw insertion. Failure to do so may result in failure of the screw or screwdriver due to excess torque during insertion.
 - Damage to the implant or screwdriver may result from failure to fully seat the driver in the screw socket or misalignment between the screwdriver and screw socket.

- Bio Only: Do not over-tighten the screw post onto the washer.

H. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and Customer Service info is intact.
- Contact your Arthrex Service if Package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbols/glossary.

I. VALIDATION

The recommended cleaning, disinfection, and sterilization methods in this DFU have been validated in compliance with federal and international guidance/standards. In accordance with ISO 17665, the "overall" approach was used for sterilization validation, and demonstrates a sterility assurance level (SAL) of 10⁻⁶. Cleaning, disinfecting, and sterilizing equipment and materials vary in performance characteristics. Therefore, it is the responsibility of the facility/end user to perform the appropriate validation testing for any use beyond recommended performance characteristics.

In accordance with EN ISO 17664 and AAMI TR30, limit values and a range of mechanical residues following cleaning have been established for the product. In assessing the level of cleaning residues following the manual cleaning and disinfection process or the machine (automated) cleaning and disinfection process, a clinically relevant method was utilized to test the safety of residuals as part of the validation protocol. Deionized (critical) water was utilized as the terminal rinse water quality to ensure that residuals will not interfere with subsequent processing steps.

Repeated processing has minimal effect on these devices. End of life is normally determined by wear and damage due to the intended use. The user assumes liability and is responsible for the use of a damaged and dirty device.

A device labeled as a Single Use device must **never** be reused. Use refers to those single use devices that have come into contact with blood, bone, tissue, or other body fluids. Any unused single use device that has been exposed to blood, bone, tissue, or body fluids must **not be reprocessed and must be discarded**.

The instructions in this DFU were developed using the guidance given in the following standards:

- ANSI/AAMI ST79, "Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities"
- ISO 17664, "Cleaning and Disinfection of Instruments to be provided by the manufacturer for the processing of resterilizable medical devices"
- ISO 17665-1: Sterilization of health care products – Moist heat – Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices
- AAMI TR30:2011: A compendium of processes, materials, test methods, and acceptance criteria for cleaning, rinsing, and disinfection of medical devices.
- AAMI S777: Containment devices for reusable medical device sterilization

J. CLEANING AND DISINFECTION

Certain Arthrex devices that may be used during this procedure are not non-sterile and must be adequately cleaned and sterilized prior to use or re-use. All devices are to be cleaned, disinfected, and sterilized prior to each application; this is required as well for the first use after delivery of the unsterile device. Arthrex cleaning and disinfection information is to be performed following cleaning, disinfection, and sterile packaging prior to use, and may be re-sterilized (if unused) following cleaning, disinfection, and sterile packaging prior to use.

Devices that are provided in a terminally sterilized configuration should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided in non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

I. STERILE PACKAGING
Singly: Single devices should be packed as to ensure that the pack is large enough to contain the device without stressing the seals. Packaging should be completed utilizing a pouch or wrap which conforms to the recommended specifications for steam sterilization as outlined below. If a wrap is utilized, it should be completed following AAMI double-wrap or equivalent guidelines with an appropriate wrap. An appropriate wrap is one that, for example, is cleared by the FDA or the local governing body at the point of use. Devices may also be placed into an approved reusable rigid sterilization container. Aesculap SteriContainer™ rigid containers with perforated bottom lids are approved for use with Arthrex, Inc. devices.

Sets: Where appropriate, cleaned, disinfected, and inspected devices should be placed into trays/cases as provided or in general-purpose sterilization trays. The total weight of trays/cases should not exceed 11.4kg/25 lbs. (Other local limits below 11.4kg/25 lbs. may apply). Trays/cases should be double wrapped following AAMI or equivalent guidelines with an appropriate wrap. An appropriate wrap is one that, for example, is cleared by the FDA or the local governing body at the point of use. Sets may also be placed into an approved reusable rigid sterilization container. Aesculap SteriContainer™ rigid containers with perforated bottom lids and lids are approved for use with Arthrex, Inc. sets.

Areas designated for specific devices shall contain only devices intended for those areas. Devices should not be stacked or placed in close contact. Only Arthrex devices should be included in the trays or cases. These validated reprocessing instructions are not applicable to trays or cases that include devices not intended to be used with Arthrex products or sets.

II. STERILIZATION
Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in the table below. Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions.

- features that might retain liquid can drain (for example, hinges should be opened and automated/sterile workflow positioned to drain).
- Using alkaline cleaning agents, a neutralization step should be applied as appropriate.
 - Run an automated wash cycle with fundamentally approved efficiency of the washer-disinfector (for example, CE marking according to EN ISO 15883 or FDA approval/clearance/registration). The following minimum recommended automated wash cycle parameters were utilized by Arthrex during the validation of these instructions.

Phase	Recirculation Time	Temperature	Detergent
Pre-Wash	5 Minutes	Cold Water	N/A
Cleaning Wash	10 Minutes	Follow detergent manufacturer's recommendation	Enzymatic or alkaline detergent
Neutralization Rise (optional)	2 Minutes	Follow detergent manufacturer's recommendation	Neutralizing agent (as needed)
Rinse	5 Minutes	Cold Water	N/A
Thermal Disinfection Rise	5 Minutes	90°C (194°F)	N/A
Drying	Minimum 6 Minutes or until visibly dry	Minimum 90°C (212°F)	N/A

- Remove the devices from the washer-disinfector following the completion of the program and check devices for visible soil. Repeat cleaning if soil is visible and re-inspect; otherwise, proceed to Sterilization according to the instructions.

IV. MANUAL CLEANING AND DISINFECTION

Following preliminary cleaning, the instructions for Manual Cleaning and Disinfection may be followed as an alternative cleaning method to Machine (Automated) Cleaning and Thermal Disinfection if an automated procedure is not available.

- After preliminary cleaning is complete, repeat steps 1-5 provided in the Preliminary Cleaning section within this DFU, including rinsing, immersion and sonication, and post-rinsing. Final rinsing should be completed with purified (critical, e.g. RO or DI) water.
- Check devices for visible soil. Repeat cleaning if soil is visible and re-inspect.
- Soak the devices for the given soaking time (provided by the disinfectant manufacturer) in disinfectant solution so that the devices are sufficiently wetted. Make sure that there is no contact between the devices. Ensure that the device is in the open position during soaking. Actuate movable parts at least five times during soaking, as applicable.
- Remove the devices from the disinfectant solution and rinse per disinfectant manufacturer's instructions.
- Dry devices thoroughly utilizing filtered medical grade air or a soft, low-linting cloth. Proceed to Sterilization section.

K. STERILIZATION

This device may be provided either sterile or non-sterile. Check the package labeling for more information. For devices that are not provided in a terminally sterilized configuration, sterilization is to be performed following cleaning, disinfection, and sterile packaging prior to use, and may be re-sterilized (if unused) following cleaning, disinfection, and sterile packaging prior to use.

Devices that are provided in a terminally sterilized configuration should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided in non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

I. STERILE PACKAGING
Singly: Single devices should be packed as to ensure that the pack is large enough to contain the device without stressing the seals. Packaging should be completed utilizing a pouch or wrap which conforms to the recommended specifications for steam sterilization as outlined below. If a wrap is utilized, it should be completed following AAMI double-wrap or equivalent guidelines with an appropriate wrap. An appropriate wrap is one that, for example, is cleared by the FDA or the local governing body at the point of use. Devices may also be placed into an approved reusable rigid sterilization container. Aesculap SteriContainer™ rigid containers with perforated bottom lids and lids are approved for use with Arthrex, Inc. devices.

Sets: Where appropriate, cleaned, disinfected, and inspected devices should be placed into trays/cases as provided or in general-purpose sterilization trays. The total weight of trays/cases should not exceed 11.4kg/25 lbs. (Other local limits below 11.4kg/25 lbs. may apply). Trays/cases should be double wrapped following AAMI or equivalent guidelines with an appropriate wrap. An appropriate wrap is one that, for example, is cleared by the FDA or the local governing body at the point of use. Sets may also be placed into an approved reusable rigid sterilization container. Aesculap SteriContainer™ rigid containers with perforated bottom lids and lids are approved for use with Arthrex, Inc. sets.

Areas designated for specific devices shall contain only devices intended for those areas. Devices should not be stacked or placed in close contact. Only Arthrex devices should be included in the trays or cases. These validated reprocessing instructions are not applicable to trays or cases that include devices not intended to be used with Arthrex products or sets.

II. STERILIZATION
Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in the table below. Sterilizers vary in design and performance characteristics. Cycle parameters and the load configuration should always be verified against the sterilizer manufacturer's instructions.

	RECOMMENDED STEAM STERILIZATION PARAMETERS			
Cycle Type	Minimum Exposure Temperature	Minimum Exposure Time	Minimum Drying Time	Minimum Cooling Time*
US Prevacum Cycle	132°C (270°F)	4 Minutes	30 Minutes	30 Minutes
UK Prevacum Cycle	134°C (273°F)	3 Minutes	30 Minutes	30 Minutes
Prevacum Cycle (Prion Cycle)	134°C (273°F)	18 Minutes	30 Minutes	30 Minutes

*Drying times vary according to load size and should be increased for larger loads.

*Cooling times vary according to the sterilizer used, device design, temperature and humidity of ambient environment, and type of packaging used. Cooling process should comply with ANSI/AAMI ST79.

Reprocessing parameters recommended by the World Health Organization (WHO) where there is concern regarding TSE/CJD contamination.

L. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

The device is made of titanium alloy or Poly (L-Lactide).

M. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

Non-sterile metal devices should be stored in a clean, dry environment. The shelf life of non-sterile devices is not limited; the devices are manufactured from non-degradable material, which does not raise any question of device stability when stored under normal conditions.

It is the responsibility of the end-user to ensure devices, once sterilized, are stored in such a way as to maintain the sterility of the device until use. Sterile, packaged devices should be stored

in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, and temperature/humidity extremes. Sterile device packages should be carefully inspected prior to opening to ensure that package integrity has not been compromised. Maintenance of sterile package integrity is generally event related. If a sterile wrap is torn, shows any evidence of tampering, or has been exposed to moisture, the device or set must be cleaned, repackaged, and sterilized.

N. INFORMATION

- In EU Only:** Procedures carried out using these devices may be used on the general population.
- In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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