English

A. DEVICE DESCRIPTION

Pencils are available in Push Button, Rocker, or Foot switching configurations. Each pencil is supplied with a coated electrode. Electrosurgical pencils are available with or without a holster.

This is a sterile, single-use product and is not designed for re-sterilization or reprocessing. Reuse and reprocessing pose risks of patient or user infection, patient injury, and may lead to device failure.

For use with maximum peak voltage of 4600 V

B. INDICATIONS

Disposable Electrosurgical Pencils are intended to remove tissue and control bleeding by use of high-frequency electrical current. The coated electrodes are intended as alternative to uncoated stainless steel electrodes or non-stick coated electrodes for use in conventional monopolar electrosurgical accessories. The coated electrodes are intended for use in situations where monopolar electrosurgical cutting and coagulation (coag) are normally used.

C. CONTRAINDICATIONS

- Use during a procedure without a properly placed ground pad.
- Use in patients with heart pacemakers, other externally attached or implanted electronic devices.
- Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Do not use for surgeries other than those indicated.

D. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- DO NOT touch the probe tip if the device has electrical current flowing to it.
- The probe tip can remain hot and cause burns after the device is deactivated.
- Fire Hazard:** Both oxygen (O₂) and nitrous oxide (N₂O) support combustion. Avoid O₂ and N₂O enriched atmospheres. Enriched atmospheres may result in fires and burns to patients or surgical personnel.
- The metal shaft of the electrode must fit completely and securely into the pencil. If metal is visible, do not use the electrode. Exposed metal may result in electric shock or burn to the patient or surgical personnel.
- Confirm proper electrosurgical settings prior to and during a procedure. Use the lowest power settings to achieve the desired effect. If increased power settings are requested, check the patient return electrode and all accessory connections before major power settings adjustments.
- Fire/Explosion Hazard:** The following substances contribute to increased fire and explosion hazards in the operating room:
 - Oxygen enriched environments.
 - Oxidizing agents such as nitrous oxide (N2O) atmospheres.
 - Verify all anesthesia circuit connections are leak free before and during use of electrosurgery.
 - Verify endotracheal tubes are leak free and that the cuff seals properly to prevent oxygen leaks.
 - If an uncouffed tube is in use, pack the throat with sponges around the uncouffed tube.
 - If possible, stop supplemental oxygen at least one minute before and during cases of electrosurgery.
 - Alcohol-based skin prepping agents and tinctures.
 - Activate the electrosurgical unit only after vapors from skin prep solutions and tinctures have dissipated.
 - Naturally occurring flammable gases (such as methane) that may accumulate in body cavities.
- Fire Hazard:** Always place the active electrode in a clean, dry, insulated safety holster when not in use.

- Electrosurgical accessories that are activated or hot from use can cause unintended burns to the patient or surgical personnel.
- Electrosurgical accessories may cause fire or burn if placed close to or in contact with flammable materials.
- Fire Hazard:** The sparking and heating associated with electrosurgery can provide an ignition source.

- Observe fire precautions at all times:

- When using electrosurgery in the presence of gases or flammable substances, prevent pooling of fluids and the accumulation of gases under surgical drapes. Tent drapes to allow vapors to mix with room air before using electrosurgery.
- Tissue buildup (eschar) on the tip of an active electrode poses a fire hazard, especially in oxygen enriched environments such as throat or mouth procedures. With sufficient heating, eschar can become a glowing ember and pose a fire hazard both as an ignition source and as a fuel. Keep the electrode clean and free of all debris.

- Facial and other body hair is flammable.
 - Water soluble surgical lubricating jelly may be used to cover hair close to the surgical site to decrease flammability.
- Conductive fluids (e.g., blood or saline) in direct contact with an active electrode or in close proximity to any active accessory may carry electrical current and cause unintended burns to the patient. This can happen as a result of either direct coupling with the active electrode or capacitive coupling between the active electrode and the external surface of the electrode insulation. Therefore, to prevent unintended burns in the presence of conductive fluids:
 - Always keep the external surface of the active electrode away from adjacent tissue while activating the electrosurgical generator.

- Buzzing the Hemostat or other instrument may cause severe burns to the surgeon or patient, including burns to the surgeon’s hands. This practice is not recommended. Some surgeons may disregard this warning and elect to buzz the hemostat during surgical procedures. The hazards of such a practice cannot be eliminated. The following precautions may slightly reduce the risks of burn:
 - Do not “buzz the hemostat” with a Needle electrode.
 - Do not lean on the patient, the table, or the retractors while buzzing the hemostat.
 - Activate cut rather than coag. Cut has a lower voltage than coag.
 - Use the lowest power setting possible for the minimum time necessary to achieve hemostasis.
 - Activate the generator after the accessory makes contact with the hemostat. Do not arc to the hemostat.
 - Firmly grasp as much of the hemostat as possible before activating the generator. This disperses the current over a larger area and minimizes current concentration at the finger tips.
 - “Buzz the hemostat” below hand level (as close as possible to the patient) to reduce the opportunity for current to follow alternate paths through the surgeon’s hands.

- When using a coated or nonstick blade electrode, place the edge of the electrode against the hemostat or other metal instrument. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device.

- A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, compromised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device. As with other electrosurgical units, electrodes and cables can provide paths for high frequency current. Position the cables to avoid contact with the patient or other leads. Other electrical equipment may experience interference when positioned near the system.
- DO NOT bend or reshape the probe; insulation damage can result.
- DO NOT use a device that has visible insulation damage. Insulation damage can cause iatrogenic injury to tissue.
- Always keep active devices away from the patient and user when not in use. Devices should always be placed in a clean, dry, non-conductive, and highly visible area within the sterile field when not in use.

- Arthrex does not assume responsibility for problems that may occur when using an Arthrex ablation probe with a third party generator. Certain generators from other manufacturers may accept the Arthrex ablation probe pin connector. It is the responsibility of the user facility to evaluate third party generators to determine the correct operating parameters and potential hazards of using Arthrex ablation probes with these generators.
- Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-company representative, and to the health authority where the incident occurred.

E. PRECAUTIONS

- Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Always use the lowest power setting that achieves the desired surgical effect. Use the active electrode for the minimum time necessary in order to reduce the possibility of unintended burn injury.
- Do not activate electrodes while in contact with or near other instruments, including cannulas. Localized burns to the patient or physician may result.
- Do not modify or add to the insulation of active electrodes.
- Activate the electrosurgical unit only when you are ready to deliver electrosurgical current and the active tip is in view and near target tissue (especially if looking through an endoscope).
- Deactivate the electrosurgical unit before the tip leaves the surgical site.
- DO NOT RESTERILIZE.
- The electrode has a coating to reduce sticking of eschar. Cleaning the electrode with an abrasive pad may damage the electrode. If the coating or electrode is damaged, discard it.
- Needle electrodes are fragile. Handle them with care to avoid damage to the needle and injury to hospital personnel.
- Needle electrodes are designed for precise low power use during surgery. Using a needle at high power setting for long periods of time may result in damage to the needle.
- Do not bend the needle electrode.
- Securely attach all components prior to use, and check the connection between the electrode and the accessory to ensure proper fit and compatibility. Improper electrode installation may result in injury to the patient or operating room personnel by arcing at the connection.
- Do not use electrode if any nicks or cuts are present on the insulation.
- If the instrument is monopolar, a dispersive electrode should be used with the generator to prevent burns/injury to the patient.

F. PACKAGING AND LABELING

- Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbols glossary.

G. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions. Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

H. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

I. INFORMATION

- In CE accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE accepting Countries:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE accepting Countries:** There are no unacceptably residual risks or uncertainties associated with the clinical use of these devices.

J. DIRECTIONS FOR USE

Instal the Electrode

- Ensure the pencil is not connected to the generator, or the generator is Off or in the Standby mode, if available.
- Grasp the insulating sleeve on the electrode. Insert the electrode into the pencil.
- Ensure the electrode is fully inserted into the pencil. The shank and insulating sleeve should fit securely into the active accessory. If the shank and/or insulating sleeve does not fit, or the insulation will not insert, do not use this electrode/accessory combination.
- Refer to the generator instruction manual for proper procedures for connecting the active accessory to the generator.
- A tip protector covers the end of some electrodes. If a tip protector is present, remove and discard it before use. Do not use it for insulation purposes.

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