English

A. DEVICE DESCRIPTION

The NanoScope™ Handpiece Kits are designed specifically for use with the NanoScope Console and together form the NanoScope System. The disposable NanoScope Handpiece Kits are not compatible with any other camera system.

The NanoScope Handpiece provides distal LED illumination to the surgical site using a fiber optic bundle surrounding a high-resolution camera sensor. The disposable NanoScope Handpiece Kit is packaged sterile and includes the NanoScope Handpiece and endoscopic instrumentation such as cannulas, obturators, trocars, and stopcocks. The camera has buttons that can be programmed using the NanoScope Console to serve different functions such as image and video capture.

The distal tip of the camera is packaged with an attached sterile cap that is used to calibrate the endoscope prior to use. The sterile cap should not be removed until the camera is inserted into the reusable NanoScope Console and calibration is successful.

B. INDICATIONS

The NanoScope System is indicated for use in diagnostic and operative arthroscopic and endoscopic procedures to provide illumination and visualization of an interior cavity of the body through either a natural or surgical opening.

C. CONTRAINDICATIONS

- Do not use for surgeries other than those indicated.
- Only the physician can evaluate the clinical factors involved with each patient and determine if the use of this device is indicated. The physician must determine the specific technique and procedure that will accomplish the desired clinical effect.
- Do not use the NanoScope Handpiece for endoscopic procedures where cannula length may limit the effectiveness of visualization.
- Patient status or conditions to be avoided include:
 - Any active infection or blood supply limitations.
 - Blood supply limitations and previous infections, which may retard healing.
 - Conditions that tend to limit the patient’s ability or willingness to restrict activities or follow directions during the healing period.
 - Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to NanoScope handpiece insertion.

D. ADVERSE EFFECTS

- As a consequence of endoscopy, damage to surrounding tissue through insertion and manipulation for visualization of internal structures can result in injury.
- Additional damage to tissue from excessive force during obturator or trocar dissection is possible.
- The camera can be damaged or degraded by vigorous use against bony or other hard surfaces.

E. WARNINGS

- Caution: Federal law restricts this device to sale by or on the order of a physician.
- This device is intended to be used by a trained medical professional.
- The NanoScope System is designed for use by medical professionals completely familiar with the required techniques and instructions for use of the equipment. Any surgeon or medical staff member involved in the receipt, storage, preparation, use, cleaning, or care of the NanoScope System or its components, should review the full documentation prior to contact with the devices.
- Do not re-sterilize this device.
- These devices should only be used in compliance with its indications for use.
- Failure to use this device in accordance with the directions for use below may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
- Surgeons should review the product-specific surgical technique prior to performing any surgery. Failure to follow these instructions can lead to:
 - Serious injuries to the patient
 - Potential injuries to the surgical team
 - Damage or malfunction of the device or accessories
- This device shall only be used with Arthrex NanoScope Console. Use of other parts or materials may degrade safety.

- Inserting an incompatible handpiece into the camera receptacle can result in damage to the NanoScope Console.
- The handpiece tube can be bent or broken if leveraged against bone or other anatomy. Always use a cannula with the camera handpiece to minimize the risk of a bent or broken camera tube.
- Before use, the outer surfaces of the NanoScope Handpiece, cannula and user accessories should be checked to ensure there are no unintended rough edges, sharp edges or protrusions which may cause injury.
- If the NanoScope Handpiece is dropped, or cable is cut or connector or camera is damaged in any way, do not use the handpiece. Open a new, undamaged, handpiece.
- Sharp trocars and cannulas must be handled carefully.
- Do not use in an Oxygen Rich Environment.
- No modifications of this equipment are allowed.
- Detailed instructions on the use and limitations of this device should be given to the patient.
- Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
- Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. PRECAUTIONS

- Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
- Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.
- The NanoScope Handpiece is not designed to be submerged (IPX4 rated).
- Surface temperatures of the NanoScope Handpiece may rise during use.
- Do not excessively bend or kink the NanoScope Handpiece cable.
- If the NanoScope Handpiece is dropped, or the cable is cut or the connector or the camera is damaged in any way, do not use the handpiece. Notify your Arthrex Sales Representative.

G. PACKAGING AND LABELING

- NanoScope Handpiece Kits should be accepted only if the factory packaging and labeling arrive intact.
- Contact Customer Service if the package has been opened or altered.
- All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

H. STERILIZATION

- This device is provided sterile. Refer to the package label for the sterilization method.
- The NanoScope Handpiece Kit is classified as a single-use device. It is delivered in a sterile package to support a single procedure with one patient.
- Reuse may pose health and/or safety risks that can include, but are not limited to cross-infection, breakage, compromised mechanical or video image/optical performance, lack of or no function, no guarantee of proper cleaning or sterilization of the device.

I. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

J. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

K. INFORMATION

- In CE accepting Countries:** Procedures carried out using these devices may be used on the general population.
- In CE accepting Countries:** There are no clinical benefits associated with the use of these devices outweigh the known clinical risks.
- In CE accepting Countries:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

L. DIRECTIONS FOR USE

Read and understand the included quick start guide (DFU-0338-SUB-XX) and the NanoScope System Manual (DFU-0298-XX) prior to using any NanoScope Handpiece Kits. The Console and Handpiece Kit are designed to be operated as a complete system.

Users of the NanoScope System are encouraged to contact their Arthrex representatives if, in their professional judgment, they require a more comprehensive surgical technique.

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