A. DEVICE DESCRIPTION

The Arthrex Autologous Conditioned Plasma (ACP) System (PCT and Cromia Edition) is comprised of an inner syringe inside an outer syringe. The inner syringe is specifically designed to fit within and function in concert with the specifically designed outer syringe. The two syringes are connected via a Luer lock connector. A threaded screw cap connected to the Luer tip of the outer syringe can be removed prior to use. The Arthrex ACP® System (PCT and Cromia Edition) was designed to be used with benchtop centrifuges with swing-out rotors that securely accommodate ACP Double Syringes (e.g., Rotofix 32 A Benchtop Centrifuge).

B. INDICATIONS

The devices comprised of the Arthrex ACP System (PCT Edition) and Arthrex ACP System (Cromia Edition) are indicated for sterile separation of non-homogenous liquids such as autologous blood or lipoaspirate for the purpose of intra-/subdermal/subcutaneous reapplication thereof for facial skin rejuvenation (i.e., skin texture, elasticity, wrinkles, pigmentation, sagging, acne and dermal scarring); aesthetic surgery; scalp disorders; and hair loss treatment (i.e., Alopecia Areata or Androgenetic Alopecia).

This device:

- Will allow the withdraw of fluids from the human body using a commercially available cannula with Luer lock connection;
- May be centrifuged up to RCF 1100g thereafter, to transfer the supernatant into the smaller syringe by pressure balance.

C. CONTRAINDICATIONS

1. Blood supply limitations and previous infections, which may retard healing.
2. Any active infection or blood supply limitations.

D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Allergies and other reactions to device materials.
3. Hematoma.
4. Damage to blood vessels and nerve damage resulting in pain or numbness from autologous sampling.

E. WARNINGS

1. Caution: Federal law restricts this device to sale by or on the order of a physician
2. This device is intended to be used by a trained medical professional.
3. Do not re-sterilize this device.
4. Follow the manufacturer's instructions when using centrifuge. Use only benchtop centrifuges with swing-out rotors that securely accommodate ACP Double Syringes (e.g., Rotofix 32 A Benchtop Centrifuge). Outcomes using other types of centrifuges are unknown.

5. Follow the manufacturer's directions for use if using an anticoagulant (e.g., Anticoagulant Citrate Dextrose Solution A (ACD-A).
6. A device labeled as a Single Use device must never be reused. Reuse may pose health and/or safety risks to the patient that can include, but are not limited to cross-infection, breakage resulting in irretrievable fragments, comprised mechanical performance due to wear, lack of or no function, no guarantee of proper cleaning or sterilization of the device.
7. Failure to use this device in accordance with the directions for use may result in device failure, render the device unsuitable for its intended use, or compromise the procedure.
8. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
9. Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative or authorized distributor for an onsite demonstration.

G. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

H. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

I. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.
Syringe and Piston: Polypropylene, TPE (Thermoplastic Elastomer), silicone fluid
Cap: Polyethylene HD (High Density)
Small Syringe: Polypropylene, Polyisoprene rubber, silicone fluid

J. STORAGE CONDITIONS

Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

K. INFORMATION

1. **In EU Only:** Procedures carried out using these devices may be used on the general population.
2. **In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
3. **In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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