A. DEVICE DESCRIPTION

The Thrombinator™ System for use with Arthrex PRP Systems (Arthrex Angel® concentrated Platelet Rich Plasma (cPRP) System or Arthrex Double Syringe (ACP) System) includes a sterile, single use device and a filter that aids in the preparation of autologous thrombin serum. Additional equipment (i.e., winged infusion withdraw needle, syringes, syringe cap(s), anticoagulant, and wet ice) used for blood drawing and processing is available and sold separately and to be supplied by the end-user in order to operate the device as intended.

B. INDICATIONS

The Thrombinator System for use with Arthrex PRP systems (Arthrex Angel concentrated Platelet Rich Plasma (cPRP) System or Arthrex Double Syringe (ACP) System) is designed for the preparation of autologous serum from anticoagulated or non-anticoagulated peripheral blood, platelet poor plasma, or platelet rich plasma (PRP) that is to be mixed with PRP and autograft or allograft bone prior to application to a bony defect for improving handling characteristics.

C. CONTRAINDICATIONS

1. Blood supply limitations and previous infections, which may retard healing.
2. Any active infection or blood supply limitations.
3. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
4. Do not use for surgeries other than those indicated.
5. Use as a dialyzer or for dialysis with a dialysate.
6. Direct connection to patient's vascular system of circulating blood volume.

D. ADVERSE EFFECTS

1. Infections, both deep and superficial.
2. Allergies and other reactions to device materials.
3. Hematoma.
4. Damage to blood vessels and nerve damage resulting in pain or numbness from autologous sampling.
5. Delayed wound healing.

E. WARNINGS

1. Caution: Federal law restricts this device to sale by or on the order of a physician
2. This device is intended to be used by a trained medical professional.
3. The autologous thrombin serum exhibits thrombin activity and therefore must not be injected into or otherwise allowed to enter the patient's vascular system.
4. Do not re-sterilize this device.
5. This is a single use device. Reuse of this device could result in failure of the device to perform as intended

and could cause harm to the patient and/or user.

6. The safety and effectiveness of the device for bone healing and hemostasis has not been established.
7. Biohazard waste, such as explanted devices, needles and contaminated surgical equipment, should be safely disposed of in accordance with the institutions policy.
8. Serious incidents should be reported to Arthrex Inc., or an in-country representative, and to the health authority where the incident occurred.

F. MRI SAFETY INFORMATION

This device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. This device has not been tested for heating, migration or image artifact in the MR environment. The safety of the device in the MR environment is unknown.

G. PRECAUTIONS

1. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Arthrex provides detailed surgical techniques in print, video, and electronic formats. The Arthrex website also provides detailed surgical technique information and demonstrations. Or, contact your Arthrex representative for an onsite demonstration.
2. The device should be operated according to the surgical technique provided with the device.
3. Surgeons must apply their professional judgment when determining the appropriate sized device based on the specific indication, preferred surgical technique, and patient history.

H. PACKAGING AND LABELING

1. Arthrex devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description and standard designation number may be found on our website at www.arthrex.com/symbolsglossary.

I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions.

Certain Arthrex instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials.

Thrombinator Construct: Polycarbonate, silicone

Filter: Polyethylene, stainless steel, styrene-butadiene copolymer, and polycarbonate

Vent: Acrylic copolymer and nylon

Beads: Borosilicate Glass

K. STORAGE CONDITIONS
Sterile devices must be stored in the original unopened packaging, away from moisture and should not be used after the expiration date.

L. INFORMATION

1. **In EU Only:** Procedures carried out using these devices may be used on the general population.
2. **In EU Only:** The clinical benefits associated with the use of these devices outweigh the known clinical risks.
3. **In EU Only:** There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

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Thrombinator™ System for use with Arthrex PRP Systems (Arthrex Angel® concentrated Platelet Rich Plasma (cPRP) System or Arthrex Double Syringe (ACP) System)

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This is not a warranty document. For all warranty information, including disclaimers, exclusions, terms, conditions and related provisions, refer to the "Arthrex U.S. Product Warranty" section of the Arthrex, Inc. website, found at www.arthrex.com whose provisions are incorporated herein by reference.

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