

# Cartiform®

Viable Osteochondral Allograft


## Preparation Guide



1 Remove package insert, patient labels and Cartiform pouch from box.



2 Peel back chevron.



3 Using aseptic technique, transfer jar into sterile field.



4 Place the sterile jar into a sterile basin.




5 Using aseptic technique, add sterile saline until volume is just below the lid. Incubate for ~7 minutes until no ice crystals are visible. **Temperature of thawing solution should not exceed 39°C (102°F). Do not thaw for longer than 30 minutes.**



6 Take the thawed Cartiform jar from the basin and unscrew the lid.



7 Using sterile forceps, remove Cartiform from the jar and place in a sterile rinse basin at room temperature containing sterile saline for 1 minute. It can be kept in sterile saline for up to 2 hours at room temperature prior to implantation.



8 Once rinsed, Cartiform is ready to use. The side with the score mark is the bottom of the graft.

*NOTE: Graft color may vary as human articular cartilage color varies.*

**Osiris®**  
THERAPEUTICS, INC.  
Smart Medicine.® Right Now.®

Please refer to the product's Instructions For Use for a full listing of warnings and directions. Cartiform is a registered trademark of Osiris Therapeutics, Inc.

**Arthrex®**

© 2015, Arthrex Inc. All rights reserved. LS1-0860-EN\_C  
U.S. Patents Pending.