


# Minimally Invasive Surgery

Cheilectomy Surgical Technique



**Arthrex**® 

## Minimally Invasive Foot Surgery

---

The Arthrex MIS product portfolio provides surgeons with all the tools necessary to perform minimally invasive or percutaneous surgery of the foot. First, Arthrex offers a dedicated, high-quality power unit with the ideal performance parameters for MIS surgery. Next, Arthrex offers an array of disposable burrs designed specifically for the various osteotomies performed during MIS procedures. Furthermore, Arthrex offers surgeons a small, reusable instrument set complete with a reusable blade handle and other reusable rasps and elevators. Finally, in instances where bony fixation is needed, Arthrex offers a comprehensive line of cannulated, headless, fully threaded Compression FT screws that allow for stable fixation of any osteotomy.



### DrillSaw Power™ System

This ergonomic, low-speed and high-torque pencil grip driver allows for safe and effective creation of osteotomies during minimally invasive or percutaneous procedures.




### Percutaneous Burrs

The wide selection of single-use burrs are designed for various procedures including cheilectomies, exostectomies, and osteotomies.



### Minimally Invasive Instrument set

This small yet functional, reusable set cuts down on cost, offering an array of the more commonly used instruments during minimally invasive or percutaneous procedures.



### Fully Threaded Compression Screws

With available diameters ranging from 2.5 mm to 7 mm and lengths from 8 mm to 140 mm, the Compression FT screws offer surgeons a variety of options for osteotomy fixation.

## Minimally Invasive Cheilectomy




Make an incision at the dorsal one-third of the metatarsal shaft approximately 15 mm proximal to the MTP joint. Using the periosteal elevator, free the dorsal capsule and soft tissue from the bone spur at the metatarsal.



Introduce the burr through the incision and begin to resect the bone spur. Rotate the burr from lateral to medial, making sure to avoid superior migration over the spur.

## Minimally Invasive Cheilectomy



**Do vs Don't:** It is important to rotate the burr into the spur as depicted in the image on the left. Do not let the burr rotate over the top of the spur (image on right).



Upon completion of the bony debridement, remove the bony debris via manipulation and compression toward the incision site. A large syringe filled with saline may be used to remove any remaining bony debris.



Finally, use the small or large rasps to further remove any bony debris from the surrounding soft tissue. It is important to note that the rasps are not used to smooth out the bone spur, but rather to gently remove any remaining bony debris.



Final confirmation of an adequate resection should be confirmed using fluoroscopy. Dorsiflexion improvement should also be noted.



Midfoot Cheilectomy: Alternatively, bone spurring at the midfoot can also be addressed and removed through a minimally invasive approach. Use the same procedural steps as previously outlined.

## Ordering Information



### MIS Instrument Set (AR-8880S)

Product Description	Item Number
Scalpel Handle, 13 cm	<b>3KL</b>
Rasp and Blunt Elevator, small	<b>AR-8880-01</b>
Rasp and Blunt Elevator, medium	<b>AR-8880-02</b>
Combination Elevator, straight and curved ends, sharp	<b>AR-8880-03</b>
MIS Instrument Case	<b>AR-8880C</b>

### DrillSaw Highspeed 200™ Set (AR-200)

Product Description	Item Number
<b>Instruments</b>	
DrillSaw Highspeed 200 set console	<b>AR-200C</b>
Motor w/ Cable 0-15,000 rpm	<b>AR-200M</b>
Irrigation Clip	<b>AR-200SP</b>
Foot Pedal	<b>OEM06202400</b>
IV Stand	<b>OEM04005900</b>
Motor Support	<b>OEM06177800</b>

### Disposables, sterile

Product Description	Item Number
Mini Scalpel Blades, sterile, qty. 10	<b>64/ST</b>
Irrigation Tubing Set, qty. 6	<b>OEM04364100</b>
<b>Osteotomies for Lesser Toe Deformity Correction</b>	
Burr, straight, sterile, 8 mm × 2 mm	<b>AR-300-B002</b>
Burr, straight, sterile, 12 mm × 2.2 mm	<b>AR-300-B003</b>
<b>Osteotomies for Hallux Valgus Correction</b>	
Burr, straight, sterile, 13 mm × 2 mm	<b>AR-300-B001</b>
Burr, straight, sterile, 19.5 × 2 mm	<b>AR-300-B201</b>
<b>Bone Resection for Hallux Valgus/Hallux Rigidus Correction</b>	
Burr, conical, sterile, 13 mm × 4.3 mm	<b>AR-300-B101</b>
Burr, straight, sterile, 13 mm × 2.9 mm	<b>AR-300-B102</b>
Burr, oval, sterile, 15 mm × 5 mm	<b>AR-300-B103</b>
<b>Chevron Osteotomy for Calcaneal Displacement</b>	
Burr, straight, sterile, 20 mm × 3.1 mm	<b>AR-300-B202</b>

### Accessories

Product Description	Item Number
Motor w/ Cable, 3.5 m	<b>AR-200M</b>
MIS Burr Adapter, 2.35 mm	<b>AR-300B</b>
Spray Clip	<b>AR-200SP</b>

### Optional

Product Description	Item Number
Sayre Elevator	<b>AR-8954-05</b>

### Compression FT Screws

Product Description	Item Number
<b>2.5 Micro Compression FT™ Screws</b>	
8 mm-14 mm (1 mm increments)	<b>AR-8725-08H – 14H</b>
16 mm-50 mm (2 mm increments)	<b>AR-8725-16H – 50H</b>
<b>3.5 Mini Compression FT™ Screws</b>	
12 mm-60 mm (2 mm increments)	<b>AR-8730-12H – 60H</b>
<b>4.0 Standard Compression FT Screws</b>	
16 mm-60 mm (2 mm increments)	<b>AR-8740-16H – 60H</b>
<b>5.0 mm Large Compression FT Screws</b>	
20 mm-50 mm (2 mm increments)	<b>AR-8750-20H – 50H</b>
55 mm-90 mm (5 mm increments)	<b>AR-8750-55H – 90H</b>
<b>7.0 mm X-Large Compression FT Screws</b>	
35 mm-120 mm (5 mm increments)	<b>AR-8770-35H – 120H</b>
125 mm-140 mm (5 mm increments)	<b>AR-8770-125HS – 140HS</b>

### Multimedia

Product Description	Item Number
Minimally Invasive Akin Osteotomy Surgical Technique, Presented by Jorge Acevedo, MD, video	<b>VID1-01430-EN</b>
Minimally Invasive Cheilectomy Sawbone Demonstration, Presented by Jorge Acevedo, MD, video	<b>VID1-01431-EN</b>
Minimally Invasive Chevron Osteotomy Sawbone Demonstration, Presented by Jorge I. Acevedo, MD, and James McWilliam, MD, video	<b>VID1-01433-EN</b>
Minimally Invasive Calcaneal Osteotomy, Presented by Jorge I. Acevedo, MD, and James McWilliam, MD, video	<b>VID1-01434-EN</b>
Minimally Invasive Calcaneal Osteotomy Surgical Technique, Presented by Jorge I. Acevedo, MD, and James McWilliam, MD, video	<b>VID1-01405-EN</b>
Minimally Invasive Chevron and Akin Osteotomy, Presented by Jorge I. Acevedo, MD, video	<b>VID1-01406-EN</b>
Minimally Invasive Cheilectomy Surgical Technique, Presented by Jorge I. Acevedo, MD, and James McWilliam, MD, video	<b>VID1-01407-EN</b>
Minimally Invasive Bunionette, video	<b>VID1-01409-EN</b>

Products advertised in this brochure/surgical technique guide may not be available in all countries. For information on availability, please contact Arthrex Customer Service or your local Arthrex representative.



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex® products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience, and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level and/or outcomes.

View U.S. patent information at [www.arthrex.com/corporate/virtual-patent-marking](http://www.arthrex.com/corporate/virtual-patent-marking)

© 2020 Arthrex, inc. All rights reserved. | [www.arthrex.com](http://www.arthrex.com) | LT1-000012-en-US\_B