

# Distal Fibular Nonunion Case

Contributed by Brian Burgess, DPM

## Background

- 30-year-old male
- Smoker, otherwise healthy
- Weber B fracture (Fig. 1)
- Treated nonsurgically
- Patient presented with pain 5 months later

## Examination and Pre-op Observations

Five months postinjury, the fractured distal fibula showed no signs of healing (Fig. 2.1). The CT scan also showed no sign of healing with minimal bone callus at the fracture site (Fig. 2.2).

## Diagnosis/Assessment/Treatment

The diagnosis of a nonunion of the left distal fibula fracture was made after X-ray and CT scan examinations. The persistent pain the patient was experiencing along with the diagnosis of a nonunion of the distal fibula fracture warranted surgery.

## Surgical Notes/Pearls

A FibuLock® fibular nail was used to provide internal fixation of the nonunion. While preparing the intramedullary canal for the implant, the reduction was stabilized with percutaneous reduction clamps. The implant talons were activated proximally and the fracture was compressed using the compression screw feature. An additional lateral/medial locking screw was placed proximally for fixation on both sides of the fracture site.

## Six Month Follow-up

At 6-month follow-up, the patient was asymptomatic and both anterior and lateral X-ray images showed healing (Fig. 3.1 and 3.2). A 1-year follow-up was scheduled.

Figure 1. X-ray of nondisplaced Weber B fracture in the distal fibula



Figure 2.1. Five-month postinjury X-ray



Figure 2.2. Five-month postinjury CT Scan



Figure 3.1. A/P 6-month follow-up X-ray



Figure 3.2. Lateral 6-month follow-up X-ray



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's Directions For Use. Postoperative management is patient specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.



**SONOMA**<sup>®</sup>  
Orthopedic Products, Inc

Sonoma Orthopedic Products, Inc.  
1388 Busch Parkway  
Buffalo Grove, IL 60089  
847-807-4378 Phone  
847-947-8082 Fax  
www.sonomaorthopedics.com

TM Trademarks and <sup>®</sup> Registered Marks of Sonoma Orthopedic Products, Inc.  
© 2018 Sonoma Orthopedic Products, Inc. All Rights Reserved USA Patents 7,846,162; 7,909,825; 7,914,533 and 7,942,875  
USA and International Patents Pending



MediTech Strategic Consultants B.V.  
Maastrichterlaan 127-129  
6291 EN Vaals