


Univers Revers™ Total Shoulder System



Celebrating Over 10 Years of
Positive Outcomes



- Evolutions to the system over the years include options for screw or post fixation, augments, and multiple lateralization and glenosphere options
- 135° neck-shaft angle optimizes impingement-free range of motion and reduces scapular notching¹
- Inlay humeral design facilitates ideal soft-tissue tensioning²
- Fracture repair features include press-fit stem geometry, 135° inclination angle, and soft-tissue fixation points for good functional outcomes and reliable tuberosity healing³

References

1. Erickson BJ, et al. *J Shoulder Elbow Surg.* 2015;24(6):988-993.
2. Haidamous G, et al. *J Shoulder Elbow Surg.* 2020;29(12):2556-2563.
3. Schmalzl J, et al. *BMC Musculoskelet Disord.* 2020;21(1):35.



Learn more about the Univers Revers
Total Shoulder System



© 2024-03 Arthrex, Inc. All rights reserved.
FL1-000698-en-US_B