




Arthrex Small External Fixator

Small Bone Fixation



Small Pins and Rods for Your Toughest Small Bone Fractures (9933)

- 2.0 mm and 3.0 mm threaded pins
- 5.0 mm carbon fiber rods
- Single and multipin clamps

Arthrex® 

Arthrex Small External Fixator

Comprehensive Pin and Rod Options




2.0 mm Pins

- 10 mm, 15 mm, and 20 mm thread lengths with varying overall lengths of 55 mm to 75 mm




3.0 mm Pins

- 10 mm, 15 mm, and 20 mm thread lengths with varying overall lengths of 75 mm to 140 mm



5.0 mm Carbon Fiber Rods

- 5 mm × 50 mm to 300 mm length



Single Pin to Rod Fixation Clamp



Double Pin to Rod Fixation Clamp



Multipin to Rod Fixation Clamp

Ordering Information

Product Description	Item Number
Carbon Fiber Rods	
Carbon Fiber Rod, 5 mm × 50 mm to 300 mm	1181-050 – 300
Clamps	
5.0 mm Small External Fixation Clamp	1400-000
2-Pin Small Multi-Pin Clamp	1406-000
Multi-Pin Clamp Extension	1407-000
3-Pin Small Multi-Pin Clamp	1409-000
Threaded Half Pins	
2.0 mm × 10 mm × 55 mm and 75 mm	1410-055 – 075
2.0 mm × 15 mm × 55 mm and 75 mm	1411-055 – 075
2.0 mm × 20 mm × 75 mm	1412-075
3.0 mm × 10 mm × 75 mm	1414-075
3.0 mm × 15 mm × 75 mm and 100 mm	1415-075 – 100
3.0 mm × 20 mm × 75 mm and 140 mm	1416-075 – 140
3.0 mm × 25 mm × 75 mm and 140 mm	1417-075 – 140

Accessories	
Bent Post, Multi-Pin Clamp	1421-000
Straight Post, Multi-Pin Clamp	1422-000
Bolt, Multi-Pin Clamp	1423-000
Small External Fixation System	1425-000
Parallel Fixation Nut	1428-000
1.5 mm Guide Pin	0109-150
Pin Caps (4 pack)	1426-000
Calibrated Drill, AO Style, 1.5 mm × 110 mm	4001-000
Calibrated Drill, AO Style, 2.0 mm × 110 mm	4002-000
Drill Guide, 1.5 mm to 2 mm	0343-000
Wrench, AO, 3/8 in - 3/16 in	5007-000
T-Handle Wrench, 3/16 in	5009-000
T-Handle Wrench, 3/8 in	5013-000
Multi-Pin Clamp Drill Guide	0357-000
Thread-In Sleeve Multi-Pin Clamp Drill Guide	0358-000

arthrex.com

© 2022 Arthrex, Inc. All rights reserved. FL1-000516-en-US_A

