


Arthrex Calcium Sulfate BioBeads

Synthetic Bone Void Filler




Overview

Arthrex Calcium Sulfate (CS) BioBeads provide a high-purity, medical-grade calcium sulfate solution that is fully resorbed and replaced with bone during the repair process. Arthrex CS BioBeads are biocompatible and biodegradable beads for filling bony voids and they can be used in an infected site.

Features and Benefits

- Bead sizes of 3 mm, 4.5 mm, 6 mm, and as well as a paste option, for case-by-case flexibility
- High-purity, medical-grade calcium sulfate powder
- Easy-to-use mixing and delivery kit
- Fast setting time of 5 minutes
- Radiopaque for placement visibility



Arthrex Calcium Sulfate BioBeads Kit

Pic.	Qty.	Description
1	1	Sterile Mixing Bowl
2	1	Bead Mold Mat
3	1	Spatula
4	1	Calcium Sulfate Powder Vial
5	1	Liquid Vial

Ordering Information

Product Description	Item Number
Arthrex Calcium Sulfate BioBeads, 5 cc	ABS-3000-05
Arthrex Calcium Sulfate BioBeads, 10 cc	ABS-3000-10
Arthrex Calcium Sulfate BioBeads, 20 cc	ABS-3000-20

arthrex.com

© 2023-08 Arthrex, Inc. All rights reserved. FL1-000692-en-US_A

