


DiamondCut™ Burr

The Diamond Standard in Arthroscopic
Femoroacetabular Impingement (FAI) Surgery



- > 19 cm (hip length)
- > 4 mm and 5 mm diameters



Learn more about
DiamondCut burr

Arthrex® 

Hip-length DiamondCut™ burrs feature a diamond-coated head designed for precise control during resection of pincer and cam deformities associated with FAI.

- > Diamond-coated burr head for controlled bone resection
- > Designed to minimize soft-tissue wrapping on the burr head
- > Engineered to minimize skipping and provide exceptional control
- > Consistent cutting in both forward and reverse

Ordering Information

| | |
|-------------------------------|------------|
| DiamondCut burr, 4 mm × 19 cm | AR-6400RBD |
| DiamondCut burr, 5 mm × 19 cm | AR-6500RBD |

