

Synergy Vision™ Control Unit

Camera Control at Your Fingertips



Take Command of the Vision Experience

The Synergy Vision™ control unit delivers a sleek, dedicated touchscreen interface that puts effortless control of the Synergy Vision camera console at the fingertips of OR staff. Fast, intuitive, and purpose-built, it streamlines how teams interact with the Synergy Vision platform.

Effortless Touchscreen Control

A modern, responsive touchscreen allows staff to quickly access and control the camera console without breaking workflow. The Synergy Vision control unit is designed to simplify actions, reduce steps, and keep attention where it matters most—on the procedure.

Control Vision Camera Functionality

Using the Web app, the OR staff can enter patient information, select a case from their worklist, start and end a case, capture images and videos, and review the case media after the procedure. They can then print and/or export their case media to Synergy.net™ management software and their facility's EHR.

Permanent Pairing Between Console and Touchscreen

The Synergy Vision control unit supports permanent pairing between the Synergy Vision camera console and the touchscreen. This ensures that the connection between these two devices persists, even through network disruptions and power loss.

Built to Fit Any OR Setup

With flexible installation options, the Synergy Vision control unit can be configured to match any environment. The touchscreen and camera console can be connected through the facility network or connected directly using a network cable.

Powered by Arthrex Software

Running on a dedicated touchscreen PC, the Synergy Vision control unit leverages Arthrex-designed software built specifically for surgical workflows. The result: clear navigation, faster adjustments, and a more streamlined vision management experience.

A Smart Alternative Interface

The Synergy Vision control unit adds choice and flexibility to the OR by serving as a dedicated alternative to the system's wired tablet for controlling the Synergy Vision camera console.

Ordering Information

Synergy Vision control unit SW license

AR-3200-1120

