


Arthroplasty and Shoulder Fracture

Next Generation in Repair and Reconstruction



Arthrex® 

Helping Surgeons Treat Their Patients Better[®]

Since its inception, Arthrex has been committed to one mission: Helping Surgeons Treat Their Patients Better. We are strategically focused on constant product innovation through scientific research, surgeon collaboration, and medical education to make less invasive surgical procedures simple, safer and more reproducible. Each year, we develop more than 1,000 new innovative products and procedures to advance minimally invasive orthopedics worldwide.

Arthrex has always remained a privately held company, which allows for the rapid evaluation of new technologies and ideas and the freedom to develop products and techniques that truly make a difference. Our experienced team of dedicated professionals represents a shared passion and commitment to delivering uncompromising quality to the health care providers who use our products and the millions of patients whose lives we impact.

The medical significance of our contributions serves as our primary benchmark of success and will continue into the future as the legacy of Arthrex.

Table of Contents

Shoulder Arthroplasty	05
Shoulder Fracture	69
Knee Arthroplasty	89

Shoulder Arthroplasty

Solutions for Partial and Total Shoulder Replacement

Eclipse™ Total Shoulder System	06
Eclipse™ Total Shoulder Instrument Set	08
Eclipse™ Condensed Instrument Set	10
Eclipse™ SpeedScap™ Implant System (AR-9400-SBK)	11
Univers™ II/Univers Apex OptiFit™ Total Shoulder Systems	12
Univers™ II Humeral Preparation Set	14
Univers™ Apex Humeral Preparation Set	15
Univers™ II/Apex Combined Humeral Preparation Set	16
Univers VaultLock® Glenoid System	18
Keeled Glenoid	19
VaultLock® and Keel Glenoid Instrument Set	20
Univers VaultLock® Augmented Glenoid System	22
VaultLock® Augmented Glenoid Instrument Set	24
Universal Glenoid™ Convertible Baseplate	27
Universal Glenoid™ Instrument Set	28
CeMend™ Shoulder Spacer Molds	31
CeMend™ Shoulder Spacer Molds Instrument Set	32
Univers Revers™ Total Shoulder System	34
Univers Revers™ Apex Humeral Instrumentation Set	36
Univers Revers™ Humeral Broach Instrumentation Set	37
Univers Revers™ CA Head and Adapter Trials	41
Univers Revers™ Modular Glenoid System	43
Modular Glenoid System, Instrument Set 1	46
Modular Glenoid System, Instrument Set 2	48
Bone Graft Preparation for the Univers Revers™ Modular Glenoid System (AR-9665S)	50
Univers Revers™ Augmented Modular Glenoid System	53
Augmented Modular Glenoid System Instrument Set	55
Univers Revers™ Consolidated Instrument Set	57
Shoulder Arthroplasty Retractor Set	61
Virtual Implant Positioning™ (VIP™) System	64
Glenoid Targeter Instrument Set	66
VIP™ Glenoid Reaming Instrument Set	67



Stemless Anatomic Shoulder Prosthesis

- Individual anatomic reconstruction of the humeral head—the prosthesis is adjusted to the cortical rim of the humeral resection at the anatomical neck
- Humeral head can be positioned independently from humeral shaft axis—particularly important in post-traumatic arthritis cases
- Variable humeral head heights for anatomic reconstruction
- Preserved bone in the humeral shaft enables primary-like setting in the event of revision arthroplasty
- Potential for less invasive exposure
- More than a decade of long-term survivorship data¹
- Proprietary cage-screw fixation provides dynamic compression at the bone-prosthesis interface and reduces medial calcar osteolysis²
- Streamlined 3-step implantation

Humeral Heads	37 mm-55 mm (2 mm increments) Humeral heads offered in cobalt-chrome and titanium (special order)
Heights	Anatomic and Extended (+2 mm) Note: Only one head height is available for 37 mm
Cage Screw	Small: 30 mm Medium: 35 mm Large: 40 mm X-Large: 45 mm
Trunnion	37 mm-55 mm (sized with humeral head)

Eclipse™ Total Shoulder System (Cont.)


Implants

Product Description	Item Number
Trunnions	
Eclipse Trunnion, slotted, 37 mm, TPS- and CaP-coated	AR-9301-37CPC
Eclipse Trunnion, slotted, 39 mm, TPS- and CaP-coated	AR-9301-39CPC
Eclipse Trunnion, slotted, 41 mm, TPS- and CaP-coated	AR-9301-41CPC
Eclipse Trunnion, slotted, 43 mm, TPS- and CaP-coated	AR-9301-43CPC
Eclipse Trunnion, slotted, 45 mm, TPS- and CaP-coated	AR-9301-45CPC
Eclipse Trunnion, slotted, 47 mm, TPS- and CaP-coated	AR-9301-47CPC
Eclipse Trunnion, slotted, 49 mm, TPS- and CaP-coated	AR-9301-49CPC
Eclipse Trunnion, slotted, 51 mm, TPS- and CaP-coated	AR-9301-51CPC
Eclipse Trunnion, slotted, 53 mm, TPS- and CaP-coated	AR-9301-53CPC
Eclipse Trunnion, slotted, 55 mm, TPS- and CaP-coated	AR-9301-55CPC
Cage Screws	
Eclipse Cage Screw, 30 mm, small	AR-9301-01
Eclipse Cage Screw, 35 mm, medium	AR-9301-02
Eclipse Cage Screw, 40 mm, large	AR-9301-03
Eclipse Cage Screw, 45 mm, x-large	AR-9301-04
CoCr Humeral Heads	
Eclipse Humeral Head, 37 mm/16 mm	AR-9337-16
Eclipse Humeral Head, 39 mm/16 mm	AR-9339-16
Eclipse Humeral Head, 39 mm/18 mm	AR-9339-18
Eclipse Humeral Head, 41 mm/18 mm	AR-9341-16
Eclipse Humeral Head, 41 mm/16 mm	AR-9341-18
Eclipse Humeral Head, 43 mm/16 mm	AR-9343-16
Eclipse Humeral Head, 43 mm/18 mm	AR-9343-18
Eclipse Humeral Head, 45 mm/17 mm	AR-9345-17
Eclipse Humeral Head, 45 mm/19 mm	AR-9345-19
Eclipse Humeral Head, 47 mm/18 mm	AR-9347-18
Eclipse Humeral Head, 47 mm/20 mm	AR-9347-20
Eclipse Humeral Head, 49 mm/18 mm	AR-9349-18
Eclipse Humeral Head, 49 mm/20 mm	AR-9349-20
Eclipse Humeral Head, 51 mm/19 mm	AR-9351-19
Eclipse Humeral Head, 51 mm/21 mm	AR-9351-21
Eclipse Humeral Head, 53 mm/20 mm	AR-9353-20
Eclipse Humeral Head, 53 mm/22 mm	AR-9353-22
Eclipse Humeral Head, 55 mm/21 mm	AR-9355-21
Eclipse Humeral Head, 55 mm/23 mm	AR-9355-23
Titanium Humeral Heads	
Eclipse Titanium Humeral Head, 37 mm/16 mm	AR-9337-16T
Eclipse Titanium Humeral Head, 39 mm/16 mm	AR-9339-16T
Eclipse Titanium Humeral Head, 41 mm/16 mm	AR-9341-16T
Eclipse Titanium Humeral Head, 43 mm/16 mm	AR-9343-16T
Eclipse Titanium Humeral Head, 45 mm/17 mm	AR-9345-17T
Eclipse Titanium Humeral Head, 47 mm/18 mm	AR-9347-18T
Eclipse Titanium Humeral Head, 49 mm/18 mm	AR-9349-18T
Eclipse Titanium Humeral Head, 51 mm/19 mm	AR-9351-19T
Eclipse Titanium Humeral Head, 53 mm/20 mm	AR-9353-20T
Eclipse Titanium Humeral Head, 55 mm/21 mm	AR-9355-21T

Disposables

Product Description	Item Number
Eclipse and Univers™ II Pin Set	AR-9207S
Eclipse Cage Screw Sizer, sterile	AR-9401-08S
Eclipse Cage Screw Sizer, sterile, condensed set	AR-9402-08S


Eclipse™ Total Shoulder Instrument Set



Eclipse Total Shoulder Instrument Case - Top Tray (AR-9400S)

Label	Product Description	Item Number	Label	Product Description	Item Number
1	Eclipse Coring Template, 37 mm	AR-9400-37	8	Eclipse Trial Head, 37 mm/16 mm	AR-9437-16
	Eclipse Coring Template, 39 mm	AR-9400-39		Eclipse Trial Head, 39 mm/16 mm	AR-9439-16
	Eclipse Coring Template, 41 mm	AR-9400-41		Eclipse Trial Head, 39 mm/18 mm	AR-9439-18
	Eclipse Coring Template, 43 mm	AR-9400-43		Eclipse Trial Head, 41 mm/16 mm	AR-9441-16
	Eclipse Coring Template, 45 mm	AR-9400-45		Eclipse Trial Head, 41 mm/18 mm	AR-9441-18
	Eclipse Coring Template, 47 mm	AR-9400-47		Eclipse Trial Head, 43 mm/16 mm	AR-9443-16
	Eclipse Coring Template, 49 mm	AR-9400-49		Eclipse Trial Head, 43 mm/18 mm	AR-9443-18
	Eclipse Coring Template, 51 mm	AR-9400-51		Eclipse Trial Head, 45 mm/17 mm	AR-9445-17
Eclipse Coring Template, 55 mm	AR-9400-55	Eclipse Trial Head, 45 mm/19 mm		AR-9445-19	
2	Cut Guide Handle	AR-9203-10		Eclipse Trial Head, 47 mm/18 mm	AR-9447-18
3	Eclipse Cage Screw Sizer	AR-9401-08		Eclipse Trial Head, 47 mm/20 mm	AR-9447-20
4	Eclipse Resection Protector, 40 mm, small	AR-9401-40		Eclipse Trial Head, 49 mm/18 mm	AR-9449-18
	Eclipse Resection Protector, 46 mm, medium	AR-9401-46		Eclipse Trial Head, 49 mm/20 mm	AR-9449-20
	Eclipse Resection Protector, 50 mm, large	AR-9401-50		Eclipse Trial Head, 51 mm/19 mm	AR-9451-19
	Eclipse Resection Protector, 55 mm, x-large	AR-9401-55		Eclipse Trial Head, 51 mm/21 mm	AR-9451-21
5	Eclipse Trial Trunnion	AR-9401-18		Eclipse Trial Head, 53 mm/20 mm	AR-9453-20
6	Eclipse Coring Reamer	AR-9401-02		Eclipse Trial Head, 53 mm/22 mm	AR-9453-22
7	Eclipse Centering Device	AR-9401-09		Eclipse Trial Head, 55 mm/21 mm	AR-9455-21
				Eclipse Trial Head, 55 mm/23 mm	AR-9455-23


Eclipse™ Total Shoulder Instrument Set (Cont.)



Eclipse Total Shoulder Instrument Case - Bottom Tray (AR-9400S)

Label	Product Description	Item Number
1	Eclipse Screwdriver	AR-9401-03
2	Eclipse Trunnion Impactor	AR-9401-05
3	Univers™/Eclipse Orientation Pin for Resection Guide	AR-9202
4	Handle for Resection Templates	AR-9215-1-02
5	Head Impactor	AR-9203-13
6	Univers II and Eclipse Humeral Resection Block	AR-9205
7	Humeral Resection Template, left Humeral Resection Template, right	AR-9200-01L AR-9200-01R
8	Eclipse Coring Reamer, secondary	AR-9401-02L
9	Eclipse and Univers II Resection Guide, small, left Eclipse and Univers II Resection Guide, small, right Eclipse and Univers II Resection Guide, large, left Eclipse and Univers II Resection Guide, large, right	AR-9401-10 AR-9401-11 AR-9401-12 AR-9401-13
10	Bone Graft Tamp	AR-13317
11	Eclipse Humeral Head Extractor	AR-9401-17
12	Eclipse Targeting Adapter	AR-9402

Eclipse™ Condensed Instrument Set



Eclipse Condensed Instrument Set (AR-9400CES)

Label	Product Description	Item Number
1	Eclipse Sizing Pin	AR-9402-08
2	Eclipse Sizer Protector, 37 mm	AR-9402-37
	Eclipse Sizer Protector, 39 mm	AR-9402-39
	Eclipse Sizer Protector, 41 mm	AR-9402-41
	Eclipse Sizer Protector, 43 mm	AR-9402-43
	Eclipse Sizer Protector, 45 mm	AR-9402-45
	Eclipse Sizer Protector, 47 mm	AR-9402-47
	Eclipse Sizer Protector, 49 mm	AR-9402-49
	Eclipse Sizer Protector, 51 mm	AR-9402-51
	Eclipse Sizer Protector, 53 mm	AR-9402-53
	Eclipse Sizer Protector, 55 mm	AR-9402-55
3	Eclipse Protector Inserter	AR-9402-10
4	Eclipse Coring Reamer	AR-9402-02
5	Eclipse Targeting Adapter	AR-9402
6	Handle for Resection Templates	AR-9215-1-02
7	Humeral Resection Template, left	AR-9200-01R
	Humeral Resection Template, right	AR-9200-01L
8	Eclipse Trunnion Impactor	AR-9402-05
9	Eclipse Humeral Head Extractor	AR-9401-17

Label	Product Description	Item Number
10	Eclipse Screwdriver	AR-9402-03
11	Eclipse Calcar Planers, small	AR-9402-P40
	Eclipse Calcar Planers, medium	AR-9402-P48
	Eclipse Calcar Planers, large	AR-9402-P56
	Short Glove Protector	AR-9402-GP
12	Eclipse Centering Device	AR-9501-09
13	Eclipse Trial Trunion	AR-9401-18
14	Head Impactor	AR-9203-13
15	Eclipse Trial Head Spacer	AR-9402-T2
	Eclipse Trial Head, 37 mm/16 mm	AR-9437-16
	Eclipse Trial Head, 39 mm/16 mm	AR-9439-16
	Eclipse Trial Head, 41 mm/16 mm	AR-9441-16
	Eclipse Trial Head, 43 mm/16 mm	AR-9443-16
	Eclipse Trial Head, 45 mm/17 mm	AR-9445-17
	Eclipse Trial Head, 47 mm/18 mm	AR-9447-18
	Eclipse Trial Head, 49 mm/18 mm	AR-9449-18
	Eclipse Trial Head, 51 mm/19 mm	AR-9451-19
	Eclipse Trial Head, 53 mm/20 mm	AR-9453-20
	Eclipse Trial Head, 55 mm/21 mm	AR-9455-21

Eclipse™ SpeedScap™ Implant System (AR-9400-SBK)




- Contains all components needed to perform a knotless double-row subscapularis repair during stemless anatomic total shoulder arthroplasty
- FiberTak® DR and 3.9 mm SwiveLock® anchors, along with their appropriate drills, are designed for open shoulder surgery by incorporating a shorter driving and shaft length compared to those designed for arthroscopic shoulder surgery
- SutureTape link and free needle are included for easy passage of repair sutures through the subscapularis
- Two 1.3 mm SutureTapes are included for tensioning sutures during subscapularis release (peel) and rotator interval closure at the end of the operation

Product Description	Item Number
Eclipse SpeedScap Implant System	AR-9400-SBK
Eclipse Targeting Adapter	AR-9402

Univers™ II/Univers Apex OptiFit™ Total Shoulder Systems

Anatomic Shoulder Prosthesis: Individually Adaptable to Anatomical Variations of the Proximal Humerus


- Variable inclination (125°-140°) and version (±10°) for anatomic reconstruction
- Multiple head diameters and heights
- Eccentric humeral head options (3.5 mm offset)
- Humeral heads offered in cobalt-chrome and titanium (special order)
- OptiFit stem offers titanium plasma spray (TPS) coating for ingrowth and a midsize length for additional stability compared to short stems and improved alignment during insertion³



Inclination



Version



Offset

Univers™ II/Univers Apex OptiFit™ Total Shoulder Systems (Cont.)


Implants

Product Description	Item Number
Trunnion Replacement Kit	AR-9100TK
Stems	
Univers II Humeral Stem, 5 mm/105 mm	AR-9100-05P
Univers II Humeral Stem, 6 mm/115 mm	AR-9100-06P
Univers II Humeral Stem, 7 mm/140 mm	AR-9100-07P
Univers II Humeral Stem, 8 mm/143 mm	AR-9100-08P
Univers II Humeral Stem, 9 mm/146 mm	AR-9100-09P
Univers II Humeral Stem, 10 mm/148 mm	AR-9100-10P
Univers II Humeral Stem, 11 mm/151 mm	AR-9100-11P
Univers II Humeral Stem, 12 mm/151 mm	AR-9100-12P
Univers II Humeral Stem, 13 mm/151 mm	AR-9100-13P
Univers II Humeral Stem, 14 mm/151 mm	AR-9100-14P
Univers II Humeral Stem, 15 mm/151 mm	AR-9100-15P
Univers Apex OptiFit Humeral Stem, 5 mm/85 mm	AR-9100-05M
Univers Apex OptiFit Humeral Stem, 6 mm/85 mm	AR-9100-06M
Univers Apex OptiFit Humeral Stem, 7 mm/85 mm	AR-9100-07M
Univers Apex OptiFit Humeral Stem, 8 mm/85 mm	AR-9100-08M
Univers Apex OptiFit Humeral Stem, 9 mm/85 mm	AR-9100-09M
Univers Apex OptiFit Humeral Stem, 10 mm/85 mm	AR-9100-10M
Univers Apex OptiFit Humeral Stem, 11 mm/85 mm	AR-9100-11M
Univers Apex OptiFit Humeral Stem, 12 mm/85 mm	AR-9100-12M
Univers Apex OptiFit Humeral Stem, 13 mm/85 mm	AR-9100-13M
CoCr Humeral Heads	
Univers II Humeral Head, 40 mm/17 mm	AR-9140-17P
Univers II Humeral Head, 42 mm/17 mm	AR-9142-17P
Univers II Humeral Head, 44 mm/17 mm	AR-9144-17P
Univers II Humeral Head, 44 mm/19 mm	AR-9144-19P
Univers II Humeral Head, 46 mm/18 mm	AR-9146-18P
Univers II Humeral Head, 46 mm/20 mm	AR-9146-20P
Univers II Humeral Head, 48 mm/19 mm	AR-9148-19P
Univers II Humeral Head, 48 mm/21 mm	AR-9148-21P
Univers II Humeral Head, 50 mm/19 mm	AR-9150-19P
Univers II Humeral Head, 50 mm/21 mm	AR-9150-21P
Univers II Humeral Head, 52 mm/20 mm	AR-9152-20P
Univers II Humeral Head, 52 mm/22 mm	AR-9152-22P
Univers II Humeral Head, 54 mm/21 mm	AR-9154-21P
Univers II Humeral Head, 54 mm/23 mm	AR-9154-23P
Univers II Humeral Head, 56 mm/22 mm	AR-9156-22P
Univers II Humeral Head, 56 mm/24 mm	AR-9156-24P
Titanium Humeral Heads	
Univers II Titanium Humeral Head, 40 mm/17 mm	AR-9140-17T
Univers II Titanium Humeral Head, 42 mm/17 mm	AR-9142-17T
Univers II Titanium Humeral Head, 44 mm/17 mm	AR-9144-17T
Univers II Titanium Humeral Head, 46 mm/18 mm	AR-9146-18P
Univers II Titanium Humeral Head, 48 mm/19 mm	AR-9148-19T
Univers II Titanium Humeral Head, 50 mm/19 mm	AR-9150-19T
Univers II Titanium Humeral Head, 52 mm/20 mm	AR-9152-20T
Univers II Titanium Humeral Head, 54 mm/21 mm	AR-9154-21T
Univers II Titanium Humeral Head, 56 mm/22 mm	AR-9156-22T

Disposables

Product Description	Item Number
#2 FiberWire® Suture (blue), 97 cm, w/ tapered needle, 26.5 mm, ½ circle	AR-7200
#2 FiberWire Suture (blue), 97 cm, w/ reverse cutting needle, 36.6 mm, ½ circle	AR-7202
#5 FiberWire Suture (blue), 97 cm	AR-7210
#5 FiberWire Suture (blue), 97 cm, w/ conventional cutting needle, 48 mm, ½ circle	AR-7211
FiberWire Suture Kit	AR-7219
FiberTape® Tendon Compression Bridge Kit	AR-7297
Univers Apex Suture Kit	AR-7298
Eclipse™ and Univers II Pin Set	AR-9207S


Univers™ II Humeral Preparation Set



Univers II Humeral Preparation Set (AR-9226UBS)

Label	Product Description	Item Number
1	Univers II Broach, 5 mm	AR-9231-05
	Univers II Broach, 6 mm	AR-9231-06
	Univers II Broach, 7 mm	AR-9231-07
	Univers II Broach, 8 mm	AR-9231-08
	Univers II Broach, 9 mm	AR-9231-09
	Univers II Broach, 10 mm	AR-9231-10
	Univers II Broach, 11 mm	AR-9231-11
	Univers II Broach, 12 mm	AR-9231-12
	Univers II Broach, 13 mm	AR-9231-13
2	Univers II Reamer, 5 mm	AR-9202-25H
	Univers II Reamer, 6 mm	AR-9202-02H
	Univers II Reamer, 7 mm	AR-9202-01H
3	Univers II Humeral Broach Alignment Guide	AR-9232
4	Univers/Eclipse™ Orientation Pin for Resection Guide	AR-9202
5	Univers II Version Rod Bracket for Broach	AR-9231-20


Univers™ Apex Humeral Preparation Set



Univers Apex Humeral Preparation Set (AR-9226ABS)

Label	Product Description	Item Number
1	Univers Apex Broach, 5 mm	AR-9231-05S
	Univers Apex Broach, 6 mm	AR-9231-06S
	Univers Apex Broach, 7 mm	AR-9231-07S
	Univers Apex Broach, 8 mm	AR-9231-08S
	Univers Apex Broach, 9 mm	AR-9231-09S
	Univers Apex Broach, 10 mm	AR-9231-10S
	Univers Apex Broach, 11 mm	AR-9231-11S
	Univers Apex Broach, 12 mm	AR-9231-12S
	Univers Apex Broach, 13 mm	AR-9231-13S
	Univers Apex Broach, 14 mm	AR-9231-14S
	Univers Apex Broach, 15 mm	AR-9231-15S
2	Univers Apex Reamer, 5 mm	AR-9202-25HS
	Univers Apex Reamer, 6 mm	AR-9202-01HS
	Univers Apex Reamer, 7 mm	AR-9202-02HS
3	Univers Apex Humeral Broach Alignment Guide	AR-9506-07
4	Univers/Eclipse™ Orientation Pin for Resection Guide	AR-9202


Univers™ II/Apex Combined Humeral Preparation Set



Uniers II/Apex Combined Humeral Preparation - Top Tray (AR-9226CC)

Label	Product Description	Item Number
1	Handle for Resection Templates	AR-9215-1-02
2	Humeral Resection Template, left Humeral Resection Template, right	AR-9200-01L AR-9200-01R
3	Intramedullary Humeral Cutting Guide	AR-9401-14
4	Uniers II/Eclipse™ Humeral Resection Block	AR-9205
5	Uniers II Pointed Stem Impactor	AR-9202-09
6	Uniers II Morse Taper Stem Impactor	AR-9202-09P
7	Uniers Screwdriver	AR-9202-10
8	Uniers II Torque Driver	AR-9224
9	Slotted Mallet	AR-9231-21
10	Uniers Reamer T-Handle, Hudson connect	AR-9202-15H

Univers™ II/Apex Combined Humeral Set



Univers II/Apex Combined Humeral Set - Bottom Tray (AR-9226CS)

Label	Product Description	Item Number
1	Univers Apex Protector Post, 5 mm	AR-9200-15S
	Univers Apex Protector Post, 6 mm	AR-9200-16S
	Univers Apex Protector Post, 7 mm	AR-9200-17S
	Univers Apex Protector Post, 8 mm	AR-9200-18S
	Univers Apex Protector Post, 9 mm	AR-9200-19S
	Univers Apex Protector Post, 10 mm	AR-9200-20S
	Univers Apex Protector Post, 11 mm	AR-9200-21S
	Univers Apex Protector Post, 12 mm	AR-9200-22S
	Univers Apex Protector Post, 13 mm	AR-9200-23S
	Univers Apex Protector Post, 14 mm	AR-9200-24S
	Univers Apex Protector Post, 15 mm	AR-9200-25S
2	Univers II Trial Head, 40/17	AR-9240-17P
	Univers II Trial Head, 42/17	AR-9242-17P
	Univers II Trial Head, 44/17	AR-9244-19P
	Univers II Trial Head, 44/19	AR-9244-19P
	Univers II Trial Head, 46/18	AR-9246-18P
	Univers II Trial Head, 46/20	AR-9246-20P
	Univers II Trial Head, 48/19	AR-9248-19P
	Univers II Trial Head, 48/21	AR-9248-21P
	Univers II Trial Head, 50/19	AR-9250-19P
	Univers II Trial Head, 50/21	AR-9250-21P
	Univers II Trial Head, 52/20	AR-9250-20P
	Univers II Trial Head, 52/22	AR-9252-22P
	Univers II Trial Head, 54/21	AR-9254-21P
	Univers II Trial Head, 54/23	AR-9254-23P
	Univers II Trial Head, 56/22	AR-9256-22P
	Univers II Trial Head, 56/24	AR-9256-24P
3	Univers II Driver for Trial Head	AR-9202-091P
4	Humeral Head Impactor	AR-9202-13

Label	Product Description	Item Number
5	Univers II Resection Protector TRL SFT, 40 mm	AR-9202-40SP
	Univers II Resection Protector, 40 mm	AR-9202-40TP
	Univers II Resection Protector TRL SFT, 45 mm	AR-9202-45SP
	Univers II Resection Protector, 45 mm	AR-9202-45TP
6	Univers II Trial Trunnion	AR-9202-27
7	Humeral Head Extractor	AR-9401-17
8	Univers II Trunnion Extractor	AR-9202-38P
9	Univers II Stem Extractor	AR-9202-41P
10	Univers Slap Hammer	AR-9202-14

Optional

Product Description	Item Number
Version Rod Bracket	AR-9231-20
Intramedullary Humeral Cutting Guide	AR-9401-1

Univers VaultLock® Glenoid System




Fluted Central Peg

- Immediate fixation
- OR efficiency
- Bone graft compatibility

Anatomic Backside Radius of Curvature (ROC)

- Matches glenoid poly to glenoid anatomy
- Subchondral bone-sparing reaming
- Simplified decision-making



Inferior Keel

- Decreased risk of cortical penetration compared to inferior pegs
- Multiple fixation features, including reverse barbs, flutes, and central cement fenestration

Superior Peg

- Enhanced immediate fixation
- Self-pressurizing design

Inline Configuration

- Combines all advantages of pegged and keeled implants, including stability and ease of preparation


Implants

Product Description	Item Number
Univers VaultLock Glenoid, small	AR-9106-01
Univers VaultLock Glenoid, medium	AR-9106-02
Univers VaultLock Glenoid, large	AR-9106-03
Univers VaultLock Glenoid, x-large	AR-9106-04

Disposables

Product Description	Item Number
DynaNite® Virtual Implant Positioning™ Glenoid Pin, nitinol, 2.8 mm	AR-5400N-400NS
Eclipse™ and Univers™ II Pin Set	AR-9207S

Keeled Glenoid




Cemented Glenoid Reconstruction

- Keel with dual fenestration for improved cement acceptance
- Reverse barbs on the keel for enhanced anchoring

Implants

Product Description	Item Number
Keeled Glenoid, cemented, small	AR-9104-01
Keeled Glenoid, cemented, medium	AR-9104-02
Keeled Glenoid, cemented, large	AR-9104-03
Keeled Glenoid, cemented, x-large	AR-9104-04


VaultLock® and Keel Glenoid Instrument Set



Univers VaultLock and Keel Glenoid Instrument Set - Top Tray (AR-9217VKS)

Label	Product Description	Item Number
1	Univers VaultLock/Keel Glenoid Pin Guide, quick connect, small	AR-9215-VK-1
	Univers VaultLock/Keel Glenoid Pin Guide, quick connect, medium	AR-9215-VK-2
	Univers VaultLock/Keel Glenoid Pin Guide, quick connect, large	AR-9215-VK-3
	Univers VaultLock/Keel Glenoid Pin Guide, quick connect, x-large	AR-9215-VK-4
2	Steinmann Pin, 2.8 mm	AR-9207
3	Quick-Connect Handle	AR-9215-1-03
4	Modular Univers VaultLock Reamer, small	AR-9228AG
	Modular Univers VaultLock Reamer, medium	AR-9229AG
	Modular Univers VaultLock Reamer, large	AR-9230AG
	Modular Univers VaultLock Reamer, x-large	AR-9231AG
5	Modular Bone Graft Reamer	AR-9229-VLG
6	Quick-Connect Drive Shaft	AR-9617
7	Glove Protector	AR-9216-4
8	Modular Cannulated Central Drill, 6 mm	AR-9216AG
9	Modular Superior Drill, 6 mm	AR-9221AG
10	Modular Inferior Drill, 4.5 mm	AR-9239AG
11	Glenoid Broach	AR-9233
12	Glenoid Punch	AR-9234
13	Univers VaultLock Drill Guide, quick connect, small	AR-9231-VL-01
	Univers VaultLock Drill Guide, quick connect, medium	AR-9231-VL-02
	Univers VaultLock Drill Guide, quick connect, large	AR-9231-VL-03
	Univers VaultLock Drill Guide, quick connect, x-large	AR-9231-VL-04
14	Univers VaultLock Glenoid Trial, small	AR-9236-01PP
	Univers VaultLock Glenoid Trial, medium	AR-9236-02PP
	Univers VaultLock Glenoid Trial, large	AR-9236-03PP
	Univers VaultLock Glenoid Trial, x-large	AR-9236-04PP
15	Glenoid Inserter/Impactor Handle	AR-9241-2
16	Secondary Glenoid Impactor	AR-9241-00
17	Glenoid Inserter/Impactor, small	AR-9241-01
	Glenoid Inserter/Impactor, medium	AR-9241-02
	Glenoid Inserter/Impactor, large	AR-9241-03
	Glenoid Inserter/Impactor, x-large	AR-9241-04

Univers VaultLock® and Keel Glenoid Instrument Set (Cont.)



Univers VaultLock and Keel Glenoid Instrument Set - Bottom Tray (AR-9217VKS)

Label	Product Description	Item Number
1	Modular Nautilus Reamer, small Modular Nautilus Reamer, medium Modular Nautilus Reamer, large Modular Nautilus Reamer, x-large	AR-9228NRM AR-9229NRM AR-9230NRM AR-9231NRM
2	Cement Pressurizer	AR-9235
3	Keel Glenoid Drill Guide, quick connect, small Keel Glenoid Drill Guide, quick connect, medium Keel Glenoid Drill Guide, quick connect, large Keel Glenoid Drill Guide, quick connect, x-large	AR-9215-K-3 AR-9215-K-5 AR-9215-K-7 AR-9215-K-9
4	Graft Compression Tool	AR-9236GT
5	Glenoid Trial Forceps	AR-9238
6	Keel Glenoid Trial, small Keel Glenoid Trial, medium Keel Glenoid Trial, large Keel Glenoid Trial, x-large	AR-9217 AR-9218 AR-9219 AR-9219-1
7	Keel Glenoid Punch	AR-9214

Uniers VaultLock® Augmented Glenoid System




Fluted Central Peg

- Immediate fixation
- OR efficiency
- Bone graft compatibility

Anatomic Backside Radius of Curvature (ROC)

- Matches glenoid poly to glenoid anatomy
- Subchondral bone-sparing reaming
- Simplified decision-making



Superior Peg

- Enhanced immediate fixation
- Self-pressurizing design

Inline Configuration

- Combines all advantages of pegged and keeled implants including stability and preparation ease

Inferior Keel

- Decreased cortical penetration compared to inferior pegs
- Multiple fixation features, including reverse barbs, flutes, and central cement fenestration

Augmented Sizes

- 15° and 25° half-wedge options available for bone preservation



Univers VaultLock® Augmented Glenoid System (Cont.)


Implants

Product Description	Item Number
Augmented Univers VaultLock Glenoid, 15°, left, small	AR-9107-01-15L
Augmented Univers VaultLock Glenoid, 15°, left, medium	AR-9107-02-15L
Augmented Univers VaultLock Glenoid, 15°, left, large	AR-9107-03-15L
Augmented Univers VaultLock Glenoid, 15°, left, x-large	AR-9107-04-15L
Augmented Univers VaultLock Glenoid, 15°, right, small	AR-9107-01-15R
Augmented Univers VaultLock Glenoid, 15°, right, medium	AR-9107-02-15R
Augmented Univers VaultLock Glenoid, 15°, right, large	AR-9107-03-15R
Augmented Univers VaultLock Glenoid, 15°, right, x-large	AR-9107-04-15R
Augmented Univers VaultLock Glenoid, 25°, left, small	AR-9107-01-25L
Augmented Univers VaultLock Glenoid, 25°, left, medium	AR-9107-02-25L
Augmented Univers VaultLock Glenoid, 25°, left, large	AR-9107-03-25L
Augmented Univers VaultLock Glenoid, 25°, left, x-large	AR-9107-04-25L
Augmented Univers VaultLock Glenoid, 25°, right, small	AR-9107-01-25R
Augmented Univers VaultLock Glenoid, 25°, right, medium	AR-9107-02-25R
Augmented Univers VaultLock Glenoid, 25°, right, large	AR-9107-03-25R
Augmented Univers VaultLock Glenoid, 25°, right, x-large	AR-9107-04-25R

Disposables

Product Description	Item Number
Angled Reamer, small	AR-9275-S
Angled Reamer, medium	AR-9275-M
Angled Reamer, large	AR-9275-L
Angled Reamer, x-large	AR-9275-XL
Augmented VaultLock Reamer, locking, small	AR-9275T-S
Augmented VaultLock Reamer, locking, medium	AR-9275T-M
Augmented VaultLock Reamer, locking, large	AR-9275T-L
Augmented VaultLock Reamer, locking, x-large	AR-9275T-XL
DynaNite® Virtual Implant Positioning™ Glenoid Pin, nitinol, 2.8 mm	AR-5400-400NS
Eclipse™ and Univers™ II Pin Set	AR-9207S


VaultLock® Augmented Glenoid Instrument Set



Univers VaultLock Augmented Glenoid Instrument Set - Top Tray (AR-9217AVS)

Label	Product Description	Item Number
1	2.8 mm Steinmann Pin	AR-9207
2	Quick-Connect Handle	AR-9215-1-03
3	Glove Protector	AR-9216-4
4	Quick-Connect Drive Shaft	AR-9617
5	Angled Reamer Sleeve, 15° Angled Reamer Sleeve, 25°	AR-9297-15 AR-9297-25
6	6 mm Central Drill 4.5 mm Inferior Drill 6 mm Superior Drill	AR-9216AG AR-9239AG AR-9221AG
7	Modular Reamer, small Modular Reamer, medium Modular Reamer, large	AR-9228AG AR-9229AG AR-9230AG
8	Angled Reamer Drive Shaft	AR-9676
9	Modular Handle	AR-9595
10	Angled Reamer Orientation Handle	AR-9678
11	Angled Reamer Depth Stop, small Angled Reamer Depth Stop, medium Angled Reamer Depth Stop, large	AR-9298-S AR-9298-M AR-9298-L
12	Univers VaultLock Augmented Broach	AR-9233AG
13	Pegged Glenoid Punch	AR-9234
14	Angled Reamer Orientation Sleeve	AR-9679
15	Pegged Glenoid Pressurizer	AR-9235
16	Glenoid Inserter/Impactor Handle	AR-9241-2
17	Glenoid Inserter/Impactor, small Glenoid Inserter/Impactor, medium Glenoid Inserter/Impactor, large	AR-9241-01 AR-9241-02 AR-9241-03


Univers VaultLock® Augmented Glenoid Instrument Set (Cont.)



Univers VaultLock Augmented Glenoid Instrument Set - Bottom Tray (AR-9217AVS)

Label	Product Description	Item Number
1	Augmented Univers VaultLock Pin Guide, 15° left (small, medium, large) Augmented Univers VaultLock Pin Guide, 25° left (small, medium, large) Augmented Univers VaultLock Pin Guide, 15° right (small, medium, large) Augmented Univers VaultLock Pin Guide, 25° right (small, medium, large)	AR-9215-2AG/4AG/6AG-15L AR-9215-2AG/4AG/6AG-25L AR-9215-2AG/4AG/6AG-15R AR-9215-2AG/4AG/6AG-25R
2	Augmented Univers VaultLock Drill Guide, 15° left (small, medium, large) Augmented Univers VaultLock Drill Guide, 25° left (small, medium, large) Augmented Univers VaultLock Drill Guide, 15° right (small, medium, large) Augmented Univers VaultLock Drill Guide, 25° right (small, medium, large)	AR-9231-01AG/02AG/03AG-15L AR-9231-01AG/02AG/03AG-25L AR-9231-01AG/02AG/03AG-15R AR-9231-01AG/02AG/03AG-25R
3	Augmented Univers VaultLock Trial, 15° left (small, medium, large) Augmented Univers VaultLock Trial, 25° left (small, medium, large) Augmented Univers VaultLock Trial, 15° right (small, medium, large) Augmented Univers VaultLock Trial, 25° right (small, medium, large)	AR-9236-01AG/02AG/03AG-15L AR-9236-01AG/02AG/03AG-25L AR-9236-01AG/02AG/03AG-15R AR-9236-01AG/02AG/03AG-25R
4	Augmented Univers VaultLock Broach Base, 15° left Augmented Univers VaultLock Broach Base, 25° left Augmented Univers VaultLock Broach Base, 15° right Augmented Univers VaultLock Broach Base, 25° right	AR-9233AG-15L AR-9233AG-25L AR-9233AG-15R AR-9233AG-25R
5	Augmented Univers VaultLock Bone Prep Trial, small (15°, 25°) Augmented Univers VaultLock Bone Prep Trial, medium (15°, 25°) Augmented Univers VaultLock Bone Prep Trial, large (15°, 25°)	AR-9236-01AG-15BP/25BP AR-9236-02AG-15BP/25BP AR-9236-03AG-15BP/25BP
6	Graft Compression Tool	AR-9236GT
7	Glenoid Trial Forceps	AR-9238


Univers VaultLock® Augmented Glenoid XL Instrument Set



Univers VaultLock Augmented Glenoid XL Instrument Set (AR-9217AVXLS)

Label	Product Description	Item Number
1	Augmented Univers VaultLock Pin Guide, x-large, 15° left	AR-9215-8AG-15L
	Augmented Univers VaultLock Pin Guide, x-large, 15° right	AR-9215-8AG-15R
	Augmented Univers VaultLock Pin Guide, x-large, 25° left	AR-9215-8AG-25L
	Augmented Univers VaultLock Pin Guide, x-large, 25° right	AR-9215-8AG-25R
2	Augmented Univers VaultLock Drill Guide, x-large, 15° left	AR-9231-04AG-15L
	Augmented Univers VaultLock Drill Guide, x-large, 15° right	AR-9231-04AG-15R
	Augmented Univers VaultLock Drill Guide, x-large, 25° left	AR-9231-04AG-25L
	Augmented Univers VaultLock Drill Guide, x-large, 25° right	AR-9231-04AG-25R
3	Augmented Univers VaultLock Trial, x-large, 15° left	AR-9236-04AG-15L
	Augmented Univers VaultLock Trial, x-large, 15° right	AR-9236-04AG-15R
	Augmented Univers VaultLock Trial, x-large, 25° left	AR-9236-04AG-25L
	Augmented Univers VaultLock Trial, x-large, 25° right	AR-9236-04AG-25R
4	Bone Prep Trial, x-large, 15°	AR-9236-04AG-15BP
	Bone Prep Trial, x-large, 25°	AR-9236-04AG-25BP
5	Angled Reamer Depth Stop, x-large	AR-9298-XL
6	Modular Reamer, x-large	AR-9231AG
7	Glenoid Inserter/Impactor, x-large	AR-9241-04

Universal Glenoid™ Convertible Baseplate



Convertible Cementless Glenoid Reconstruction


- Baseplate may be combined with an anatomic polyethylene inlay for the advantages of anatomic reconstruction with the stability of screw fixation
- With the addition of a glenosphere to an existing baseplate, anatomic reconstructions can be conveniently converted to a reverse prosthesis
- Provides an anatomic option for patients with more severe glenoid wear patterns

Implants

Product Description	Item Number
Universal Glenoid Baseplate, small	AR-9120-01PC
Universal Glenoid Baseplate, medium	AR-9120-02PC
Universal Glenoid Baseplate, large	AR-9120-03PC
Universal Glenoid Inlay, small	AR-9121-01
Universal Glenoid Inlay, medium	AR-9121-02
Universal Glenoid Inlay, large	AR-9121-03
Universal Glenoid Inlay, small PLUS	AR-9121-04
Universal Glenoid Inlay, medium PLUS	AR-9121-05
Universal Glenoid Inlay, large PLUS	AR-9121-06
Universal Glenoid Glenosphere, 42	AR-9504L
Universal Glenoid Glenosphere, 42+4 LAT	AR-9504L-04
Universal Glenoid Glenosphere, 42+2.5 INF	AR-9504M-02
Universal Glenoid Glenosphere, 39	AR-9504M
Universal Glenoid Glenosphere, 39+4 LAT	AR-9504M-04
Universal Glenoid Glenosphere, 39+2.5 INF	AR-9504S-02
Universal Glenoid Glenosphere, 36	AR-9504S
Universal Glenoid Glenosphere, 36+4 LAT	AR-9504S-04
Universal Glenoid Glenosphere, 36+2.5 INF	AR-9504S-INF
Universal Glenoid Peripheral Locking Screw, 24 mm	AR-9145-24
Universal Glenoid Peripheral Locking Screw, 30 mm	AR-9145-30
Universal Glenoid Peripheral Locking Screw, 36 mm	AR-9145-36
Universal Glenoid Peripheral Locking Screw, 42 mm	AR-9145-42
Universal Glenoid Peripheral Locking Screw, 48 mm	AR-9145-48
Universal Glenoid Central Screw, 15 mm	AR-9165-15
Universal Glenoid Central Screw, 20 mm	AR-9165-20
Universal Glenoid Central Screw, 25 mm	AR-9165-25

Product Description	Item Number
Univers Revers™ Peripheral Screw, nonlocking, 4.5 mm × 24 mm	AR-9145-24NL
Univers Revers Peripheral Screw, nonlocking, 4.5 mm × 30 mm	AR-9145-30NL
Univers Revers Peripheral Screw, nonlocking, 4.5 mm × 36 mm	AR-9145-36NL
Univers Revers Peripheral Screw, nonlocking, 4.5 mm × 42 mm	AR-9145-42NL
Univers Revers Peripheral Screw, nonlocking, 4.5 mm × 48 mm	AR-9145-48NL
Univers Revers Central Screw, nonlocking, 6.5 mm × 15 mm	AR-9165-15NL
Univers Revers Central Screw, nonlocking, 6.5 mm × 20 mm	AR-9165-20NL
Univers Revers Central Screw, nonlocking, 6.5 mm × 25 mm	AR-9165-25NL


Universal Glenoid™ Instrument Set



Universal Glenoid Instrument Set - Top Tray (AR-9501GS)

Label	Product Description	Item Number
1	Central Screw Depth Gauge	AR-9165G
2	Primary Post Reamer	AR-9126RP
3	Final Reamer, small Final Reamer, medium Final Reamer, large	AR-9126-01 AR-9126-02 AR-9126-03
4	Correction Post Reamer	AR-9126RC
5	Universal Glenoid Baseplate Impactor (Trinity Guide)	AR-9165CDG
6	Torque Indicating Adapter	AR-9545-T15H
7	Glenoid Drill Guides, small Glenoid Drill Guides, medium Glenoid Drill Guides, large	AR-9125-1/AR-9125-10 AR-9125-2/AR-9125-20 AR-9125-3/AR-9125-30
8	2.8 mm Guide Pin, qty. 2	AR-9165K
9	Glenoid Drill Guide Handle	AR-9125H
10	Ratcheting Modular Handle	AR-1999HH
11	Baseplate Impactor Drill, 2.8 mm	AR-9165DDG
12	Baseplate Impactor Central Screw Tap	AR-9165TDG
13	T15 Driver Shaft, long	AR-9545-T15-03


Universal Glenoid™ Instrument Set (Cont.)



Universal Glenoid Instrument Set - Middle Tray (AR-9501GS)

Label	Product Description	Item Number
1	Peripheral Screw Drill Guide	AR-9145DG
2	T15 Driver Shaft, short	AR-9545-T15-01
3	T15 Driver Shaft, medium	AR-9545-T15-02
4	Peripheral Screw Drill, 2.5 mm	AR-9145K
5	Baseplate Coring Reamer, small Baseplate Coring Reamer, medium Baseplate Coring Reamer, large	AR-9127-01 AR-9127-02 AR-9127-03
6	Glenosphere Inserters	AR-9542
7	Glenoid Forceps	AR-9240
8	Glenosphere Forceps	AR-9544
9	Trial Glenospheres, 36 mm	AR-9540SM, AR-9540SM-INF, AR-9540SM-04
10	Trial Glenospheres, 39 mm	AR-9540MD, AR-9540SM-02, AR-9540MD-04
11	Trial Glenospheres, 42 mm	AR-9540LG, AR-9540MD-02, AR-9540LG-04


Universal Glenoid™ Instrument Set (Cont.)



Universal Glenoid Instrument Set - Bottom Tray (AR-9501GS)

Label	Product Description	Item Number
1	Baseplate Extractor	AR-9120E
2	Slotted Mallet	AR-9231-21
3	Glenosphere Extractor/Chisel	AR-9123GE
4	Cup Alignment Guides	AR-9532-36/39/42
5	Revers Humeral Assembly Station	AR-9532
6	Universal Glenoid Inlay Extractor	AR-9121E

CeMend™ Shoulder Spacer Molds



Arthrex Bone Cement

- Available in 4 options, including high viscosity and medium viscosity, both with and without gentamicin

Configurable Molds for Patient Optimization


- Allows for complete customization of a monoblock spacer with antibiotic dilution for complex two-stage revision cases
- 8 mm, 10 mm, 12 mm, and 14 mm diameters
- 42 mm and 48 mm head sizes
- 15 mm and 18 mm head heights
- Metal core for stability with extensions available to accommodate stem lengths from 90 mm to 170 mm

Disposables

Product Description	Item Number
Arthrex Bone Cement HV	AR-901-HV
Arthrex Bone Cement HV+G	AR-901-HVG
Arthrex Bone Cement MV	AR-901-MV
Arthrex Bone Cement MV+G	AR-901-MVG
Arthrex Bone Cement Open Mixing Bowl	AR-901-1
Arthrex Bone Cement Closed Mixing Bowl	AR-901-2
Arthrex Bone Cement Mixing Syringe	AR-901-3
Arthrex Bone Cement Syringe Gun, disposable	AR-901-4
CeMend Shoulder Spacer Mold 8-42	AR-902-0842M
CeMend Shoulder Spacer Mold 10-42	AR-902-1042M
CeMend Shoulder Spacer Mold 10-48	AR-902-1048M
CeMend Shoulder Spacer Mold 12-42	AR-902-1242M
CeMend Shoulder Spacer Mold 12-48	AR-902-1248M
CeMend Shoulder Spacer Mold 14-48	AR-902-1448M
CeMend Core Extension	AR-902-EXT




CeMend™ Shoulder Spacer Molds Instrument Set



CeMend Spacer Mold Set – Top Tray (AR-902-S)

Label	Product Description	Item Number
1	CeMend Mold Base	AR-902-01
2	CeMend Base Compression Bolt	AR-902-08
3	CeMend Head Insert, 42 mm CeMend Head Insert, 48 mm	AR-902-10 AR-902-09
4	Pliers	AR-902-12
5	Rasp	AR-902-11

CeMend™ Shoulder Spacer Molds Instrument Set (Cont.)



CeMend Spacer Mold Set – Bottom Tray (AR-902-S)

Label	Product Description	Item Number
1	CeMend Mold Insert Press, 8-42	AR-902-02
2	CeMend Mold Insert Press, 10-42	AR-902-03
3	CeMend Mold Insert Press, 10-48	AR-902-04
4	CeMend Mold Insert Press, 12-42	AR-902-05
5	CeMend Mold Insert Press, 12-48	AR-902-06
6	CeMend Mold Insert Press, 14-48	AR-902-07

Univers Revers™ Total Shoulder System

Humeral Stems

- Rectangular proximal filling stem for rotational stability
- Calcium phosphate (HA or CaP)-coated proximal stem
- 1 mm size increments for optimized press-fit
- Traditional-length stems from 95 mm to 147 mm
- Univers Revers Apex stems from 60 mm to 65 mm
- Revision length stems available at 180 mm
- 135° and 155° inclination angles in one universal stem body



SutureCups, Liners, and Spacers

- 33 mm, 36 mm, 39 mm, and 42 mm options
- Suture holes for tuberosity and rotator cuff attachment
- Calcium phosphate (CaP)-coated cup
- Centered and offset options for bone preservation and anatomic placement
- 2 polyethylene liner and 4 titanium spacer thicknesses for optional soft-tissue balancing
- Standard and constrained liner options
- Combination liners allow upsizing the glenosphere by one size with respect to the humeral cup
- Adapters available for revision to anatomic hemi-shoulder replacement with cuff arthropathy (CA) humeral head




Univers Revers™ Total Shoulder System (Cont.)

Implants

Product Description	Item Number
Stems	
Univers Revers Stem, size 4/33, monoblock 135°	AR-9501-04-135P
Univers Revers Stem, size 5	AR-9501-05P
Univers Revers Stem, size 6	AR-9501-06P
Univers Revers Stem, size 7	AR-9501-07P
Univers Revers Stem, size 8	AR-9501-08P
Univers Revers Stem, size 9	AR-9501-09P
Univers Revers Stem, size 10	AR-9501-10P
Univers Revers Stem, size 11	AR-9501-11P
Univers Revers Stem, size 12	AR-9501-12P
Univers Revers Stem, size 13	AR-9501-13P
Univers Revers Stem, size 14	AR-9501-14P
Univers Revers Stem, size 15	AR-9501-15P
Univers Revers Stem, size 6, revision (180 mm)	AR-9501-06RCPC
Univers Revers Stem, size 9, revision (180 mm)	AR-9501-09RCPC
Univers Revers Stem, size 12, revision (180 mm)	AR-9501-12RCPC
Univers Revers Apex Stem, 6	AR-9501-06S
Univers Revers Apex Stem, 7	AR-9501-07S
Univers Revers Apex Stem, 8	AR-9501-08S
Univers Revers Apex Stem, 9	AR-9501-09S
Univers Revers Apex Stem, 10	AR-9501-10S
Univers Revers Apex Stem, 11	AR-9501-11S
Univers Revers Apex Stem, 12	AR-9501-12S
Univers Revers Apex Stem, 13	AR-9501-13S
Univers Revers Apex Stem, 14	AR-9501-14S
Univers Revers Apex Stem, 15	AR-9501-15S
SutureCups	
Univers Revers SutureCup, 33 (neutral)	AR-9502F-33CPC
Univers Revers SutureCup, 33 (+2 mm left)	AR-9502F-33LCPC
Univers Revers SutureCup, 33 (+2 mm right)	AR-9502F-33RCPC
Univers Revers SutureCup, 36 (neutral)	AR-9502F-36CPC
Univers Revers SutureCup, 36 (+2 mm left)	AR-9502F-36LCPC
Univers Revers SutureCup, 36 (+2 mm right)	AR-9502F-36RCPC
Univers Revers SutureCup, 39 (neutral)	AR-9502F-39CPC
Univers Revers SutureCup, 39 (+2 mm left)	AR-9502F-39LCPC
Univers Revers SutureCup, 39 (+2 mm right)	AR-9502F-39RCPC
Univers Revers SutureCup, 42 (neutral)	AR-9502F-42CPC
Univers Revers SutureCup, 42 (+2 mm left)	AR-9502F-42LCPC
Univers Revers SutureCup, 42 (+2 mm right)	AR-9502F-42RCPC
Humeral Inserts	
Humeral Insert, 33+3 mm	AR-9503XS-03
Humeral Insert, 33+6 mm	AR-9503XS-06
Humeral Insert, 33+3 mm, constrained	AR-9503XS-03C
Humeral Insert, 33+6 mm, constrained	AR-9503XS-06C
Humeral Insert, 39+3 mm	AR-9503M-03
Humeral Insert, 39+3 mm, constrained	AR-9503M-03C
Humeral Insert, 39+6 mm	AR-9503M-06
Humeral Insert, 39+6 mm, constrained	AR-9503M-06C
Humeral Insert, 42+3 mm	AR-9503L-03
Humeral Insert, 42+3 mm, constrained	AR-9503L-03C
Humeral Insert, 42+6 mm	AR-9503L-06
Humeral Insert, 42+6 mm, constrained	AR-9503L-06C

Product Description	Item Number
Humeral Insert, 36+3 mm	AR-9503S-03
Humeral Insert, 36+3 mm, constrained	AR-9503S-03C
Humeral Insert, 36+6 mm	AR-9503S-06
Humeral Insert, 36+6 mm, constrained	AR-9503S-06C
Humeral Spacers	
Humeral Spacer, 33+6 mm	AR-9504-06
Humeral Spacer, 33+9 mm	AR-9504-09
Humeral Spacer, 33+12 mm	AR-9504-12
Humeral Spacer, 33+15 mm	AR-9504-15
Humeral Spacer, 36+6 mm	AR-9555-06
Humeral Spacer, 36+9 mm	AR-9555-09
Humeral Spacer, 36+12 mm	AR-9555-12
Humeral Spacer, 36+15 mm	AR-9555-15
Humeral Spacer, 39+6 mm	AR-9505-06
Humeral Spacer, 39+9 mm	AR-9505-09
Humeral Spacer, 39+12 mm	AR-9505-12
Humeral Spacer, 39+15 mm	AR-9505-15
Humeral Spacer, 42+6 mm	AR-9550-06
Humeral Spacer, 42+9 mm	AR-9550-09
Humeral Spacer, 42+12 mm	AR-9550-12
Humeral Spacer, 42+15 mm	AR-9550-15
Combination Humeral Inserts	
Combination Humeral Insert, 33+3/36	AR-9503-3336-3
Combination Humeral Insert, 33+6/36	AR-9503-3336-6
Combination Humeral Insert, 33+3C/36, constrained	AR-9503-3336-3C
Combination Humeral Insert, 33+6C/36, constrained	AR-9503-3336-6C
Combination Humeral Insert, 36+3/33	AR-9503-3633-3
Combination Humeral Insert, 36+6/33	AR-9503-3633-6
Combination Humeral Insert, 36+3/33, constrained	AR-9503-3633-3C
Combination Humeral Insert, 36+6/33, constrained	AR-9503-3633-6C
Combination Humeral Insert, 36+3/39	AR-9503-3639-3
Combination Humeral Insert, 36+6/39	AR-9503-3639-6
Combination Humeral Insert, 36+3/39, constrained	AR-9503-3639-3C
Combination Humeral Insert, 36+6/39, constrained	AR-9503-3639-6C
Combination Humeral Insert, 39+3/42	AR-9503-3942-3
Combination Humeral Insert, 39+6/42	AR-9503-3942-6
Combination Humeral Insert, 39+3/42, constrained	AR-9503-3942-3C
Combination Humeral Insert, 39+6/42, constrained	AR-9503-3942-6C
Combination Humeral Insert, 42+3/45	AR-9503-4245-3
Combination Humeral Insert, 42+6/45	AR-9503-4245-6
Combination Humeral Insert, 42+3/45, constrained	AR-9503-4245-3C
Combination Humeral Insert, 42+6/45, constrained	AR-9503-4245-6C


Univers Revers™ Apex Humeral Instrumentation Set



Univers Revers Apex Humeral Instrumentation Set (AR-9501-RAB)

Label	Product Description	Item Number
1	Revers Apex Broach/Trial, size 6	AR-9510-06SM
	Revers Apex Broach/Trial, size 7	AR-9510-07SM
	Revers Apex Broach/Trial, size 8	AR-9510-08SM
	Revers Apex Broach/Trial, size 9	AR-9510-09SM
	Revers Apex Broach/Trial, size 10	AR-9510-10SM
	Revers Apex Broach/Trial, size 11	AR-9510-11SM
	Revers Apex Broach/Trial, size 12	AR-9510-12SM
	Revers Apex Broach/Trial, size 13	AR-9510-13SM
	Revers Apex Broach/Trial, size 14	AR-9510-14SM
	Revers Apex Broach/Trial, size 15	AR-9510-15SM
2	Univers Revers Humeral Alignment Guide Boom	AR-9512-AG-1
3	Univers Revers Humeral Alignment Guide Body	AR-9512-AG-2
4	Univers Revers Apex Impactor/Extractor Adapter	AR-9511-2


Univers Revers™ Humeral Broach Instrumentation Set



Univers Revers Humeral Broach Instrumentation Set - Bottom Tray (AR-9501HBS)

Label	Product Description	Item Number
1	Lever-Lock Broach Handle	AR-9510-2
2	Version Rod	AR-9510-01
3	Humeral Broach Alignment Guide	AR-9232
4	Monoblock Broach, size 4/33	AR-9510-04M
5	Revers Humeral Broach, size 5 Revers Humeral Broach, size 6 Revers Humeral Broach, size 7 Revers Humeral Broach, size 8 Revers Humeral Broach, size 9 Revers Humeral Broach, size 10 Revers Humeral Broach, size 11 Revers Humeral Broach, size 12 Revers Humeral Broach, size 13	AR-9510-05MDM AR-9510-06M AR-9510-07M AR-9510-08M AR-9510-09M AR-9510-10M AR-9510-11M AR-9510-12M AR-9510-13M
6	Monoblock Stem Trial, size 4/33	AR-9523-04-135


Univers Revers™ Humeral Instrument Set 1



Univers Revers Humeral Instrument Set 1 (AR-9501HS-1)

Label	Product Description	Item Number
1	Osteotomy Guide Pin	AR-13303-2.4
2	Humeral Resection Guide Assembly	AR-9507RGDP
3	Universal T-Handle	AR-9202-15H
4	Humeral Resection Protector, extra small Humeral Resection Protector, small Humeral Resection Protector, medium Humeral Resection Protector, large	AR-9509-XS AR-9509-S AR-9509-M AR-9509-L
5	2.4 mm Guidewires	AR-1250L
6	Cannulated Drill, 6 mm	AR-1206L
7	IM Reamer, 5 mm IM Reamer, 6 mm IM Reamer, 7 mm IM Reamer, 8 mm IM Reamer, 9 mm IM Reamer, 10 mm IM Reamer, 11 mm IM Reamer, 12 mm IM Reamer, 13 mm	AR-9202-25H AR-9202-02H AR-9202-01H AR-9202-28H AR-9202-29H AR-9202-30H AR-9202-31H AR-9202-32H AR-9202-33H


Univers Revers™ Humeral Instrument Set 2



Univers Revers Humeral Instrument Set 2 - Top Tray (AR-9501HS-2)

Label	Product Description	Item Number
1	Stem/Cup Impactor/Extractor Handle	AR-9512
2	Univers Modular Reamer Pins, 135° Univers Modular Reamer Pins, 135°, right Univers Modular Reamer Pins, 135°, left	AR-9510RG-C AR-9510RG-R AR-9510RG-L
3	Reamer/Guide, size 4/33, 135°	AR-9510-04-135
4	Humeral Cup Reamer, 33 mm Humeral Cup Reamer, 36 mm Humeral Cup Reamer, 39 mm Humeral Cup Reamer, 42 mm	AR-9508-33 AR-9508-36 AR-9508-39 AR-9508-42
5	Liner/Glenosphere Impactor	AR-9531
6	Humeral Impactor/Extractor Adapter	AR-9511-2
7	Assembly Stand	AR-9532
8	Cup Screw Alignment Guide, 33 mm Cup Screw Alignment Guide, 36 mm Cup Screw Alignment Guide, 39 mm Cup Screw Alignment Guide, 42 mm	AR-9532-33 AR-9532-36 AR-9532-39 AR-9532-42
9	Pointed Impactor	AR-9202-09
10	Univers Revers Reamer Guides, size 5 Univers Revers Reamer Guides, size 5, left Univers Revers Reamer Guides, size 5, right	AR-9510RG-5C AR-9510RG-5R AR-9510RG-5L
11	T15 Driver Shaft, short	AR-9545-T15-01


Univers Revers™ Humeral Instrument Set 2 (Cont.)



Univers Revers Humeral Instrument Set 2 - Bottom Tray (AR-9501HS-2)

Label	Product Description	Item Number
1	135° Trial Cups, 33 mm	AR-9523-135/135R/135L
2	135° Trial Cups, 36 mm	AR-9522-135/135R/135L
3	135° Trial Cups, 39 mm	AR-9520-135/135R/135L
4	135° Trial Cups, 42 mm	AR-9521-135/135R/135L
5	Trial Spacers, 33 mm	AR-9554-06/09/12/15
6	Trial Spacers, 36 mm	AR-9551-06/09/12/15
7	Trial Spacers, 39 mm	AR-9552-06/09/12/15
8	Trial Spacers, 42 mm	AR-9553-06/09/12/15
9	Trial Liners, 33 mm	AR-9530XS-03/06
10	Trial Liners, 36 mm	AR-9530S-03/06
11	Trial Liners, 39 mm	AR-9530M-03/06
12	Trial Liners, 42 mm	AR-9530L-03/06
13	Constrained Trial Liners, 33 mm	AR-9530XS-03C/06C
14	Constrained Trial Liners, 36 mm	AR-9530S-03C/06C
15	Constrained Trial Liners, 39 mm	AR-9530M-03C/06C
16	Constrained Trial Liners, 42 mm	AR-9530L-03C/06C
17	Combo Trial Liners, 33 mm/36 mm	AR-9603-3336-3/6
18	Combo Trial Liners, 36 mm/39 mm	AR-9603-3639-3/6
19	Combo Trial Liners, 39 mm/42 mm	AR-9603-3942-3/6
20	Constrained Combo Trial Liners, 33 mm/36 mm	AR-9603-3336-3C/6C
21	Constrained Combo Trial Liners, 36 mm/39 mm	AR-9603-3639-3C/6C
22	Constrained Combo Trial Liners, 39 mm/42 mm	AR-9603-3942-3C/6C

Univers Revers™ CA Head and Adapter Trials



Univers Revers CA Head and Adapter Trials (AR-9501H-CA155S)

Label	Product Description	Item Number
1	155° Trial Cups, 33 mm	AR-9523-155/155R/155L
2	155° Trial Cups, 36 mm	AR-9522-155/155R/155L
3	155° Trial Cups, 39 mm	AR-9520-155/155R/155L
4	155° Trial Cups, 42 mm	AR-9521-155/155R/155L
5	Univers Modular Reamer Pins, 155°	AR-9511RG-C/R/L
6	Univers Revers CA Trial Adapter, 33 mm	AR-9522-33ARCAT
7	Univers Revers CA Trial Adapter, 36 mm	AR-9522-36ARCAT
8	Univers Revers CA Trial Adapter, 39 mm	AR-9522-39ARCAT
9	Univers Revers CA Trial Adapter, 42 mm	AR-9522-42ARCAT
10	Univers Revers CA Trial Head, 40/17 mm	AR-9540-17RCAT
11	Univers Revers CA Trial Head, 44/17 mm	AR-9544-17RCAT
12	Univers Revers CA Trial Head, 50/19 mm	AR-9550-19RCAT
13	Univers Revers CA Trial Head, 56/22 mm	AR-9556-22RCAT

Univers Revers™ CA Head and Adapter Trials (Cont.)

Disposables

Product Description	Item Number
Univers Revers Sterile Pin Set	AR-9507S

Optional Instruments

Product Description	Item Number
Humeral Resection Block, 135°/155° (superolateral)	AR-9507RGSL-1
Superolateral Version Guide Adapter	AR-9507RGSL-2
Univers Revers Reduction Tool	AR-9545
Superolateral Humeral Resection Instrument Set	AR-9507RGSL
Humeral Resection Block, 135°/155° (superolateral)	AR-9507RGSL-1
Superolateral Version Guide Adapter	AR-9507RGSL-2
Superolateral Version Downrod	AR-9507RGSL-3
Univers Revers Reduction Tool	AR-9545
Large Humeral (14/15 mm) Instrument Set	AR-9501LHS
IM Reamer, 14 mm	AR-9202-34H
IM Reamer, 15 mm	AR-9202-35H
Broach/Trial, 14 mm	AR-9510-14
Broach/Trial, 15 mm	AR-9510-15

Product Description	Item Number
Revision Humeral Instrument Set	AR-9501RHS
IM Reamer, 6 mm, revision	AR-9506-06R
IM Reamer, 9 mm, revision	AR-9506-09R
IM Reamer, 12 mm, revision	AR-9506-12R
Broach/Trial, 6 mm, revision	AR-9510-06R
Broach/Trial, 9 mm, revision	AR-9510-09R
Broach/Trial, 12 mm, revision	AR-9510-12R

Univers Revers™ Modular Glenoid System



Baseplates

- 24 mm and 28 mm diameters
- 0 mm, +2 mm, +4 mm lateral augmentation options
- BioSync® matrix with interconnected porosity for biological fixation⁴
- 4 peripheral screw holes accommodate both 5.5 mm locking and 4.5 mm nonlocking screws



Univers Revers™ Modular Glenoid System (Cont.)

Central Fixation

- Hybrid central screw for compression and ingrowth
- Central post for press-fit and ingrowth
- Numerous lengths to address various glenoid vault depths (15 mm to 35 mm)
- BioSync® matrix with interconnected porosity for biological fixation⁴



Glenospheres

- 33 mm, 36 mm, 39 mm, 42 mm, and 45 mm diameters
- +0 mm, +4 mm lateralized, and +2.5 mm eccentric offset options
- Titanium options available for patients with cobalt alloy intolerance (special order)



Univers Revers™ Modular Glenoid System (Cont.)

Implants


Product Description	Item Number
Glenoid Baseplates	
Baseplate, monoblock screw, 24 mm	AR-9560-24S
Baseplate, monoblock post, 24 mm	AR-9560-24P
Baseplate, modular, 24 mm	AR-9560-24
Lateralized Baseplate, modular, 24 mm+2	AR-9560-24-2
Lateralized Baseplate, modular, 24 mm+4	AR-9560-24-4
Baseplate, modular, 28 mm	AR-9560-28
Lateralized Baseplate, modular, 28 mm+2	AR-9560-28-2
Lateralized Baseplate, modular, 28 mm+4	AR-9560-28-4
Modular Central Fixation	
Modular Central Screw, 20 mm	AR-9561-20S
Modular Central Screw, 25 mm	AR-9561-25S
Modular Central Screw, 30 mm	AR-9561-30S
Modular Central Screw, 35 mm	AR-9561-35S
Modular Central Post, 20 mm	AR-9561-20P
Modular Central Post, 25 mm	AR-9561-25P
Modular Central Post, 30 mm	AR-9561-30P
Modular Central Post, 35 mm	AR-9561-35P
4.5 mm Nonlocking Peripheral Bone Screws	
Nonlocking Screws, 4.5 mm × 16 mm	AR-9562-16NL
Nonlocking Screws, 4.5 mm × 20 mm	AR-9562-20NL
Nonlocking Screws, 4.5 mm × 24 mm	AR-9562-24NL
Nonlocking Screws, 4.5 mm × 28 mm	AR-9562-28NL
Nonlocking Screws, 4.5 mm × 32 mm	AR-9562-32NL
Nonlocking Screws, 4.5 mm × 36 mm	AR-9562-36NL
Nonlocking Screws, 4.5 mm × 40 mm	AR-9562-40NL
Nonlocking Screws, 4.5 mm × 44 mm	AR-9562-44NL
Nonlocking Screws, 4.5 mm × 48 mm	AR-9562-48NL
5.5 mm Locking Peripheral Bone Screws	
Locking Screws, 5.5 mm × 16 mm	AR-9563-16
Locking Screws, 5.5 mm × 20 mm	AR-9563-20
Locking Screws, 5.5 mm × 24 mm	AR-9563-24
Locking Screws, 5.5 mm × 28 mm	AR-9563-28
Locking Screws, 5.5 mm × 32 mm	AR-9563-32
Locking Screws, 5.5 mm × 36 mm	AR-9563-36
Locking Screws, 5.5 mm × 40 mm	AR-9563-40
Locking Screws, 5.5 mm × 44 mm	AR-9563-44
Locking Screws, 5.5 mm × 48 mm	AR-9563-48
CoCr Glenospheres	
Glenosphere, 33/24	AR-9564-2433
Glenosphere, 33+4 LAT/24	AR-9564-2433-LAT
Glenosphere, 36/24	AR-9564-2436
Glenosphere, 36+4 LAT/24	AR-9564-2436-LAT
Glenosphere, 36+2.5 INF/24	AR-9564-2436-INF
Glenosphere, 39/24	AR-9564-2439
Glenosphere, 39+4 LAT/24	AR-9564-2439-LAT
Glenosphere, 39+2.5 INF/24	AR-9564-2439-INF

Product Description	Item Number
CoCr Glenospheres	
Glenosphere, 42/24	AR-9564-2442
Glenosphere, 42+4 LAT/24	AR-9564-2442-LAT
Glenosphere, 42+2.5 INF/24	AR-9564-2442-INF
Glenosphere, 45/24	AR-9564-2445
Glenosphere, 45+4 LAT/24	AR-9564-2445-LAT
Glenosphere, 45+2.5 INF/24	AR-9564-2445-INF
Glenosphere, 36/28	AR-9564-2836
Glenosphere, 36+4, LAT/28	AR-9564-2836-LAT
Glenosphere, 39/28	AR-9564-2839
Glenosphere, 39+4 LAT/28	AR-9564-2839-LAT
Glenosphere, 39+2.5 INF/28	AR-9564-2839-INF
Glenosphere, 42/28	AR-9564-2842
Glenosphere, 42 LAT/28	AR-9564-2842-LAT
Glenosphere, 42+2.5 INF/28	AR-9564-2842-INF
Glenosphere, 45/28	AR-9564-2845
Glenosphere, 45+4 LAT/28	AR-9564-2845-LAT
Glenosphere, 45+2.5 INF/28	AR-9564-2845-INF
Titanium Glenospheres	
33/24 Ti Glenosphere	AR-9564-T2433
33+4 LAT/24 Ti Glenosphere	AR-9564-T2433-LAT
36/24 Ti Glenosphere	AR-9564-T2436
36+4 LAT/24 Ti Glenosphere	AR-9564-T2436-LAT
36+2.5 INF/24 Ti Glenosphere	AR-9564-T2436-INF
39/24 Ti Glenosphere	AR-9564-T2439
39+4 LAT/24 Ti Glenosphere	AR-9564-T2439-LAT
39+2.5 INF/24 Ti Glenosphere	AR-9564-T2439-INF
42/24 Ti Glenosphere	AR-9564-T2442
42+4 LAT/24 Ti Glenosphere	AR-9564-T2442-LAT
42+2.5 INF/24 Ti Glenosphere	AR-9564-T2442-INF

Disposables

Product Description	Item Number
Glenosphere Screw, extra	AR-9565
MGS Sterile Pin Kit	AR-9607S
3 mm MGS Sterile Drill Bit	AR-9628S

Modular Glenoid System, Instrument Set 1




Modular Glenoid System, Instrument Set 1 - Top Tray (AR-9615C-1)

Label	Product Description	Item Number
1	Glenoid Scrapette	AR-9601
2	Pin Guide, 24 mm Pin Guide, 28 mm	AR-9616-24 AR-9616-28
3	Glove Protector Sleeve	AR-9216-4
4	Modular Reamer Shaft	AR-9617
5	Guide Pin, 2.8 mm	AR-9616
6	Guide Pin Depth Gauge	AR-9616-DG
7	Manual Driver Shaft	AR-9631
8	Torque-Indicating Adapter	AR-9545-T15H
9	Glenosphere Sizing Guide	AR-9664-SG
10	Primary Glenoid Reamer, 24 mm Primary Glenoid Reamer, 28 mm	AR-9618-24 AR-9618-28
11	Peripheral Reamer, 33 mm Peripheral Reamer, 36 mm Peripheral Reamer, 39 mm Peripheral Reamer, 42 mm	AR-9619-33 AR-9619-36 AR-9619-39 AR-9619-42

Label	Product Description	Item Number
12	Inferior Glenoid Reamer	AR-9632-INF
13	Ratcheting Handle, hybrid Hudson	AR-1999HH
14	Tap Caddy Drill, 10 mm Tap, 15 mm Tap, 20 mm Tap, 25 mm Tap, 30 mm Tap, 35 mm	AR-9615C-TC AR-9620-10D AR-9621-15T AR-9621-20T AR-9621-25T AR-9621-30T AR-9621-35T
15	Drill Caddy Drill, 15 mm Drill, 20 mm Drill, 25 mm Drill, 30 mm Drill, 35 mm	AR-9615C-DC AR-9620-15D AR-9620-20D AR-9620-25D AR-9620-30D AR-9620-35D


Modular Glenoid System, Instrument Set 1 (Cont.)



Modular Glenoid System, Instrument Set 1 - Bottom Tray (AR-9615C-1)

Label	Product Description	Item Number
1	Taper Assembly Press	AR-9622
2	Baseplate Inserter, threaded	AR-9623
3	Baseplate Inserter, cannulated	AR-9623-C
4	Drill, 3 mm	AR-9628
5	Drill Guide, locking screw	AR-9626
6	Drill Guide, nonlocking screw	AR-9627
7	Peripheral Screw Depth Gauge	AR-9629
8	Speed Guide for 24 mm Baseplate Speed Guide for 28 mm Baseplate	AR-9630-24NL AR-9630-28NL
9	Glenosphere Inserter	AR-9624
10	Driver, Hex, 3 mm	AR-9625
11	Combination Humeral Trial, 36+3/33 Combination Humeral Trial, 36+6/33	AR-9603-3633-3 AR-9603-3633-6
12	Constrained Combination Humeral Trial, 36+3/33 Constrained Combination Humeral Trial, 36+6/33	AR-9603-3633-3C AR-9603-3633-6C

Modular Glenoid System, Instrument Set 2



Modular Glenoid System, Instrument Set 2 - Top Tray (AR-9615C-2)

Label	Product Description	Item Number
1	Mallet, slotted	AR-9231-21
2	Baseplate Extractor and Extractor Shaft	AR-9659-R
3	Glenosphere Distractor	AR-9658-R
4	Central Post Extractor	AR-9660-R
5	Central Post/Screw Trephine Adapter	AR-9661-R
6	Central Post/Screw Trephine	AR-9662-R
7	Driver for Central Screw, 4 mm	AR-9663-R
8	Baseplate Extractor Wrench	AR-9659-RW

Optional: MGS Size 45 Instrument Set (AR-9615-45S)

Label	Product Description	Item Number
9	Trial Glenosphere, 45/24	AR-9664-2445
	Trial Glenosphere, 45+4 LAT/24	AR-9664-2445-LAT
	Trial Glenosphere, 45+2.5 INF/24	AR-9664-2445-INF
10	Trial Glenosphere, 45/28	AR-9664-2845
	Trial Glenosphere, 45+4 LAT/28	AR-9664-2845-LAT
	Trial Glenosphere, 45+2.5 INF/28	AR-9664-2845-INF
11	Over Baseplate Reamer, 45/24	AR-9632-2445
	Over Baseplate Reamer, 45/28	AR-9632-2845
12	45 mm Peripheral Reamer	AR-9619-45
13	Inferior Glenoid Reamer, 45	AR-9632-45INF
14	Combination Humeral Trial, 42+3/45	AR-9603-4245-3
	Combination Humeral Trial, 42+6/45	AR-9603-4245-6

Modular Glenoid System, Instrument Set 2 (Cont.)



Modular Glenoid System, Instrument Set 2 - Bottom Tray (AR-9615C-2)

Label	Product Description	Item Number
1	Trial Glenosphere, 33/24	AR-9664-2433
	Trial Glenosphere, 33+4 LAT/24	AR-9664-2433-LAT
2	Trial Glenosphere, 36/24	AR-9664-2436
	Trial Glenosphere, 36+4 LAT/24	AR-9664-2436-LAT
	Trial Glenosphere, 36+2.5 INF/24	AR-9664-2436-INF
3	Trial Glenosphere, 39/24	AR-9664-2439
	Trial Glenosphere, 39+4 LAT/24	AR-9664-2439-LAT
	Trial Glenosphere 39+2.5 INF/24	AR-9664-2439-INF
4	Trial Glenosphere, 42/24	AR-9664-2442
	Trial Glenosphere, 42+4 LAT/24	AR-9664-2442-LAT
	Trial Glenosphere, 42+2.5 INF/24	AR-9664-2442-INF
5	Trial Glenosphere, 36/28	AR-9664-2836
	Trial Glenosphere, 36+4 LAT/28	AR-9664-2836-LAT
6	Trial Glenosphere, 39/28	AR-9664-2839
	Trial Glenosphere, 39+4 LAT/28	AR-9664-2839-LAT
	Trial Glenosphere, 39+2.5 INF/28	AR-9664-2839-INF
7	Trial Glenosphere, 42/28	AR-9664-2842
	Trial Glenosphere, 42+4 LAT/28	AR-9664-2842-LAT
	Trial Glenosphere, 42+2.5 INF/28	AR-9664-2842-INF
8	Over Baseplate Reamer, 36/24	AR-9632-2436
	Over Baseplate Reamer, 39/24	AR-9632-2439
	Over Baseplate Reamer, 42/24	AR-9632-2442
	Over Baseplate Reamer, 36/28	AR-9632-2836
	Over Baseplate Reamer, 39/28	AR-9632-2839
	Over Baseplate Reamer, 42/28	AR-9632-2842

Bone Graft Preparation for the Univers Revers™ Modular Glenoid System (AR-9665S)




Technique Precision in Harvesting, Shaping, and Preparing Bony Augments

- Humeral sizing instruments available in both 135° and 155°
- Glenoid defect gauge to determine the angle/ thickness of the graft required
- Adjustable variable drill for precision in central vault drilling length

Disposables

Product Description	Item Number
Angled Reamer, small	AR-9675-S
Angled Reamer, medium	AR-9675-M
Angled Reamer, large	AR-9675-L
Angled Reamer, x-large	AR-9675-XL


Bone Graft Instrument Set



Bone Graft Instrument Set - Top Tray (AR-9665S)

Label	Product Description	Item Number
1	Pin Guide, 135°, 2.8 mm, small	AR-9672-135S
	Pin Guide, 135°, 2.8 mm, large	AR-9672-135L
2	Pin Guide, 155°, 2.8 mm, small	AR-9672-155S
	Pin Guide, 155°, 2.8 mm, large	AR-9672-155L
3	Pin Guide, handle, 2.8 mm	AR-9672
4	Univers Revers Version Rod	AR-9510-01
5	Guide Pin, 2.8 mm	AR-9165K
6	Univers Revers Flat Reamer, cannulated, 33 mm diameter, small	AR-9128RF-01
7	Hole Saw, cannulated, Ø25 mm, w/ Ø9.5 mm, central bit	AR-9669-25
	Hole Saw, cannulated, Ø30 mm, w/ Ø9.5 mm, central bit	AR-9669-30

Bone Graft Instrument Set (Cont.)



Bone Graft Instrument Set - Bottom Tray (AR-9665S)

Label	Product Description	Item Number
1	Glenoid Defect Gauge, housing	AR-9673
2	Graft Preparation Station, cart	AR-9670-2
3	Glenoid Defect Gauge, Ø25 mm graft angle leg Glenoid Defect Gauge, Ø25 mm graft thickness leg	AR-9673-25A AR-9673-25T
4	Glenoid Defect Gauge, Ø30 mm graft angle leg Glenoid Defect Gauge, Ø30 mm graft thickness leg	AR-9673-30A AR-9673-30T
5	Graft Preparation Station, base	AR-9670-1
6	Reamer Adapter, 10° Reamer Adapter, 12.5° Reamer Adapter, 15° Reamer Adapter, 17.5° Reamer Adapter, 20° Reamer Adapter, 22.5° Reamer Adapter, 25°	AR-9674-10 AR-9674-125 AR-9674-15 AR-9674-175 AR-9674-20 AR-9674-225 AR-9674-25
7	Modular Reamer Sleeve	AR-9674
8	Angled Reamer, orientation sleeve	AR-9679
9	Angled Reamer, orientation handle	AR-9678
10	Angled Reamer, drive shaft	AR-9676
11	Variable Depth Central Post Drill, shaft	AR-9677
12	Variable Depth Central Post Drill, adjustment collar	AR-9677-1
13	Graft Preparation Raise Plate, Ø25 mm graft Graft Preparation Raise Plate, Ø30 mm graft	AR-9670-3 AR-9670-4

Univers Revers™ Augmented Modular Glenoid System



Baseplates

- 10° and 20° full-wedge
- 15°, 25°, and 35° half-wedge
- 24 mm and 28 mm diameters
- Available with 2 mm of additional lateral offset
- Standard or oblique augment orientation
- BioSync® matrix with interconnected porosity for biological fixation⁴
- 4 peripheral screw holes accommodate both 5.5 mm locking and 4.5 mm nonlocking screws

Central Fixation

- Central post for press-fit and ingrowth
- Numerous lengths to address various glenoid vault depths (20 mm to 40 mm)
- BioSync matrix with interconnected porosity for biological fixation⁴

Glenospheres

- Compatible with all MGS glenospheres



Standard augment orientation



Oblique augment orientation

Univers Revers™ Augmented Modular Glenoid System (Cont.)

Implants


Product Description	Item Number
Full-Wedge	
24 mm Baseplate, 10° full-wedge, oblique	AR-9580-2410
24 mm Baseplate, 20° full-wedge, oblique	AR-9580-2420
24 mm Baseplate, 10° full-wedge, +2 lateralized, oblique	AR-9580-2410-2
24 mm Baseplate, 20° full-wedge, +2 lateralized, oblique	AR-9580-2420-2
24 mm Baseplate, 10° full-wedge	AR-9580-2410S
24 mm Baseplate, 20° full-wedge	AR-9580-2420S
24 mm Baseplate, 10° full-wedge, +2 lateralized	AR-9580-2410-2S
24 mm Baseplate, 20° full-wedge, +2 lateralized	AR-9580-2420-2S
28 mm Baseplate, 10° full-wedge, oblique	AR-9580-2810
28 mm Baseplate, 20° full-wedge, oblique	AR-9580-2820
28 mm Baseplate, 10° full-wedge, +2 lateralized, oblique	AR-9580-2810-2
28 mm Baseplate, 20° full-wedge, +2 lateralized, oblique	AR-9580-2820-2
28 mm Baseplate, 10° full-wedge	AR-9580-2810S
28 mm Baseplate, 20° full-wedge	AR-9580-2820S
28 mm Baseplate, 10° full-wedge, +2 lateralized	AR-9580-2810-2S
28 mm Baseplate, 20° full-wedge, +2 lateralized	AR-9580-2820-2S
Half-Wedge	
24 mm Baseplate, 15° half-wedge, oblique	AR-9581-2415
24 mm Baseplate, 25° half-wedge, oblique	AR-9581-2425
24 mm Baseplate, 35° half-wedge, oblique	AR-9581-2435
24 mm Baseplate, 15° half-wedge, +2 lateralized, oblique	AR-9581-2415-2
24 mm Baseplate, 25° half-wedge, +2 lateralized, oblique	AR-9581-2425-2
24 mm Baseplate, 35° half-wedge, +2 lateralized, oblique	AR-9581-2435-2
24 mm Baseplate, 15° half-wedge	AR-9581-2415S
24 mm Baseplate, 25° half-wedge	AR-9581-2425S
24 mm Baseplate, 35° half-wedge	AR-9581-2435S
24 mm Baseplate, 15° half-wedge, +2 lateralized	AR-9581-2415-2S
24 mm Baseplate, 25° half-wedge, +2 lateralized	AR-9581-2425-2S
24 mm Baseplate, 35° half-wedge, +2 lateralized	AR-9581-2435-2S
28 mm Baseplate, 15° half-wedge, oblique	AR-9581-2815
28 mm Baseplate, 25° half-wedge, oblique	AR-9581-2825
28 mm Baseplate, 35° half-wedge, oblique	AR-9581-2835
28 mm Baseplate, 15° half-wedge, +2 lateralized, oblique	AR-9581-2815-2
28 mm Baseplate, 25° half-wedge, +2 lateralized, oblique	AR-9581-2825-2
28 mm Baseplate, 35° half-wedge, +2 lateralized, oblique	AR-9581-2835-2
28 mm Baseplate, 15° half-wedge	AR-9581-2815S
28 mm Baseplate, 25° half-wedge	AR-9581-2825S
28 mm Baseplate, 35° half-wedge	AR-9581-2835S
28 mm Baseplate, 15° half-wedge, +2 lateralized	AR-9581-2815-2S
28 mm Baseplate, 25° half-wedge, +2 lateralized	AR-9581-2825-2S
28 mm Baseplate, 35° half-wedge, +2 lateralized	AR-9581-2835-2S

Product Description	Item Number
Central Posts	
Modular Post, 20 mm	AR-9582-20
Modular Post, 25 mm	AR-9582-25
Modular Post, 30 mm	AR-9582-30
Modular Post, 35 mm	AR-9582-35
Modular Post, 40 mm	AR-9582-40

Disposables

Product Description	Item Number
Angled Reamer, small	AR-9675-S
Angled Reamer, medium	AR-9675-M
Angled Reamer, large	AR-9675-L
Angled Reamer, x-large	AR-9675-XL
Augmented MGS Reamer, locking, small	AR-9675T-S
Augmented MGS Reamer, locking, medium	AR-9675T-M
Augmented MGS Reamer, locking, large	AR-9675T-L
Augmented MGS Reamer, locking, x-large	AR-9675T-XL
3 mm Drill	AR-9628S


Augmented Modular Glenoid System Instrument Set



Augmented Modular Glenoid System Instrument Set - Top Tray (AR-9579S)

Label	Product Description	Item Number
1	Modular Handle	AR-9595
2	Pilot Drill	AR-9596-P
3	Angled Reamer, orientation handle	AR-9678
4	Augmented Reamer Drive Shaft, locking	AR-9676T
5	Angled Reamer, orientation sleeve	AR-9679
6	Baseplate Inserter	AR-9623-A
7	Full-Wedge Augment Depth Gauge	AR-9599F
8	Full-Wedge Augment Sizer and Trial Caddy 24 mm Baseplate, 10° full-wedge augment sizer 24 mm Baseplate, 20° full-wedge augment sizer 28 mm Baseplate, 10° full-wedge augment sizer 28 mm Baseplate, 20° full-wedge augment sizer	AR-9579C-FAC AR-9592-2410 AR-9592-2420 AR-9592-2810 AR-9592-2820
9	24 mm Baseplate, 10° full-wedge augment trial 24 mm Baseplate, 20° full-wedge augment trial 28 mm Baseplate, 10° full-wedge augment trial 28 mm Baseplate, 20° full-wedge augment trial	AR-9594-2410 AR-9594-2420 AR-9594-2810 AR-9594-2820
10	Angled Reamer Sleeve, 10°	AR-9597-10
11	Angled Reamer Sleeve, 20°	AR-9597-20
12	Modular Drills, 10° full, 20 mm Modular Drills, 10° full, 25 mm Modular Drills, 10° full, 30 mm Modular Drills, 10° full, 35 mm Modular Drills, 10° full, 40 mm Modular Drills, 20° full, 20 mm Modular Drills, 20° full, 25 mm Modular Drills, 20° full, 30 mm Modular Drills, 20° full, 35 mm Modular Drills, 20° full, 40 mm	AR-9596-1020F AR-9596-1025F AR-9596-1030F AR-9596-1035F AR-9596-1040F AR-9596-2020F AR-9596-2025F AR-9596-2030F AR-9596-2035F AR-9596-2040F

Augmented Modular Glenoid System Instrument Set (Cont.)



Augmented Modular Glenoid System Instrument Set - Bottom Tray (AR-9579S)


Label	Product Description	Item Number
1	Paleo Reamer	AR-9596
2	Trephine Adapter	AR-9661-A
3	Half-Wedge Depth Gauge	AR-9599H
4	24 mm Baseplate, 15° half-wedge augment sizer 24 mm Baseplate, 25° half-wedge augment sizer 24 mm Baseplate, 35° half-wedge augment sizer 28 mm Baseplate, 15° half-wedge augment sizer 28 mm Baseplate, 25° half-wedge augment sizer 28 mm Baseplate, 35° half-wedge augment sizer	AR-9591-2415 AR-9591-2425 AR-9591-2435 AR-9591-2815 AR-9591-2825 AR-9591-2835
5	24 mm Baseplate, 15° half-wedge augment trial 24 mm Baseplate, 25° half-wedge augment trial 24 mm Baseplate, 35° half-wedge augment trial 28 mm Baseplate, 15° half-wedge augment trial 28 mm Baseplate, 25° half-wedge augment trial 28 mm Baseplate, 35° half-wedge augment trial	AR-9593-2415 AR-9593-2425 AR-9593-2435 AR-9593-2815 AR-9593-2825 AR-9593-2835
6	20 mm Modular Drill, half 25 mm Modular Drill, half 30 mm Modular Drill, half 35 mm Modular Drill, half 40 mm Modular Drill, half	AR-9596-20H AR-9596-25H AR-9596-30H AR-9596-35H AR-9596-40H
7	Angled Reamer Sleeve, 15°	AR-9597-15
8	Angled Reamer Sleeve, 25°	AR-9597-25
9	Retroverted Pin Guide	AR-9680
10	24/28 mm Depth Stop, 15° 24 mm Depth Stop, 25° 28 mm Depth Stop, 25° 24 mm Depth Stop, 35° 28 mm Depth Stop, 35°	AR-9598-15 AR-9598-2425 AR-9598-2825 AR-9598-2435 AR-9598-2835

Univers Revers™ Consolidated Instrument Set



- Streamlined technique and instrumentation
- Reduced back-table burden
- Reduced sterilization cost
- Maintain patient-centric implant options with the modularity of the Univers Revers system
- 2 trays for standard reverse, one glenoid and one humeral
- Additional tray for augmented glenoid


Univers Revers™ Consolidated Instrument Set (Cont.)



Consolidated Humeral Set – Top Tray (AR-9501CHS)

Label	Product Description	Item Number
1	Trial Humeral Spacers, 33 mm	AR-9554-06/09/12/15
2	Trial Humeral Spacers, 36 mm	AR-9551-06/09/12/15
3	Trial Humeral Spacers, 39 mm	AR-9552-06/09/12/15
4	Trial Humeral Liners, 33 mm	AR-9530XS-03/06/03C/06C
5	Trial Humeral Liners, 36 mm	AR-9530S-03/06/03C/06C
6	Trial Humeral Liners, 39 mm	AR-9530M-03/06/03C/06C
7	Trial Combo Liners, 33/36	AR-9603-3336-3/6/3C/6C
8	Trial Combo Liners, 36/39	AR-9603-3639-3/6/3C/6C
9	Trial Combo Liners, 39/42	AR-9603-3942-3/6/3C/6C
10	Humeral Cup Reamer, 33 mm Humeral Cup Reamer, 36 mm Humeral Cup Reamer, 39 mm	AR-9508-33 AR-9508-36 AR-9508-39
11	Pointed Stem Impactor	AR-9202-09
12	Stem/Cup Inserter/Extractor Adapter Handle	AR-9512
13	Driver Shaft, T-15, Short	AR-9545-T15-01


Univers Revers™ Consolidated Instrument Set (Cont.)



Consolidated Humeral Set – Top Tray (AR-9501CHS)

Label	Product Description	Item Number
1	Univers Revers Broach, size 5	AR-9510-05MDM
	Univers Revers Broach, size 6	AR-9510-06M
	Univers Revers Broach, size 7	AR-9510-07M
	Univers Revers Broach, size 8	AR-9510-08M
	Univers Revers Broach, size 9	AR-9510-09M
	Univers Revers Broach, size 10	AR-9510-10M
2	Revers Apex Broach, size 6	AR-9510-06SM
	Revers Apex Broach, size 7	AR-9510-07SM
	Revers Apex Broach, size 8	AR-9510-08SM
	Revers Apex Broach, size 9	AR-9510-09SM
	Revers Apex Broach, size 10	AR-9510-10SM
3	Univers Modular Reamer Pins	AR-9510RG-C/R/L
	Univers Revers Reamer Guides, size 5	AR-9510RG-5C/5R/5L
4	IM Reamer, 5 mm	AR-9202-25H
	IM Reamer, 6 mm	AR-9202-02H
5	Monoblock Broach, 4/33	AR-9510-04M
6	Reamer/Guide, 4 mm, 135°	AR-9510-04-135
7	Univers Revers Apex Impactor/Extractor	AR-9511-2
8	T-Handle – Hudson Connect	AR-9202-15H
9	Humeral Resection Protector, extra small	AR-9509-XS
	Humeral Resection Protector, small	AR-9509-S
	Humeral Resection Protector, medium	AR-9509-M
	Humeral Resection Protector, large	AR-9509-L
10	Revers Lever-Locking Broach Handle	AR-9510-2
11	Version Rod	AR-9510-01

Univers Revers™ Consolidated Instrument Set (Cont.)



Consolidated Glenoid Set – Top Tray (AR-9615CGS)


Label	Product Description	Item Number
1	Glenoid Scrapette	AR-9601
2	Pin Guide, 24 mm	AR-9616-24
3	Glove Protector Sleeve	AR-9216-4
4	Modular Reamer Shaft	AR-9617
5	Guide Pin, 2.8 mm	AR-9616
6	Guide Pin Depth Gauge	AR-9616-DG
7	Manual Driver Shaft	AR-9631
8	Glenosphere Distractor	AR-9658-R
9	Torque-Indicating Adapter	AR-9545-T15H
10	Tap Caddy Drill, 10 mm Tap, 15 mm Tap, 20 mm Tap, 25 mm Tap, 30 mm Tap, 35 mm	AR-9615C-TC AR-9620-10D AR-9621-15T AR-9621-20T AR-9621-25T AR-9621-30T AR-9621-35T

Label	Product Description	Item Number
11	Drill Caddy Drill, 15 mm Drill, 20 mm Drill, 25 mm Drill, 30 mm Drill, 35 mm	AR-9615C-DC AR-9621-15D AR-9621-20D AR-9621-25D AR-9621-30D AR-9621-35D
12	Glenosphere Sizing Guide	AR-9664-SG
13	24 mm Primary Reamer	AR-9618-24
14	Over Baseplate Reamer Caddy, consolidated set Over Baseplate Reamer, 36/24 Over Baseplate Reamer, 39/24 Over Baseplate Reamer, 42/24	AR-9615CGC-01 AR-9632-2436 AR-9632-2439 AR-9632-2442
15	Ratcheting Handle, Hybrid Hudson	AR-1999HH

Shoulder Arthroplasty Retractor Set




- More than 25 standard and specialty retractors
- Most extensive shoulder arthroplasty-specific retractor grouping on the market



Shoulder Arthroplasty Retractor Set - Top Tray (AR-9260S-53)

Label	Product Description	Item Number
1	Darrach Elevator, 14 in (L) × 1 in (W)	AR-9260-37
2	Hohmann Retractor, wide	AR-9260-34
3	Hohmann Retractor, narrow	AR-9260-33
4	Darrach Elevator, 10 in (L) × 5/8 in (W)	AR-9260-36
5	Darrach Elevator, 10 in (L) × 3/8 in (W)	AR-9260-51
6	Volkmann Bone Hook	AR-9260-38
7	Volkmann "Rake" Retractor	AR-9260-39


Shoulder Arthroplasty Retractor Set (Cont.)



Shoulder Arthroplasty Retractor Set - Middle Tray (AR-9260S-53)

Label	Product Description	Item Number
1	Thin Glenoid Retractor, narrow tip	AR-9260-30
2	Thin Glenoid Retractor, wide tip	AR-9260-31
3	Hohmann Retractor, curved	AR-9260-35
4	Cobra Retractor	AR-9260-49
5	Richardson Retractor, deep	AR-9260-42
6	Richardson Retractor, wide	AR-9260-41
7	Richardson Retractor, narrow	AR-9260-40

Shoulder Arthroplasty Retractor Set (Cont.)



Shoulder Arthroplasty Retractor Set - Bottom Tray (AR-9260S-53)

Label	Product Description	Item Number
1	Hohmann Retractor, bent	AR-9260-48
2	Fukuda Retractor, small	AR-9260-52
3	Fukuda Retractor, large	AR-9260-32
4	Gelpi Retractor, sharp tip	AR-9260-44
5	Kolbel Retractor, ring handle	AR-9260-45
6	Kolbel Retractor Blades, large	AR-9260-47
7	Kolbel Retractor Blades, small	AR-9260-46
8	Gelpi Retractor, blunt tip	AR-9260-43
9	Browne Deltoid Retractor	AR-9260-50

Virtual Implant Positioning™ (VIP™) System



Optimized for Surgeon Control and Capability

- Team of dedicated surgical planners use the patient's uploaded CT data to develop a preoperative plan that is subsequently reviewed and adjusted, as necessary, by the surgeon
- Plan and easily switch between anatomic and reverse total shoulder arthroplasty cases
- Capability to plan both glenoid and humeral components
- Metal segmentation allows for planning in complex revision cases
- Optimize implant selection through osteophyte removal and range of motion simulation
- Intraoperative flexibility with multi-plan approval
- 2-day turnaround of plan PDF and reusable transfer technology settings once plan is approved
- Web-based portal available on desktop, mobile, and iOS tablet devices

MyVIP™ Texting Platform

- Delivers timely, automated text messages to guide patients through their presurgical journey
- Educate patients and ease anxiety-ridden waiting time
- Based on patients' MyVIP experience alone⁵:
 - 91% would refer their surgeon to a friend or family member
 - 97% understand why a CT scan was needed
 - 94% were satisfied with communication and information


Virtual Implant Positioning™ (VIP™) System (Cont.)

Reusable, Adjustable, Patient-Specific Instrumentation

- Facilitates accurate intraoperative reproduction of the surgeon-approved preoperative plan
- The glenoid targeter (a) is used to transfer the planned glenoid guide pin position and trajectory from the software to the patient's anatomy
- The VIP glenoid reamer (b) can be used to replicate the desired depth of ream from the VIP plan to the patient's anatomy when used in conjunction with the modular VIP glenoid reamers
- A 3D glenoid model (a) of the patient's anatomy can be ordered



Glenoid Targeter Instrument Set



Glenoid Targeter Instrument Set (AR-5400S)

Label	Product Description	Item Number
1	Glenoid Targeter Leg, size 10	AR-5400-10
	Glenoid Targeter Leg, size 11	AR-5400-11
	Glenoid Targeter Leg, size 12	AR-5400-12
	Glenoid Targeter Leg, size 13	AR-5400-13
	Glenoid Targeter Leg, size 14	AR-5400-14
	Glenoid Targeter Leg, size 15	AR-5400-15
	Glenoid Targeter Leg, size 16	AR-5400-16
	Glenoid Targeter Leg, size 17	AR-5400-17
	Glenoid Targeter Leg, size 18	AR-5400-18
	Glenoid Targeter Leg, size 19	AR-5400-19
	Glenoid Targeter Leg, size 20	AR-5400-20
	Glenoid Targeter Leg, size 21	AR-5400-21
	Glenoid Targeter Leg, size 22	AR-5400-22
	Glenoid Targeter Leg, size 23	AR-5400-23
	Glenoid Targeter Leg, size 24	AR-5400-24
	Glenoid Targeter Leg, size 25	AR-5400-25
2	VIP™ Glenoid 3D Model Base	AR-VIPBASE
3	2.8 mm Pin, stainless steel	AR-5400-400


Glenoid Targeter Instrument Set - Case Insert (DE-VIPTARG2-IN)

Label	Product Description	Item Number
4	Glenoid Targeter Shaft	AR-5400-01
5	Glenoid Targeter Handle	AR-5400-02
6	Glenoid Targeter Capture Nut	AR-5400-03
7	Glenoid Targeter Locking Nut	AR-5400-04

Disposables

Product Description	Item Number
DynaNite® VIP Glenoid Pin, nitinol, 2.8 mm	AR-5400-400NS

VIP™ Glenoid Reaming Instrument Set



VIP Glenoid Reaming Set - Reusable Instruments (AR-5410S)

Label	Product Description	Item Number
1	Virtual Implant Positioning™ Glenoid Reamer, pilot	AR-5410-01
2	VIP Glenoid Reamer, Univers VaultLock® glenoid, small VIP Glenoid Reamer, Univers VaultLock glenoid, medium VIP Glenoid Reamer, Univers VaultLock glenoid, large VIP Glenoid Reamer, Univers VaultLock glenoid, X-large	AR-5410-VLS AR-5410-VLM AR-5410-VLL AR-5410-VLXL
3	VIP Glenoid Reamer, Univers Revers™ MGS, 24 mm VIP Glenoid Reamer, Univers Revers MGS, 28 mm	AR-5410-24 AR-5410-28
4	VIP Glenoid Reamer, paleo	AR-5410-P
5	Glove Protector Sleeve	AR-9216-4
6	Modular Reamer Shaft	AR-9617

Disposables

Product Description	Item Number
VIP Glenoid Reamer, Univers Revers Augmented MGS, small VIP Glenoid Reamer, Univers Revers Augmented MGS, medium VIP Glenoid Reamer, Univers Revers Augmented MGS, large	AR-5410-AMGS-S AR-5410-AMGS-M AR-5410-AMGS-L
VIP Augmented MGS Reamer, locking, small VIP Augmented MGS Reamer, locking, medium VIP Augmented MGS Reamer, locking, large	AR-5410T-AMGS-S AR-5410T-AMGS-M AR-5410T-AMGS-L

Shoulder Fracture

Solutions for Proximal Humerus Fracture Repair

AOS Humeral Plating System Options	70
AOS Humeral Plating Instrument Set.....	71
Humeral Nail System	74
Humeral Nail Instrument Set	75
SuturePlate™ Humeral Plate	78
SuturePlate™ Instrument Set	79
Arthrex Proximal Humeral Plating System	80
Arthrex Proximal Humeral Plating Instrument Set	81
FxBridge™ Tuberosity Repair System (AR-9517).....	86
Uniers Revers™ Fracture Adapter.....	87

AOS Humeral Plating System Options



ALPHA Plate

The ALPHA plate (a) is a unique anatomical, side-specific plate with a proximal contour that allows for sparing of the deltoid insertion and avoidance of dissection/soft-tissue stripping at the subdeltoid space.

- Designed to help avoid postoperative adhesions/stiffness and promote preservation of the vascular supply to the bone
- Preserves native deltoid insertion for potential future conversion to a reverse shoulder prosthesis
- The ALPHA plate's proprietary dual curvature design creates a near 90° opposition from proximal cluster and anterior shaft, providing for superior torsional stability in fractures with shaft extension, while avoiding the superior flare of the lateral epicondyle of the humerus.
- Convergent proximal screw pattern allows for increased spread in the humeral head while creating an intramedullary strut to resist excessive settling or subsidence of the humeral head fragment
- Shaft suture holes facilitate biceps tenodesis


130° and 95° Proximal Plates

The proximal humeral plate is available in a distal fitting, limited impingement 130° configuration (b), or a more proximal fracture buttressing 95° construct (c).

- Five plate lengths are offered in the 130° configuration
- 3 to 14-hole constructs
- 4 mm cancellous locking screws in the proximal cluster and 3.5 mm locking or nonlocking cortical screws for distal fixation in the shaft.

Note: 95° proximal humeral plate is available in 3- and 5-hole constructs


AOS Humeral Plating Instrument Set



Proximal Humeral Plate Main Tray – Top Tray (9912-100I)

Label	Product Description	Item Number
Instruments		
1	Cancellous Tap, 3.5 mm, AO style	0255-200
2	Proximal Screw Depth Gauge	0522-100
3	Calibrated Step Drill, AO style, 3.5/2.5 mm × 130 mm	0254-000
4	Calibrated Drill, long, 2.5 mm × 130 mm	0250-000
5	Alignment Guide, proximal humeral plate Bolt, T15, alignment guide, proximal humeral plate	1250-200 1251-200
6	Cortical Tap, 3.5 mm, AO style	0226-100
7	Hook Tip Depth Gauge, locking/nonlocking	0523-100
8	Drill Guide, nonlocking/locking, 2.5 mm	0318-100
9	Calibrated Drill, AO style, 2.5 mm × 95 mm	0237-100
10	Screw Sheath Drill Guide, 2.5 mm Drill Guide, 3.5 mm	0316-100 0327-000 0331-000
11	Drill Guide, 2.5 mm	0324-200


AOS Humeral Plating Instrument Set (Cont.)



Proximal Humeral Plate Main Tray – Bottom Tray (9912-100I)

Label	Product Description	Item Number
1	Trocar Tip Guidewire, 2 mm × 150 mm	0102-150
2	T15 Driver, AO style T15 Driver, large Hudson	0432-400 0452-200
3	Hex Driver, 5 mm	5001-000
4	Screw Gripper, 3.5 mm/4 mm screws	0433-000
5	Ratcheting Small Axial Handle, AO quick connect, cannulated	5026-100
6	Solid T15 Driver	0448-300

AOS Humeral Plating Instrument Set (Cont.)



Proximal Humeral Plate Expansion Kit (9941-000I)

Label	Product Description	Item Number
Instruments		
1	Alignment Guide, 95° PHP	1288-000
2	Guide Pin, 1.5 mm	0109-150
3	Alignment Guide Locking Bolt, proximal humeral plate	1251-200
4	Alignment Guide, ALPHA plate, left	1296-000
5	Alignment Guide, ALPHA plate, right	1297-000

Implants

Product Description	Item Number
95° Plates	
95° Proximal Humeral Plate, 3H	3022-003
95° Proximal Humeral Plate, 5H	3022-005
ALPHA Plates	
ALPHA Plate, left, 4H	3024-004
ALPHA Plate, left, 5H	3024-005
ALPHA Plate, left, 7H	3024-007
ALPHA Plate, left, 11H	3024-011
ALPHA Plate, left, 14H	3024-014
ALPHA Plate, left, 16H	3024-016
ALPHA Plate, left, 18H	3024-018
ALPHA Plate, right, 4H	3025-004
ALPHA Plate, right, 5H	3025-005
ALPHA Plate, right, 7H	3025-007
ALPHA Plate, right, 11H	3025-011
ALPHA Plate, right, 14H	3025-014
ALPHA Plate, right, 16H	3025-016
ALPHA Plate, right, 18H	3025-018

Product Description	Item Number
135° Plates	
Proximal Humeral Plate, 3H	3012-003
Proximal Humeral Plate, 5H	3012-005
Proximal Humeral Plate, 8H	3012-008
Proximal Humeral Plate, 11H	3012-011
Proximal Humeral Plate, 14H	3012-014
Screws	
Cortical Screw, 3.5 mm × 20 mm-32 mm (2 mm increments), 35 mm, 40 mm, 45 mm	8010-200-450
Cortical Locking Screw, 3.5 mm × 20 mm-36 mm (2 mm increments)	8014-200-360
Cortical Bone Screw, 3.5 mm × 20 mm-32 mm (2 mm increments), 35 mm, 40 mm, 45 mm	8110-020-045
Cortical Screw, double lead lock, 3.5 mm × 20 mm-36 mm (2 mm increments)	8114-020-036
Cancellous Locking Screw, fully threaded, 4 mm × 22 mm-60 mm (2 mm increments), 65 mm, 70 mm	8124-022-070
Cancellous Locking Screw, partially threaded, 4 mm × 22 mm-60 mm (2 mm increments), 65 mm, 70 mm	8128-022-070


Humeral Nail System



The Arthrex Humeral Nail System consists of both short and long intramedullary nails. The 15 cm short nail is best suited for simple 2-part fractures of the proximal humerus involving the surgical neck of the humerus and for proximal third humeral fractures without comminution. The long humeral nail is primarily intended for humeral shaft fractures, which are inherently not prone to varus malposition.

- 6° proximal bend facilitates nail insertion with an incision just medial to the rotator cuff insertion
- Threaded proximal holes in every nail eliminate screw back out
- Self-reaming tip allows for less bone disruption at the entry site and decreases OR time (short nail)
- 5 mm of compression is accomplished intraoperatively through the jig, allowing for exact and reproducible compression (long nail)


Humeral Nail Instrument Set



Humeral Nail Instrument Set – Top Tray (9903-200S)

Label	Product Description	Item Number
1	Entry Reamer, cannulated, 10.5 mm	0221-100
2	Guidewire Depth Gauge	0515-000
3	Flexible Reamer Shaft	0222-200
4	Reamer Head, 6 mm Reamer Head, 6.5 mm Reamer Head, 7 mm Reamer Head, 7.5 mm Reamer Head, 8 mm Reamer Head, 8.5 mm Reamer Head, 9 mm Reamer Head, 9.5 mm Reamer Head, 10 mm	0223-060 0223-065 0223-070 0223-075 0223-080 0223-085 0223-090 0223-095 0223-100
5	Guidewire Gripper, 2.0 mm/3.0 mm	0419-200
6	Guidewire, 3.2 mm × 330 mm Trocar Tip Guidewire, 2.0 mm × 508 mm	AR-9095-40 0102-500
7	Reduction Tool	0807-000
8	Straight Awl, cannulated	0218-200
9	Soft Tissue Protector	0649-000


Humeral Nail Instrument Set (Cont.)



Humeral Nail Instrument Set – Middle Tray (9903-200S)

Label	Product Description	Item Number
1	Humeral Nail Drill Guide, 2.9 mm	0604-000
2	Wrench, 3/4-1/2 Combination	0813-000
3	Cortical Tap, 3.5 mm	0226-200
4	Drill, AO Style, short, 2.9 mm Drill, calibrated, 2.5 mm Drill, calibrated, AO style, 3.5 mm	0224-100 0217-100 0220-100
5	Hook Tip Depth Gauge	0513-000
6	Compression Nut Driver, 3.5 mm	0458-100
7	Hex Driver, 2.5 mm Hex Driver, 3.5 mm	AR-9903-01 AR-9903-02
8	Hook Tip Depth Gauge Sheath	0602-000
9	Distal Drill Guide, 2.5 mm Proximal Drill Guide, 3.5 mm	0307-000 0309-000
10	Screw Sheath	0601-100
11	Screw Trocar	0607-000
12	T-Handle, Quick Connect	AR-9505-26
13	Screw Gripper	0413-000
14	Impactor Pad	0808-200
15	Humeral Nail Extractor	0803-000
16	Hex Driver, short, 9/32	0407-100
17	Jig Locking Bolt	1231-100
18	Targeting Module	1228-200
19	Insertion Guide	1230-100

Humeral Nail Instrument Set (Cont.)



Implants

Product Description	Item Number
Nails	
Proximal Humeral Nail, 8 mm × 15 cm	1170-150
Humeral Nail, 7 mm × 20 cm	1171-200
Humeral Nail, 7 mm × 22.5 cm	1171-225
Humeral Nail, 7 mm × 25 cm	1171-250
Humeral Nail, 7 mm × 27.5 cm	1171-275
Humeral Nail, 7 mm × 30 cm	1171-300
Humeral Nail, 8 mm × 20 cm	1172-200
Humeral Nail, 8 mm × 22.5 cm	1172-225
Humeral Nail, 8 mm × 25 cm	1172-250
Humeral Nail, 8 mm × 27.5 cm	1172-275
Humeral Nail, 8 mm × 30 cm	1172-300
Humeral Nail, 9 mm × 20 cm	1173-200
Humeral Nail, 9 mm × 22.5 cm	1173-225
Humeral Nail, 9 mm × 25 cm	1173-250
Humeral Nail, 9 mm × 27.5 cm	1173-275
Humeral Nail, 9 mm × 30 cm	1173-300
Screws	
Cancellous Screw, 5 mm × 25 mm-60 mm (5 mm increments)	8022-025-060
Cortical Screw, 3.5 mm × 17.5 mm-40 mm (2.5 mm increments)	8010-175-400
End Cap, 2 mm	1074-020
End Cap, 5 mm	1074-050
Fracture Reduction Screw	1175-000

Disposables

Product Description	Item Number
Guide Pin, 3.2 mm	0100
Guidewire w/ Trocar Tip, 2 mm	0102-500
Guidewire, ball nose, 2 mm	0103-800
Drill Bit, calibrated, green, 2.5 mm	0217-100
Drill Bit, calibrated, gold, 3.5 mm	0220-100
Drill Bit, short, 2.9 mm	0224-100
Tap, cortical, 3.5 mm	0226
Drill Bit, calibrated, blue, 3 mm	0230

SuturePlate™ Humeral Plate



The SuturePlate humeral plate is an anatomically designed, low-profile, titanium polyaxial locking plate and screw system.

- Multiple chamfered suture eyelets along the margin of the plate allow the surgeon to pass FiberWire® suture after the plate is fixed to bone
- Polyaxial bushings offer surgeon the freedom to direct screws within the anatomical template, based on the fracture pattern and bone quality for better fixation
- Compatible with variable angle locking and nonlocking screws


Implants

Product Description	Item Number
Humeral SuturePlate, 3H	AR-14003
Humeral SuturePlate, 5H	AR-14005
Cortical Locking Screws, 3.5 mm × 16 mm-40 mm (2 mm increments)	AR-14116-40
Cancellous Locking Screws, 4 mm × 24 mm-56 mm (2mm increments)	AR-14224-56
Cortical Screws, 3.5 mm × 16 mm-40 mm (2 mm increments), 45 mm, 50 mm	AR-14116NL-40NL

Disposables

Product Description	Item Number
Compatible FiberWire Suture w/ Needles	
FiberWire Suture Kit for SuturePlate	AR-7214
#5 FiberWire Suture, 38 in w/ cutting and tapered needles	AR-7213
#2 FiberWire Suture (blue), 38 in w/ tapered needle	AR-7200
#2 FiberWire Suture (blue), 38 in w/ two tapered needles	AR-7205
#2 FiberWire Suture (white/black), 38 in w/ two tapered needles	AR-7205T
#2 FiberWire Suture and TigerWire® Suture (one blue, one white/black), 38 in w/ tapered needle	AR-7208


SuturePlate™ Instrument Set



SuturePlate Humeral Plate Instrument Set (AR-14003S)


Label	Product Description	Item Number
1	Terminally Threaded Pin, 2.8 mm	AR-2521
2	Ball Wire, 1.6 mm × 15 cm K-Wire, 1.6 mm × 15 cm	AR-14016B AR-14016
3	Hex Screw Driver Shaft	AR-14025
4	Easy-Out	AR-1994-25
5	Temporary Compression Device	AR-14023TCD
6	Drill Bit, 2.5 mm	AR-8943-13
7	Proximal Screw Depth Gauge	AR-14013
8	Pin Extractor	AR-14016PE
9	Drill Guide	AR-14029
10	Threaded Depth Gauge, qty. 2	AR-14012
11	SuturePlate Handle	AR-14024
12	Tear Drop Handle w/ AO Connection	AR-2001AO
13	Torque Driver	AR-14025AO

Arthrex Proximal Humeral Plating System



- Limited-contact option for fixation of proximal humerus fractures
- Divergent proximal screw cluster for optimal articular reconstruction
- Proximal screw holes can accept 4.0 mm cancellous or 3.5 mm reinforced VAL/KreuLock™ locking screws with the option to use 3.5 mm cortical or locking screws distally
- FiberTape® cerclage buttons may be used in the shaft to facilitate fracture reduction
- Low-profile plate and screw interaction
- Proximal holes for suture engagement with the plate, with cutouts to prevent plate lifting

Arthrex Proximal Humeral Plating Instrument Set



Titanium VA Proximal Humeral Plating System, Top Tray

Label	Product Description	Item Number
Proximal humeral plate, screw caddy, KreuLock™ reinforced		
1	Proximal Humeral Plate, screw caddy, KreuLock	AR-9953C-SC-VCL
	Proximal Humeral Plate, screw caddy, KreuLock lid	AR-9953C-SC-VCL01
	VAL KreuLock Comp. Screw, reinforced, Ti, 3.5 mm × 10-70 mm	AR-8935RVCL-XX
	Cortical Screw, Ti, 3.5 mm × 10-70 mm	AR-8935-XX
	Cancellous Screw, Ti, 4 mm × 10-60 mm	AR-8940-XX
	Low Profile Cortical Screw, Ti, 4.5 mm × 18-40 mm	AR-8545-XX
	Low Profile Locking Screw, Ti, 4.5 mm × 18-40 mm	AR-8545L-XX
	Washer, Ti, 7 mm	AR-8740W
	Cerclage Button, 3.5 mm	AR-9953CB-T35*
	Large BB-Tak	AR-8970-09
	Large BB-Tak, threaded	AR-8970-09T
	Screw Forceps	AR-8941F


Label	Product Description	Item Number
Proximal Humeral Plate Instruments		
2, 3, 4	Bending Iron, proximal humeral plate	AR-9953-08
	Drill Guide, locking, 4.5 mm	AR-8970-01
	Aiming Block Assembly, proximal humerus	AR-9953-01
	Aiming Block Captured Bolt, T15 (not pictured)	AR-9953-01B
	Reduction Nut, proximal humeral plate	AR-9953-05
	Drill Bit, calibrated, reusable 2.8 mm	AR-9953DB-28RU**,**
	Drill Bit, calibrated, reusable 2.5 mm	AR-9953DB-25RU**,**
	Drill Bit, calibrated, 3 mm	AR-8970-30
	4.5 mm Drill Bit	AR-8970-45
	4.5 mm Profile Drill	AR-18900-11
	T20 Hexalobe Driver (AO)	AR-8970-03
	Drill Guide, 3 mm/4.5 mm	AR-8970-02
	Depth Measuring Device, 4.5 mm/5.5 mm	AR-8970-07
	Screwdriver T20 Hexalobe, straight	AR-8970-04
	Depth Probe, 3.5 mm/4.5 mm	AR-9953-06

Note: 4.5 mm instruments and implants intentionally left out of set

*Available sterile, with "S" designation, eg, AR-9953DB-28S.

** Available nonsterile, with no designation, eg, AR-9953DB-28.

Arthrex Proximal Humeral Plating Instrument Set (Cont.)




Titanium VA Proximal Humeral Plating System, Middle Tray

Label	Product Description	Item Number
Proximal Humeral Plate, screw caddy, VAL reinforced		
1	Proximal Humeral Plate, screw caddy, VAL Proximal Humeral Plate, screw caddy, VAL lid VAL Screw, reinforced, Ti, 3.5 × 10-70 mm Cortical Screw, Ti, 3.5 × 10-70 mm Cancellous Screw, Ti, 4.0 × 10-60 mm Low Profile Cortical Screw, Ti, 4.5 × 18-40 mm Low Profile Locking Screw, Ti, 4.5 × 18-40 mm Washer, Ti, 7 mm Cerclage Button, 3.5 mm Large BB-Tak Large BB-Tak, threaded Screw Forceps	AR-9953C-SC-VAL AR-9953C-SC-VAL01 AR-8935RV-XX AR-8935-XX AR-8940-XX AR-8545-XX AR-8545L-XX AR-8740W AR-9953CB-T35* AR-8970-09 AR-8970-09T AR-8941F
Proximal Humeral Plate, lateral tray		
2, 3	Lateral Proximal Humeral Plate, 3h, 91 Mm Lateral Proximal Humeral Plate, 5h, 114 Mm Lateral Proximal Humeral Plate, 7h, 153 Mm Lateral Proximal Humeral Plate, 10h, 200 Mm Lateral Proximal Humeral Plate, 12h, 231 Mm Lateral Proximal Humeral Plate, 14h, 263 Mm Lateral Proximal Humeral Plate, 17h, 309 Mm Lateral Proximal Humeral Template, 5h Lateral Proximal Humeral Template, 17h	AR-9953UH-03* AR-9953UH-05* AR-9953UH-07* AR-9953UH-10* AR-9953UH-12* AR-9953UH-14* AR-9953UH-17* AR-9953UH-05T AR-9953UH-17T

Note: 4.5 mm implants intentionally left out of set

*Available sterile, with "S" designation, eg, AR-9953UH-03S.


Arthrex Proximal Humeral Plating Instrument Set (Cont.)



Titanium VA Proximal Humeral Plating System, Bottom Tray

Label	Product Description	Item Number
1	Proximal Humeral Plate, anatomic case	AR-9953C-02

Arthrex Proximal Humeral Plating Instrument Set (Cont.)




General Instruments 2.7 mm/3.5 mm, Top Tray

Label	Product Description	Item Number
2.7 mm Instruments and Disposables		
1	Torque Limiting Adapter, AO, 0.8 Nm	AR-18700-39
	Locking Drill Guide, 2.7 mm	AR-9954-06
	Holding Sleeve, 2.7 mm	AR-18800-32
	Drill Guide, 2 mm/2.7 mm	AR-18800-22
	VAL Drill Guide, double-ended, 2.7 mm	AR-9954-07
	Driver Shaft, T10	AR-18800-24
	Driver Shaft, T10, long	AR-18800-26
	Solid Screwdriver, T10	AR-18800-27
	Tap, 2.7 mm	AR-18800-34
	Countersink, 2.7 mm	AR-18800-23
	Drill Bit, 2 mm, short, AO	AR-18800-18
Drill Bit, 2.7 mm	AR-18800-20	
Drill Bit, 2 mm, calibrated, AO	AR-9954-20	
General Instruments and Disposables		
2	Depth Device, low profile, 2.7 mm/3.5 mm/4 mm	AR-8943-15
	Ratcheting Handle, AO QC	AR-8700RH
	Sharp Hook	AR-8943-21
	Freer Elevator	AR-8943-19
	Hohmann Retractor, double ended (6 mm/10 mm)	AR-18700-53
	Periosteal Elevator, curved blade, 6 mm	AR-8943-20
	BB-Tak	AR-13226
	BB-Tak, threaded	AR-13226T
	Guidewire w/ Trocar Tip, 0.053 in (1.35 mm)	AR-8943-01
	Guidewire w/ Trocar Tip, 0.062 in (1.6 mm)	AR-8941K
	Guidewire w/ Trocar Tip, 0.078 in (2 mm)	AR-8945K


Label	Product Description	Item Number
3.5 mm Instruments and Disposables		
3	Solid Screwdriver, T15	AR-8943-10
	Driver Shaft, T15, long	AR-9953-T15
	T15 Hexalobe Driver	AR-8941DH
	VAL Drill Guide, double ended, 3.5 mm	AR-8940VDG
	Lag Drill Guide, 3.5 mm/2.5 mm	AR-8963-21
	Drill Guide, 2.5 mm/3.5 mm	AR-8943-14
	Torque Limiting Attachment, 1.5 Nm, QC	AR-8963TL-01
	Countersink, 3.5 mm/4 mm	AR-8950-03
	Tap, 3.5 mm	AR-8935T
	Drill Bit, 2.5, short	AR-8943-30
	Drill Bit, 3.5, short	AR-4160-35
	Drill Bit, 4 mm	AR-4160-40
	Drill Bit, 2.5, calibrated, AO	AR-9954-25
	Drill Bit, 2.8, calibrated, AO	AR-9954-28
	Drill Bit, 3.5, long	AR-8964-10

Arthrex Proximal Humeral Plating Instrument Set (Cont.)



General Instruments 2.7 mm/3.5 mm, Bottom Tray

Label	Product Description	Item Number
3.5 mm Instruments and Disposables		
1	Drill Bit, calibrated, AO, long, 2.8 mm Drill Bit, calibrated, AO, long, 2.5 mm Drill Tip K-Wire, calibrated, 2.8 mm Drill Tip K-Wire, calibrated, 2.5 mm Drill Tip K-Wire, calibrated, 2 mm Drill Tip K-Wire, 1.6 mm MIS Tissue Sleeve Depth Gauge Profile Drill, reinforced, AO, long, 3.5 mm BB-Tak, threaded, long MIS Tissue Sleeve Trocar, long MIS Drill Guide, threaded, long, 3.5 mm MIS Drill Guide, non-threaded, long, 3.5 mm MIS Tissue Sleeve	AR-9954-28L AR-9954-25L AR-9954-28GW AR-9954-25GW AR-9954-20GW AR-18800-36 AR-9954-05 AR-18900-10 AR-14016B AR-9954-04 AR-9954-02 AR-9954-03 AR-9954-01
General Instruments		
2	Locking Bending Tower, 3.5 mm Bending Iron Bending Iron, double-ended	AR-9954-09 AR-8943-18 AR-18800-44
Hand Instruments		
3	Needlenose Pliers Weber Clamp, small Weber Clamp Submuscular Tissue Elevator, distal tibia Lobster Claw	AR-8916-24 AR-5050-05 AR-8943-24 AR-8963-10 AR-8943-23



- All-inclusive kit streamlines repairs during a case
- Three distinct suture colors aid in suture management for a simple, reproducible repair
- Provides an environment for tuberosity healing, which improves postoperative range of motion (ROM)⁶
- FiberTape® suture is 27% more resistant to tissue pull-through compared to round #2 suture⁷
- 135° humeral inclination increases tuberosity healing rate and stability⁸

Univers Revers™ Fracture Adapter

- Compatible with all sizes of Univers Revers humeral stems
- Suture holes for tuberosity fixation
- Machine-textured, grit-blasted, and hydroxyapatite (HA)-coated for bony ongrowth
- Narrow and wide sizes available:
 - Size 5: +11 mm
 - Sizes 6-9: +9.5 mm
 - Sizes 10-15: +7.5 mm
- Heights available as standard and +5 mm

Implants

Product Description	Item Number
Fracture Adapter, sizes 6-9	AR-9501-FX6N
Fracture Adapter, sizes 6-9, wide	AR-9501-FX6W
Fracture Adapter, sizes 6-9, tall	AR-9501-FX6NT
Fracture Adapter, sizes 6-9, wide, tall	AR-9501-FX6WT
Fracture Adapter, sizes 10-15	AR-9501-FX10N
Fracture Adapter, sizes 10-15, wide	AR-9501-FX10W
Fracture Adapter, sizes 10-15, tall	AR-9501-FX10NT
Fracture Adapter, sizes 10-15, wide, tall	AR-9501-FX10WT
Fracture Adapter, size 5	AR-9501-FX5N
Fracture Adapter, size 5, wide	AR-9501-FX5W
Fracture Adapter Trunnion	AR-9501FTK

Fracture Adapter Case (AR-9510FX)

Product Description	Item Number
Trial Fracture Adapter, size 5, narrow	AR-9510-FX5N
Trial Fracture Adapter, size 5, wide	AR-9510-FX5W
Trial Fracture Adapter, sizes 6-9, narrow	AR-9510-FX6N
Trial Fracture Adapter, sizes 6-9, wide	AR-9510-FX6W
Trial Fracture Adapter, sizes 6-9, narrow, tall	AR-9510-FX6NT
Trial Fracture Adapter, sizes 6-9, wide, tall	AR-9510-FX6WT
Trial Fracture Adapter, sizes 10-15, narrow	AR-9510-FX10N
Trial Fracture Adapter, sizes 10-15, wide	AR-9510-FX10W
Trial Fracture Adapter, sizes 10-15, narrow, tall	AR-9510-FX10NT
Trial Fracture Adapter, sizes 10-15, wide, tall	AR-9510-FX10NW
Fracture Adapter Impactor	AR-9510-F1
Fracture Adapter Extractor	AR-9510-F2
Fracture Adapter Extractor Shaft	AR-9510-F3

Disposables

Product Description	Item Number
FiberWire® Suture Kit	AR-7219
#2 FiberWire Suture (blue), 97 cm, w/ tapered needle, 26.5 mm, ½ circle, qty. 12	AR-7200
#2 FiberWire Suture, 97 cm, w/ reverse cutting needle, 36.6 mm, ½ circle, qty. 12	AR-7202
#5 FiberWire Suture (blue), 97 cm	AR-7210
#5 FiberWire Suture, 97 cm, w/ conventional cutting needle, 48 mm, ½ circle, qty. 12	AR-7211



Knee Arthroplasty

Solutions for Partial Knee Replacement

iBalance® UKA System.....	90
UKA Instrumentation Set 1	93
UKA Instrumentation Set 2	97
iBalance® PFJ System.....	99
iBalance® PFJ Instrument Set.....	100
Patellar Resection System.....	102



Unicondylar Knee Arthroplasty

- Only UKA implant that allows independent flexion and extension gap balancing
- Patented keel design resists anterior shear forces and anchors the tibial implants
- Tibial cortical support is preserved through all aspects of instrument fixation
- Independent control of flexion and extension space provide uncompromising joint line control, ligament balancing, bone preservation, and kinematic restoration
- Femoral thickness of 7 mm distal, posterior, and at midflexion
- Tibial thickness options range from 8 mm to 12 mm in 1 mm increments; 14 mm option also available
- Overall composite thickness between 15 mm and 21 mm

iBalance® UKA System (Cont.)

Implants

Product Description	Item Number
Femoral Components	
Femoral Component, size 1, LM	AR-501-UFLA
Femoral Component, size 2, LM	AR-501-UFLB
Femoral Component, size 3, LM	AR-501-UFLC
Femoral Component, size 4, LM	AR-501-UFLD
Femoral Component, size 5, LM	AR-501-UFLE
Femoral Component, size 6, LM	AR-501-UFLF
Femoral Component, size 1, RM	AR-501-UFRA
Femoral Component, size 2, RM	AR-501-UFRB
Femoral Component, size 3, RM	AR-501-UFRC
Femoral Component, size 4, RM	AR-501-UFRD
Femoral Component, size 5, RM	AR-501-UFRE
Femoral Component, size 6, RM	AR-501-UFRF
Tibial Components	
Tibial Tray Component, size 1, LM	AR-511-T1L
Tibial Tray Component, size 2, LM	AR-511-T2L
Tibial Tray Component, size 3, LM	AR-511-T3L
Tibial Tray Component, size 4, LM	AR-511-T4L
Tibial Tray Component, size 5, LM	AR-511-T5L
Tibial Tray Component, size 6, LM	AR-511-T6L
Tibial Tray Component, size 1, RM	AR-511-T1R
Tibial Tray Component, size 2, RM	AR-511-T2R
Tibial Tray Component, size 3, RM	AR-511-T3R
Tibial Tray Component, size 4, RM	AR-511-T4R
Tibial Tray Component, size 5, RM	AR-511-T5R
Tibial Tray Component, size 6, RM	AR-511-T6R

Product Description	Item Number
Polyethylene Components, vit E	
Tibial Bearing, size 1, 8 mm, vit E	AR-521-TBA8
Tibial Bearing, size 1, 9 mm, vit E	AR-521-TBA9
Tibial Bearing, size 1, 10 mm, vit E	AR-521-TBA0
Tibial Bearing, size 1, 11 mm, vit E	AR-521-TBA1
Tibial Bearing, size 1, 12 mm, vit E	AR-521-TBA2
Tibial Bearing, size 2, 8 mm, vit E	AR-521-TBB8
Tibial Bearing, size 2, 9 mm, vit E	AR-521-TBB9
Tibial Bearing, size 2, 10 mm, vit E	AR-521-TBB0
Tibial Bearing, size 2, 11 mm, vit E	AR-521-TBB1
Tibial Bearing, size 2, 12 mm, vit E	AR-521-TBB2
Tibial Bearing, size 3, 8 mm, vit E	AR-521-TBC8
Tibial Bearing, size 3, 9 mm, vit E	AR-521-TBC9
Tibial Bearing, size 3, 10 mm, vit E	AR-521-TBC0
Tibial Bearing, size 3, 11 mm, vit E	AR-521-TBC1
Tibial Bearing, size 3, 12 mm, vit E	AR-521-TBC2
Tibial Bearing, size 3, 14 mm, vit E	AR-521-Tb(A-F)4
Tibial Bearing, size 4, 8 mm, vit E	AR-521-TBD8
Tibial Bearing, size 4, 9 mm, vit E	AR-521-TBD9
Tibial Bearing, size 4, 10 mm, vit E	AR-521-TBD0
Tibial Bearing, size 4, 11 mm, vit E	AR-521-TBD1
Tibial Bearing, size 4, 12 mm, vit E	AR-521-TBD2
Tibial Bearing, size 5, 8 mm, vit E	AR-521-TBE8
Tibial Bearing, size 5, 9 mm, vit E	AR-521-TBE9
Tibial Bearing, size 5, 10 mm, vit E	AR-521-TBE0
Tibial Bearing, size 5, 11 mm, vit E	AR-521-TBE1
Tibial Bearing, size 5, 12 mm, vit E	AR-521-TBE2
Tibial Bearing, size 6, 8 mm, vit E	AR-521-TBF8
Tibial Bearing, size 6, 9 mm, vit E	AR-521-TBF9
Tibial Bearing, size 6, 10 mm, vit E	AR-521-TBF0
Tibial Bearing, size 6, 11 mm, vit E	AR-521-TBF1
Tibial Bearing, size 6, 12 mm, vit E	AR-521-TBF2

Optional Instruments


Product Description	Item Number
Disposable Tibial Cutting Guide Left	AR-621-TL2
Disposable Tibial Cutting Guide Right	AR-621-TR2
iBalance UKA Tibial Fixation Buttress	AR-611-5
iBalance UKA Tibial Trialing Shim, 2/3 mm	AR-611-7
iBalance UKA Tibia Trialing Shim, 4/5 mm	AR-611-13
iBalance UKA Femoral Peg Drill, ø 0.312 mm	AR-601-FPDO
iBalance UKA Tibial Peg Drill, ø 0.312 mm	AR-601-TPDO
iBalance UKA Tibial Trial Handle	AR-611-18
EM Tibial Guide, long distal body	AR-623-32
Tibial Recut Guide	AR-601-TCCG-2

iBalance® UKA System (Cont.)

Optional Spacer Block Organization (AR-611-SC1 and AR-611-SC2)

Product Description	Item Number	AR-611-SC1	AR-611-SC2	AR-611-SCA
iBalance UKA Space Block Caddy, 11 spaces	AR-611-C3	✓		
iBalance UKA Space Block Caddy, 12 spaces	AR-611-C4		✓	
iBalance UKA Spacer Block, 6 mm	AR-611-SB06	✓		
iBalance UKA Spacer Block, 7 mm	AR-611-SB07	✓		
iBalance UKA Spacer Block, 8 mm	AR-611-SB08	✓	✓	
iBalance UKA Spacer Block, 9 mm	AR-611-SB09	✓	✓	
iBalance UKA Spacer Block, 10 mm	AR-611-SB10	✓	✓	
iBalance UKA Spacer Block, 11 mm	AR-611-SB11		✓	
iBalance UKA Spacer Block, 12 mm	AR-611-SB12		✓	
iBalance UKA Spacer Block, 14 mm	AR-611-SB14		✓	
iBalance UKA Spacer Block, 15 mm	AR-611-SB15	✓	✓	✓
iBalance UKA Spacer Block, 16 mm	AR-611-SB16	✓	✓	✓
iBalance UKA Spacer Block, 17 mm	AR-611-SB17	✓	✓	✓
iBalance UKA Spacer Block, 18 mm	AR-611-SB18	✓	✓	✓
iBalance UKA Spacer Block, 19 mm	AR-611-SB19		✓	
iBalance UKA Spacer Block, 21 mm	AR-611-SB21		✓	
iBalance UKA Spacer Block Caddy, adjustable	AR-611-C5			✓
Adjustable Spacer Block	AR-611-20			✓
Lamina Spreader	AR-1340T			✓

UKA Instrumentation Set 1




UKA Instrumentation Set 1 (Cont.)

UKA Instrumentation Set 1, Top Tray (AR-611-C1)

Label	Product Description	Item Number
1	iBalance® UKA Spacer Block, 6 mm	AR-611-SB06
	iBalance UKA Spacer Block, 7 mm	AR-611-SB07
	iBalance UKA Spacer Block, 8 mm	AR-611-SB08
	iBalance UKA Spacer Block, 9 mm	AR-611-SB09
	iBalance UKA Spacer Block, 10 mm	AR-611-SB10
	iBalance UKA Spacer Block, 15 mm	AR-611-SB15
	iBalance UKA Spacer Block, 16 mm	AR-611-SB16
	iBalance UKA Spacer Block, 17 mm	AR-611-SB17
	iBalance UKA Spacer Block, 18 mm	AR-611-SB18
2	iBalance UKA Modular Spacer Block, 5 mm	AR-611-SM05
	iBalance UKA Modular Spacer Block, 10 mm	AR-611-SM10
3	iBalance Patella Calipers	AR-602-48
4	iBalance UKA Quick-Connect Handle, short	AR-611-8
5	Short Distal Body, EM tibial guide	AR-623-31
6	UKA Tibial Stylus, adjustable	AR-611-17
7	Ankle Clamp EM Tibial Guide	AR-623-33
8	Proximal Body, EM tibial guide	AR-623-30
9	Tibia Measuring Guide	AR-611-14
10	Angel Wing, narrow body	AR-623-76
11	iBalance UKA Tibia Trialing Shim, $\frac{2}{3}$ mm	AR-611-7
12	iBalance UKA Shim, $\frac{1}{5}$ mm	AR-611-13
13	iBalance UKA Cement Removal Tool	AR-611-3
14	iBalance UKA Tibial Lug Hole Stabilizer	AR-611-12
15	iBalance UKA Offset Attachment Alignment, Drop Rod	AR-611-1
16	Threaded Pin, headed, short	AR-611-15
17	iBalance UKA Trocar Pin Driver, smooth	AR-613-50
18	iBalance UKA Headless Pin Driver	AR-613-91
19	iBalance UKA Femoral Peg Drill, Ø0.260	AR-611-11
20	iBalance UKA Tibia Peg Drill	AR-611-10
21	iBalance UKA Quick-Connect Handle, short	AR-611-8

UKA Instrumentation Set 1 (Cont.)



UKA Instrumentation Set 1, Bottom Tray (AR-611-S1)


Label	Product Description	Item Number
1	iBalance® UKA, Trial Bearing Puller	AR-601-TBP0
2	iBalance UKA Tibial Bearing Trial, size 1, 10 mm iBalance UKA Tibial Bearing Trial, size 1, 11 mm iBalance UKA Tibial Bearing Trial, size 1, 12 mm iBalance UKA Tibial Bearing Trial, size 1, 14 mm iBalance UKA Tibial Bearing Trial, size 1, 8 mm iBalance UKA Tibial Bearing Trial, size 1, 9 mm	AR-601-TBA0 AR-601-TBA1 AR-601-TBA2 AR-601-TBA4 AR-601-TBA8 AR-601-TBA9
3	iBalance UKA Tibial Bearing Trial, size 2, 10 mm iBalance UKA Tibial Bearing Trial, size 2, 11 mm iBalance UKA Tibial Bearing Trial, size 2, 12 mm iBalance UKA Tibial Bearing Trial, size 2, 14 mm iBalance UKA Tibial Bearing Trial, size 2, 8 mm iBalance UKA Tibial Bearing Trial, size 2, 9 mm	AR-601-TBB0 AR-601-TBB1 AR-601-TBB2 AR-601-TBB4 AR-601-TBB8 AR-601-TBB9
4	iBalance UKA Tibial Bearing Trial, size 3, 10 mm iBalance UKA Tibial Bearing Trial, size 3, 11 mm iBalance UKA Tibial Bearing Trial, size 3, 12 mm iBalance UKA Tibial Bearing Trial, size 3, 14 mm iBalance UKA Tibial Bearing Trial, size 3, 8 mm iBalance UKA Tibial Bearing Trial, size 3, 9 mm	AR-601-TBC0 AR-601-TBC1 AR-601-TBC2 AR-601-TBC4 AR-601-TBC8 AR-601-TBC9
5	iBalance UKA Tibial Bearing Trial, size 4, 10 mm iBalance UKA Tibial Bearing Trial, size 4, 11 mm iBalance UKA Tibial Bearing Trial, size 4, 12 mm iBalance UKA Tibial Bearing Trial, size 4, 14 mm iBalance UKA Tibial Bearing Trial, size 4, 8 mm iBalance UKA Tibial Bearing Trial, size 4, 9 mm	AR-601-TBD0 AR-601-TBD1 AR-601-TBD2 AR-601-TBD4 AR-601-TBD8 AR-601-TBD9
6	iBalance UKA Tibial Bearing Trial, size 5, 10 mm iBalance UKA Tibial Bearing Trial, size 5, 11 mm iBalance UKA Tibial Bearing Trial, size 5, 12 mm iBalance UKA Tibial Bearing Trial, size 5, 14 mm iBalance UKA Tibial Bearing Trial, size 5, 8 mm iBalance UKA Tibial Bearing Trial, size 5, 9 mm	AR-601-TBE0 AR-601-TBE1 AR-601-TBE2 AR-601-TBE4 AR-601-TBE8 AR-601-TBE9

UKA Instrumentation Set 1 (Cont.)

UKA Instrumentation Set 1, Bottom Tray (AR-611-S1)

Label	Product Description	Item Number
7	iBalance® UKA Tibial Bearing Trial, size 6, 10 mm	AR-601-TBF0
	iBalance UKA Tibial Bearing Trial, size 6, 11 mm	AR-601-TBF1
	iBalance UKA Tibial Bearing Trial, size 6, 12 mm	AR-601-TBF2
	iBalance UKA Tibial Bearing Trial, size 6, 14 mm	AR-601-TBF4
	iBalance UKA Tibial Bearing Trial, size 6, 8 mm	AR-601-TBF8
	iBalance UKA Tibial Bearing Trial, size 6, 9 mm	AR-601-TBF9
8	Dilating Rasp, slot, meniscal allograft	AR-2963BR
9	iBalance UKA, femoral impactor	AR-611-6
10	iBalance UKA SportsPlasty™ Tibial Keel Punch	AR-611-9
11	iBalance UKA Tibial Impactor	AR-611-4
12	iBalance UKA Alignment Drop Rod w / Sleeve, 3/16 in × 18 in	AR-611-2
13	iBalance UKA, Alignment Rod, 3/16 in × 18 in	AR-601-AR00
14	iBalance UKA Base Plate, size 1	AR-601-TBP1
	iBalance UKA Base Plate, size 2	AR-601-TBP2
	iBalance UKA Base Plate, size 3	AR-601-TBP3
	iBalance UKA Base Plate, size 4	AR-601-TBP4
	iBalance UKA Base Plate, size 5	AR-601-TBP5
	iBalance UKA Base Plate, size 6	AR-601-TBP6


UKA Instrumentation Set 2



UKA Instrumentation Set 2, Right Medial Inner Tray (AR-611-C2R)

Label	Product Description	Item Number
1	iBalance® UKA RM/LL Distal Cut 4 mm Block	AR-611-DR04
	iBalance UKA RM/LL Distal Cut 5 mm Block	AR-611-DR05
	iBalance UKA RM/LL Distal Cut 6 mm Block	AR-611-DR06
	iBalance UKA RM/LL Distal Cut 7 mm Block	AR-611-DR07
	iBalance UKA RM/LL Distal Cut 8 mm Block	AR-611-DR08
	iBalance UKA RM/LL Distal Cut 9 mm Block	AR-611-DR09
	iBalance UKARM/LL Distal Cut 10 mm Block	AR-611-DR10
2	iBalance UKA RM/LL Post Cut 3 mm Block	AR-611-PR03
	iBalance UKA RM/LL Post Cut 4 mm Block	AR-611-PR04
	iBalance UKA RM/LL Post Cut 5 mm Block	AR-611-PR05
	iBalance UKA RM/LL Post Cut 6 mm Block	AR-611-PR06
	iBalance UKA RM/LL Post Cut 7 mm Block	AR-611-PR07
	iBalance UKA RM/LL Post Cut 8 mm Block	AR-611-PR08
	iBalance UKA RM/LL Post Cut 9 mm Block	AR-611-PR09
	iBalance UKA RM/LL Post Cut 10 mm Block	AR-611-PR10
3	iBalance UKA Femoral Finish Guide, size 1, RM/LL	AR-611-CR1
	iBalance UKA Femoral Finish Guide, size 2, RM/LL	AR-611-CR2
	iBalance UKA Femoral Finish Guide, size 3, RM/LL	AR-611-CR3
	iBalance UKA Femoral Finish Guide, size 4, RM/LL	AR-611-CR4
	iBalance UKA Femoral Finish Guide, size 5, RM/LL	AR-611-CR5
	iBalance UKA Femoral Finish Guide, size 6, RM/LL	AR-611-CR6
4	iBalance UKA RM/LL Tibial Finish Guide 1	AR-611-TR1
	iBalance UKA RM/LL Tibial Finish Guide 2	AR-611-TR2
	iBalance UKA RM/LL Tibial Finish Guide 3	AR-611-TR3
	iBalance UKA RM/LL Tibial Finish Guide 4	AR-611-TR4
	iBalance UKA RM/LL Tibial Finish Guide 5	AR-611-TR5
	iBalance UKA RM/LL Tibial Finish Guide 6	AR-611-TR6
5	iBalance UKA Femoral Trial, size 1, RM/LL	AR-601-FTRA
	iBalance UKA Femoral Trial, size 2, RM/LL	AR-601-FTRB
	iBalance UKA Femoral Trial, size 3, RM/LL	AR-601-FTRC
	iBalance UKA Femoral Trial, size 4, RM/LL	AR-601-FTRD
	iBalance UKA Femoral Trial, size 5, RM/LL	AR-601-FTRE
	iBalance UKA Femoral Trial, size 6, RM/LL	AR-601-FTRF

UKA Instrumentation Set 2 (Cont.)



UKA Instrumentation Set 2, Left Medial Inner Tray (AR-611-C2L)

Label	Product Description	Item Number
1	iBalance® UKA LM/RL Distal Cut 4 mm Block iBalance UKA LM/RL Distal Cut 5 mm Block iBalance UKA LM/RL Distal Cut 6 mm Block iBalance UKA LM/RL Distal Cut 7 mm Block iBalance UKA LM/RL Distal Cut 8 mm Block iBalance UKA LM/RL Distal Cut 9 mm Block iBalance UKA LM/RL Distal Cut 10 mm Block	AR-611-DL04 AR-611-DL05 AR-611-DL06 AR-611-DL07 AR-611-DL08 AR-611-DL09 AR-611-DL10
2	iBalance UKA LM/RL Post Cut 3 mm Block iBalance UKA LM/RL Post Cut 4 mm Block iBalance UKA LM/RL Post Cut 5 mm Block iBalance UKA LM/RL Post Cut 6 mm Block iBalance UKA LM/RL Post Cut 7 mm Block iBalance UKA LM/RL Post Cut 8 mm Block iBalance UKA LM/RL Post Cut 9 mm Block iBalance UKA LM/RL Post Cut 10 mm Block	AR-611-PL03 AR-611-PL04 AR-611-PL05 AR-611-PL06 AR-611-PL07 AR-611-PL08 AR-611-PL09 AR-611-PL10
3	iBalance UKA Femoral Finish Guide, size 1 LM/RL iBalance UKA Femoral Finish Guide, size 2 LM/RL iBalance UKA Femoral Finish Guide, size 3 LM/RL iBalance UKA Femoral Finish Guide, size 4 LM/RL iBalance UKA Femoral Finish Guide, size 5 LM/RL iBalance UKA Femoral Finish Guide, size 6 LM/RL	AR-611-CL1 AR-611-CL2 AR-611-CL3 AR-611-CL4 AR-611-CL5 AR-611-CL6
4	iBalance UKA LM/RL Tibial Finish Guide 1 iBalance UKA LM/RL Tibial Finish Guide 2 iBalance UKA LM/RL Tibial Finish Guide 3 iBalance UKA LM/RL Tibial Finish Guide 4 iBalance UKA LM/RL Tibial Finish Guide 5 iBalance UKA LM/RL Tibial Finish Guide 6	AR-611-TL1 AR-611-TL2 AR-611-TL3 AR-611-TL4 AR-611-TL5 AR-611-TL6
5	iBalance UKA Femoral Trial, size 1, LM/RL iBalance UKA Femoral Trial, size 2, LM/RL iBalance UKA Femoral Trial, size 3, LM/RL iBalance UKA Femoral Trial, size 4, LM/RL iBalance UKA Femoral Trial, size 5, LM/RL iBalance UKA Femoral Trial, size 6, LM/RL	AR-601-FTLA AR-601-FTLB AR-601-FTLC AR-601-FTLD AR-601-FTLE AR-601-FTLF



Patellofemoral Joint Arthroplasty

- Simplified onlay system with the benefits of trochlear inlay techniques
- Efficient, single-tray system designed for use in the ASC minimizes sterilization costs and turnover times and helps ensure predictable, reproducible patellofemoral arthroplasty procedures
- Precisely prepare the trochlea to match the underside geometry of the implant with the anterior femoral resection guide and captured jig
- Low-profile, anatomic trochlear implant—available in 4 sizes—articulates with a round dome, 3-peg patellar component


Implants

Product Description	Item Number
Trochlear Components	
iBalance Patellofemoral Implant, left, size 1	AR-502-1L
iBalance Patellofemoral Implant, left, size 2	AR-502-2L
iBalance Patellofemoral Implant, left, size 3	AR-502-3L
iBalance Patellofemoral Implant, left, size 4	AR-502-4L
iBalance Patellofemoral Implant, right, size 1	AR-502-1R
iBalance Patellofemoral Implant, right, size 2	AR-502-2R
iBalance Patellofemoral Implant, right, size 3	AR-502-3R
iBalance Patellofemoral Implant, right, size 4	AR-502-4R
Patellar Components	
iBalance Patella Implant Dome, 27 mm/8 mm, vit E	AR-524-PSA8
iBalance Patella Implant Dome, 30 mm/8 mm, vit E	AR-524-PSB8
iBalance Patella Implant Dome, 32 mm/8 mm, vit E	AR-524-PS328
iBalance Patella Implant Dome, 34 mm/9 mm, vit E	AR-524-PSC9
iBalance Patella Implant Dome, 37 mm/10 mm, vit E	AR-524-PSD0
iBalance Patella Implant Dome, 40 mm/10 mm, vit E	AR-524-PSE0

iBalance PFJ Finishing Blades (Each kit includes 3 blades [distal, medial, and proximal])

Product Description	Item Number
Finishing Blades, size 1 mm	AR-602-26
Finishing Blades, size 2 mm	AR-602-27
Finishing Blades, size 3 mm	AR-602-28
Finishing Blades, size 4 mm	AR-602-29

iBalance® PFJ Instrument Set




iBalance PFJ Instrument Set, Top Tray (AR-602-S)

Label	Product Description	Item Number
1	Anterior Cortex Alignment Guide	AR-602-5
2	Anterior Cut Guide Base, left Anterior Cut Guide Base, right	AR-602-8L AR-602-8R
3	Anterior Cut Guide	AR-602-9
4	Distal/Proximal Stylus	AR-602-15
5	Anterior Cut Stylus	AR-602-10
6	IM Drill	AR-602-6
7	Femoral Peg Drill, Ø 0.312 in	AR-602-42
8	Angel Wing	AR-613-76
9	Finishing Guide, size 1 Finishing Guide, size 2 Finishing Guide, size 3 Finishing Guide, size 4	AR-602-11 AR-602-12 AR-602-13 AR-602-14
10	Trochlear Reamer, size 1 Trochlear Reamer, size 2 Trochlear Reamer, size 3 Trochlear Reamer, size 4 Blade Stabilizer, size 1 Blade Stabilizer, size 2 Blade Stabilizer, size 3 Blade Stabilizer, size 4	AR-602-16 AR-602-17 AR-602-18 AR-602-19 AR-602-20 AR-602-21 AR-602-22 AR-602-23

Label	Product Description	Item Number
11	iBalance Femoral Trial, size 1, left	AR-602-1L
	iBalance Femoral Trial, size 1, right	AR-602-1R
	iBalance Femoral Trial, size 2, left	AR-602-2L
	iBalance Femoral Trial, size 2, right	AR-602-2R
	iBalance Femoral Trial, size 3, left	AR-602-3L
	iBalance Femoral Trial, size 3, right	AR-602-3R
	iBalance Femoral Trial, size 4, left	AR-602-4L
	iBalance Femoral Trial, size 4, right	AR-602-4R
	iBalance PFJ Pegged Femoral Trial, size 1, left	AR-602-50
	iBalance PFJ Pegged Femoral Trial, size 2, left	AR-602-51
	iBalance PFJ Pegged Femoral Trial, size 3, left	AR-602-52
	iBalance PFJ Pegged Femoral Trial, size 4, left	AR-602-53
	iBalance PFJ Pegged Femoral Trial, size 1, right	AR-602-54
	iBalance PFJ Pegged Femoral Trial, size 2, right	AR-602-55
	iBalance PFJ Pegged Femoral Trial, size 3, right	AR-602-56
	iBalance PFJ Pegged Femoral Trial, size 4, right	AR-602-57
12	iBalance PFJ Pegged Drill Stabilizer, size 1	AR-602-38
	iBalance PFJ Pegged Drill Stabilizer, size 2	AR-602-39
	iBalance PFJ Pegged Drill Stabilizer, size 3	AR-602-40
	iBalance PFJ Pegged Drill Stabilizer, size 4	AR-602-41

iBalance® PFJ Instrument Set (Cont.)




iBalance PFJ Instrument Set, Bottom Tray (AR-602-S)

Label	Product Description	Item Number
1	iBalance UKA Alignment Rod, 3/16 in × 18 in	AR-601-AR00
2	iBalance Handle, quick connect	AR-613-8
3	iBalance PFJ Blade Handle	AR-602-25
4	iBalance PFJ Femoral Impactor Head	AR-602-43
5	iBalance PFJ Universal Impactor Handle	AR-603-31
6	iBalance UKA Fixation Pin, 1/8 in × 3 in	AR-601-FP30
7	iBalance PFJ Threaded Pin, headed	AR-603-72
8	iBalance Patella Drill Guide, 27 mm iBalance Patella Drill Guide, 30 mm iBalance Patella Drill Guide, 34 mm iBalance Patella Drill Guide, 37 mm	AR-602-44 AR-602-45 AR-602-46 AR-602-47
9	iBalance PFJ Dome Patella Trial, size 27, 8 mm thick iBalance PFJ Dome Patella Trial, size 30, 8 mm thick iBalance PFJ Dome Patella Trial, size 34, 9 mm thick iBalance PFJ Dome Patella Trial, size 37, 10 mm thick	AR-603-77 AR-603-79 AR-603-81 AR-603-83
10	Patella Peg Drill	AR-603-92
11	iBalance PFJ Pin Driver	AR-613-91
12	iBalance PFJ Patella Cement Clamp	AR-603-85
13	iBalance PFJ Blade Stabilizer Handler	AR-602-24
14	iBalance Patella Calipers	AR-602-48

Not Pictured

Product Description	Item Number
iBalance PFJ Flexion Anterior Cut Guide	AR-602-58
iBalance PFJ Dome Patella Trial, size 40, 10 mm thick	AR-603-84
iBalance Patella Drill Guide Combo	AR-602-49

Patellar Resection System



Patellar Resection System (AR-612-S)

Label	Product Description	Item Number
1	Patella Resection Instrument Case	AR-612-C1
2	Quick-Connect Handle	AR-613-8
3	Patella Trial, 34 mm Patella Trial, 34 mm Patella Trial, 34 mm Patella Trial, 37 mm Patella Trial, 40 mm	AR-603-77 AR-603-79 AR-603-81 AR-603-83 AR-603-84
4	Patella Drill Guide, 27 mm Patella Drill Guide, 30 mm Patella Drill Guide, 34 mm Patella Drill Guide, 37 mm Patella Drill Guide, 40 mm	AR-602-44 AR-602-44 AR-602-46 AR-602-47 AR-602-49
5	O-ring Attachment	AR-613-95
6	Patella Calipers	AR-602-48
7	Patella Peg Drill	AR-603-92
8	Patella Multitool	AR-613-93
9	Patella Resection Guide	AR-612-19

References

1. Magosch P, Lichtenberg S, Habermeyer P. Survival of stemless humeral head replacement in anatomic shoulder arthroplasty: a prospective study. *J Shoulder Elbow Surg.* 2021;30(7):e343-e355. doi:10.1016/j.jse.2020.09.034
2. Alikhah A, Imiolczyk J, Krukenberg A, Scheibel M. Screw fixation in stemless shoulder arthroplasty for the treatment of primary osteoarthritis leads to less osteolysis when compared to impaction fixation. *BMC Musculoskelet Disord.* 2020;21(1):295. doi:10.1186/s12891-020-03277-3
3. Bents EJ, Werner BC, Griffin JW, Raiss P, Denard PJ. A radiographic analysis of proximal humeral anatomy in patients with primary glenohumeral arthritis and implications for press-fit stem length. *J Clin Med.* 2022;11(10):2867. doi:10.3390/jcm11102867
4. Sites Medical Research and Development. Bone ingrowth performance of OsteoSync™ Ti. Available at: <https://www.sitesmedical.com/wp-content/uploads/2024/08/The-Bone-Ingrowth-Performance-of-OsteoSync-Ti.pdf>
5. Arthrex, Inc. Data on file (MyVIP pilot study). Naples, FL; 2021.
6. Erickson BJ, Shishani Y, Bishop ME, Romeo AA, Lederman E, Gobeze R. Tuberosity repair in reverse total shoulder arthroplasty for fracture using a stem-based double-row repair: a cadaveric biomechanical study. *J Am Acad Orthop Surg.* 2020;28(23):e1059-1065. doi:10.5435/JAAOS-D-19-0066
7. Arthrex, Inc. Data on file (LA0239). Naples, FL; 2008.
8. Schmalzl J, Jessen M, Holschen M, et al. Tuberosity healing improves functional outcome following primary reverse shoulder arthroplasty for proximal humeral fractures with a 135° prosthesis. *Eur J Orthop Surg Traumatol.* 2020;30(5):909-916. Dot:10.1007/s00590-020-02649-8



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.

Products may not be available in all markets because product availability is subject to the regulatory approvals and medical practices in individual markets. Please contact Arthrex if you have questions about the availability of products in your area.

arthrex.com

© 2026-01 Arthrex, Inc. All rights reserved. LB1-000312-en-US_F



Arthrex manufacturer,
authorized representative,
and importer information
(Arthrex eIFUs)



US patent information

