

Synergy.net


*A Breakthrough in
OR Management Visibility*



Synergy.net


Synergy.net is a browser-based software application designed to work in conjunction with Synergy^{HD3} and expand the capabilities of Synergy^{HD3} beyond a single operating room.

OR Dashboard



Manage all ORs from one central location

Preferences



Automatically update surgeon preferences for Arthrex devices

OR Command



Live Case View



Watch live surgery without having to be in the OR

OR Scheduler



Improve awareness and time management

Functionalities of Synergy.net

- Increases nurse-management capabilities with real-time views into each OR. Ability to schedule OR rooms and input surgeon and device settings once and have the settings automatically pushed throughout the facility to all Synergy^{HD3} devices.
- Interfaces with EHRs and other HL7-based systems
- Provides one single database for all cases captured by all Synergy^{HD3} devices
- Improves DICOM capabilities, which includes transfer of video to a PACS
- Allows staff to mine procedural data to help improve efficiencies
- Offers network printing from any Synergy^{HD3} device

Part Numbers

AR-3200-1030
AR-3200-1034

Synergy.net Management Software
Synergy.net Management Software HL7



©2014, Arthrex Inc. All rights reserved.
LB1-0695-EN_A

<http://www.arthrex.com/imaging-resection/imaging-consoles>