

QuadPro™ Tendon Harvester

Product Highlight




Arthrex pioneered quadriceps tendon (QT) autograft harvesting in 2014 with the minimally invasive QT harvest instrument set, allowing a much easier and more reproducible harvest of the quad tendon.

The QuadPro tendon harvester is the next evolution in graft harvesting for ACL reconstruction (ACLR). It was specifically engineered to allow for efficient, safe graft harvesting while reducing the morbidity and challenges associated with traditional quad harvesting techniques.

Features and Benefits

- **Reproducible Harvesting:** Available in multiple sizes to accommodate surgeon preferences and each patient's needs. The transparent handle enables direct visualization of the graft during harvesting to enable accurate harvest length.
- **Minimally Invasive Harvesting:** The sharp cylindrical tip safely and reliably cores out a smooth cylindrical graft. The graft is easily amputated through the cutting window. New harvesting technique reduces graft site morbidity and overall procedure time.
- **Versatility:** Compatible with either all-soft-tissue or bone block (BQT) harvesting techniques, and harvesting can be performed using a small transverse incision or a traditional longitudinal incision.



QuadPro™ Tendon Harvester

Ordering Information


QuadPro Tendon Harvester

Product Description	Item Number
QuadPro Tendon Harvester, 8 mm	AR-2386-08
QuadPro Tendon Harvester, 9 mm	AR-2386-09
QuadPro Tendon Harvester, 10 mm	AR-2386-10
QuadPro Tendon Harvester, 11 mm	AR-2386-11

Additional ACL Reconstruction Products

Product Description	Item Number
FiberTag® TightRope® and FlipCutter® III Implant System	AR-1288RTT-FC3
ACL FiberTag TightRope Implant	AR-1588RTT
ACL FiberTag TightRope ABS Implant	AR-1588TNT
Concave TightRope ABS Button 11 mm, round	AR-1588TB-3
Concave TightRope ABS Button 14 mm, round	AR-1588TB-4
Concave TightRope ABS Button 20 mm, round	AR-1588TB-5
ACL Backup Fixation Kit w/ BioComposite SwiveLock® Anchor	AR-1593-BC
ACL Backup Fixation Kit w/ PEEK SwiveLock Anchor	AR-1593-P
FlipCutter III Drill	AR-1204FF
ACL ToolBox Instrument Set	AR-1900S
PCL ToolBox Set	AR-1269S
RetroConstruction™ Drill Guide System Instrument Set	AR-1510S
GraftClamp, Graft Preparation Instrument	AR-2386T

Products may not be available in all markets because product availability is subject to the regulatory approvals and medical practices in individual markets. Please contact your Arthrex representative if you have questions about the availability of products in your area.



arthrex.com

© 2021 Arthrex, Inc. All rights reserved. DOC1-000528-en-US_A

