



FlexiGRAFT®

Cancellous Plugs

Clinical Overview

FlexiGraft Cancellous Plugs are sterile plugs of 100% cancellous bone for use in backfilling OATS® procedures

Applications

- Used to back-fill donor sites created during an osteochondral autograft transfer system (OATS) procedures
- Trephine arthrodesis of tarsal-metatarsal joint
- Other procedures/bone voids that require a small diameter bone

Why Use

- **Versatile:** FlexiGraft Cancellous Plugs are available in diameters 6 - 11 mm and length can be trimmed to properly back-fill donor recovery site
- **Biocompatible:** FlexiGraft Cancellous Plugs are comprised of 100% human bone, which is remodeled naturally during the bone healing process
- **Compatible:** When used in conjunction with the Arthrex OATS Harvester, the graft should be 1 mm larger than the harvester used to provide a line-to-line fit in the donor site. (8 mm Harvester = 9 mm allograft plug) Refer to OATS Backfill Plugs Brochure from Arthrex
- **Convenience:** Created to precisely fit the recovery site created when using the Arthrex OATS system
- **Fully-hydrated, Ambient Storage Available:** FlexiGraft Cancellous Plugs are also offered in Preservon®, a proprietary glycerol-based preservation technology that allows allograft bio-implants to be stored in a fully hydrated state at ambient temperature. Preservon eliminates lengthy thawing and rehydrating times and does not require freezer storage
- **Sterile:** FlexiGraft Cancellous Plugs are sterilized using proprietary and patented Allowash XG® technology. This provides a sterility assurance level of 10⁻⁶ without compromising the graft's inherent biomechanical properties





FlexiGraft Cancellous Plugs with non-viable cartilage

Set length of 16 mm ▶ Frozen Storage/5 year shelf life

Frozen	Diameter
FCPD7	7 mm
FCPD8	8 mm
FCPD9	9 mm
FCPD10	10 mm
FCPD11	11 mm

FlexiGraft Cancellous Plugs without cartilage

Set length of 16 mm ▶ Ambient Storage/5 year shelf life

Preservon	Diameter
PCPD6	6 mm
PCPD7	7 mm
PCPD8	8 mm
PCPD9	9 mm
PCPD10	10 mm
PCPD11	11 mm

68-60-097-02 .01

The LifeNet Health logo, FlexiGraft, Preservon, and Allowash XG are registered trademarks of LifeNet Health. The Arthrex logo and OATS are registered trademarks of Arthrex, Inc.

©2018 LifeNet Health, Virginia Beach, VA. All rights reserved.



North America
1.888.847.7831
orders@lifenethealth.org

Europe
+ 43 1 375002710
eu_orders@lifenethealth.eu

Latin America □ **Asia** □ **Middle East**
1.757.464.4761 ext. 2000
internat.orders@lifenethealth.org

www.LifeNetHealth.org
www.LifeNetHealth.eu