

# AlloSync™ Expand

Demineralized Cortical Fibers



**Arthrex®** 

## Features and Benefits

The unique geometry of AlloSync™ Expand 100% demineralized bone is ideal for intraoperative handling and controlled expansion into bone voids. AlloSync Expand fibers come preloaded in a syringe that allows for consistent hydration of the graft with biologic fluids, such as bone marrow aspirate.

### 100% Demineralized Bone Fibers

- No added fillers for maximum demineralized bone content and osteoinductive potential
- Specific fiber geometry provides exceptional handling and controlled expansion
- Lyophilized fibers extend shelf life while preserving the osteoinductive potential

### Expands to Fill Gaps

- Wicks blood, bone marrow, and other physiological fluids that allow the graft to expand and improve fill

### Cellular Highways

- Fibers have demonstrated superior bone-forming capacity compared to standard particulate demineralized bone marrow<sup>1</sup>
- Entangled fibers create a 3D interconnected matrix to promote cell migration and fusion

### Simplicity of Hydration

- Luer lock portal delivers a simple yet thorough hydration process
- Flexibility to select various hydration fluids

## Ordering Information

Product Description	Item Number
AlloSync Expand, 1 cc	ABS-2017-01
AlloSync Expand, 2.5 cc	ABS-2017-02
AlloSync Expand, 5 cc	ABS-2017-05
AlloSync Expand, 10 cc	ABS-2017-10



Refer to the instructions for use for indications and important safety information.

Warning: Applicable laws restrict these products to sale by or on the order of a physician.

## Reference

1. Martin GJ, Boden SD, Titus L, Scarborough NL. New formulations of demineralized bone matrix as a more effective graft alternative in experimental posterolateral lumbar spine arthrodesis. *Spine*. 1999;24(7):637-645. doi:10.1097/00007632-199904010-00005.

To order, please call Arthrex Customer Service at 1-800-934-4404.

Products may not be available in all markets because product availability is subject to the regulatory approvals and medical practices in individual markets. Please contact Arthrex if you have questions about the availability of products in your area.

View U.S. patent information at [www.arthrex.com/corporate/virtual-patent-marking](http://www.arthrex.com/corporate/virtual-patent-marking)

**[www.arthrex.com](http://www.arthrex.com)**

© 2020 Arthrex, Inc. All rights reserved. FL1-000125-en-US\_B

