

# Lift-Assist® Beach Chair Positioner System

Effortless Patient Positioning

## Easy-to-Use and Protective Head Positioning

- Universal™ Head Positioner (UHP), single-use, simple to apply disposable kit
- Comfortable and secure
- Free-sliding head positioner protects the patient's neck when raising and lowering the table

## Unobstructed Posterior and Lateral Access

- Completely removable fold-away shoulder panels
- Clear access to the operative shoulder

## Lightweight with Lift-Assist

- Quick and easy patient repositioning
- Lift-assist piston supports up to 80 lb (36.3 kg) of the patient's upper body weight, enabling nearly effortless positioning
- Class leading 500 lb (227 kg) patient weight limit
- Lateral slide feature enables the surgeon to move the patient closer to the edge of the OR table for increased surgical site access
- Lightweight, under 30 lb (13.6 kg) without panels



Easy to adjust and position with lift-support pistons.



Safe and supportive. Simple adjustments can be made using lift-support handle.



Lies flat for anesthesia access.



Beach Chair with Reusable Pad Set



Rail Clamp



Intubation Pad and Plate



Reusable Security Strap



LPS Arm Support



Counter Traction Strap



Easy Lock Socket



Universal Head Positioner (UHP) Disposable



UHP w/ Handles



UHP w/o Handles



Reusable Pad Set



UHP Replacement Strap



Cart



Beach Chair w/ Components

## Ordering Information

Product Description	Item Number
<b>Included:</b>	
Beach Chair with Reusable Pad Set	AR-1627
Intubation Pad and Plate	AR-1627-10
Reusable Security Strap	AR-1627-15
Rail Clamps (2 ea)	AR-1627-11
<b>Required:</b>	
LPS Arm Support	AR-1627-01
Counter Traction Strap	AR-1627-03
Easy Lock Socket	AR-1627-12
UHP, w/ handles	AR-1627-05
UHP, w/o handles	AR-1627-20
UHP Disposable (Box of 10)	AR-1627-06
<b>Recommended:</b>	
Reusable Pad Set	AR-1627-09
Cart	AR-1627-13
UHP Replacement Strap	AR-1627-17



Scan to view more product information

**arthrex.com**

© 2021 Arthrex, Inc. All rights reserved. FL1-000368-en-US\_B

