

# Torpedo™ Shaver Blades

Synergy<sup>Resection</sup>™ System




Arthrex® 

# Portfolio of Torpedo™ Shaver Blades

---

## A Size and an Angle for Any Joint Space Application


The best-in-class Torpedo shaver blades have a unique tapered tip and dual intercutting windows for precise and efficient resection across all joint spaces. Explore the unique portfolio offerings available to meet your specific needs.




**4.2 mm**  
Hip-length curved  
Torpedo™ blade

**4.2 mm**  
Hip-length  
Torpedo blade

**5.0 mm**  
Torpedo blade





## Shoulder Procedure Application

---

### **Labral Bony Preparation**

- Tapered tip for improved access

### **Labral Debridement**


- Curved option to reach inferior labrum
- Smooth window for precise cut

### **Biceps Debridement**

- 2× inner cutting windows to decrease clogging

### **Subacromial Decompression**

- 2× inner cutting windows for increased efficiency while resecting bursal and frayed rotator cuff tissue



## Knee Procedure Application

---

### **Chondroplasty**


- Smooth window configuration for creating a stable border

### **Cartilage Repair**

- Clinically tested for high cell viability when harvesting chondrocytes with GraftNet™ tissue collector<sup>1</sup>

### **Meniscectomy**

- Razor sharp blade for clean resection
- Curved option and tapered tip for reaching posterior meniscus in tight knees



## Hip Procedure Application

---

### **Labral Debridement**


- Curved option to support two portal techniques

### **Peritrochanteric Space Procedures**

- Straight option available to fit down a 4.5 mm metal cannula

### **Capsular Debridement for Improved Repair**

- Efficient, precise window leaves a smooth edge to help improve capsular closure



## Ankle Procedure Application

---

### **Lateral Ankle Instability**

- Operates on Nano handpiece (small hub) with 7 cm length for improved control and precision


### **Osteochondral Lesions**

- Low-profile tapered tip for better navigation in tight joint spaces

**Note: Multiple hub options are available to support all Arthrex shaver handpieces.**

#### Reference

1. Arthrex, Inc. LA1-000143-en\_US-A. Naples, FL; 2021.



**Arthrex**<sup>®</sup> 

©2022-10 Arthrex, Inc. All rights reserved.  
LB1-000391-en-US\_A