

JumpStart®

ANTIMICROBIAL WOUND DRESSING

Infection Control and Cosmetic Closure—All in One Dressing

- › Polyester substrate with embedded elemental silver and zinc microcell batteries
- › JumpStart dressing combats biofilms and bacteria to lower infection risk¹
- › Bioelectric dressings can deliver 35% faster healing with visibly better scars²
- › Up to 7-day wear time, resulting in fewer dressing changes and lower cost of care³
- › May be used in conjunction with negative-pressure wound therapy (NPWT) systems

Combining wireless electrical dressing with 5-day NPWT




cut cost of care by

20%

vs NPWT alone³



The microcurrent technology of JumpStart dressing can reduce wound infection risk and promote healing²



Learn more about
JumpStart dressing

Arthrex®

Ordering Information

JumpStart® Contact Layer Dressing

Dressing Dimensions	Qty/Box	Item Number
1.5 in × 8 in	10	ABS-4005
1.5 in × 10 in	10	ABS-4006
2 in × 5 in	10	ABS-4025

JumpStart FlexEFIT® Dressing

Dressing Dimensions	Adhesive Dimensions	Qty/Box	Item Number
1.5 in × 4.5 in	2.5 in × 6.3 in	5	ABS-4060-05
1.5 in × 4.5 in	2.5 in × 6.3 in	10	ABS-4060-10

JumpStart ClearFit® Dressing

Dressing Dimensions	Qty/Box	Item Number
4 in × 4 in	10	ABS-4061
6 in × 6 in	10	ABS-4062
3.5 in × 2.25 in	10	ABS-4063

Dressing Dimensions	Qty/Box	Item Number
JumpStart pin-site dressing kit	5	ABS-4059

Biologic CT Kits

Product Description	Qty/Box	Item Number
FlexEFit Sternal Set		ABS-4060L-T
FlexEFit dressing, 1.5 in × 4.5 in	10	
Energel wound hydrogel	10	
FlexEFit – Angel PRP Sternal Set		ABS-4060L-61T
FlexEFit dressing, 1.5 in × 4.5 in	10	
Energel wound hydrogel	10	
Angel PRP kit	1	
JumpStart Sternal Set		ABS-4006L-T
JumpStart dressing, 1.5 in × 10 in	10	
Energel wound hydrogel	10	
JumpStart – Angel PRP Sternal Set		ABS-4006L-61T
JumpStart dressing, 1.5 in × 10 in	5	
Energel wound hydrogel	10	
Angel PRP kit	1	

Product Description	Qty/Box	Item Number
Energel® wound hydrogel	10	AGL-L075-10

Products advertised in this brochure / surgical technique guide may not be available in all countries. For information on availability, please contact Arthrex Customer Service or your local Arthrex representative.

References

- Kim H, Park S, Housler G, Marcel V, Cross S, Izadjoo M. *Mil Med.* 2016;181(5 Suppl):184-190. doi:10.7205/MILMED-D-15-00157
- Blount AL, Foster S, Rapp DA, Wilcox R. *J Burn Care Res.* 2012;33(3):354-357. doi:10.1097/BCR.0b013e31823356e4
- Ghatak PD, Schlanger R, Ganesh K, et al. *Adv Wound Care (New Rochelle).* 2015;4(5):302-311. doi:10.1089/wound.2014.0615