


Bone Graft Delivery

Brochure



Arthrex® 

BioXpress™ Graft Delivery Device



- Efficient delivery of bone graft through 10 cm or 15 cm cannula
- Second internal plunger clears cannula contents
- Luer-to-luer connector allows for delivery into smaller syringes
- Funnel ensures clean bone graft transfer
- **Device can be used with:** AlloSync™ Pure demineralized bone matrix (DBM), AlloSync Expand demineralized cortical fibers, AlloSync gel, AlloSync putty, AlloSync CB paste, AlloSync CB putty, and AlloSync demineralized cancellous chips

Product Description	Item Number
Blunt Tip Cannula, 10 cm	ABS-10053-10
Angled Tip Cannula, 10 cm	ABS-10053-10-45
Blunt Tip Cannula, 15 cm	ABS-10053-15
Angled Tip Cannula, 15 cm	ABS-10053-15-45





Blunt Tip Cannula



Angled Tip Cannula




Ankle Arthroscopy Instrument Set—Bone Graft Delivery Cannula



- Included with ankle arthroscopy instrument set
- Tapers from 5.5 mm to 7.7 mm
- Associated tamp for graft delivery (available separately)
- **Set can be used with:** AlloSync™ Pure DBM, AlloSync Expand demineralized cortical fibers, AlloSync gel, AlloSync putty, AlloSync CB paste, AlloSync CB putty, and AlloSync chips

Product Description	Item Number
Ankle Arthroscopy Instrument Set	AR-8656S
Cannula, bone graft	AR-8656-71
Tamp, bone graft	AR-8656-70

IntraOsseous BioPlasty® (IOBP) Kit



- Multiple sizes (8 ga - 13 ga) of open tip delivery cannulas allow for percutaneous graft delivery
- Includes mixing syringe with luer-to-luer connector
- ABS-2020-OT includes 1 cc delivery syringes
- **Kit can be used with:** AlloSync™ Pure DBM and AlloSync gel

Product Description	Delivery Cannula	Item Number
IOBP® Core Decompression and Delivery Kit, open tip	Open Tip, 8 ga × 4.5 in	ABS-2000-OT
IOBP Core Decompression and Delivery Kit, foot and ankle	Open Tip, 13 ga × 4.5 in	ABS-2020-OT

Mixing Syringe



- Mixing syringe allows for thorough mixing of bone graft with autologous fluids
- Male luer connection delivers bone graft or transfers to delivery cannulas
- 14 cc syringe included in IOBP® kits
- 3 cc syringe included in mixing and delivery kits
- 14 cc syringe is also available individually packaged
- **Mixing syringe can be used with: AlloSync™ Pure DBM and AlloSync gel**

Product Description	Delivery Cannula	Item Number
IOBP Core Decompression and Delivery Kit, open tip	Open Tip, 8 ga × 4.5 in	ABS-2000-OT
IOBP Core Decompression and Delivery Kit, foot and ankle	Open Tip, 13 ga × 4.5 in	ABS-2020-OT
BioCartilage® Mixing and Delivery Kit, large joint	Tuohy Needle w/ ArthroPaddle™ Device	ABS-1000-L
BioCartilage Mixing and Delivery Kit, small joint	Tuohy Needle w/ ArthroPaddle Device	ABS-1000-S
Mixing Syringe w/ Luer Cap	N/A	ABS-2000





This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience, and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level and/or outcomes.

View U.S. patent information at www.arthrex.com/corporate/virtual-patent-marking

arthrex.com

© 2022 Arthrex, Inc. All rights reserved. LB1-000352-en-US_A